

# FABRICATION ENTERPRISES

physical and occupational therapy products since 1974







exercise ADL evaluation electrotherapy

hot/cold DME

### 4 easy ways to order









1-800-431-2830 1-914-345-9300



mail to
Fabrication Enterprises
PO Box 1500
White Plains, NY 10602 USA

Fabrication Enterprises (FEI) is a manufacturer, importer and master distributor of products for physical and occupational therapy, rehabilitation, athletic training, sports medicine and homecare. Since 1974 our products have been sold to hospitals, clinics, fitness centers, professionals and individuals by a network of dealers throughout the world. With our *new packaging*, some products can now be sold through retail distribution as well!

Many of our products are manufactured in the USA and are "CE marked" for sale in Europe and beyond. Most products you see in our catalog are in stock and available for immediate shipment. FEI can manufacture and supply many products with your private-label and can even supply products not in our catalog. Please call Jason (ext.1215) to discuss.

Fabrication is the manufacturer of the following product lines:

Cando® exercise band and tubing (look for the new packaging!)
Cando® Digi-Flex®, Digi-Extend®, Digi-Squeeze® and Digi-Extend n'Squeeze®
Cando® inflatable balls, foam rollers, and other balance products
Cando® Theraputty™ exercise putty
Cuff®, Pouch™ and SoftGrip™ exercise weight
Cando® and Armasport® exercise mats
Baseline® measurement instruments
Relief Pak® hot and cold products
Point-Relief™ topical analgesics and massage products
WaxWel™ paraffin bath and refills
Skillbuilders® positioning products
Dipsters® patientwear

The prices in this catalog are suggested "list" prices. Call for our most current dealer price list.

Please contact me directly (ext.1200) or Jason (ext.1215) regarding additional pricing for products you may want to include in your catalog or on your website. Photography for all products is available on CD-Rom and on our FTP site. Contact Andrew (ext.1225) for further information.

Sincerely.

Elliott Goldberg, President



Elliott Goldberg President ext.1210 elliott@Fab-Ent.com



Linda Goldberg
Customer Service ext. 1205
linda@Fab-Ent.com



Andrew Goldberg Marketing ext. 1225 andrew@Fab-Ent.com



Jason Drucker Sales ext. 1215 jason@Fab-Ent.com

# Fabrication Enterprises ~ 2011-2012

strength training and rehabilitation, balance and coordination, exercise/therapy mats, range-of-motion, flexibility	exercise	4-60
strength, range-of-motion, posture, edema, sensory, body fat,	evaluation	61-91
vital sign, timers and manual dexterity	evaluation	01-91
anatomical model, anatomical chart	models & charts	92-95
cold therapy: packs, chilling unit and wraps, hot therapy: moist heat	hot / cold	96-109
packs, paraffin units/wax, diathermy, topical analgesics: spray/gel	not / colu	30-103
electrotherapy, iontophoresis, compression therapy, continuous passive motion, traction	treatment modalities	110-131
whirlpools, disposable patient wear for hydrotherapy,	hydrotherapy	132-135
aquatic strength and endurance training	пустотнегару	102 100
splinting pan, splinting sheet material, kinesiology tape	splint material and taping	136-137
manual/electric massagers, massage lotion/warmers, massage tables	massage	138-141
treatment tables, mat tables, positioning rolls/wedges, supports, stools, lamps, glass/glassless mirrors, tilt tables, work tables, ambulation training	treatment furniture	142-153
pediatric walkers, canes and crutches, adult walkers, rollators, wheelchairs and transport chairs, wheelchair cushions	DME mobility	154-159
	·	
commodes, potty seats, raised toilet seats, bathing chairs, bath benches, bathroom aids	DME bathroom	160-163
splinting and supports for various body parts, drop foot braces	splinting & supports	164-165
swing seats, helmets, rolls and wedges, positioning shapes/pillows,		100 100
modular seating systems, standers, car seats, strollers, crawlers	pediatric positioning	166-180
ball pit, mat shapes, balance systems, balance beams, inflatable balance balls	pediatric balance	181-183
non-slip material, hook/loop, dressing aids, reachers, eating/drinking	45	194 10F
aids, bedroom aids, medication aids, hygienic aids, transfer/gait belts	ADL	184-195
patient alarms, patient sensors	patient alarms	196
alphabetical listing of brand/product names with page numbers	indov	197-208
aipriabetical ilsuing of brand/product frames with page numbers	index	131-200



### **NEW PRODUCTS**



Cold n'Hot™ treatment compresses for the convenience of cold and hot therapy in one pack

page 106



The best value in exercise bands: the new **Perf 100™** in the economical 100 yard bulk dispenser with ready-tear perforations

page 106



The Cando® Plus antimicrobial foam rollers protect against microbes and odors

page 39

### **NEW PRODUCTS**



Slide the new **Absolute+Axis™** attachment onto your 12 inch goniometer to add vertical and horizontal accuracy

page 77



For temporary pain relief try Point Relief® ColdSpot® and HotSpot® topical analgesic gels

page 98



The Baseline® Sit and Reach Box enables testing against four different standards; deluxe model has variable start points for different arm/leg lengths page 89

### **NEW PRODUCTS**



Use **Digi-Extend n' Squeeze**® for complete flexion and extension exercise and to develop isolated finger strength page 20



Tough enough for clinics, school and fitness centers, ArmaSport®
Mats have sanitized coating, are useable in water, and have hanging grommets for convenient storage

page 50



For grip strength testing, the **Baseline® HD™ dynamometer** features heavy duty construction and a 5 year warranty

page 63





# **Easier Than Ever!**

- Bulk Sized Roll
  - 100-yard Roll Dispenser Box
- Perfect Size for all Exercises
  - 60 5-foot Segments per Roll
- Easy to Dispense
  - Just Pull and Separate
- Save Time and Money

Low Powder Latex\* Band

x-light	yellow	10-5191
light	red	10-5192
medium	green	10-5193
heavy	blue	10-5194
x-heavy	black	10-5195
	·	·

\*CAUTION latex may cause allergic reactions

### **No Latex Formula Band**

x-light	yellow	10-5691
light	red	10-5692
medium	green	10-5693
heavy	blue	10-5694
x-heavy	black	10-5695



Bulk roll is perforated every 5 feet (total of 60 exercisers per roll)



this



not this







# low powder exercise band

### for small clinics

- an economic solution for upper and lower body exercise
- lightweight, compact and portable
- durable
- effective when used alone, or with handles and anchors
- color-coded according to the exercise band's resistance
- Cando<sup>®</sup> is licensed to use the Thera-Band<sup>®</sup> color sequence so resistance level is easily identified
- available in 8 levels of resistance for progressive exercise
- uniform width each band is 5" wide
- time-tested and user-tested in clinics since 1987

Cando® low powder exercise band is made of natural rubber latex, which may cause allergic reactions. For those with latex sensitivities, we have formulated a synthetic latex-free band (page 6).

Thera-Band® colors licensed by Goldberg from Hygenic Corp. without endorsement.



tan	xx-light	10-5210	
yellow	x-light	10-5211	
red	light	10-5212	
green	medium	10-5213	
blue	heavy	10-5214	
black	x-heavy	10-5215	
silver	xx-heavy	10-5216	
gold	xxx-heavy	10-5217	
yel red g	re blu blk set	10-5218	

### single exerciser



tan	xx-light	10-5200	
yellow	x-light	10-5201	
red	light	10-5202	
green	medium	10-5203	
blue	heavy	10-5204	
black	x-heavy	10-5205	
silver	xx-heavy	10-5206	
gold	xxx-heavy	10-5207	
yel red g	gre blu blk set	10-5208	

### 40-piece dispenser



- dispense individually wrapped bands
- 40 bands per box

tan	xx-light	10-5240	
yellow	x-light	10-5241	
red	light	10-5242	
green	medium	10-5243	
blue	heavy	10-5244	
black	x-heavy	10-5245	
silver	xx-heavy	10-5246	
gold	xxx-heavy	10-5247	
yel red	gre blu blk set	10-5248	

### PEP™ pack



- contains 3 different 4-foot elastic exercise bands
- contains an exercise chart

### contents

easy	yellow	red	green	10-5280
moderate	green	blue	black	10-5282
challenging	black	silver	gold	10-5284






# PERF 100

### 100 yard roll with perforations in dispenser box

- economic bulk size roll is best value
- easy to dispense; just pull and separate
- 60 5-foot perfect-size pieces for all exercises







# no latex exercise band

- used in clinics for rehabilitation, conditioning and training
- specially formulated to mimic the stretch characteristics of latex band
- available in 8 strengths to allow for progressive exercise
- color-coded according to the exercise band's resistance level
- identifiable as "no latex" by checkered pattern on packaging







- effective exerciser when used alone or with handles, anchors and stirrups
- each band is 5" wide
- lightweight, compact and portable

### continuous roll

# NEW roll with perforations



### 50-yard tan xx-light 10-5620 x-light light medium heavy

yellow 10-5621 10-5622 10-5623 10-5624 x-heavy 10-5625 xx-heavy 10-5626 xxx-heavy 10-5627 yel red gre blu blk set 10-5628



6 yards

	o-yar	a
tan	xx-light	10-5610
yellow	x-light	10-5611
red	light	10-5612
green	medium	10-5613
blue	heavy	10-5614
black	x-heavy	10-5615
silver	xx-heavy	10-5616
gold	xxx-heavy	10-5617
yel red g	gre blu blk set	10-5618







- easy to dispense
- sized for all exercises



no scissors necessary

### 100-yard roll with perforations

yellow	x-light	10-5691
red	light	10-5692
green	medium	10-5693
blue	heavy	10-5694
black	x-heavy	10-5695
yel red	gre blu blk	set 10-5698

# pre-cut 4 foot exerciser



### 40-exerciser box

tan	xx-light	10-5640
yellow	x-light	10-5641
red	light	10-5642
green	medium	10-5643
blue	heavy	10-5644
black	x-heavy	10-5645
silver	xx-heavy	10-5646
gold	xxx-heavy	10-5647
vel red	are blu blk set	10-5648



### single exerciser

	_	
tan	xx-light	10-5600
yellow	x-light	10-5601
red	light	10-5602
green	medium	10-5603
blue	heavy	10-5604
black	x-heavy	10-5605
silver	xx-heavy	10-5606
gold	xxx-heavy	10-5607
vel red	re blu blk set	10-5608





- thread band through strap; secure 10-5330 a handle pair
- loop at one end, nub or disc at other end
- place in doorway
- close door to secure

10-5320 nub anchor ea disc anchor ea 10-5325

- double loop: place band through one loop and other loop around object or foot
  - d stirrup anchor ea 10-5310

Thera-Band® colors licensed by Goldberg from Hygenic Corp. without endorsement.



# AccuForce<sup>™</sup> band

 specially designed exercise band allows user to know how much force is used when exercising


- simple to use:
  - 1. stretch Cando® AccuForce™ exercise band
  - 2. when a rectangle forms a square, the force indicated on the square has been exerted
- used for upper and lower body exercises
- AccuForce™ band is 5" wide
- available in 5 progressive levels of resistance
- made of high quality latex



### continuous roll pre-cut piece Cando<sup>®</sup> **AccuForce**<sup>™</sup> **Band** Color / Resistance 40 4-Foot Pieces **4-Foot Pieces** Chart 50-Yard 6-Yard 10-5901 10-5921 10-5911 10-5941 yellow x-light 10-5902 10-5912 10-5942 light 10-5922 10-5903 10-5913 10-5943 green medium 10-5923 10-5904 10-5914 10-5944 10-5924 heavy 10-5905 10-5915 10-5945 black x-heavy 10-5925 yel red gre blu blk set 10-5948 10-5908 10-5928 10-5918

### Cando® loops

- made of high quality latex exercise band pre-formed into a loop
- can be used to perform upper and lower body exercises
- available in 7 progressive levels of resistance and 3 sizes
- lengths represent loop when flat:
- loops are 3" wide
- FULL BODY SET contains 10", 15" and 30" loop of specified level for full-body workout



0" loop is equivalent to 20" of flat band cut and tied in a knot

		10" loop each	15" loop each	30" loop each	full body set each
yellow	x-light	10-5251	10-5261	10-5291	10-5841
red	light	10-5252	10-5262	10-5292	10-5842
green	medium	10-5253	10-5263	10-5293	10-5843
blue	heavy	10-5254	10-5264	10-5294	10-5844
black	x-heavy	10-5255	10-5265	10-5295	10-5845
silver	xx-heavy	10-5256	10-5266	10-5296	10-5846
gold	xxx-heavy	10-5257	10-5267	10-5297	10-5847
		bulk loon	packages		
		·		30" loop 10-pack	Contract of the Contract of th
yellow	x-light	10-5251-10	10-5261-10	10-5291-10	Care College
red	light	10-5252-10	10-5262-10	10-5292-10	
green	medium	10-5253-10	10-5263-10	10-5293-10	14
blue	heavy	10-5254-10	10-5264-10	10-5294-10	

10-5265-10

10-5266-10

10-5267-10

10-5295-10

10-5296-10

10-5297-10

individual exercise loops



10-inch loop

(upper body)



x-heavy

xx-heavy

xxx-heavy


10-5255-10

10-5256-10

10-5257-10

30-inch loop

(full body)



# exercise tubing

- · an economic solution for rehabilitation and strength training
- use for upper and lower body exercises
- lightweight resistance tubing is compact and portable
- time-tested and user-tested in clinics since the late 1980s
- effective when used alone, or with handles and anchors
- · color-coded to indicate resistance level
- Cando<sup>®</sup> is licensed to use the Thera-Band<sup>™</sup> color sequence, so resistance level is easily identified
- available in 8 levels of resistance
- available in both low powder (latex) and latex-free formulas

For those with latex sensitivities we have formulated a synthetic latex-free tubing (bottom of page).

Thera-Band® colors licensed by Goldberg from Hygenic Corp. without endorsement.

### Cando® Low Powder exercise tubing

### PEP™ packs



### 100-foot box

tan	xx-light	10-5520
yellow	x-light	10-5521
red	light	10-5522
green	medium	10-5523
blue	heavy	10-5524
black	x-heavy	10-5525
silver	xx-heavy	10-5526
gold	xxx-heavy	10-5527
yel red o	gre blu blk set	10-5529

### 25-foot bag

tan	xx-light	10-5510
yellow	x-light	10-5511
red	light	10-5512
green	medium	10-5513
blue	heavy	10-5514
black	x-heavy	10-5515
silver	xx-heavy	10-5516
gold	xxx-heavy	10-5517
yel red	10-5519	

\*Cando® low powder exercise tubing is made of natural rubber latex, which may cause allergic reactions.

### **Cando**® **PEP**<sup>TM</sup> tubing packs



- with 6-foot latex tubing as indicated
- illustrated instruction/exercise chart

### **PEP™ Tubing**

		_		
easy	yellow	red	green	10-5380
moderate	green	blue	black	10-5382
challenging	black	silver	gold	10-5384

### Cando® No Latex exercise tubing



# 100-foot box

tan	xx-light	10-5720
yellow	x-light	10-5721
red	light	10-5722
green	medium	10-5723
blue	heavy	10-5724
black	x-heavy	10-5725
silver	xx-heavy	10-5726
gold	xxx-heavy	10-5727
vel red	gre blu blk set	10-5729

# 25-foot bag

tan	xx-light		10-5/10
yellow	x-light		10-5711
red	light	10-5712	
green	medium		10-5713
blue	heavy		10-5714
black	x-heavy		10-5715
silver	xx-heavy		10-5716
gold	xxx-heavy		10-5717
yel red g	10-5719		

### accessories



- · thread band or tubing through strap and secure
  - a handle pair
- 10-5330
- loop at one end, nub or disc at other end b nub anchor ea 10-5320
- place in doorway, close door to secure c disc anchor ea 10-5325
- double loop: thread band through one loop; place other loop around object or
  - d stirrup anchor ea 10-5310



# exercise tubing



### Cando® Tubing-with-Handles™ exerciser

- preassembled Cando® Tubing-with-Handles™ exercisers available in three sizes to accommodate upper and lower body workouts
- lightweight, compact and portable
- perform strengthening workouts on-the-go, even in your hotel room!
- two durable foam covered plastic handles for comfort and a firm grip
- made with high quality latex Cando<sup>®</sup> tubing
- FULL BODY SET contains 18", 36" and 48" tubing; for total body workout





		18 inch*	36 inch*	48 inch*	full body set	
tan	xx-light	10-5540	10-5550	10-5560	10-5820	
yellow	x-light	10-5541	10-5551	10-5561	10-5821	
red	light	10-5542	10-5552	10-5562	10-5822	
green	medium	10-5543	10-5553	10-5563	10-5823	0
blue	heavy	10-5544	10-5554	10-5564	10-5824	
black	x-heavy	10-5545	10-5555	10-5565	10-5825	
silver	xx-heavy	10-5546	10-5556	10-5566	10-5826	
yel red g	gre blu blk set	10-5549	10-5559	10-5569	10-5829	

\*length of tubing from end to end when unit is laid flat and is taut

### Cando® BowTie™ exerciser

- preassembled Cando® BowTie™ exercisers are lightweight, take up minimal room and are portable
- color-coded to indicate resistance level
- two durable foam covered handles for comfort and a firm grip
- made with high quality latex Cando® tubing
- used for strength training, conditioning and
- FULL BODY SET contains 14", 22" and 30" tubing; for total body workout



		14 inch*	22 inch*	30 inch*	full body set
tan	xx-light	10-5570	10-5580	10-5590	10-5830
yellow	x-light	10-5571	10-5581	10-5591	10-5831
red	light	10-5572	10-5582	10-5592	10-5832
green	medium	10-5573	10-5583	10-5593	10-5833
blue	heavy	10-5574	10-5584	10-5594	10-5834
black	x-heavy	10-5575	10-5585	10-5595	10-5835
silver	xx-heavy	10-5576	10-5586	10-5596	10-5836
<mark>yel</mark> red gre	blu blk set	10-5579	10-5589	10-5599	10-5839

<sup>\*</sup>length of tubing from end to end when unit is laid flat and is taut





# bungee cords and exercisers

# bungee exercise cords and accessories

- color-coded bungee cords used for upper and lower body exercises
- each cord has snap hook on one end AND door jamb nub and loop anchor on other end
- attach snap hook to D-ring handles, straps or mounts
- grip loop, use as a foot stirrup, or place nub in a door jamb





bungee cord not included

(d) door mount: vertical strap is 16 ft. long to attach around door

(e) wall mount: vertical strap is 5.5 ft. long, horizontal strap is 16" long with 5 ring attachments

### a foam padded handle webbing strap 10-3221 a foam padded handle 10-3220 b adjustable large strap 12-0422 c adjustable medium strap 12-0420 adjustable small strap 12-0424 adjustable shoulder harness 12-0427 adjustable head harness 12-0429 7-foot bungee cord door mount 10-5150 e wall mount 10-5152 b adjustable large strap adjustable medium strap

Color	Strength	4 foot cord	7 foot cord
tan yellow red green	xx-light	10-5810	10-5800
	x-light	10-5811	10-5801
	light	10-5812	10-5802
	medium	10-5813	10-5803
blue	heavy	10-5814	10-5804
black	x-heavy	10-5815	10-5805
silver	xx-heavy	10-5816	10-5806
gold	xxx-heavy	10-5817	10-5807
yel red	gre blu blk set	10-5819	10-5809

### tubing with cuffs exercisers

- continuous tubing loop with two padded ankle/wrist cuffs
- preassembled exerciser has a circumference of 35"
- each exerciser is made with high quality latex Cando<sup>®</sup> tubing
- used for strength training, conditioning and rehabilitation
- use on arms or legs
- available in 8 levels of resistance, color-coded



blue	neavy	10-5764
black	x-heavy	10-5765
silver	xx-heavy	10-5766
gold	xxx-heavy	10-5767

<mark>yel red gre</mark> set 10-5769 blu blk





for lower body workout



### adjustable Cando® bands

The resistance of each band can be increased or decreased by adjusting length of band used. "Quick-attachment" hooks attach to the "D-ring" on the handle and strap accessories.



3 band kit includes: carrying bag, instruction booklet, 1 yellow, 1 red, 1 green adjustable band, 2 handles, and 1 anchor strap

# adjustable Cando $^{\scriptsize \odot}$ kits (bands, 2 handles, 1 anchor strap, 1 exercise booklet)

10-3234 2 band moderate (green, blue)
10-3233 2 band easy (yellow, red)
10-3232 3 band (red, green, blue)
10-3231 4 band (yellow, red, green blue)
10-3230 5 band (yellow, red, green, blue, black)

### adjustable Cando® bands with adjustment clips and hooks

		— 1 band —	10 bands
		i baila	10 bando
yellow	x-light	10-3201	10-3201-10
red	light	10-3202	10-3202-10
green	medium	10-3203	10-3203-10
blue	heavy	10-3204	10-3204-10
black	x-heavy	10-3205	10-3205-10

### adjustable Cando® band accessories with D-ring

10-3220 handle (pair) 10-3221 leg attachment strap (25 in.), each 10-3222 anchor strap (16 in.), each



# Point-Relief™ Be-Better™ kits



### combine exercise therapy with pain relief in convenient kits



general rehab & pain relief kits

### shoulder rehab & pain kits

- for home exercise programs and individual patient use
- includes the combination of Cando® exercise products to increase strength, and ColdSpot™ pain relieving spray to increase comfort

### kit includes:

- Cando<sup>®</sup> exercisers
- Cold Spot<sup>™</sup> relief spray
- · retail packaging



shoulder

48" yellow Cando® band Cando® shoulder pulley 2 oz. Cold Spot™ spray 10-5160

### targeted rehab kits



- for home exercise programs and individual patient use
- includes the combination of Cando® exercise products to increase strength and Cold Spot™ pain relieving spray

### kits include:

Cando<sup>®</sup> exercisers

to increase comfort

- Cold Spot<sup>™</sup> relief spray
- · retail packaging



light 48" yellow & red Cando® bands 10-5161 2 oz. Cold Spot™ spray

48" red & green Cando® bands 3 oz. Cold Spot™ roll-on moderate

10-5162

- · designed to target specific areas
- includes the combination of Cando® exercise products to increase strength

and ColdSpot™ pain relieving spray to increase comfort

### lower back

topical relief

· retail packaging

48" red tubing w/ handles, door anchor, 10-5165 2 oz. Cold Spot™ spray

### shoulder

48" red tubing w/ handles, door anchor, 10-5166 3 oz. Cold Spot™ roll-on

48" red tubing w/ handles, 15" red loop, 10-5167 door anchor, 2 oz. Cold Spot™ spray

### knee/hip

10" red loop, 3 oz. Cold Spot™ roll-on 10-5168

48" red tubing w/ handles, door anchor, 10-5169 red twist-n-bend bar, 3 oz. Cold Spot™ roll-on

### lower body

48" blue exercise tube, loop stirrups, door anchor 10-5163

For additional information on Point Relief Pain Relieving Topical Analgesic products see page 98 in the 2011 Fabrication Enterprises catalog, or visit www.Point-Relief.com. Cando® band and tubing offered in the Point Relief Be Better kits contains natural rubber latex which may cause allergic reactions.



# standard adjustable handle

- makes exercising with band and tubing easier!
- exercise band or tubing attaches to foam handles for comfort-
- each handle has 6" long foam covered rigid bar
- · heavy-duty webbing for durability
- easy to load and unload band or tubing into handles



10-5330

1 pair standard

10-5330-10 10-5330-50

10 pair standard 50 pair standard

### MarV<sup>™</sup> handles

- MarV™ handle has 3 grips: channel grip, vertical grip and stirrup grip
- use with exercise tubing or as an overdoor pulley



channel

10-5347

10-5347-10

10-5347-50



1 pair 10 pair 50 pair



10-5300 1 pair 10-5300-10 10 pair 10-5300-50 50 pair

### foam handles


- thread band or tubing through foam covered rigid handle
- tie band/tubing into a loop
- handle uniformly distributes the pressure across hand

# band and tubing handles

### accessory kit

includes accessories to enhance a band or tubing workout

standard adjustable foam handles - pair



accessory kit Hold-Rite™ handles

• with plastic or foam padded hand grips

10-5334

· self-locking "adjuster" holds band or tubing in place

 use self-locking "adjuster" to increase or decrease resistance by shortening or lengthening band/tubing

### plastic HoldRite™ handles

10-5340 1 pair 10-5340-10 10 pair 10-5340-50 50 pair

### foam covered HoldRite™ handles

10-5450 1 pair 10-5450-10 10 pair 10-5450-50 50 pair

### **D-ring handles**

• webbing exercise handle can be used by threading the band or tubing

through the ring and knotting



10-3220 1 pair 10-3220-10 10 pair 10-3220-50

### sports handle

- 12" foam covered rigid handle
- 3 attachment positions
- use with band/tubing anchored to object or foot



both ends

center 10-5336 1 each 10-5336-10 10 each 10-5336-50 50 each

either end




# band and tubing anchors



### loop stirrup








create "luggage tag" loop around sturdy anchor

- loop stirrup has two webbing loops: one large, one small
- place band or tubing through small webbing loop and put other loop around any stationary object such as your foot, a door anchor, or a wall anchor
- can be used in conjunction with handles to facilitate a wide range of exercises



10-5310 1 each 10-5310-10 10 each 10-5310-50 50 each

### band Klips™

• form the exercise band into a loop and secure with clamp



Band Klips™ 5½" Klip™ (each) 10-5360 1 band klip 10-5360-10 10 band klips

Band Klips™ 3½" Klip™ (each) 10-5370 1 band klip 10-5370-10 10 band klips

### tubing Klips™

• can be used to clamp tubing loop and to prevent one end of tubing from being pulled through an anchor



Klips™ for tan-red (pair) 10-5365 1 tube klip 10-5365-10 10 tube klips

Klips™ for green-black pair 10-5375 1 tube klip 10-5375-10 10 tube klips \*for silver and gold tubing use 3.5" band Klip™. Colors refer to Cando®

### door disc anchor strap



- · place flexible disc in door jamb at any height, even on top of door!
- · close door to anchor disc
- easy to load and unload band or tubing into anchor strap
- heavy-duty webbing for durability

10-5325 1 each 10-5325-10 10 each 10-5325-50 50 each



### anchor stirrups with nub



- · facilitates anchoring band or tubing in door jamb or around anchor
- webbing strap has a loop on one end and a nub on the other
- place nub in door jamb at any height, even on top of door!
- close door to workout!

10-5320 1 each 10-5320-10 10 each 10-5320-50 50 each

# D-ring velcro strap

- position strap around ankle, wrist or thigh and secure with hook and loop
- D-ring can be used as an anchor for the band or tubing to perform various exercises



10-3222 16 inch, 1 each 10-3221 25 inch, 1 each



### extremity cuff strap

- secure band or tubing into fastening system
- loop cuff around body part



16" diameter webbing 10-5356 1 each 10-5356-10 10 each 10-5356-50 50 each

25" diameter webbing 10-5346 1 each 10-5346-10 10 each 10-5346-50 50 each





# WalSlide™ exercise station

### Cando<sup>®</sup> WalSlide<sup>™</sup> exercise station original 1.0



- easy assembly, simple to use, takes up very little space
- anchor on track allows for unlimited height positions
- heavy-duty aluminum frame perfect for the professional setting
- has very small wall footprint
- all Cando WalSlide units ship UPS
- band, tubing and handles sold separately

### **Optional Features:**

- overhead exercises can be performed using horizontal overhead section; positions at any height along WalSlide™ and up to 2 feet from the wall
- use optional *pulley attachment system* to perform shoulder stretching and range-of-motion exercises
- pulley system attaches to overhead piece of WalSlide™ and consists of a pulley with heavy duty rope and comfortable handle grips

### Complete WalSlide™ exercise station accessories pulley attachment · contains top overhead section, wall piece, sliding vertical station with sliding band/tubing anchor 10-5088 anchors for both top and wall pieces 10-5091 3 feet 10-5092 4 feet · sliding attachment mechanism allows overhead additional sliding anchor 10-5093 5 feet section to vary in height 10-5094 6 feet 10-5089 horizontal overhead section 10-5096 complete, 3 feet • with band/tubing anchor and sliding attach-10-5097 complete, 4 feet ment mechanism for height adjustable overcomplete, 5 feet 10-5098 head section 10-5099 complete, 6 feet 10-5095 overhead section

# Cando<sup>®</sup> WalSlide 2.0<sup>™</sup> exercise station



- anchor (included) moves along center channel and locks in place by tightening with screwdriver
- adjust anchor to any height along 36" channel
- works with any elastic band or tubing
- ships UPS

10-5070 WalSlide 2.0 10-5071 WalSlide 2.0 add'l anchor


### Cando® shoulder pulley set



- helps increase range of shoulder motion
- can be used for prevention and rehabilitation of shoulder injuries
- · mount in a doorway
- with tubing, handles, pulley, door anchor

	ye	llow		10-	5040	) x-lig	ıht, ea.
red		10-	504	ligh	t, ea.		
	gr	een	n 10-504		5042	2 med	dium, ea.
	b	lue		10-	5043	3 hea	vy, ea.
	bl	ack		10-	5044	1 x-he	eavy, ea.
	yel	red	gre	blu	blk	10-5048	5 pc. set

# Cando<sup>®</sup> over door exercise bar



- use standing, sitting, or balancing
- exercise superficial muscles and the deeper lying muscles
- exercise bar can be used isometrically and dynamically

yellow	10-5060	x-light, ea.
red	10-5061	light, ea.
green	10-5062	medium, ea.
blue	10-5063	heavy, ea.
black	10-5064	x-heavy, ea.
yel red g	re blu blk 10-	5068 5 pc. set



# exercise stations

### Web-Slide™ wall rail

- Web-Slide™ exercise rail is easy to install, use and maintain
- each 28" section has 5 slots to accommodate anchor straps
- perform upper and/or lower band/tubing exercises
- purchase three exercise rail sections (as pictured) for range of motion
- accessories sold separately

### Web Slide exercise rail

10-5090 web slide section 28" 10-5109 full web slide 84"

### anchor straps

10-5320 strap (1 each) 10-5320-10 strap (10 each) 10-5320-50 strap (50 each)

### Lifeline™ wall pulley





Wall Gym 4000


- pulley mechanism positioned to allow for exercises to be performed at floor level, shoulder height and anywhere in between
- includes three graded resistance latex tubes and two handles
- · attach included jogging belt harness and the wall pulley exerciser becomes a treadmill

single channel (2000)	10-2940
double channel (4000)	10-2941
neck strap	10-2942
foot strap	10-2943
ankle/wrist strap	10-2944

### Thera-Band® exercise station



Pro Exercise Station consists of the base [a]. Wall Station consists of the wall piece [b]. Rehab station consists of [a] and [b]



0-1570	pro exercise station [a]
0-1571	wall station [b]
0-1572	rehab station [a]+[b]
0-1573	rehab and wellness station
ccessory	tubing with connectors

color	12" tubing	18" tubing	24" tubing
yel	10-1577	10-1578	10-1579
red	10-1580	10-1581	10-1582
green	10-1583	10-1584	10-1585
blue	10-1586	10-1587	10-1588

assist strap with D-ring 10-1590 handles (2) with D-ring 10-1591 10-1592 padded 3' bar with D-ring

# storage racks

- 50 yard exercise band is made with a core or a cradle (universal)
- · universal racks can accommodate both the core and coreless 50 yard bands and all 6
- · core band rack dispenser can accommodate 50 yard bands that are rolled on a core
- universal tube rack comes with an easy-touse loading mechanism to facilitate the wrapping of tubing onto the dispensing spindles

### band racks (bands sold separately) vertical dispenser with shelf

10-5101 universal, 6 roll 10-5100 core only, 6 roll [a]

### horizontal dispenser with shelf

10-5103 universal, 5 roll 10-5102 core only, 5 roll [b]

### multi-use weight rack

mobile MaxRac™ rack/storage unit

10-5120 without mirror [c] 10-5121 with mirror

# [b]







- 50-yard rolls of resistive exercise band
- rod-in-tube cradle facilitates smooth dispensing

10-5104 single rack (holds 5 rolls) 10-5105 duplex rack (holds 10 rolls)

band rack with 50-yard Cando® band roll set

(1 each: yel red gre blu blk )
10-5229 rack w/ 5 rolls low powder band 10-5629 rack w/ 5 rolls no latex band

10-5949 rack w/ 5 rolls AccuForce® band 10-5106 scissors

### Cando® band wood dispenser rack

- holds five 100-yard Perf-100<sup>™</sup> dispenser boxes or any 50 yard dispenser boxes
- · easiest, most convenient rack for storage and dispensing of band



band rack (5 roll capacity, ships k/d)

10-5111 deluxe rack

10-5480 rack w/ 5 Perf 100 low powder rolls 10-5481 rack w/ 5 Perf 100 no latex rolls





# exercise bands, loops, and kits

### Thera-Band® latex exercise band



 5" band can be used independently or as part of a progressive exercise program

Color	Strength	6 yards	50 yards	30 5-foot pieces
tan	xx-light	10-1016	10-1019	
yellow	x-light	10-1000	10-1006	10-1551
red	light	10-1001	10-1007	10-1552
green	medium	10-1002	10-1008	10-1553
blue	heavy	10-1003	10-1009	10-1554
black	x-heavy	10-1004	10-1010	10-1555
silver	xx-heavy	10-1017	10-1018	
gold	xxx-heavy	10-1080	10-1081* *	*25-yard band
	•			-

# Thera-Band® rehab kits

- rehab kits contain everything you need to enhance recovery with effective home treatment
- each kit includes Thera-Band® resistance products and accessories, BioFreeze® pain reliever, and a comprehensive guide with DVD



<b>lower body:</b> 6' blue exercise tube, 1 loop stirrup, 2 exercise handles, 1 door anchor	10-2303
lower back: 5' red tube with handles, door anchor, BioFreeze® 2 oz. spray, DVD and care guide	10-2305
<b>shoulder</b> : 5' red tube with handles, door anchor, BioFreeze® 3 oz. roll-on, DVD and care guide	10-2306
<b>neck:</b> 5' red tube with handles, 15" red loop, door anchor, BioFreeze® 2 oz. spray, DVD and care guide	10-2307
<b>knee / hip:</b> 12" red loop, BioFreeze® 3 oz. roll-on, DVD and care guide	10-2308
hand / wrist: 5' red tube with handles, door anchor, red FlexBar®, BioFreeze® 3 oz. roll-on, DVD & care guide	10-2309

### Thera-Band® no latex exercise band



- Thera-Band® bands are also available in a latex-free formulation
- latex-free Thera-Band® is 4" wide and comes in 25 yard dispenser packages

Color	Strength	25 yar
yellow	x-light	10-1236
red	light	10-1237
green	medium	10-1238
blue	heavy	10-1239
black	x-heavy	10-1240

### Thera-Band® active recovery kit

- home exercise kit
- packaged for resale, these convenient, professional kits are ideal for a wide range of home exercise programs



### Thera-Band® exercise band sets

- each set is in retail display packaging
- each piece is 5 ft. in length
- set contains one each of the colors indicated below



### exercise band sets

10-1036lightyellowredgreen10-1037heavyblueblack

latex-free exercise band sets

10-1041 light yellow red green 10-1042 heavy blue black

### general rehab & pain relief (Beginner)

 2 Thera-Band<sup>®</sup> resistance bands – yellow & red, 5 ft. pre-cut, BioFreeze<sup>®</sup> pain relieving spray, 2 oz.

10-2301 beginner

### general rehab & pain relief (Intermediate)

• 2 Thera-Band® resistance bands – red & green, 5 ft. pre-cut, BioFreeze® pain relieving spray, 2 oz.

10-2302 intermediate

### shoulder rehab & pain relief

 Thera-Band® shoulder pulley features an innovative design for improved patient feedback, Thera-Band® resistance band – yellow, 5 ft. pre-cut, BioFreeze® pain relieving spray, 2 oz.

10-2300 shoulder



# exercise bands, tubing and loops



### latex exercise tubing

REP®

cord tubing



- tubing is made of natural rubber latex
- resistance levels are obtained by varying the thickness of the Thera-Band® tubing
- color-coded in Thera-Band® colors increasing in difficulty to seven levels from tan through silver
- color-coding system provides positive reinforcement and feedback for gauging results
- available in 100-foot and 25-foot rolls, and retail packs containing one 5-foot piece of each color indicated





1			
NEP Cords	T-man	SEP Cords	
	0:	M	1

- REP® cord tubing is latex-free
- level 1 is lightest, 5 is heaviest

	dispenser boxes			
color/le	vel	25 feet	100 feet	
peach	1	10-5670	10-5675	
orange	2	10-5671	10-5676	
green	3	10-5672	10-5677	
blue	4	10-5673	10-5678	
plum	5	10-5674	10-5679	

### retail packs

- each pack contains 5 foot pieces of tubing
- light <u>yellow</u> red green

heavy blue black

10-1327 light 10-1328 heavy

			~~9	
	color	strength	25 feet	100 feet
	tan	xx-light	10-1310	10-1320
	yellow	x-light	10-1311	10-1321
	red	light	10-1312	10-1322
ĺ	green	medium	10-1313	10-1323
	blue	heavy	10-1314	10-1324
	black	x-heavy	10-1315	10-1325
	silver	xx-heavy	10-1316	10-1326

### Thera-Band<sup>®</sup> latex exercise loops

# Thera-Band

- "lay flat lengths" of 8", 12", and 18" can be used to support a broad range of upper and lower body exercises
- looped bands have same resistance level per color as a single length of Thera-Band<sup>®</sup> exercise band tied in a loop
- each loop is 3" wide

color	strength	8" loop	12" loop	18" loop
yellow	x-light	10-1931	10-1941	10-1951
red	light	10-1932	10-1942	10-1952
green	medium	10-1933	10-1943	10-1953
blue	heavy	10-1934	10-1944	10-1954

### REP® band no latex exercise band



- all REP® bands are latex-free
- $\bullet\,$  REP® bands are 4" in width and come in 6 yard and 50 yard dispenser boxes
- 5 color-coded resistance levels from peach (easiest) to plum (hardest)

color	level	6 yard	50 yard
peach	1: x-light	10-1074	10-1089
orange	2: light	10-1075	10-1090
green	3: medium	10-1076	10-1091
blue	4: heavy	10-1077	10-1092
plum	5: x-heavy	10-1078	10-1093





# hand exercisers

flexion for individual fingers and the entire hand!

### Cando<sup>®</sup> Digi-Flex<sup>®</sup> exerciser

### great for clinical use

- develops isolated finger strength, flexibility, and coordination
- for rehabilitation of carpal tunnel, arthritis, stroke, fractures, tendon injury, nerve lacerations, tennis elbow and more
- · comes with exercise pamphlet







with plastic stand with metal stand

### force to compress finger / hand

yellow	x-light	1.5 lbs. / 5.0 lbs.	10-0740
red	light	3.0 lbs. / 10.0 lbs.	10-0741
green	medium	25.0 lbs. / 16.0 lbs.	10-0742
blue	heavy	7.0 lbs. / 23.0 lbs.	10-0743
black	x-heavy	9.0 lbs. / 31.0 lbs.	10-0744
-	gre blu blk	set of 5	10-0745
yel red	gre blu blk	set with plastic stand [a]	10-0746
yel red	gre blu blk	set with metal stand [b]	10-0749

10-0747 plastic stand unit 10-0750 metal stand unit

### Cando® Digi-Squeeze® exerciser

- develops isolated finger strength, flexibility, and coordination
- for rehabilitation of carpal tunnel, arthritis, stroke, fractures, tendon injury, nerve lacerations, tennis elbow and more
- comes with exercise pamphlet



3 different sizes for small, medium and large hands



different resistances for varied strengthening requirements

stand with product, see below 10-1920 stand only

	small	medium	large
yellow x-light	10-1970	10-1980	10-1990
red light	10-1971	10-1981	10-1991
green medium	10-1972	10-1982	10-1992
blue heavy	10-1973	10-1983	10-1993
black x-heavy	10-1974	10-1984	10-1994
vol rod gro blu blk			

 yel
 red
 gre
 blu
 blk
 set of 5
 10-1975
 10-1985
 10-1995

 yel
 red
 gre
 blu
 blk
 with stand
 10-1976
 10-1986
 10-1996

### Cando® Via® exerciser

### perfect for patient give-aways

- soft finger pads and palm rest
- ergonomic design

yellow

 non-intimidating design appeals to children and seniors



x-light

light

medium

heavy



10-0738.

### force to compress

5.5 lbs. / 22.0 lbs


ilalia	
4.0 lbs.	10-0735
10.0 lbs.	10-0736
16.0 lbs.	10-0737
	4.0 lbs. 10.0 lbs. 16.0 lbs.

black	x-heavy	7.0 lbs. / 28.0 lbs.	10-0739
yel red g	gre blu blk	set of 5	10-0734
yel red	gre blu blk	set of 5 with plastic stand	10-0733


# some exercises for the DigiFlex®, Pro®, Via® and VariGrip® exercisers















# hand exercisers

flexion for individual fingers and the entire hand!



### Cando® VariGrip® exerciser

### individual adjustable finger pads

- each finger individually adjusts
- use "thumb screw" to adjust tension
- palm rest doubles as "full-hand" topper
- soft finger pads and palm rest
- great for post-surgery rehab





adjusts easily by turning thumb screw

min strength	max strength
finger / hand	finger / hand (lbs.)

yellow	x-light	1.0 / 4.0	1.5 / 6.0	10-0540
red	light	2.5 / 10.0	3.0 / 12.0	10-0541
green	medium	4.0 / 16.0	5.0 / 20.0	10-0542
blue	heavy	5.5 / 22.0	7.0 / 28.0	10-0543
black	x-heavy	7.5 / 30.0	9.0 / 36.0	10-0544
yel red	gre blu blk	set of 5		10-0545
yel red	gre blu blk	set of 5 with p	lastic stand	10-0546

### Cando® rubber band hand exerciser

- ideal for finger flexion exercises
- vary resistance by increasing or decreasing number of rubber bands, or by changing to a different color-coded set of bands
- range-of-motion may be decreased by using flexion limiter
- comes with 5 red bands or with 5 sets of color-coded bands (25 bands - 5 of each color)

### exerciser w/ 5 red (light) bands

	latex	latex-free
1 set	10-0800	10-1860
10 sets	10-0800-10	10-1860-10
50 sets	10-0800-50	10-1860-50

### exerciser w/ 25 bands (5 each color)

	latex	latex-free
1 set	10-1800	10-1865
10 sets	10-1800-10	10-1865-10
50 sets	10-1800-50	10-1865-50

### 25 additional bands

		latex	latex-free
tan	xx-light	10-1820	10-1850
yellow	x-light	10-1821	10-1851
red	light	10-1822	10-1852
green	medium	10-1823	10-1853
blue	heavy	10-1824	10-1854
25 band	l set*	10-1825	10-1855
and the second	anne India	Lui e	









### Cando® Pro® exerciser

for larger hands for sports medicine for athletic training

- larger body
- heavier tensions
- 5 color-coded resistances
- soft finger pads
- contoured palm rest



### force to compress

		illiger / Haria	
yellow	x-light	3.0 lbs. / 12.0 lbs.	10-0835
red	light	5.0 lbs. / 20.0 lbs.	10-0836
green	medium	7.0 lbs. / 28.0 lbs.	10-0837
blue	heavy	9.0 lbs. / 36.0 lbs.	10-0838
black	x-heavy	11.0 lbs. / 44.0 lbs.	10-0839
yel red	gre blu blk	set of 5	10-0834
vel red	are blu blk	set of 5 with plastic stand	10-0833

### Cando® fixed resistance grips





sold in pairs

- resistance sequence on grips covers needs from earliest stage of rehabilitation after a stroke (yellow) to advanced athletic training (black)
- force required to squeeze grip doubles with each level

		ibs	
yellow	x-easy	3	10-1804
red	easy	6	10-1805
green	medium	12	10-1806
blue black	hard x-hard	24 48	10-1807 10-1808
yel red	gre blu blk	5 pair set	10-1809





### Digi-Extend® finger extension exerciser





- facilitates extension and flexion exercises
- · isolate extension of all fingers and individual finger joints
- develop isolated finger strength, flexibility, and coordination
- vary resistance by changing to a different color-coded rubber band
- 4 color-coded band resistance levels or joint to work at a comfortable level and build as increased strength is achieved
- for rehabilitation of carpal tunnel, arthritis, stroke, fractures, tendon injury, nerve lacerations, tennis elbow and more...
- · exercise pamphlet included

### Digi-Extend® exerciser

Digi-Extend® unit with color-coded bands and instructions

10-0775 Digi-Extend® unit

### Digi-Extend® clinic set

set contains 4 Digi-Extend® units with color-coded bands and metal stand

10-0757 clinic set w/stand 10-0756 stand only

### additional latex-free bands

 tan
 xx-light
 10-1850

 yellow
 x-light
 10-1851

 red
 light
 10-1852

 green
 medium
 10-1853


 blue
 heavy
 10-1854

25 latex-free bands - 5 each color tan yel red gre blu set 10-1855

# Digi-Extend n' Squeeze®

- combines best features of Digi-Extend® and Digi-Squeeze®
- helps to develop isolated finger strength, flexibility and coordination
- provides complete reciprocal exercise program (flexion and extension)
- ergonomic ball-shape positions fingers and hand for ideal pressure distribution
- 5 color-coded flexion balls accommodate lightest to the heaviest requirements
- 4 color-coded extension bands accommodate lightest to the heaviest requirements
- balls available in 3 sizes to accommodate all hand sizes
- metal stand holds 5 balls
- comes with exercise pamphlet

		small	medium	large
yellow	x-light	10-2270	10-2280	10-2290
red	light	10-2271	10-2281	10-2291
green	medium	10-2272	10-2282	10-2292
blue	heavy	10-2273	10-2283	10-2293
black	x-heavy	10-2274	10-2284	10-2294
yel red g	re blu blk	10-2275	10-2285	10-2295
-	re blu blk	10-2276	10-2286	10-2296
VVILIT	stariu	10-1920	stand only	





# finger exercisers

# Xtensor™ finger extension exerciser



- Xtensor<sup>™</sup> can promote recovery from tennis elbow, tendonitis, effects of a stroke, and hand trauma rehabilitation
- Xtensor™ applies tension when opening the hand using an optimal line of motion replicating the opposite movement to closing the hand
- instructions included

10-0960B blue Xtensor 10-0960Y yellow Xtensor

10-0962 foam cushion pad for Xtensor 2 comfort pads are included with each Xtensor

### HandMaster Plus™



- strengthens hand muscles
- HandMaster Plus<sup>™</sup> can be used to treat: carpal tunnel syndrome, tennis elbow, golfer's elbow, RSI's (repetitive stress injuries), tendonitis, DeQuervain's syndrome, Dupuytren's contracture, osteoarthritis, osteoporosis, stroke, neuropathy, post-surgical rehabilitation, sprain/dislocation/fracture rehabilitation, subluxation and adhesion, and circulatory concerns
- combines flexion and extension
- instructions included

10-0784 soft (purple ball, yellow cord) 10-0785 med. (red ball, black cord) 10-0786 hard (orange ball, gray cord)

10-0787 set of 3 progressive resistances



## hand exercise / stress balls



### Cando® gel hand exercise balls

- hand exerciser for grip strength, dexterity, mobility, and fine and gross motor skills
- heat or chill for use in hot or cold therapy

### standard (circular)

,			
tan	xx-light	10-1490	
yellow	x-light	10-149	
red	light	10-1492	
green	medium	10-1493	
blue	heavy	10-149	
black	x-heavy	10-149	

tan yel red gre blu blk 10-1496

### large ergonomic (egg-shaped)

, ,	, , ,	. ,
tan	xx-light	10-1890
yellow	x-light	10-1891
red	light	10-1892
green	medium	10-1893
blue	heavy	10-1894
black	x-heavy	10-1895
tan vol	od gro blu blk	40 4000





### Cando® memory foam ball exerciser

- memory foam balls are perfect for early stages of hand therapy
- color-coded balls available in 3 levels

		dia	— each—	dozen
yellow	x-easy	1.0"	10-0776	10-0776-12
red	easy	1.5"	10-0777	10-0777-12
green	medium	2.0"	10-0778	10-0778-12
yel red	gre 3 b	all set	: 10-0779	10-0779-12



### Thera-Band® gel balls

balls exercise hands, fingers and forearms

### standard size

otania		
yellow	x-light	10-1330
red	light	10-1331
green	medium	10-1332
blue	heavy	10-1333
black	x-heavy	10-1334

### large size

9		
red	light	10-1341
green	medium	10-1342
blue	heavy	10-1343
black	v hoavy	10.13//



### Eggsercizer® balls

• for exercise and progressive rehabilitation of hand and arm muscles

### individual packaging

10-1290	x-light	orange
10-1291	light	green
10-1292	medium	blue
10-1293	heavy	purple



### foam ball exerciser

 soft foam balls are perfect for the early stages of hand rehabilitation

10-0821	3" each
10-0821-12	3" dozen
10-0823	4" each
10-0823-12	



### foam hand grip

- contoured to fit comfortably in palm
- moderate resistance provides hand, finger and forearm exercise
- sold in pairs

10-0812 foam grip, pair



### Exerball™ exerciser

• has ring to fit around finger

10-0796 soft 10-0797 medium 10-0798 firm





### power hand grip

• adjustable grippers **pro grip** has 25 spring settings from 5-150 lbs.

**deluxe grip** has 15 settings from 10-100 lbs.

**adjustable grip** has 10 settings from 5-50 lbs.

10-0815 pro 10-0813 deluxe 10-0816 adjustable



### big grip

• exercise grip requires flexion to fully compress

10-0791	10 lbs.
10-0792	15 lbs.
10-0793	30 lbs.
10-0794	60 lbs.
10-0790	set of 4



### 6 spring hand grip

- vary tension by using up to six springs
- use for finger, hand, and forearm exercises

10-0801 grip, each



# color-coded resistive hand exercise material

tan yellow red green blue black

### Cando<sup>®</sup> Theraputty™

- standard in resistive hand exercise material
- each color-coded putty has a different consistency ranging from xx-soft for strengthening the weakest grasp to extra-firm
- available in convenient easy-to-open plastic containers (2, 3, 4 and 6 ounce) for patients to use at home and a choice of larger containers (1, 5 and 50 pound) for the hospital or clinic
- Cando® Theraputty™ material is gluten, casein and latex free











standard putty - each										
	tan		yellow		red	green		blue	black	
	xx-soft		x-soft		soft	medium		firm	x-firm	
2 oz.	10-0956		10-0900		10-0901	10-0902		10-0903	10-1466	
3 oz.	10-0957		10-0967		10-0968	10-0969		10-0970	10-1467	
4 oz.	10-0958		10-0905		10-0906	10-0907		10-0908	10-1468	
6 oz.	10-0959		10-0910		10-0911	10-0912		10-0913	10-1469	
1 lb.	10-0994		10-0918		10-0919	10-0920		10-0921	10-1470	
5 lb.	10-0995		10-0923		10-0924	10-0925		10-0926	10-1471	
50 lb.	10-1460		10-1461		10-1462	10-1463		10-1464	10-1473	

putty - set							
complete set							
1 each level (6 total)							
tan yel red gre blu blk							
<b>2 oz.</b> 10-1480							
<b>3 oz.</b> 10-1481							
<b>4 oz</b> . 10-1482							
<b>6 oz</b> . 10-1483							
<b>1 lb.</b> 10-1484							
<b>5 lb.</b> 10-1485							

ех	tra	containers
	emp	ty containers
		and lids
2 oz.	25	10-0940
3 oz.	25	10-0948
4 oz.	25	10-0941
6 oz.	25	10-0942
1 lb.	10	10-0943
5 lb.	10	10-0944

### Cando Plus<sup>®</sup> Theraputty<sup>™</sup> - anti-microbial formula

- same resistance levels and packaging choices as standard Cando® Theraputty™ exercise material
- non-toxic anti-microbial additive resists growth of odor-causing microbes
- anti-microbial putty retains the clean, non-greasy feel that has made Cando® Theraputty™ so popular
- Cando<sup>®</sup> anti-microbial Theraputty<sup>™</sup> material is gluten, casein and latex free









	putty - set								
(	complete anti-microbial set								
	1	each level (6 total)							
	tan	yel red gre blu blk							
2	oz.	10-2606							
3	oz.	10-2616							
4	oz.	10-2626							
6	oz.	10-2636							
1	lb.	10-2646							
5	lb.	10-2656							

Theraputty™ and Theraputty™ stain removal suggestions

remove from hair: Use baby oil or a generous amount of hair conditioner with warm water. Massage around the scalp and area containing the Theraputty™ and then gently comb it through. After Theraputty™ is removed, wash hair with shampoo.

remove from fabric and carpeting: Scrape off excess Theraputty™ with a dull knife or metal spoon. Spray with WD-40 and let stand a few minutes. Respray with WD-40 and wipe off stain with cotton balls. Wipe remaining residue or stain with damp sponge or cloth moistened with liquid dish soap.



		anti-	microbia	al putty -	each	
	tan	yellow	red	green	blue	black
	xx-soft	x-soft	soft	medium	firm	x-firm
2 oz.	10-2600	10-2601	10-2602	10-2603	10-2604	10-2605
3 oz.	10-2610	10-2611	10-2612	10-2613	10-2614	10-2615
4 oz.	10-2620	10-2621	10-2622	10-2623	10-2624	10-2625
6 oz.	10-2630	10-2631	10-2632	10-2633	10-2634	10-2635
1 lb.	10-2640	10-2641	10-2642	10-2643	10-2644	10-2645
5 lb.	10-2650	10-2651	10-2652	10-2653	10-2654	10-2655
50 lb.	10-2660	10-2661	10-2663	10-2663	10-2664	10-2665



# takes putty beyond flexion<sup>sm</sup>

Puttycise® tools can be pushed, pulled or turned through putty to produce specified exercises. Changing the putty resistance will increase or decrease the difficulty level. The puttycise® family consists of 5 tools which all function by inserting the tool base into a resistant putty.

- simulates functional activity
- teaches joint protection techniques
- easy to use, easy to teach



KNOB TURN opening a door, tightening a iar. etc.









**CAP TURN** opening bottle, turning on stove, etc.

### Puttycise® tools

10-2820 5-piece tool set 10-2821 5-piece tool set with carry bag and cards

10-2810 knob turn tool 10-2811 L-bar tool

10-2812 peg turn tool 10-2813 key turn tool 10-2814 cap turn tool

10-2819 carry bag 10-2818 instruction cards

10-2840 putty board

### Puttycise® Theraputty™ kits

each kit includes:

 4 Theraputty<sup>™</sup> resistances (1 each as indicated below)

set of 5 Puttycise® tools (1 each)

· set of instructions

easy	mea	nard
10-2822	10-2828	10-2834
10-2824	10-2830	10-2836
10-2826	10-2832	10-2838
	10-2822 10-2824	10-2822 10-2828 10-2824 10-2830

1 each putty tan yellow red green easy: medium: 1 each putty yellow red green 1 each putty red green hard:

### Cando<sup>®</sup> variable strength Theraputty™

- work pliable pellets into base putty to increase resistance
- pack contains 4 pellets

### base putty with pellet packs

10-2670 21/2 oz. base putty, 1 pack 10-2671 1 lb. base putty, 6 packs 10-2672 5 lb. base putty, 32 packs

### base putty only

10-2680 21/2 oz. base putty 10-2681 1 lb. base putty 10-2682 5 lb. base putty

pellet packs (4 pellets per pack)

10-2690 1 pellet pack

### Cando® scented Theraputty

- · makes therapy more pleasant for youngsters
- · fun way for children to develop grip strength and stay attentive and participate in physical therapy
- tan (vanilla) is xx-soft, yellow (banana) is x-soft, red (cherry) is soft, green (apple) is medium

	tan	yellow	red	green
	vanilla	banana	cherry	apple
	xx-soft	x-soft	soft	medium
2 oz.	10-2760	10-2761	10-2762	10-2763
1 lb.	10-2770	10-2771	10-2772	10-2773
5 lb.	10-2780	10-2781	10-2782	10-2783



10-2797

5 lb.

Cando<sup>®</sup> microwavable Theraputty™

• new formula allows you to warm exercise putty in a microwave


### Cando<sup>®</sup> glitter Theraputty™

- glitter Theraputty<sup>™</sup> sparkles!
- hand therapy becomes visually stimulating

yellow	red	
x-soft	soft	
10-2764	10-2765	2 oz.
10-2774	10-2775	1 lb.
10-2784	10-2785	5 lb.



yellow	blue	
soft	firm	
10-2710	10-2720	2 oz.
10-2711	10-2721	3 oz.
10-2712	10-2722	4 oz.
10.0710	40.0700	•
10-2713	10-2723	6 oz.
10-2714	10-2724	1 lb.
10-2715	10-2725	5 lb.







## hand and wrist exercisers

# Cando® web - finger flexion, extension and wrist exerciser

• perform flexion, extension, opposition and supination exercises

• resistance may be modified by adjusting the hand position, depth of finger insertion, or by switching to a different resistance (color) web

• web is available in both the large 14" diameter size and the small 7" diameter size

· color-coded resistance levels

sets include 1 of each color

large 14" diameter		diameter	low-powder	no-latex
	tan	xx-light	10-0850	10-0870
	yellow	x-light	10-0851	10-0871
	red	light	10-0852	10-0872
	green	medium	10-0853	10-0873
	blue	heavy	10-0854	10-0874
	black	x-heavy	10-0855	10-0875
	tan yel r	ed gre blu blk set	10-0856	10-0876

small 7"	diameter	low-powder	no-latex
tan	xx-light	10-0860	10-0880
yellow	x-light	10-0861	10-0881
red	light	10-0862	10-0882
green	medium	10-0863	10-0883
blue	heavy	10-0864	10-0884
black	x-heavy	10-0865	10-0885
tan yel r	red gre blu blk set	10-0866	10-0886

### multi-resistance 14" large diameter, low powder

yellow green x-light/medium 10-0857 10-0859 light/heavy

### FingerWeights™ exerciser





- · weight belts for each finger
- each weight is adjustable with up to 3 10-gram inserts

10-0470 5pc set colored 10-0471 5pc set white 10-0472 10pc set colored 10-0473 10pc set white

### graded pinch exercisers



- · improve pinch strength and coordination
- includes 5 sets of colorcoded graded resistance pinchpins, three horizontal rods, and a vertical rod for shoulder ROM exercise

### Complete set

10-0840 set / 35 pins

### Replacement pinchpins (7ea.)

10-0841 1 lb 10-0842 2 lb. 10-0843 4 lb. green 10-0844 6 lb 10-0845 8 lb

### Psytech™ wrist exerciser

10-0859 multi-resistance



- wrist flexion and extension exerciser strengthens grip and forearms
- adjustable resistance to 20foot pounds

10-0819 exerciser

### Jux-a-Cisor



 ROM exerciser for finger, hand, wrist, elbow, and shoulder

10-1230 Jux-a-Cisor

### Thera-Band® Progressive Hand Trainer™



- exercise band is held in place by two locking rings
- to change resistance change the sheet color
- full set comes with rings and six elastic sheets (tan through black)

### set

10-1276 all six sheets

### sheet refills (6 each)

10-1270 xx-light tan 10-1271 yellow x-light 10-1272 light red 10-1273 medium green 10-1274 heavy blue 10-1275 x-heavy black



# hand, wrist and body exercisers




### Cando® Twist-Bend-Shake™ bar



- change resistance by varying distance between hands and/or by changing to a different resistance bar
- TWIST: with both hands on the bar, twist by flexing one wrist and extending the other
- BEND: place both hands on or near the ends of the bar. Keep your hands in position and exert force to bend the bar
- TWIST-N-BEND: combine Twist and Bend methods
- SHAKE: place hand(s) at center or end, shake the bar
- OSCILLATION: hold bar at center or end, oscillate in circular motion

		36-inch l
tan	xx-light	10-1530
yellow	x-light	10-1531
red	light	10-1532
green	medium	10-1533
blue	heavy	10-1534
black	x-heavy	10-1535
tan yel r	ed gre blu blk set	10-1536



### Cando® Twist-n-Bend™ bar








- 12" twist
- 12" bend
- lightweight, portable, used for hand, wrist and shoulder
- promotes coordination, increases ROM
- consistent bar diameters throughout resistances means you use the same muscles with each bar
- 24" lengths are easier to bend; ideal for older people because it is easier to bend and twist

		12-inch long	24-inch long
tan	xx-light	10-1510	10-1520
yellow	x-light	10-1511	10-1521
red	light	10-1512	10-1522
green	medium	10-1513	10-1523
blue	heavy	10-1514	10-1524
black	x-heavy	10-1515	10-1525
tan yel re	ed gre blu blk set	10-1516	10-1526

flexible mesh 12" bar cover 10-1519

### Bodyblade® exerciser



- improve strength, balance, coordination, posture, and core stabilization with low level vibration and inertia
- delivers progressive resistance

		length
10-1540	Pro, black	60"
10-1541	Classic, black	48"
10-1542	CXT, yellow	40"
10-1543	Lite, yellow	32"
10-1544	Cardio, yellow, 2 ea.	32"
10-1545 10-1546	Small Stand, for 6 bla Large Stand, for 26 b	

<sup>\*</sup>each unit comes with a wall chart and DVD

### Flexi-Bar® stick



- improve strength, balance, coordination, posture, and core stabilization with low level vibration and inertia
- soft grip in center of bar
- easy-to-use exerciser
- portable

10-1960	oscillator rod 1
10-1961	oscillator rod 2
10-1962	oscillator rod 3

### Thera-Band® FlexBar®





- improve strength, balance, coordination, posture, and core stabilization with low level vibration and inertia
- wide scope of application broad applications for both rehab and sports specific training
- versatile ability to work bending, twisting, and oscillation movements with one piece of equipment
- progressive resistance provides easy and immediate motivation
- portable easy to carry/use anywhere

yellow	x-light	10-1350
red	light	10-1351
green	medium	10-1352
blue	heavy	10-1353



# the Cuff®

# wrist and ankle weights

### the **DRIGINAL** Cuff® weight

- flexible weight contours to snugly fit the ankle, wrist... and even thigh
- long closure strap for comfortable secure fit
- heavy-duty materials and double-stitching
- color-coded weight levels
- marked with both pounds and kilograms
- grommeted for hanging on any weight rack

### original Cuff® individual weights

1/4 lb.	white	10-0200
1/2 lb.	walnut	10-0201
3/4 lb.	orange	10-0202
1 lb.	blue	10-0203
11/2 lb.	olive	10-0204
2 lb.	white	10-0205
2½ lb.	red	10-0206
3 lb.	gold	10-0207
4 lb.	turquoise	10-0208
5 lb.	black	10-0209
6 lb.	beige	10-0210
7 lb.	lemon	10-0211
7½ lb.	orange	10-0212
8 lb.	red	10-0213
9 lb.	parchment	10-0214
10 lb. 12½ lb. 15 lb. 20 lb. 25 lb.		10-0215 10-0216 10-0217 10-0218 10-0219

### original Cuff® weight sets

### Functional set (7 pieces)

1 each: 1, 2, 3, 4, 5, 7½, 10 functional set 10-0250 set with rack 10-0251

### Economy set (20 pieces)

2 each: ¼, ½, ¾, 1, 1½, 2, 2½, 3, 4, 5 economy set 10-0252 set with rack 10-0253

### Standard set (24 pieces)

2 ea.: ½, ½, ¾, 1, 1½, 2, 2½, 3, 4, 5, 7½, 10 standard set 10-0254 set with rack 10-0255

### Deluxe set (32 pieces)

2 ea.: 1/4, 1/2, 3/4, 1, 11/2, 2, 21/2, 3, 4, 5, 6, 7, 71/2, 8, 9, 10


deluxe set 10-0256 set with rack 10-0257


### Knee-Hip-Leg set (8 pieces)

2 each: 10, 12½, 15, 20

knee-hip-leg 10-0258







### mobile weight rack



- heavy-duty steel rack helps keep weights and accessories organized
- holds up to 500 pounds on the 2-sided pegboard and tray
- 3-inch wheels allow rack to be easily moved (even over carpet)
- assembly required; assembled rack measures 24" x 35" x 48"
- ships UPS

10-0150 mobile weight cart





# wrist and ankle weights



### the **DELUXE** Cuff® weight

- the Deluxe Cuff® weight has all the features of the original Cuff® weight PLUS an extra-long loop strap with extra hook tabs
- extra long loop strap allows for firm closure on the largest wrist and ankle and facilitates use on larger areas such as the thigh
- extra hook tabs add additional fastening points and prevent extra-long strap from dangling

### Deluxe Cuff® individual weights

⅓ lb.	white	10-2500
½ lb.	walnut	10-2501
¾ lb.	orange	10-2502
1 lb.	blue	10-2503
1½ lb.	olive	10-2504
2 lb.	white	10-2505
2½ lb.	red	10-2506
3 lb.	gold	10-2507
4 lb.	turquoise	10-2508
5 lb.	black	10-2509
6 lb.	beige	10-2510
7 lb.	lemon	10-2511
7½ lb.	orange	10-2512
8 lb.	red	10-2513
9 lb.	parchment	10-2514
10 lb.	brown	10-2515
12½ lb.	olive	10-2516
15 lb.	tan	10-2517
20 lb.	blue	10-2518
25 lb	green	10-2519

### Deluxe Cuff® weight sets

### Functional set (7 pieces)

1 each: 1, 2, 3, 4, 5, 71/2, and 10 10-2550 functional set 10-2551 set with rack

### Economy set (20 pieces)

2 each: 1/4, 1/2, 3/4, 1, 11/2, 2, 21/2, 3, 4, and 5 10-2552 economy set set with rack 10-2553

### Standard set (24 pieces)

2 each: 1/4, 1/2, 3/4, 1, 11/2, 2, 21/2, 3, 4, 5, 71/2, and 10 standard set 10-2554

set with rack 10-2555

### Deluxe set (32 pieces)

2 ea: 1/4, 1/2, 3/4, 1, 11/2, 2, 21/2, 3, 4, 5, 6, 7, 71/2, 8, 9,10

10-2556 deluxe set

set with rack 10-2557

Knee-Hip-Leg set (8 pcs) 2 each: 10, 121/2, 15, 20

knee-hip-leg set 10-2558

### what makes the **DELUXE** Cuff® weight different from the **DRIGINAL** Cuff® weight?





### adjustable weight



- weight is varied by inserting or removing the individual weight packets
- consists of 5 weight packets and weight
- for example, a 5-lb. Pouch™ weight includes five 1-lb. weight packets



change number of weight packets to adjust the total overall weight

### The Pouch™ weight sets

2½ lb.	red	10-0300
5 lb.	black	10-0301
10 lb.	brown	10-0302
15 lb.	tan	10-0303
20 lb.	blue	10-0304
25 lb.	green	10-0305







- easy-to-grip, easy-to-store color-coded Cando® exercise Wate™ bars come in weights ranging from 1 to 25 lbs.
- · use in supine, sitting or standing position to increase strength and ROM
- use to improve endurance, balance and motor skills
- use in therapy, group classes, aerobics, yoga and pilates

# weighted exercise bars

- anti-roll end cap to prevent Wate™ bar from rolling around and becoming a safety hazard









stand	ard wat	e™ ba	r - eacl
tan	10-1600	1 lb.	.45 kg
tan stripe	10-1601	1½ lb.	.68 kg
yellow	10-1602	2 lb.	.91 kg
yellow stripe	10-1603	2½ lb.	1.14 kg
red	10-1604	3 lb.	1.36 kg
red stripe	10-1605	3½ lb.	1.59 kg
green	10-1606	4 lb.	1.82 kg
green stripe	10-1607	4½ lb.	2.05 kg
blue	10-1608	5 lb.	2.27 kg
blue stripe	10-1609	6 lb.	2.73 kg
black	10-1610	7 lb.	3.18 kg
black stripe	10-1623	7½ lb.	3.41 kg
silver	10-1611	8 lb.	3.64 kg
silver stripe	10-1612	9 lb.	4.09 kg
gold	10-1613	10 lb.	4.55 kg

### heavy wate™ bar - each

10-1614	blue	121/2 lb. 5.68 kg
10-1615	blue	15 lb. 6.82 kg
10-1616	blue	17½ lb. 7.95 kg
10-1617	blue	20 lb. 9.09 kg
10-1618	blue	221/2 lb.10.23 kg
10-1619	blue	25 lb.11.36 kg

# wate™ bars - sets iped: half pounds 틸 solid: Introductory set (3 pieces) 1 each 1, 2, 3 10-1624 introductory set Lite™ set (5 pieces)

1 each 1, 2, 3, 4, 5

10-1620 Lite™ set Standard set (10 pcs)

1 ea. 1, 2, 3, 4, 5, 6, 7, 8, 9, 10 10-1621 standard set

Deluxe set (16 pcs)

1 ea. 1, 2, 3, 4, 5, 6, 7, 8, 9, 10,12½,15,17½,20,22½,25 10-1622 deluxe set

mini wate™ bars Mini-bars • 14" long, dumbbell size • same features as the larger bars • sold in pairs

yellow	10-1650	1 lb. each (pair)
blue	10-1651	11/2 lb. each (pair)
green	10-1652	2 lb. each (pair)
blue	10-1653	21/2 lb. each (pair)
black	10-1654	3 lb. each (pair)

### exercise bar and tubing combo




- foam covered exercise bar
- · tubing attached at each side of bar
- wind up tubing to increase resistance
- · yellow tubing included

10-2979 combo exercise bar

### weighted rol-up™ exercise bar





- use to exercise upper and lower body
- resistance is increased by rotating the bar ("winding-up" the tubing)
- standard with 8 ft. strip of mediumweight green Cando® resistive exercise tubing and 2 tube Klips™

10-2981	1 lb.
10-2982	5 lb.
10-2983	10 lb.
10-2984	15 lb.
10-2985	20 lb.
10-2986	25 lb

### wate™ bar racks and carts









[o]		[5]
10-1640 [A 10-1641	-	k wall rack ok wall rack
10-1645 10-1643 [B		vertical wall rack vertical wall rack
10-1642 [C	[] 8 bar h	norizontal wall rack
10-1630 [[	)] 15 slot	t cart (16x13x12")




# hand-held weights

# SoftGrip<sup>™</sup> hand weights

- comfortable SoftGrip<sup>™</sup> hand weight is an alternative to standard dumbbells
- easier to grasp and less intimidating than traditional dumbbell
- · comfortable, soft, not floppy
- color-coded to indicate weight
- stores on standard Cuff® weight rack
- made in USA

_	Thade in 6071					
		lb.	individual	weight pairs		
	tan	1/2	10-0350	10-0350-2		
	yellow	1	10-0351	10-0351-2		
	red	11/2	10-0352	10-0352-2		
	green	2	10-0353	10-0353-2		
	blue	21/2	10-0354	10-0354-2		
	black	3	10-0355	10-0355-2		
	silver	4	10-0356	10-0356-2		
	gold	5	10-0357	10-0357-2		
	tan	6	10-0358	10-0358-2		
	yellow	7	10-0359	10-0359-2		
	red	$7\frac{1}{2}$	10-0360	10-0360-2		



### SoftGrip™ hand weight sets


Lite™ Set (10-piece) 2 each: 1, 2, 3, 4, 5 lbs.

10-0383

Standard Set (16-piece)

2 each: ½, 1, 1½, 2, 2½, 3, 4, 5 lbs.

functional set





# Cando® vinyl coated cast-iron dumbbells

- vinyl coated, color-coded cast iron dumbbells are ideal for upper body exercise
- · coated weight does not scratch floor
- · coating gives weight comfortable feel
- clearly marked in pounds and kilograms

lb.		individual	weight pairs
1	pink	10-0550	10-0550-2
2	purple	10-0551	10-0551-2
3	green	10-0552	10-0552-2
4	silver	10-0553	10-0553-2
5 6 7 8	blue pink purple black	10-0554 10-0555 10-0556 10-0557	10-0554-2 10-0555-2 10-0556-2 10-0557-2
9 10 15	yellow orange	10-0558 10-0559 10-0560	10-0558-2 10-0559-2 10-0560-2
20	not coated	10-0561	10-0561-2

### weight sets Lite™ set (10 pieces)

2 each: 1, 2, 3, 4, 5 lbs. 10-0563 economy 10-0564 with wall rack 10-0578 with floor rack

### standard set (20 pieces) 2 each: 1 through 10 lbs.

10-0565 standard 10-0566 with 2 wall racks 10-0579 with floor rack

### weight racks

wall rack (holds 10 each) 10-0575 wall rack

floor rack (holds 10 pair) 10-0577 floor rack

### studio rack (holds 1100 pounds)

- mobile studio rack has casters
- holds 10 columns of dumbbells

10-0581 studio rack



# Cando<sup>®</sup> Handy Grip<sup>™</sup> balls

- achieve the benefits of traditional weight training without grasping metal dumbbells
- · ergonomic shape makes for easy grasping and textured surface adds grip
- enables a variety of strength training activities and functional movements
- · secure ball to hand with adjustable strap for firm grasp
- available in 1 to 5 lb. weights

tan	1 lb.	10-0440
yellow	2 lb.	10-0441
red	3 lb	10-0442
green	4 lb.	10-0443
blue	5 lb.	10-0444













### Cando<sup>®</sup> Wate™ hand ball



- · color-coded weighted balls
- small enough and flexible enough to grasp with one hand, yet big enough to hold with two hands
- increase in weight from 1.1 lb. through 6.6 lbs. (0.5 to 3 kgs.) while maintaining a constant 5" (12.7cm) diameter
- use a bicycle pump to inflate/deflate to make ball easier to grasp

### constant diameter (5" diameter)

		(	,
tan	½ kg	1.1 lb.	10-3160
yellow	1 kg	2.2 lb.	10-3161
red	1½ kg	3.3 lb.	10-3162
green	2 kg	4.4 lb.	10-3163
blue	2½ kg	5.5 lb.	10-3164
black	3 kg	6.6 lb.	10-3165
tan vel r	ed are blu	hlk set	10-3166

### wate racks

10-3167 2-tier holds 3 balls each tier 10-3168 1-tier holds 6 balls

### Thera-Band® Soft Weights™



- Soft Weights<sup>™</sup> are designed with a soft cover and pliable filler free of any hard edges or rims
- each weight is 41/2" in diameter
- use for strengthening, plyometrics, balance training, mobilization and functional grip strength development
- supports one or two hand activities
- comfortable to hold
- · easily fits in gym bag

### constant diameter (41/2" diameter)

3150
3151
3152
3153
3154
3155

tan yel red gre blu blk set 10-3156

# medicine balls

### leather medicine ball



· available with increasing diameter (traditional medicine balls) or with constant 71/2" diameter (plyometric balls)

### increasing diameter

10-3062	2 kg	4.4 lb.	9" dia.
10-3063	3 kg	6.6 lb.	9½" dia.
10-3064	4 kg	8.8 lb.	10" dia.
10-3065	5 kg	11 lb.	10½" dia.
10-3066	6 kg	3.2 lb.	11" dia.

### constant diameter (71/2" diameter)


10-3090	2 kg	4.4 lb.	
10-3091	3 kg	6.6 lb.	
10-3092	4 kg	8.8 lb.	
10-3093	5 kg	11 lb.	
10-3094	6 kg	13.2 lb.	

### Cando® dual handle medicine ball

- for abdominal strengthening, rotational movement and upper and lower body strength training
- dual, soft grip handle is easy to grip with one or two hands
- handles enable wide ranges of movements
- each ball is 9" / 23cm diameter
- · color-coded to indicate weight

kg / lb					
yellow	2.7 / 6	10-3180			
red	5.5 / 12	10-3181			
green	8.2 / 18	10-3182			
blue	10.9 / 24	10-3183			
black	13.6 / 30	10-3184			









### Cando® vinyl-coated kettlebell





- for rehabilitation, strength training and toning
- vibrant colors for easy weight identification
- · coating reduces potential floor scratching
- available in convenient weight increments

	weig			diameter		
	(kg / l	b)	(	cm / inch)		
yellow	2.27 /	5		20.3 / 8	10	-3191
red	3.4 / 7	7.5		20.3 / 8	10	-3192
green	4.5 /	10		22.9 / 9	10	-3193
blue	6.8 /	15		22.9 / 9	10	-3194
black	9.1 / 2	20		25.4 / 10	10	-3195
silver	9.1 / 2	25		27.9 / 11	10	-3196
gold	13.6 /	30		27.9 / 11	10	-3197
yel red	gre blk gld	set of	5 (5,	7.5, 10, 20, 30	) 10	-3199
studio rack, stationary 660 pound capacity				10	-3198	


studio rack, stationary 660 pound capacity

**Fabrication Enterprises** 

# balls, rebounders and racks







- balanced, uniform weight
- color-coded for easy identification
- constructed with rubber to bounce

	weight	diameter	
	(kg / lb)	(cm / inch)	
tan	.45 / 1	20.3 / 8	10-3140
yellow	.90 / 2	20.3 / 8	10-3141
red	1.8 / 4	20.3 / 8	10-3142
green	3.2 / 7	22.9 / 9	10-3143
blue	5.0 / 11	22.9 / 9	10-3144
black	6.8 / 15	25.4 / 10	10-3145
silver	9.1 / 20	27.9 / 11	10-3147
gold	13.6 / 30	27.9 / 11	10-3148
yel red	gre blu blk set	(2, 4, 7, 11, 15 lb.)	10-3146

### Cando<sup>®</sup> P.T. pliable medicine balls



- textured tactile surface is easy to grip, grasp and catch
- increase/decrease air in ball to enhance grasp
- use for upper body rehabilitation and fitness exercises
- · designed to throw against plyometric rebounders

	weight (kg / lb)	diameter (cm / inch)	
tan	.45 / 1	12.7 / 5	10-3170
yellow	.90 / 2	12.7 / 5	10-3171
red	1.8 / 4	12.7 / 5	10-3172
green	3.2 / 7	17.8 / 7	10-3173
blue	5.0 / 11	17.8 / 7	10-3174
black	6.8 / 15	22.9 / 9	10-3175
silver	9.1 / 20	22.9 / 9	10-3177
gold	13.6 / 30	22.9 / 9	10-3178
yel red gr	e blu blk set (2,	4, 7, 11, 15 lb.)	10-3176

### adjustable angle rebounder

- returns ball with speed, force and angle consistent with throw
- angle of inclination adjusts



### Cando® rectangle rebounder

- angle adjust: 0-60 degrees (any angle)
- 39" diagonal surface (27.5" x 27.5")

**square rebounder (left)** 10-3110 square rebounder

with 5 PT balls: 2,4,7,11,15 lbs. 10-3132 with 5 balls 10-3131 with 5 balls & vertical rack

### Cando® round rebounder

- angle adjust: 0, 28, 33, 39 & 42 degrees
- with floor tray to store balls
- 27" diameter rebounder surface

round rebounder (right) 10-3112 round rebounder

with 5 PT balls: 2,4,7,11,15 lbs.

10-3136 with 5 balls

10-3135 with 5 balls & vertical rack

### 5 ball stand

- store 5 medicine balls
- made from high strength steel tubing
- small foot-print: 12" x 12"

vertical ball rack 10-3120

horizontal ball rack 10-3121

\*vertical rack shown with rubber medicine balls (sold separately)

\*minimal assembly required for all rebounders and medicine ball racks - tools included

\*ships by UPS







# weights and racks

### chest weight pulley systems

- · chest pulleys are an important part of many rehabilitation programs; these pulleys come with handles, weights, hardware indicated
- handle is mounted at chest height on single handle unit; chest and floor height on double handle unit; and chest, floor and shoulder height on triple handle unit (ceiling height of 10 feet is needed for triple handle); choose from single or dual tower set-ups
- select weight level by inserting "pin" into weight stack at increments of 2 pounds
- exercisers are wheelchair accessible
- handles can be adjusted to any starting height; adjust the rope slip to lengthen or shorten the pulley rope



single tower / single handle




single handle



dual tower / double handle



triple handle



chest pulley accessories

① chest pulley

2 floor pulley 3 overhead (arm) pulley single handle chest pulley

single weight tower 10-0660

dual weight tower 10-0667

double handle chest/floor pulley

single weight tower 10-0661

dual weight tower 10-0668

triple handle arm/chest/floor

single weight tower 10-0662

dual weight tower 10-0669

10-0673 curl bar 10-0674 twin ball

10-3220 webbing handle (pair)10-0664 metal handle (pair)10-0663 2 pound weight

each weight tower contains 10 pounds of weights (5 x 2-lb weights)

### leg/ankle exerciser



- for lower leg and ankle exercises
- use weights to control level
- use for load-resisting inversion, eversion, dorsi-flexion and plantar-flexion exercises
- optional accessories include: two 5 lb. and four 2½ lb. disk weights

10-1110 leg/ankle exerciser10-1119 with 20 lbs. of weights

### parallel dip



- parallel dip is wheelchair accessible
- designed to improve upper body strength, especially the muscles of the deltoids and triceps

10-0691 EN-304 exerciser

### iron weight plates / racks





### individual weight plates (free weights)

10-0600 1½ pound 10-0601 2½ pound 10-0602 5 pound 10-0603 7½ pound 10-0604 10 pound 10-0606 25 pound

weight plate sets

10-0601-4 10 pound (4 ea. 21/2 lb.)

weight rack

10-0640 stationary rack 10-0616 mobile rack




# bi-directional pedal exercisers



### Cando® pedal exerciser





- portable lower and upper body exerciser
- use when exercise bicycle is not available or is contra-indicated
- wide leg spread and rubber non-marring feet keep exerciser stable when patient has little muscle control
- use wrap around adjustable straps to secure feet or arms

### adjustable resistance

 tension knob turns clockwise towards "High" to increase difficulty or counter clockwise towards "Low" to reduce difficulty



### bi-directional pedal exercisers

10-0710 preassembled

10-0711 assembly required (KD)

### digital exerciser



- digital counter keeps track of amount of exercise performed
- folds for easy storage
- used for both upper and lower body exercise
- pedal action is bidirectional

10-0712 digital

### exerciser with handle



- handle provides added stability while using pedaler
- four anti-slip rubber pads prevent sliding and protect floor
- safe and gentle form of low impact exercise
- stimulates circulation
- ideal for toning leg muscles
- tension can be adjusted
- attractive chrome finish
- handle can be removed easily for storage

10-0713 with handle

### Cando® chair cycle for clinic and home use







- for individuals who cannot sit safely on a stationary bike
- resistance adjusts from 0-29 lbs. (0-13 kgs.)
- length can be adjusted to accommodate various leg lengths
- ideal for hospital, nursing home, or clinical settings
- deluxe model retains features of standard model PLUS unit has adjustable pedals

### pedal exercisers

10-0720 Cando® standard

10-0722 Cando® deluxe w/ adjustable pedals

# optional accessories for deluxe exerciser

10-0724 upper body kit 10-0725 long leg brace





10-0720

standard



### vestibular stepper



- double chambered inflatable stepper can be used for lower extremity exercises, balance training and vestibular orientation
- use to improve coordination and motor skills
- two chambers are connected by a valve which adjusts the air-flow (level of difficulty)

30-1935 stepper





# Magneciser™ ROM / strength exercise systems

### ultra quiet magnetic resistance for smooth operation

- easy to assemble
- bi-directional operation
- variable resistance (turn knob to change resistance)
- improve circulation, muscle strength, joint range-of-motion and coordination
- 5-function digital display standard on all units
- wall mountable units include mounting hardware



5-function digital fitness computer display: speed, time exercised, calories burned, distance "traveled" and scan

### upper and lower extremity pedal exerciser table top upper extremity wheelchair accessible • for muscles and joint of shoulder, arm, wrist and upper body • arm length adjusts from 13" (33 cm) to 22" (56 cm) optional non-skid table-top pad • for upper or lower body exercises 01-8030 pedal 10-0714 table-top 10-0714PAD table pad 3 in 1 exerciser vertical shoulder exerciser rotation/supination (shoulder wheel) with 3 attachments: wrist shoulder constructed with a elbow steel base for searm length cure wall mounting adjusts 13" (33 cm) to **₹** arm length adjusts 22" (56 cm) 13" (33 cm) to 22' height ad-(56 cm) justs height adjusts • adjustable height allows for sitting, standing, or wheelchair access • arm length adjusts from 13" (33 cm) to 22" (56 cm) easy access for wheelchair 10-0715 10-0716 3 in1 01-8025 shoulder upper/lower



# pedal and arm exercisers / UBE and LBE

#### pedal and arm exercisers



#### Cando® magneciser™

- quiet magnetic, bi-directional resistance system
- for upper/lower body
- comes with a fully automatic 5-function digital display showing speed, elapsed time, calories, distance and scan
- adjustable base for stability and storage
- small and portable

01-8030 magneciser



# Cando® deluxe pedal exerciser

- bi-directional exerciser can be used for upper/lower body
- resistance can be adjusted with one easy-to-use dial
- fully automated electronic monitor with LED display
- lightweight with non-skid footpads to ensure safety

01-8010 deluxe pedal



#### **Hudson® UBE**

- UBE 2000 is a one-direction, lightweight, portable exerciser that offers ultra-quiet electromagnetic resistance for smooth, maintenance-free use
- eight resistance levels
- computer provides feedback for pulse rate, speed, trip distance, total distance, elapsed time, and calories burned
- comfortable pedal handles with oval gripping surface are more ergonomic than traditional round handles

10-2870 UBE 2000



# Hudson® motorized table

- UBE-T table has motorized height adjustment
- wired remote enables push button elevation adjustment
- height adjusts from 27-39"
- adjustable brackets located on tabletop enables table to work safely with many UBEs
- table can be used with many chairs and wheelchairs

\*pictured with UBE (sold separately)

10-2880 UBE table



#### **Hudson® UBE-BDP**

- bi-directional, belt-driven electromagnetic resistance
- ergonomically designed handles
- wireless heart telemetry
- electronically selectable work-load range, displayed in watts
- 5 programmed workout templates and manual mode
- 5 heart rate control workouts

10-2873 UBE-BDP



#### **Hudson® UBE-BD**

- bi-directional UBE-BD has same sturdy build and features as the UBE 2000
- work arms and shoulders in two directions for a better, more balanced workout
- from fully welded steel frame to the onetouch computer, it is designed to withstand abuse in any clinical setting
- silent, smooth to operate
- eight workload levels for a conditioning or rehabilitation workout

10-2871 UBE-BD



#### **Hudson® UBE-PAT powered**

- bi-directional passive assist trainer facilitates the user's upper extremity range of motion while providing motorized passive assistance exercise
- three motorized speeds (15, 30, and 45 rpm), and three program run times (15, 30, and 45 minutes)
- LCD shows elapsed time, distance, speed, and pulse rate. Wireless remote control allows therapist to change time and/or speed of PAT without having to touch device or user

10-2872 UBE-PAT



# endorphin upper and lower body exercisers

All Endorphin units come standard with:

- pulse meter that displays distance, time, speed, calories and pulse
- UBE comes with comfort hand grip
- LBE standard with pedals
- UBE comes with body stabilizing chest pad:
  - adjustable to vertical or horizontal position
  - with 7" range of motion

	bi-direc-	chest	external		, , programmable	work	display screen		
	tional pedal		power 110 volt	nandles	resistance	levels WATTS	pulse meter	workout selection	screen size
e1	Y	Υ	N	comfort	8 manual		Υ	N	small
e2	Y	Y	Y	comfort	8 electronic plus manual	19-48	Υ	Υ	large
е3	Y	Υ	Y	comfort	8 electronic plus manual	6-20	Υ	Υ	large
e4	Υ	Υ	Υ	multi- position	8 electronic plus manual	19-48	Υ	Υ	large

#### **UBE 300**



- use on table or counter
- unit is 28" long by 16" wide
- multi-position grip has 3 preset hand positions. It also can be set to simulate rowing

#### comfort grip

10-3600

[e3] 10-3610 [e2] 10-3605

multi-position grip

10-3617

#### **UBE 355 on stand**



- · with height adjustable platform
- crank arm adjustable from 36"-50" from floor allowing unit to be used in seated or standing position

#### comfort grip

10-3620 [e1] 10-3630

[e2] 10-3625

multi-position grip

10-3637



- · with height adjustable platform
- semi-recumbent seat adjusts front to rear for proper user positioning; rotates 360° to provide easy access and egress
- seat assembly may be quickly removed for wheelchair accessibility

#### comfort grip

[e1] 10-3660 둭

10-3670 [e3]

[e2] 10-3665 🚐 multi-position grip

10-3677 [e4]

#### **LBE 300**



- · place unit on floor, in front of chair, bench or wheelchair

•	unit	is 28"	long by	y 16"
		pe	dal	
[e	1]	10-	3640	
[e	3]	10-	3650	
[e	2]	10-	3645	
		pla	te	
[e	1]	10-	3641	
[e	3]	10-	3651	
[e	2]	10-	3646	

#### **LBE + UBE 370**



- · comes with upper and lower ergometer
- height adjustable platform
- semi-recumbent adjustable seat

#### pedal

	poul.	
[e1]	10-3680	comfort grip 📮
[e3]	10-3690	comfort grip 📮
[e2]	10-3685	comfort grip 📮
[e4]	10-3697	multi-position grip 📮
	plate	
[e1]	10-3720	comfort grip 🕞
[e3]	10-3730	comfort grip 🕞
[e2]	10-3725	comfort grip 🚐
[e4]	10-3737	multi-position grip

#### accessories





foot pedal



multi-position handgrip



10-2925 pedal 10-2921 plate 10-2927 multi-position grip

10-2922 quad grip 10-2923 table 📮



# exercise bikes, elliptical, treadmill

# Hudson® upright or recumbent exercise bikes 2200 3300 display for 2200 and 3300 exercise bikes magnetic resistance fitness computer displays time,


- 8-level magnetic resistance; very low startup load
- fitness computer displays time, speed, distance, calories, pulse
- contact HR monitor
- adjusts for heights 4'10"-6'4"
- 300 pound weight capacity

upright 10-6000 UE 2200 recumbent 10-6020 RE 2200

- speed, distance, calories, pulse
- contact&wireless HR monitor
- manual mode, 3 hill profiles, 3 interval training styles, 6 levels of auto HR control
- 3 position seat adjustment (for upright only)
- adjusts for heights 4'10"-6'4"
- 300 pound weight capacity

upright 10-6001 unidirectional 10-6002 bidirectional recumbent 10-6021 unidirectional

10-6022 bidirectional



- 5-300 watts of work load
- fitness computer displays time, speed, distance, calories, pulse

3200 / 3500

• 350 pound weight capacity

upright 10-6003 HEC3200 recumbent 10-6023 HEC3500

- · ergonomically designed
- belt driven with self-powered generator system

400


- · fitness computer displays time, speed, distance, calories, watts, and pulse
- hand grip or wireless pulse
- 5-400 watts of work load
- 350 pound weight capacity

upright 10-6004 HEC400 recumbent 10-6024 HEC400R

#### Hudson® elliptical bike

- exercise upper and lower body, and heart
- high inertia drive system
- fitness computer displays time, speed, pedal RPM, distance, calories, watts, and pulse
- user programs include manual mode, 2 hill, mountain profiles, 3 styles of interval training, 5 levels of automatic heart-rate control
- · hand grip and wireless pulse
- 300 lbs. capacity

10-6040 EL3300 elliptical





#### treadmill

- extra long 22"x60" running surface
- 425 lb user weight capacity
- 0-12mph speed, 0-15 incline
  - contact hear rate sensors
    - accessory tray
      - safety key stop
      - blue matrix LCD display; 1/4 mile track display
      - red LED window displays speed, distance, time, calories, incline, pulse and more
      - 12 built-in programs include manual, preset, custom and heart rate control
        - light commercial warranty



additional treadmills available, see website or contact us





#### round & half-round foam rollers















- for positioning, balance, postural and muscle re-education, spinal stabilization, body awareness and coordination, and ranging and strengthening activities
- half-round roller for ankle stretching and knee rehabilitation; use half roll as rocker board by standing on flat side; diameter refers to the width of the half-round profile
- foam rollers and half rollers are strong enough to be used by large adults and active children
- available in 48", 36" and 12" lengths with either a round (circular) or half-round (semi-circular) profile
- can be used alone or with other exercisers such as Cando® exercise band, tubing, weights and balls

#### which style is right for you?

- white rollers best seller, economic choice for quality roller
- *TufCoat*<sup>™</sup> foam rollers for *washability*, *firmness* and *durability*
- blue foam rollers have polyethylene foam like white rollers
- EVA foam rollers have an attractive finish and firm feel
- black foam rollers made of molded foam for increased firmness





- · white foam rollers are economic, easy to use, and versatile
- 6" diameter roller is our most popular
- new 8" Jumbo roller is extra sturdy and extra large!

6" white			8'	8" white jumbo	
6x36	30-2100	round	8x36	30-2260	
6x36	30-2110	half-round	8x36	30-2270	
6x12	30-2101	round	8x12	30-2261	
6x12	30-2111	half-round	8x12	30-2271	

4" white			3" white	
4x36	30-2102	round	3x36	30-2106
4x36	30-2112	half-round	3x36	30-2116
4x12	30-2103	round	3x12	30-2107
4x12	30-2113	half-round	3x12	30-2117

- TufCoat™ rollers have coated exterior that is washable, durable and adds an extra layer of firmness
- interior is standard white foam

	6" Tufcoat™			4" Tufcoat™
6x36	30-2230	round	4x36	30-2235
6x36	30-2240	half-round	4x36	30-2245
6x12	30-2231	round	4x12	30-2236
6x12	30-2241	half-round	4x12	30-2246

- blue foam rollers use same industry standard polyethylene foam as white foam rollers
- attractive blue color keeps roller looking new

6" blue						
round	30-2150					
half-round	30-2152					
round	30-2151					
half-round	30-2153					
	round half-round round	round 30-2150 half-round 30-2152 round 30-2151				

- EVA rollers are made from an extra firm foam
- · hard foam maintains soft feel
- · attractive mottled design

6" EVA foam						
6x36	round	30-2200				
6x36	half-round	30-2210				
6x12	round	30-2201				
6x12	half-round	30-2211				

• black molded foam rollers are very hard / strong

6" black composite						
6x36	round	30-2280				
6x36	half-round	30-2290				
6x12	round	30-2281				
6x12	half-round	30-2291				



### round & half-round foam rollers





# Cando<sup>®</sup> Plus antimicrobial foam roller

- PE foam made with anti-microbial additive to protect against microbes and odors
- brown tint differentiates anti-microbial foam from standard foam

#### full round, roller

30-2330 6"x36", round 30-2332 6"x12", round

#### half-round, roller

30-2331 6"x36", half-round 30-2333 6"x12", half-round



- extruded as a half circle, not sliced in half; this process gives the flat side a skin
- half-round roller does not bow or warp

#### half round, extruded roller

30-2340 6"x36", half-round 30-2341 6"x12", half-round

#### star roller



- durable white polyethylene foam
- octagonal gear-shape prevents unit from rolling around the floor
- multiple units interlock like gears

30-2220 star roller 6"x36" 30-2221 star roller 6"x12"



#### 6" wide white foam roller



#### new sizes available!

- 48" x 6" for extra length applications
- 30" x 6" inexpensive to ship performs as 36" roller
- 24" x 6" performs as 12" or 36" roller
- 18" x 6" easy to store

#### full round, roller

30-2140 6"x30", round 30-2141 6"x24", round 30-2142 6"x18", round 30-2109 6"x48", round

#### half-round, roller

30-2143 6"x30", half-round 30-2144 6"x24", half-round 30-2145 6"x18", half-round 30-2108 6"x48", half-round

#### foam roller storage



- convenient wall storage for 6 foam rollers
- polished stainless steel construction
- load at the top, remove from the bottom
- ships unassembled
- 36" x 10" x 40", 17 lbs

30-2180 foam roller wall storage

- convenient storage for 12 foam rollers
- heavy duty PVC with casters ideal for fitness studio
- 24" x 24" x 30"



30-2181 upright foam roller storage

#### inflatable rollers



- ball rollers inflate like inflatable ball
- fully inflate to increase instability, and partially inflate for added stability

#### inflatable ball roller

30-2080 17" L x 7" D - green 30-2080 30" L x 7" D - blue 30-2081 28" L x 9" D - black

#### foam roller covers



- convenient covers slip over a foam roller
- covers are washable and protect foam rollers against dirt and nicks

#### full round cover

30-2090 6" L x 36" D, round 30-2091 4" L x 36" D, round

#### half-round cover

30-2094 6" L x 36"D, half-round 30-2095 4" L x 36"D, half-round



all components of the clinic MVP set hang on wall for easy storage!



MVP™ balance system

#### complete MVP™ set

10-1904 complete set

#### set includes:

10 ball wall rack, 10 balls - 2 each level, 16" diameter board, 20" diameter board, 30" diameter board

#### additional platforms and balls

#### platforms

10-1776 16" diameter, hard 10-1786 20" diameter, medium 30" diameter, easy 10-1796 20"x40" board 10-1696

- · circular boards have 3 ball-insets
- rectangular board has 5 ball-insets
- boards have slots to tie exercise band or tubing for a more dynamic vestibular challenge

#### 2 ball sets

C	hallenge	height	pair
yellow	x-easy	1"	10-1760-2
red	easy	1.5"	10-1761-2
green	moderate	2"	10-1762-2
blue	hard	2.5"	10-1763-2
black	x-hard	3"	10-1764-2

#### 10 ball sets

10 ball set

10-1766

yel red gre blu blk 2 each

10 ball set with rack yel red gre blu blk 2 each

10-1743

# 10 ball rack

10 ball rack only


10-1744





# Multi-Axial platform system





weight rack with five color-coded instability balls: yellow, red,

weight rod (for any of the five sockets on either side of board)

Professional System pictured above

green, blue and black

2.5 lb. weight plates

[c]

[ d ]

two-sided professional Multi-Axial board

motion limiter complete with clips and padding

disk weights and a wall rack. Optional motion limiter is used to limit ROM.



#### Multi-Axial Platform System - PRO

10-1730	board with balls, rods
10-1731	board with balls, rods, weights, rack
10-1735	hoard with halls rods weights tub

#### instability balls (hemispheres)

	challenge	height	each
yellow	x-easy	1"	10-1760
red	easy	1.5"	10-1761
green	moderate	2"	10-1762
blue	hard	2.5"	10-1763
black	x-hard	3"	10-1764

yel red gre blu blk	5 ball set	10-1765
yel red gre blu blk	5 ball set with rack	10-1742

Multi-Axial board only

2.5 lb. weight plate

weight rod (each)

4 each 2.5 lb. weight plate

#### instability ball wall rack

accessories

wall rack for 5 balls and board

10-1732

10-1741

10-1733

10-0601

10-1770

10-1780

10-0601-4

10-1740

#### professional instability balls model can can be used for be hung for easy storage!

both Multi-Axial Platform System and MVP™



#### Dycem foot pad for added traction Multi-Axial Platform System - Home

storage tub

motion limiter



- board is elliptical with off-centered pivot
- · select board for specific foot, and weight class

adult model (over 120 pounds) - blue home board for right foot 10-1752 10-1753 home board for left foot

child model (under 120 pounds) - yellow 10-1750 home board for right foot 10-1751 home board for left foot

non-skid accessory

10-1780 Dycem foot pad for added traction



#### Cando® heavy-duty balance boards


- · strengthen and condition lower body
- improve ankle ROM, balance, body awareness, and develop joint mobility and reaction skills
- · durable, heavy-duty and non-skid
- · lightweight, compact for portability

#### round wobble (15" diameter)

10-1103 intermediate, 2" ball

#### 10-1099 advanced 3" ball

rectangular rocker (15x18") 10-1100 basic: 2" flat bottom 10-1101 intermediate: 21/2" bottom 10-1102 advanced: 31/2" bottom



#### Cando® modular stability trainer

#### Cando<sup>®</sup> Board-on-Stone™ balance trainer







- combines the vestibular disc and non-slip board to make a stability trainer
- stack discs or use smaller board for added difficulty
- stand on disc and shift weight to perform stability exercises

#### board and stability disc (35cm disc)

10-1705 16" advanced 20" intermediate 10-1706 10-1707 30" beginner

- · use to challenge balance skills
- to increase difficulty inflate stone, use larger stone, or smaller board
- place non-slip board on top of stone
- stone can be placed flat side up or down

10-1725 7" stone and 16" board 10-1726 7" stone and 20" board 7" stone and 30" board 10-1727

10-1722 13" stone and 16" board 10-1723 13" stone and 20" board 10-1724 13" stone and 30" board


#### boards only - use with disc and balance stone

10-1700 16" circular 10-1701 20" circular 10-1702 30" circular

see page 45 for disc and stones

#### Cando® economy wobble board

balance boards



- 14" diameter wobble board makes balance training fun
- · raised circles on surface makes the board slip resistant

10-1745 wobble board

#### Chango® balance boards



- 2-in-1 board (rocker and wobble) improves balance and stimulates muscles
- 4" center ball is for aggressive training, 3" ball for beginner

rectangular board (left) 10-1716 with 4" ball 10-1717 with 3" + 4" balls

diameter circular board (right) 10-1715 with 3" + 4" balls

#### Chango® balance paws





- attach to feet, use while walking/standing for balance and proprioceptive training
- attach to hands for strengthening of shoulder and back muscles

10-1710 advanced 2" (pair) 10-1711 beginner 11/2" (pair)

#### foam eva balance boards



- · soft cushioned balance pad
- · low-impact balance training

rectangular 30-2170

rectangular beginner board

circular

30-2171 circular intermediate board



### balance trainers

#### each pad in the series of 5 progressive balance pads has a different level of "give" and offers a different vestibular challenge

- yellow is the least unstable (easiest), and black is the most unstable (most challenging) surface
- made from molded foam
- combination of foam and cut-out design yields comfort and functionality
- pads can also be used as core-strengthening vestibular seating cushions

pads measure 20" diameter patent pending



		each —	2-each
yellow x-e	asy	30-2161	30-2161-2
red eas	sy	30-2162	30-2162-2
green mo	derate	30-2163	30-2163-2
blue har	d	30-2164	30-2164-2
black x-h	ard	30-2165	30-2165-2
yel red gre blu	blk set	30-2166	30-2166-2


#### Core-Tex<sup>™</sup> balance platform



- Core-Tex<sup>™</sup> encourages movement and change of direction to maximize the integration of the body's receptors from the foot and ankle all the way up to neck, eyes and inner ears
- available with and without hand rail

Core-Tex™ standard 30-2190 30-2191 Core-Tex™ with hand rail

#### Fitter® wobble boards



- square rocker board has beveled bottom for single plane of motion
- circular boards have adjustable ball
- balance poles optional

#### square (rocker board)

10-1104 easy 20"

10-1133 2 poles

#### circle (wobble board)

medium 10-15°, 20" 10-1128 10-1129 advanced 15-20°, 16" 10-1132 set of 3 boards

#### Thera-Band® rocker boards



- · anti-slip base design with a sensory stimulation surface to offer a new performance level to balance training
- · light weight molded plastic


#### rocker board (front-back)

10-1184 rocker

wobble board (all directions)

10-1185 wobble

#### vestibular board



- rocking on board stimulates balance control and vestibular responses
- boards are covered with non-slip surface
- side-to-side or end-to-end motion

#### rocker board 30"x 60"x 12"

32-2020 side-to-side 32-2021 forward/back 32-2022 combination 🕰

rocker board 18"x 18"x 5"

32-2023 combination





#### ArmaSport® balance pad



- · free of phthalate plasticizers, latex and triclosan
- · sanitized coating protects against microbes and unpleasant odors
- resists sweat and saliva
- · closed cell foam provides excellent elasticity and energy absorption
- coarse surface for traction on bottom, top surfaces

32-1500B 1 pad 32-1500B-10 10 pads

# balance pads

#### Cando® low impact aerobic pad

- · designed to be the most effective and functional pad for low impact aerobic exercises and group classes
- both top and bottom surfaces are textured for non-slip, secure footing
- economically priced
- colorful pads make the perfect cushioned exercise surface in any setting
- available in Y (yellow), R (red), G (green), B (blue) and BLK (black)

pads measure 20"diameter, extra thick 21/2" depth



	each	case (10 pads) -
yellow	30-2120Y	30-2120Y-10
red	30-2120R	30-2120R-10
green	30-2120G	30-2120G-10
blue	30-2120B	30-2120B-10
black	30-2120BLK	30-2120BLK-10

#### Airex® balance pads



- balance pads for vestibular training and exercises
- · the difficulty level can be increased by stacking pad
- individual and group use

pad 16" x 20" 30-1910 21/4" thick 30-1910-20 20 per case

balance pad plus (w/non-slip backing) 30-1915 21/4" thick 30-1915-20 20 per case

extra-large pad 16" x 40" x 21/4" 30-1908 XL each 30-1908-10 XL case-10

#### Airex® balance wedge



balance wedge is used for balance training, coordination and reaction training, motor-skill training to maintain balance, and posture stability training

wedge 20" x 16" x 21/2" 30-1917 21/2" thick (each) 30-1917-2 21/2" thick (pair) 30-1917-40 case of 40

#### Airex® balance beam



- 64" soft foam trapezoidal beams can be positioned wide-side down (91/2") for greater stability or narrowside down (41/2") for a greater vestibular challenge
- unstable base helps improve upper/lower body stabilization and coordination
- for core conditioning, rehabilitation, gait and balance

30-1913 balance log (ea.) 30-1913-10 case of 10

#### Thera-Band® stability trainer





- used to improve balance, postural stability, proprioception, and coordination
- unsteady surface challenges body to maintain posture and balance
- softer pad provides more vestibular challenge

30-2131 green, firm 30-2131-2 green (2 ea.)

30-2130 blue, soft 30-2130-2 blue (2 ea.)

#### inflatable

black, x-soft 30-2132 30-2132-2 black (2 ea.)



# vestibular balance discs & wedges



#### Cando® vestibular disc - 60 cm / 23.6" dia.





large/standard disc



- large vestibular disc can be used as a balance board or sitting disc, depending on how much it is inflated
- inflates with standard pump
- tactile side offers user feedback and stimulation

#### 60 cm diameter (23.6 inch) disc

 30-1868R
 red

 30-1868G
 green

 30-1868B
 blue

 30-1868Y
 yellow

 30-1868BLK
 black

30-1867 60 cm, cover

#### Cando® vestibular disc - 35 cm / 13.8" diameter











flat side

nubby side

- inflatable vestibular disc mimics both the movement and shape of the inflatable ball when used on any seat
- use it on the floor as a standing disc for balance training, proprioception and strengthening of the lower extremities
- each disc has one nubby side for tactile feedback and stimulation and one flat side

#### 35 cm diameter (13.8 inch) disc

30-1870R red 30-1870G green 30-1870B blue 30-1870Y yellow 30-1870BLK black

30-1869 35 cm cover

#### Cando® balance stones



- stone can be used to practice and reinforce coordination, proprioception, and balance
- stones can be filled with water or sand to add weight for added safety and stability
- stones can be used with the curved side down as a vestibular balance "board"

#### 1 each

30-1982 7" diameter 17½ cm 30-1983 13" diameter 33 cm

set of 6

30-1982-6 7" diameter 17½ cm 30-1983-6 13" diameter 33 cm

#### Cando® vestibular wedge



- inflatable wedge requires active participation from user for correct sitting posture
- dynamic seating and static positioning
- inflate or deflate wedge to vary difficulty

30-1925 large wedge (15x15") 30-1926 small wedge (10x10")

#### Bosu™ balance exerciser

- · versatile balance trainer
- use as a board or a ball
- inflate ball to change in-



stability; 25" diameter

30-1900 w/ 1 hr video, 124 pg. manual

30-1901 w/ ½ hr video, 4 pg. manual

#### Cando® core trainer dome



- inflatable dome on rigid base can be used flatside up or down
- dome has 21" diameter
- complete with pump, resistance tubing and instructional DVD

30-1902 dome







# inflatable exercise balls

#### Cando® inflatable exercise balls - standard



- use to improve balance, coordination, flexibility and strength, or even just for fun
- used by therapists as an aid for vestibular movement and equilibrium therapy, by fitness
  professionals as an aid for stretching and strength programs, as seats, and even by
  pregnant women in birthing classes
- non-slip surface is ribbed for extra security
- under inflate to give a soft, mushy feel or fully inflate for a firm, bouncy feel
- supports 300 lbs. (136kg)

retail hox

retail box			
30-1801B	45cm	17.7in.	yellow
30-1802B	55cm	21.6in.	orange
30-1803B	65cm	25.6in.	green
30-1804B	75cm	29.5in.	red
30-1805B	85cm	33.5in.	blue



polybagge	ed		711
30-1800	30cm	11.8in.	blue V
30-1801	45cm	17.7in.	yellow
30-1802	55cm	21.6in.	orange
30-1803	65cm	25.6in.	green
30-1804	75cm	29.5in.	red
30-1805	85cm	33.5in.	blue
30-1806	95cm	37.4in.	red
30-1841	105cm	41.3in.	blue
30-1807	120cm	47.3in.	orange
30-1808	150cm	59.1in.	yellow
30-1809	100	70.9in.	

#### Cando® deluxe ABS ball



- deluxe balls are made with a specially formulated slow-deflate material
- non-slip surface ribbed for traction
- "mottled" color differentiates these from standard balls
- supports 600 lbs. (272kg)

#### retail box

30-1851B	45cm	17.7in.	yellow
30-1852B	55cm	21.6in.	orange
30-1853B	65cm	25.6in.	green
30-1854B	75cm	29.5in.	red
30-1855B	85cm	33.5in.	blue

#### polybagged

45cm	17.7in.	yellow
55cm	21.6in.	orange
65cm	25.6in.	green
75cm	29.5in.	red
85cm	33.5in.	blue
	55cm 65cm 75cm	45cm 17.7in. 55cm 21.6in. 65cm 25.6in. 75cm 29.5in. 85cm 33.5in.

#### Cando® super-duty exercise balls



- strongest, thickest, heaviest Cando® inflatable exercise ball
- made of extra heavy duty phthalate free material
- non-slip surface is ribbed for extra security
- under inflate to give a soft, mushy feel or fully inflate for a firm, bouncy feel
- supports 600 lbs. (275 kg)

cm	in.		retail box	polybagged
45	17.7	yellow	30-1961B	30-1961
55	21.6	orange	30-1962B	30-1962
65	25.6	green	30-1963B	30-1963
75	29.5	red	30-1964B	30-1964
85	33.5	blue	30-1965B	30-1965

#### Cando® sensi-ball



- sensi-ball has a specially textured "nubby" surface that adds tactile and visual elements to ball therapy
- used to improve balance, coordination, flexibility and strength, or even just for fun
- supports 300 lbs. (136kg)

cm	in.		polybagged
55	21.6	orange	30-1773
65	25.6	green	30-1776
75	29.5	red	30-1777
85	35.5	blue	30-1778
100	39.4	red	30-1779

#### Cando® body ball



- strong, light, easy to grasp and have a large inflation range
- attractive pastel colors

cm	in.	polybagged
55	21.6	30-1742
65	25.6	30-1743
75	29.5	30-1744
85	35.5	30-1745
95	38.0	30-1746



### inflatable exercise rolls & balls



#### Cando® inflatable saddle rolls



- easier to control because motion is limited to forward and backward direction
- aids in developing muscles, coordination and balance
- the "saddle" offers extra support for the patient
- supports up to 300 lbs.

	cm-HxL	ın-HxL	
30-1725	40x90	16x35	yellow
30-1726	50x100	20x39	orange
30-1727	60x110	24x43	green
30-1728	70x120	28x47	red
30-1729	80x130	32x51	blue



Cando® inflatable

- straight roll only moves in a side-to-side manner making it easier to use than the ball
- · added length of roll makes it possible for two people to be on one roll
- · 300 lb. weight capacity
- measures height x length

cm-HxL in-HxL 30-1780 40x90 16x35 yel 30-1781 50x100 20x39 ora 30-1782 60x110 24x43 gre 30-1783 70x120 28x47 red

#### Cando® jump ball



- jump ball is ideal for proprioceptive and vestibular stimulation
- ball has heavy-duty grip handles designed to give secure control at all times
- · available in four sizes
- · 300 lb. weight capacity

	cm	in	
30-1825	40	15.8	yellow
30-1826	45	17.7	red
30-1827	50	19.7	green
30-1828	55	21.7	blue

#### Cando<sup>®</sup> sensi-saddle roll



- textured "nubby" surface adds tactile and visual elements to ball therapy
- easier to control than balls
- built-in saddle gives user added sense of security
- supports up to 300 lbs. cm/in

30-1735 40/15.8 yellow 30-1736 50/19.7 orange 30-1737 60/23.6 green

#### Cando<sup>®</sup> Cushy Air™ ball



- · soft, cushy material inflates to desired firmness
- light and easy to grasp
- 7-10" diameter in yellow, red, green, blue and black
- 13-17" diameter in yellow

hand ball cm in color 30-1740\* 25 7-10 Y,R,G, B,BLK

30-1741 45 13-17 Y

\* specify color when ordering: Y-yellow, R-red, G-green, B-blue, BLK-black

#### Cando® ball with feet



- stabilizing "feet" keep ball from rolling when not in use
- weight of person makes the "feet" retract into ball so it can function like any other inflatable exercise ball
- non-slip surface is ribbed for extra security
- supports up to 300 lbs.

cm in 30-1891 45 17.7 yellow 30-1892 55 21.0 orange 30-1893 65 25.6 green

30-1894 75 29.5 red

#### Cando<sup>®</sup> economy ball sets



· economy ball set come with ball, 6" inflator pump, and chart in a retail box

30-1845 55cm ball/ pump 30-1846 65cm ball/ pump 30-1847 75cm ball/ pump

#### ball sizing chart

٠.	gested		
ball diameter		body h	eight
45cm	17.7in.	under 56in.	142cm
55cm	21.6in.	57in 63in.	143cm -160cm
65cm	25.6in.	64in 72in.	161cm -182cm
75cm	29.5in.	73in 79in.	183cm -200cm
85cm	33.5in.	over 80in.	200cm

#### inflation

- ball size reflects maximum inflation
- maximum size may be realized by finishing the inflation process after letting the ball rest 24 hours





### inflatable ball chairs & exercisers

#### exercise ball stool/trainer



- versatile exercise ball chair holds a 50cm ball (included) securely in place
- use as a chair, or use as an inflatable exercise dome
- casters have locking mechanism to make chair stationary
- lock casters when using during vestibular exercises . . . you do not want exercise ball chair to roll away
- includes ball

30-1796 mobile exercise ball stool

#### plastic ball chair

- mobile ergonomic seating system combines the benefits of ball therapy with the convenience of a chair base
- ball held securely in place by the base
- features removable back, without arms
- child and adult sizes available
- casters glide on carpet and hard surfaces
- includes ball; colors may vary
- · light assembly required



#### plastic chair without arms

30-1792 adult size w/ 55 cm ball 30-1795 child size w/ 35 cm ball

30-1797 locking casters, 2 each

#### Thera-Band® exercise balls





balls are boxed, color-coded and include an exercise chart

- regular (left) standard exercise ball
- PRO Series SCP™ formula (right)
   offers slow deflation and enhanced
   performance

cm	in.		standard
45	17.7	yellow	30-1881
55	21.6	red	30-1882
65	25.6	green	30-1883
75	29.5	blue	30-1884
85	33.5	silver	30-1885
cm	in.		PRO SCP™
<b>cm</b> 45	<b>in.</b> 17.7	yellow	<b>PRO SCP™</b> 30-1876
		yellow red	
45	17.7	•	30-1876
45 55	17.7 21.6	red	30-1876 30-1877

#### PhysioGymnic ball and roll



- premium Italian balls and rolls
- 300 pound weight capacity

#### inflatable molded vinyl BALL

30-1713	12 (30cm) pink	
30-1700	18" (45cm) yellow	
30-1701	22" (55cm) orange	
30-1702	26" (65cm) green	
30-1709	30" (75cm) red	
	30" (75cm) red 34" (85cm) blue	
	,	
30-1703	34" (85cm) blue	

#### inflatable molded vinyl ROLL

30-1720 12" (30cm) blue 30-1721 16" (40cm) red 30-1722 22" (55cm) yellow 30-1723 28" (70cm) blue 30-1724 34" (85cm) red

#### metal ball chair



- mobile ergonomic seating system combines the benefits of ball therapy with the convenience of a chair base
- ball held securely in place by the base
- features removable back
- available with or without arms
- · casters glide on carpet or hard surfaces
- includes 55cm ball
- light assembly required

#### metal ball chair

30-1791 with arms 30-1790 without arms



### exercise balls and accessories

# Cando® dynamic ball stabilizer base

#### base has 4 built-in anchors and 2 grooves for exercise band and tubing

- lightweight, portable gym ball base measures 19.5" diameter and fits 45, 55, 65 and 75 cm exercise balls
- use base to stabilize exercise ball and prevent ball from rolling around gym
- · use base with band or tubing for dynamic workout

30-1848 dynamic base



Cando® ball stabilizer bases




- plastic base provides extra protection and stability for the ball during use
- ideal for individuals who are not prepared for instability of therapy ball
- prevents ball from rolling away

30-1840 30-1842 small (for 45cm, 55cm) medium (for 65cm, 75cm)

#### PVC ball rack 30-1833

30-1831







· efficiently store balls for quick and easy access

#### wall mount

wall mount, 64"W x 18"D x 2" H (3 ball) 30-1831

stationary

floor stationary, 62"W x 20"D x 42" H (6 ball) 30-1833

mobile

30-1832 floor mobile, 62"W x 20"D x 72" H, 3 shelf rack with 3" swivel casters

please inquire about other rack sizes

30-1832



- for optimal inflation, use electric first then fully inflate with manual pump
- double piston and electric pump have inflation and deflation valve
- colors may vary

#### manual pumps

double piston (a) 30-1051 6" hand (b) 30-1047 12" hand (c) 30-1048 30-1038 small bellows (d) 30-1039 large bellows (e)

#### electric pumps

30-1054 inflator/deflator (f) 30-1829 ball plug remover

#### Cando® ball stacker



- efficient way to store balls while conserving floor space
- stack 35cm 95cm balls
- set includes 3 adjustable rings
- stack balls up to three high

30-1830 3 ring stacker set





# premium exercise mats

#### ArmaSport® mats are better...

- tough enough for clinic, school and fitness center
- extra thick and extra soft
- maximum breaking acceleration is better than permitted standards (tested with DIN EN 12503 and the TUV GS certificate)
- free of phthalate plasticizers, latex and triclosan
- sanitized treatment provides reliable, permanent protection against bacteria and fungal attack, without the use of triclosan
- sanitized coating protects against microbial attack and unpleasant odor
- resists sweat and saliva
- closed-cell structure of ArmaSport® mats do not absorb moisture
- ribbed for traction during use and ventilation when stored
- mats are tear proof, scratch resistant and do not take a compression set





#### six sizes, three colors and three thicknesses:

- 24"x48" Fit mat, 0.4" and 0.6" thick
- 24"x56" Body mat, 0.4" and 0.6" thick
- 24"x72" Star mat, 0.4" and 0.6" thick
- 32"x72" Top mat, 0.6" thick
- 40"x72" Elite mat, 0.6" thick
- 40"x78" Power mat, 0.6" and 1.0" thick

#### which thickness is right for you?

- 0.4" thick mats for pilates, yoga, home fitness
- 0.6" thick mats for fitness, studios, clubs and competitions
- 1.0" thick mats for physical therapy and school sport



#### more reasons ArmaSport® mats are better:

- hang for storage using integral flexible grommets (0.6" fit, body, star and top mats)
- closed-cell rubberized latex-free foam provides excellent elasticity, energy absorption and resistance
- use in the water during hydrotherapy





# premium exercise mats









ArmaSport®

storage rack to hang mat

adjustable width mat hanger


32-1491

wall rack, holds 10 mats

32-1490

floor rack with casters, holds 20 mats





<sup>\*</sup>with flexible grommets for hanging




# essional exercise line

# exercise, flotation, yoga mats

- Airex® mats are lightweight, comfortable and do not slide on the floor
- made from closed-cell foam with an "integrally molded skin" for extra toughness
- transverse ribbing pattern is molded into the mat to resist tearing and to give the mat enough flexibility to be "rolled-up" for carrying and storage
- available either individually boxed or in case packaging
- ideal for individual workouts
- tough enough for use in fitness centers and therapy clinics


<b>size</b> 1"x78"x39" 1"x78"x39"	name/color Hercules, green Hercules, blue	individually boxed 32-1230G 32-1230B	case pack (quanting 32-1230G-6 (6) 32-1230B-6 (6)	ty)
5/8"x78"x48"	Atlas, green	32-1234G	32-1234G-10 (10)	
5/8"x78"x48"	Atlas, red	32-1234R	32-1234R-10 (10)	
5/8"x72"x39"	Corona, green	32-1236G	32-1236G-10 (10)	
5/8"x72"x39"	Corona, red	32-1236R	32-1236R-10 (10)	
5/8"x72"x39"	Corona, blue	32-1236B	32-1236B-10 (10)	
5/8"x72"x23"	Coronella, green	32-1238G	32-1238G-10 (10)	
5/8"x72"x23"	Coronella, red	32-1238R	32-1238R-10 (10)	
5/8"x72"x23"	Coronella, blue	32-1238B	32-1238B-10 (10)	
2/5"x72"x23"	Fitline 180, aqua	32-1247A	32-1247A-15 (15)	
2/5"x56"x23"	Fitline 140, aqua	32-1248A	32-1248A-20 (20)	
3/5"x48"x23"	Fitness 120, blue	32-1240B	32-1240B-20 (20)	
3/10"x75"x23"	Piloga, purple	32-1232P	32-1232P-15 (15)	
3/10"x75"x23"	Piloga, black	32-1232BLK	32-1232BLK-15 (15)	





hercules blue

atlas red





coronella green

corona blue





fitline 180/140



Airex® Corona mat in the water!



Airex® Fitline 180 great for family exercise



Airex® Fitline 180 used for yoga/pilates



Airex® Hercules being used by trainer and athlete




#### exercise mats



- cushioned exercise mats offer superior comfort, quality and durability, yet are lightweight and sturdy enough for commercial use in health and fitness clubs
- colorful, non-leaching, soft, resilient closed-cell foam material
- waterproof and sanitary; can even be used in pool
- latex-free

26" wi	individual mat de x 72" long x .6" thick	bulk pack case of 1		
yellow	30-2310Y	30-2310Y-10		
red	30-2310R	30-2310R-10		
green	30-2310G	30-2310G-10		
blue	30-2310B	30-2310B-10		
black	30-2310BLK	30-2310BLK-10		
	<b>a</b>	. ,		

#### Cando® individual exercise/rest mats



- therapy mat is available with or without centerfold
- choose polyurethane (PU) foam or EnviroSafe™ foam
- foam mats are covered with 18oz. tough fire-retardent phthalate-free PVC coated vinyl: cover is sewn around foam

38-2412

#### center-folding mat

wxl	1" polyurethane	2" polyurethane	2" EnviroSafe™				
2x4 feet	38-0500	38-0400	38-2400				
2x5 feet	38-0501	38-0401	38-2401				
2x6 feet	38-0502	38-0402	38-2402				
non-fol	non-folding mat						
wxl	1" PU IFD100	2" PU IFD100	2" EnviroSafe™				
2x4 feet	38-0510	38-0410	38-2410				
2x5 feet	38-0511	38-0411	38-2411				

#### EnviroSafe™ foam mats are an eco-friendly alternative:

• 100% toxin free, zero VOC's or off gas, hypoallergenic, latexfree, mold and mildew proof, 100% recyclable with no leaching chemicals

38-0412

medium-firm

2x6 feet 38-0512





#### Cando® eco friendly PER premium mat



#### eco friendly premium yoga Mat


- · lightweight, yet cushioned and strong
- 0.25" thick mat offers extra padding with the convenience of being thin enough to roll-up and take anywhere
- mats roll up into compact roll to fit in gym bag
- PER (polymer environmental resin) yoga mats comply with the European Environmental Protection test standard
- no phthalates, dioxin, furons, phenol, or heavy metals
- non-skid mat surface offers traction during exercise and prevents mat from sliding on flooring
- mats measure 68" long x 24" wide x 1/4" thick

individual mat bulk pack case of 10

24" wide x 68" long x 1/4" thick 30-2401B

30-2401B-10

#### Thera-Band® personal exercise mats



- lightweight, portable, enhanced tear resistance
- new fabric-like finish offers distinct comfort and durability
- rolls up to carry or store

	size	color	each
32-1380B	24"x75"x0.6"	blue	
32-1380G	24"x75"x0.6"	green	
32-1380R	24"x75"x0.6"	red	
32-1383B	24"x75"x1"	blue	
32-1383G	24"x75"x1"	green	
32-1383R	24"x75"x1"	red	
32-1386B	40"x75"x0.6"	blue	
32-1386G	40"x75"x0.6"	green	
32-1386R	40"x75"x0.6"	red	



#### exercise mats

- mats are covered with 18oz. tough fire-retardent phthalate-free PVC coated vinyl; cover is sewn around foam
- vinyl coated nylon cover is strong, mildew resistant, anti-bacterial and wipes clean with a damp cloth
- mat is made with square side walls and is reinforced at all stress points for added durability
- PE and PU Foam Options:
  - poly-ethylene (PE) foam 1 3/8" thick PE mats have a hard feel, used for more strenuous exercises
  - poly-urethane (PU) foam 2" thick PU mats have a softer feel, used when cushioning is most important

#### Cando® accordion fold mat







- mat folds every 2 feet to minimize storage space when not in use
- velcro ends allow mats to safely connect to other mat to "expand" and form a larger mat surface
- mat folds in 2 foot sections for easy storage and transportation

WxL	number sections	color selection	1 3/8" PE foam	2" PU foam	1 3/8" EnviroSafe™	2" EnviroSafe™
4x4 feet 4x4 feet	2 2	1 color 2 color	38-0060 38-0070	38-0000 38-0010	38-2060 38-2070	38-2000 38-2010
4x6 feet	3	1 color	38-0061	38-0001	38-2061	38-2001
4x6 feet	3	2 color	38-0071	38-0011	38-2071	38-2011
4x8 feet	4	1 color	38-0062	38-0002	38-2062	38-2002
4x8 feet	4	2 color	38-0072	38-0012	38-2072	38-2012
4x8 feet	4	rainbow	38-0082	38-0022	38-2082	38-2022
5x4 feet 5x4 feet	2 2	1 color 2 color	38-0093 38-0094	38-0090 38-0091	38-2093 38-2094	38-2090 38-2091
5x10 feet	5	1 color	38-0063	38-0003	38-2063	38-2003
5x10 feet	5	2 color	38-0073	38-0013	38-2073	38-2003
5x10 feet	5	rainbow	38-0083	38-0023	38-2083	38-2023
6x12 feet	6	1 color	38-0064	38-0004	38-2064	38-2004
6x12 feet	6	2 color	38-0074	38-0014	38-2074	38-2014
6x12 feet	6	rainbow	38-0084	38-0024	38-2084	38-2024

• EnviroSafe<sup>™</sup> foam option - EnviroSafe<sup>™</sup> foam mats are the eco-friendly alternative with the following features: 100% toxin free, zero VOC's or off gas, hypoallergenic, latex-free, mold and mildew proof, 100% recyclable with no leaching chemicals. EnviroSafe<sup>™</sup> foam available in 1 3/8" thick medium firm, and 2" thick softer feel. EnviroSafe<sup>™</sup> material conforms to California Fire Retardent code 117 T.B.

#### Cando<sup>®</sup> Add-a-Mat<sup>™</sup> extension panels



- increase the length of any accordion fold mat by adding an "add-a-mat" panel
- velcro ends allow for easy attachment to any accordion fold mat

#### Add-A-Mat extension sections (panels)

WxL	1 3/8" PE foam	2" PU foar
4x2 feet	38-0120	38-0100
5x2 feet	38-0121	38-0101
6x2 feet	38-0122	38-0102

1 3/8" EnviroSafe™	
38-2120	38-2100
38-2121	38-2101
38-2122	38-2100 38-2101 38-2102



## exercise mats



#### Cando® center-fold and non-folding mat

# centerfold



- · reinforced handles can be used to carry and store mat
- see information about foam options and mat covers on page 54
- choose color from chart on page 54

#### center fold mat with handles

WxL	1 3/8" PE foam	2" PU foam	1 3/8" EnviroSafe™	2" EnviroSafe™
4x4 feet	38-0315	38-0305	38-2215	38-00
4x6 feet	38-0310	38-0300	38-2210	38-2200
4x7 feet	38-0316	38-0306	38-2216	38-2206
4x8 feet	38-0311	38-0301	38-2211	38-2201
➡ 5x7 feet	38-0312	38-0302	38-2212	38-2202
	38-0313	38-0303	38-2213	38-2203
\Rightarrow 6x8 feet	38-0317	38-0307	38-2217	38-2207
➡ 6x12 feet	38-0314	38-0304	38-2214	38-2204

#### non-folding mat with handles

WxL	1 3/8" PE foam	2" PU foam	1 3/8" EnviroSafe™	2" EnviroSafe™
4x4 feet	38-0215	38-0205	38-2315	38-2305
4x6 feet	38-0210	38-0200	38-2310	38-2300
4x7 feet	38-0216	38-0206	38-2316	38-2306
4x8 feet	38-0211	38-0201	38-2311	38-2301
⇒ 5x7 feet ⇒ 5x10 feet ⇒ 6x8 feet ⇒ 6x12 feet	38-0213 38-0217	38-0202 38-0203 38-0207 38-0204	38-2312 38-2313 38-2317 38-2314	38-2302 38-2303 38-2307 38-2304

**EnviroSafe™** foam mats are the eco-friendly alternative with the following features:

• 100% toxin free, zero VOC's or off gas, hypoallergenic, latex-free, mold and mildew proof, 100% recyclable with no leaching chemicals. Eco-friendly foam is medium-firm

#### custom projects



#### custom cut-outs

• call for your custom-sized mats and panels today!



#### custom graphics

• digital printing available for indoor and outdoor mats




#### indoor and outdoor wall padding

• call for more details!

Truck Shipments - These mats ship by truck. Unlike our other products that ship out immediately, mats take 3-5 weeks.



- lightweight and portable
- adult and child mat comes with a carry bag and has die-cut handles
- adult mat has 2 pieces that fit together in jigsaw pattern
- "air channel" construction allows mat to conform to body position



38-0550 72"x20"x1.2" (180x50x3cm) - each 38-0550-10 case of 10 adult mats

roll-up Cando<sup>®</sup> Cushy-air™ mat [c]

72"x20"x0.6" (180x50x1.5cm) - each 38-0552

38-0552-10 case of 10 roll-up mats












## rolls, wedges and mats

- positioning and vestibular shapes covered with heavy-duty vinyl-coated nylon, filled with high compression urethane foam
- inexpensive shapes are ideal for institutional use in physical therapy departments, home use for positioning rolls and wedges, special education programs, pre-school programs, physical education classes and the gym







- wedge shaped mats can be used to facilitate gross motor activities like rolling, tumbling and crawling
- safe, soft mats encourage participation
- incline mats with straps can be used for adaptive positioning

LxWxH	incline wedge without strap	positioning wedge with strap
2ft x 4ft x 14"	31-2023	
3ft x 2ft x 14"	31-2020	31-2032
6ft x 3ft x 16"	31-2021	31-2033
6ft x 4ft x 16"	31-2022	31-2034
- 7ft x 5ft x 18"	31-2024	31-2035

thick mats give added support and cushioning to increase the confidence level of the child using the mat

• side vents allow air to flow in and out

WxL	4-inch thick	8-inch thick
4x6ft (non-fold)	32-1000	32-1004
4x7ft (non-fold)	32-1008	32-1009
4x8ft (non-fold)	32-1001	32-1005
5x10ft (non-fold)	32-1002	32-1006
5x10ft (folding)	32-1011	32-1014
6x12ft (non-fold)	32-1003	32-1007
6x12ft (folding)	32-1012	32-1005
6x12ft (folding)	32-1012	

Truck Shipments - These mats ship by truck. Allow 3-5 weeks for delivery.



# exercise and rehabilitation positioning Skillbuilders® premium rolls and wedges have integral waterproof coating that is easily washed for years of use



#### Skillbuilders® wedges

- · wedges are ideal positioning shapes for exercise therapy
- · put under back, arm or leg of patient for comfortable and functional positioning
- provides necessary leverage to promote exercise regime; i.e. under the knee to enable weighted ankle lifts

wi	dth x length x h	neight	incline
30-1010	20x22x4"	small	11°
30-1011	20x22x6"	small	16°
30-1012	20x22x8"	small	20°
30-1013	20x22x10"	small	25°
30-1014	24x26x6"	medium	9°
30-1015	24x26x8"	medium	13°
30-1016	24x26x10"	medium	17°
30-1017	24x26x12"	medium	22°



#### quad boards

#### wood

- · provide support by forcing the proper amount of knee (or arm) flexion
- used on a treatment table or floor
- · board folds flat when not in use

#### steel

- · each section can be adjusted independently to match exact femoral and tibial lengths
- adjustability of this quad board allows for the knee to be fixed to any degree
- · chrome-plated steel, folds flat

10-1140	wood board
10-1141	wood board, padded
10-1142	steel board



#### Skillbuilders® rolls

- rolls, like wedges, are used to foster exercise therapy and general positioning
- roll can be placed under the knee to allow weighted ankle lifts


	length x d	iameter
30-1000	24x4"	small
30-1001	24x6"	small
30-1002	24x8"	small
30-1003	36x8"	medium
30-1004	36x10"	medium
30-1005	36x12"	medium
30-1006	48x12"	large
30-1007	48x14"	large
30-1008	48x16"	large

#### Skillbuilders® roll rack

30-1009 5-piece roll rack




Skillbuilders® knee supports



• supports knees and lower legs when user is in supine position

		WxLxH
30-1250	wedge	22x13x17"
30-1251	crescent	22x13x9"



#### Skillbuilders® half rolls

- half rolls combine stability of a wedge with the positioning versatility of a roll
- use with hook on bottom to anchor roll
- provides stable base for knee flexion and extension

	width x length	elevation
30-1230	4x18"	2"
30-1231	4x18"	3"
30-1232	6x24"	41/2"
30-1233	8x30"	6"



Pron-Pillo® face cradle



- · provides comfort when patient is lying in prone position
- · elevated base and nose cutout allow patient to breathe freely

00-4206 pron-pillo



# pulleys & range of motion (ROM) exercises

#### overdoor pulley exerciser, single



- pulley exercisers increase shoulder ROM
- · position clamp or strap in top door jamb

50-1023 pulley with strap 50-1023-25 pulley with strap (25ea)

50-1024-25

50-1024 pulley with clamp pulley with clamp (25ea)

#### overdoor pulley exerciser, double

- · heavy duty shoulder pulley is popular as a physical therapy exerciser to increase shoulder range-of-motion
- · deluxe exerciser has 2 pulleys and bracket for smooth and easy oper-
- portable exerciser can be easily transported and fitted to most doors

50-1010 deluxe pulley



#### color-coded overdoor pulleys




clamp

50-1035 pulley w/strap 50-1035-25 pulley w/strap (25ea)

50-1036 pulley w/clamp 50-1036-25 pulley w/clamp (25ea)

#### shoulder wheel



- range-of-motion arc can be varied from 10" to 39" by moving handle
- resistance varied by turning resistance knob
- height of unit can adjust 26" to accommodate seated and standing therapy

10-1150 shoulder wheel

deluxe skateboard

#### skate and powder board



· exercise skate may be used on floor, table or powder board

10-1130 6"x6" skate 10-1131 6"x12" skate

• facilitates various exercises. Portable board has a 30"x40" top and folds from its 7" height to 2"

10-1120 powder board

# finger/shoulder ladders





# standard arm skate



ergo arm skate

- figure 8 pattern for ROM
- vinyl cover, metal rim

10-1191 skate board

 resistance "pulley" system w/ 9 weights (0.4 lbs each) 10-1192 resistance sys.

10-1193 ergo arm skate

system (left)

wide hook/loop straps

• can be used with resistance

- · padded fastening straps
- removeable padded surface
- 6" x 12" surface

10-1194 std arm skate

- · offers progressive mobility of arm at shoulder, elbow and wrist joints
- wood unit has 36 steps at 11/4" intervals for shoulder, elbow, wrist, finger: 2"W x 54"L
- plastic steps at 11/4" intervals color-coded every 6 intervals: 6"W x 47"H

10-1160 wood (left) 10-1167 plastic (right)



stretching and balance products



# Cando<sup>®</sup> elastic stretch strap

the strap that stretches as you stretch

- dynamic stretch strap is made with elastic pockets that allow the user to ease into a stretch
- permits dynamic stretching
- designed for "facilitated stretching", an active form of stretching that uses isometric contractions to achieve greater flexibility gains than from static stretching
- reduces the need for a stretching partner

10-1385 strap 10-1385-25 pack of 25



- features six progressive handgrips to provide controlled stretch for improved flexibility and ROM
- increase stretch by moving to next handgrip position
- lightweight and portable

10-1380 Stretch-Rite strap



#### **CoreStretch™**

- stretches muscles, tendons and ligaments in back and legs
- lightweight and portable

10-1376 core stretch



# flexibility exercisers

- stability and support for stretching lower leg
- rocker stretches achilles tendon and plantar fascia
- can be used while standing or sitting
- motion enhancement system allows patients to safely add resistance while performing lower extremity stretching, strengthening and ROM exercises and includes light, medium and heavy cords (only)

10-1170 single leg (a) 10-1173 bilateral (b) 10-1174 single leg motion enhancement system (c)



# adjustable incline boards

- 14x18" with non-slip tread
- 5°-10°-15°-20°-25° elevation adjustment
- non-slip bottom
- 400 lb. capacity

10-1182 adjustable board



# progressive incline boards

- 161/4" wide x 60" long
- 4 units, each 16½" x15" with non slip treads
- units elevation measurement 15°, 20°, 25°, 30°

10-1183 set of 4



#### four-level incline board

- board increases extension in the lower and upper back while stretching the calves
- made of a composite plastic rated over 500 lbs.
- adjusts to 4 levels and folds flat (3-inches) to store easily
- levels:10°, 15°, 20° and 30°
- last slot on the board makes an excellent hand-hold for stretching
- non-slip top is 14"x14"

10-1179 4-level board



#### AnkleTough™ ankle exerciser

- for performing eversion, dorsiflexion and other ankle exercises
- each resistance is available in an eight-pack...or choose a set with one each of the four resistances

#### box of 8 pieces

10-1370 x-light white 10-1371 light yellow 10-1372 medium red 10-1373 heavy blue

#### box of 4 pieces

10-1374 1 ea. resistance





#### multi-reach shoulder arc

- range of motion arch featuring multiple reach capabilities
- · adjustable angle for greater flexibility and movement

10-1168 multi-reach



#### horizontal ring tree

- for people with neurological disorders, burns and orthopedic problems
- user can sit or stand while increasing ROM, strength, balance and perceptual skills. 3 ring sizes and 3 rods (12", 24" and 36")
- includes 36 plastic rings (12 ea. 11/8", 13/4" and 2")

10-1158 horizontal tree



#### single curved shoulder arc

- arc is 29" at its highest point
- includes 24 movable colored tabs for ROM exercises, crossing the midline, pinch strengthening and eye-hand coordination
- base measures 3" x 36" and also accommodates double curved arc
- base has two holes for optional tabletop attachment, or can be secured using clamps or vice grip

10-1154 single curve

# range-of-motion exercisers



#### double curved shoulder arc

- allows customizing program for people with limited ROM
- moving the colored tabs from one side of the arc to the other achieves full upper extremity ROM and improved motor planning and visual tracking skills
- use while standing for balance
- 241/2" at highest point, 71/2" at lowest point
- includes base and 24 tabs

10-1155 double curve



#### climbing board

- board and bar frame with 7 rubber-tipped extensions on each side that are 4" apart
- 38" bar has handle grips
- can be attached to table by clamps or through base
- · accessory kit combines activities to improve perceptual motor skills with activities they increase upper extremity ROM, strength and coordination

10-1152 board and bar 10-1153 accessory kit



#### incline board

- molded-plastic board can be permanently mounted or clamped on a table top
  - vary angle to provide different shoulder extension
  - 31" high, 81/2" wide, 24" base
  - bilateral or single handle
  - handle piece: 2<sup>3</sup>/<sub>4</sub>"x5<sup>1</sup>/<sub>2</sub>"x8<sup>1</sup>/<sub>2</sub>"
  - · weights can be added for increased resistance

10-1161 incline board



#### vertical ring tree

- increase ROM, eye-hand coordination, manipulation, endurance and sequencing
- base: 213/4"x141/2"x11/2" deep
- includes 6 plastic rods: 3 18" high and 3 13" rods
- diameters between 1/4" and 1"
- includes 12 rings for each rod

10-1156 vertical tree 10-1157 72 extra rings



#### semi-circular pegboard

- for people with hand injuries, neurological involvement or orthopedic injuries
- increase upper extremity ROM, strength, coordination and endurance, or improve cognitive and visual skills
- 3 tiers hold 43 pegs
- pegs come in 2 lengths and three diameters
- unit measures 26"x14"x2"

10-1159 peg board





# is better

# here's why:

- won't leak
- accurate
- repeatable
- made in USA





# here's how:

- one-piece body
- rugged double-wall bellows
- durable metal gauge



# choose your warranty









# DASELINE EVALUATION INSTRUMENTS

# Baseline® pneumatic squeeze dynamometer



- squeeze bulb dynamometer is the inexpensive way to measure hand and finger strength
- available with and without a maximum force indicator (reset) that remains at the maximum reading until reset
- · unit is calibrated in PSI

12-0291 with reset, 30 psi

12-0293 with reset, 15 psi

12-0290 30 psi 12-0292 15 psi

# grip strength dynamometers

# Baseline® pneumatic bulb hand dynamometer/pinch gauge



- measures both grip and pinch strength
- three different sized (1", 1.5" and 2") interchangeable squeeze bulbs are fitted to the gauge
- · lightweight and portable
- · maximum reading remains until reset
- unit is calibrated in PSI (0-30)
- similar to the Martin Vigorometer

12-0297 3-bulb dynamometer with reset, 30 psi

# Baseline® smedley spring dynamometer



- smedley spring-type dynamometer has been used successfully for many years to measure grip strength
- adjustable handle comfortably accommodates any hand size
- · calibrated in both pounds and kilograms
- · with protective carrying case

12-0281 adult (220 lb/100kg) 12-0282 child (110 lb/50kg)

#### Jamar® dynamometers



#### hydraulic

- hydraulic dynamometer gives accurate and repeatable grip strength readings
- adjustable 5-position handle
- maximum strength indicator remains until reset
- reads lbs. and kgs.

12-0600 200lb/90kg



#### digital Jamar+

- for routine screening of grip strength and initial and ongoing evaluation of clients with hand trauma and dysfunction
- sturdy aluminum body construction with scratch resistant UV coating
- easy-to-read LCD display can be set to display lbs. or kgs.

12-0604 200lb/90kg

#### Baseline® digital smedley dynamometer



- digital, lightweight smedley spring dynamometer offers accurate grip strength readout
- 5-preset grip sizes
- captures and displays maximum grip force
- assesses results according to age and gender
- saves/stores results for fast retrieval with up to 19 users

12-0286 200lb/90kg



## hydraulic hand dynamometers

- Baseline® hydraulic hand dynamometer gives accurate grip strength readings without the subject being able to "feel" the handle move
- the Baseline<sup>®</sup> hand dynamometer is engineered to be durable, so it will last for years and years
- · maximum reading remains until reset
- readings in lbs. and kgs. CE certified
- Baseline® dynamometers have been manufactured in the USA since 1978
- five position handle and body design yield results that are consistent with published Baseline® and Jamar® studies
- includes protective carrying case
- internationally accepted design ensures reliability, user convenience and measurement repeatability

# DASELINE® EVALUATION INSTRUMENTS

#### options

ER™ extended range 3.5" dial HiRes™ 3.5" dial

LCD 2.5" LCD gauge
Lite™ 1 year warranty
Standard 2 year warranty
HD™ 5 year warranty

#### 200 lb. capacity



12-0241 200 lb. LiTE™ with case




12-0240 200 lb. standard w/case



12-0221 200 lb. HD™ w/case



12-0243 200 lb. HiRes™ w/case






# DASELINE EVALUATION INSTRUMENTS

## pinch gauges

precision instruments for measuring finger strength







#### Baseline<sup>s</sup> mechanical pinch gauge



- measure tip, key and palmer pinch strength in both pounds and kilograms
- measurements are accurate and repeatable
- · results are consistent with published Markowitz studies
- indicator remains at the maximum reading until reset
- protective case

#### for orthopedic & sports medicine

12-0200 30 lb. with case, blue 12-0201 60 lb. with case, red

#### for weak and damaged hand

12-0203 10 lb. with case, silver 12-0202 2 lb. with case, gold

#### Psytech™ finger flexion/extension gauge





- lightweight plastic pinch gauge measures both finger flexion and extension
- · separate scales record each measurement in both pounds and kilograms
- maximum reading remains until reset
- flexion capacity is 20 pounds; extension capacity is 5 pounds

12-0230

flexion extension gauge

#### Jamar<sup>®</sup> hydraulic pinch gauge



- accurate and repeatable pinch strength readings (tip, key and palmer)
- maximum strength indicator remains after each reading until reset
- dual scale shows pounds and kilograms

12-0601

50lb./22.5kg.

#### Jamar® digital pinch gauge



- easy-to-read LCD display
- push-button console includes button to zero last reading stored and two toggle buttons used to accurately store up to 99 readings
- wrist strap and one 3V battery included
- 10" x 5" x 2" (25 x 13 x 5cm)

12-0607

50lb./22.5kg.



# hydraulic pinch gauges

- Baseline® hydraulic pinch gauge uses a hydraulic system to assure convenience, product reliability and measurement accuracy and repeatability
- therapist can support the pinch gauge during testing, yielding a more accurate result for all pinch tests
- use for tip, key and palmer tests
- Baseline® pinch gauges made in USA
- the results are consistent with published Baseline® and Jamar® studies
- · max reading remains until reset
- readings in lbs. and kgs.
- includes protective carrying case
- CE certified

#### options

ER™ extended range 3.5" dial HiRes™ 3.5" dial


LCD

2.5" LCD gauge Lite™ 1 year warranty

Standard 2 year warranty HD™ 5 year warranty

#### 50 lb. standard capacity





- 50 lb. capacity
- 21/2" diameter dial
- · attractive red finish



our most popular style

- used since the 1980's
- 50 lb. capacity
- 2¹/₂" diameter head
- mottled blue finish



**Standard** 

- extra heavy-duty parts
- 200 lb. capacity
- 21/2" diameter dial
- sleek black finish



 $HD^{\mathsf{m}}$ 



digital LCD

12-0226 LiTE™ 50 lb. with case

12-0235 50 lb. standard w/case 12-0222 HD™ 50 lb. with case

• 31/2" diameter head

- 500 lb. capacity
- mottled blue finish

50 lb. HiRes™ w/case



• 31/2" diameter head

100 lb. capacity

mottled blue finish

12-0228 100 lb. ER™ HiRes™



#### 100 lb. capacity



- 2<sup>1</sup>/<sub>2</sub>" LCD head
- electronic zero calibration system
- max clear button
- toggle maximum and current reading
- lbs. / kgs. toggle
- 100 lb. capacity
- · mottled blue finish
- 2 AAA batteries included

12-0237 100 lb. digital with case





12-0239



# VALUATION

### hand evaluation sets











standard 7-piece 12-0100

standard 3-piece 12-0103

12-0112 300/100 lb. digital



#### Baseline® 7-piece hand evaluation sets

- instruments stored in protective carrying case
- · grip strength is measured by the hydraulic hand dynamometer and pinch strength by the mechanical pinch gauge
- other measurements are made with the stainless steel goniometer (12-1010), 2-point discriminator with 3rd point (12-1481), Wartenburg pinwheel (12-1450), finger circumference gauge (12-1220) and functional finger motion gauge (12-1060)

#### choice of (dynamometer / pinch gauge)

12-0100	200/30 lb. standard	(12-0240	12-0200)
12-0102	200/30 lb. LiTE™	(12-0241	12-0200)
12-0117	200/30 lb. HD™	(12-0221	12-0200)

300/60 lb. ER™ HiRes™ (12-0246 12-0201) 12-0111 (12-0247 12-0201) 12-0113 300/60 lb. digital

#### Baseline® 3-piece hand evaluation sets

- instruments stored in protective carrying case
- includes choice of hydraulic hand dynamometer (below), or hydraulic pinch gauge (below) and 6" stainless steel goniometer (12-1010)

#### choice of (dynamometer / pinch gauge)

	200/50 lb. standard 200/50 lb. LiTE™	(12-0240 (12-0241	,
	200/50 lb. HD™	(12-0221	,
12-0110	300/100 lb. ER HiRes™	(12-0246	12-0228)



# Baseline® 8-piece evaluation sets

8-piece portable carrying case contains:

- 200 lb. hand dynamometer
- 50 lb. pinch gauge
- finger goniometer (12-1010)
- finger circumference gauge (12-1221)
- 180° 6" plastic goniometer (12-1005HR)
- measuring tape (12-1560)
- Discrim-a-Gon™ (12-1492)
- tactile™ monofilament set (12-1662)

#### choice of (dynamometer / pinch gauge)

12-0128 standard (12-0240 12-0235) 12-0126 LiTE™ (12-0241 12-0226) 12-0127 HD™ (12-0221 12-0222)



# Jamar<sup>®</sup> 3-piece <u>ha</u>nd evaluation sets

- portable carrying case contains instruments necessary for a hand evaluation
- 200 lb. hydraulic hand dynamometer, 50 lb. pinch gauge and 6" finger goniometer

#### 3-piece hand evaluation set

12-0602 Jamar® hydraulic (shown) 12-0605 Jamar® digital



(12-0247 12-0237)

### evaluation sets

- · hand dynamometer
- · mechanical pinch gauge
- small joint goniometer
- 180° 6" plastic goniometer
- flexible tape measure
- finger circumference gauge (cm)
- Touch-Test™ monofilaments (5 each)
- Discrim-a-Gon™

#### "Hands On" evaluation set

12-0603 Jamar® hydraulic (shown) 12-0606 Jamar® digital



## hydraulic wrist dynamometers

# DASELINE EVALUATION INSTRUMENTS

# wrist strength can now be measured

- simply and accurately measures the strength of the wrist muscles during flexion, extension, abduction and adduction and the forearm muscles during supination and pronation
- wrist/forearm dynamometer features the time proven hydraulic system used in the industry accepted Baseline® and Jamar® hand dynamometers
- accommodates the Baseline® dual-grip (12-0389) handle
- maximum reading remains until reset
- · measurement in pounds and kilograms
- standard with protective carrying case
- made in USA with a 1-year manufacturer's warranty
- CE certified

#### Baseline® wrist/forearm 500 lb. / 225 kg. dynamometer

includes door knob grip, table mount and protective carrying case

12-0250 analog (dial) readout 12-0251 digital (LCD) readout

#### Baseline® wrist evaluation sets



- protective carrying case contains instruments to evaluate the strength and range-of-motion of the wrist and forearm
- includes a wrist inclinometer (12-0502), wrist dynamometer (below), door knob grip (12-0261), and mounting bracket (12-0266)


analog gauge with shovel handle and mounting bracket


#### 4-piece wrist set

12-0114 12-0115 500 lb. analog 500 lb. digital














# MMT manual muscle tester grip strength dynamometer

# take manual muscle test readings with your grip strength dynamometer





### computer compatible dynamometers


computer interfaced dynamometers are now available in 2 styles:



• dual output: dial gauge with analog signal

single output: analog signal only

electronic output 300 lb. hand dynamometer



electronic output 500 lb. wrist dynamometer

electronic output 500 lb. push-pull dynamometer









Hand dynamometer contains a transducer that generates electric signal that can be read into a computer.

dial gauge with analog output 12-0021 grip [ a ]

analog output only 12-0248 grip [b] Pinch gauge contains a transducer that generates electric signal that can be read into a computer.

dial gauge with analog output 12-0023 pinch

analog output only 12-0238 pinch

Wrist dynamometer contains a transducer that generates electric signal that can be read into a computer.

dial gauge with analog output 12-0027 wrist

analog output only 12-0267 wrist Push-pull dynamometer contains a transducer that generates electric signal that can be read into a computer.

dial gauge with analog output 12-0025 push-pull

analog output only 12-0384 push-pull

#### Baseline® trade-in/upgrade program

Now is the time to trade-in or upgrade your old dynamometer, pinch gauge or force gauge (any brand, any condition) for a new Baseline® measurement instrument complete with a molded carrying case.

#### trade-in program

trade-in SKU	old product description	new product description
12-0240T	dynamometer	Baseline 200lb. hydraulic dynamometer
12-0235T	pinch gauge	Baseline 50lb. hydraulic pinch gauge
12-0200T	pinch gauge	Baseline 30lb. mechanical pinch gauge
12-0201T	pinch gauge	Baseline 60lb. mechanical pinch gauge
12-0394T	force gauge	Baseline 250lb. dial force gauge

000
CUS

cost

### upgrade program

trade-in SKU	old product description	new product description
12-0243T	dynamometer	Baseline 200lb. hydraulic dynamometer, large dial
12-0246T	dynamometer	Baseline 300lb. hydraulic dynamometer, large dial
12-0247T	dynamometer	Baseline 300lb. hydraulic digital dynamometer
12-0239T	pinch gauge	Baseline 50lb. hydraulic pinch gauge, large dial
12-0228T	pinch gauge	Baseline 95lb. hydraulic pinch gauge, large dial
12-0237T	pinch gauge	Baseline 100lb. hydraulic digital pinch gauge
12-0399T	force gauge	Baseline 250lb. digital force gauge

Call for a return authorization number. A packing slip will be faxed to you. Include it with your old measurement instrument. Send your old dynamometer or pinch gauge (prepaid) to Baseline Corporation at the address below.

Send to: Baseline trade-in/upgrade - Baseline Corporation, 3 Westchester Plaza, STE 111 Elmsford, NY 10523 (USA)



# DASELINE® ALUA ION INSTRUMENTS

A simple, easy-to-use, ergonomically designed instrument that **objectively** measures push, pull and lift forces for manual muscle testing, functional capacity evaluation and job task evaluation at a remarkably **affordable** price. Because the instrument is lightweight, small and **portable**, you can perform precise, objective evaluations in your office, at the client's location, or in the field. Ergonomically designed dynamometer is easy to grasp while testing small forces. The easy-to-attach single or dual grip handle can be used when measuring larger forces. Can be used with lift platform to perform lifting evaluation.

#### muscle strength measurement

This hand-held dynamometer lets you *objectively* measure manual muscle strength.

#### job task analysis

Measure actual push, pull and lift forces needed to perform a particular task (function).

#### functional capacity evaluation

Quantitatively evaluate an individual's push, pull or lift capacity to perform a given task (function).

### push-pull dynamometers


#### use for push or pull applications





muscle strength testing

job task evaluation



#### hydraulic



- comes with accessory pack: 3 push pads (padded curved, padded straight, and 1cm² circular), 1 pull hook and 1 snaplock hook
- cushioned carrying case
- maximum remains until rest
- optional accessories available (page 71)
- 1 year warranty; CE certified


#### analog (Dial) readout

12-0392 50 lb./22.5 kg. 12-0393 100 lb./45 kg. 12-0394 250 lb./115 kg. 12-0388 500 lb./225 kg.

#### digital (LCD) readout

12-0397 50 lb./22.5 kg. 12-0398 100 lb./45 kg. 12-0399 250 lb./115 kg. 12-0387 500 lb./225 kg.

#### electronic



- comes with accessory pack: 3 push pads, 1 pull hook, 1 snap hook, protective carrying case
- state-of-the-art "load-cell" technology to assure the highest standards of accuracy, repeatability and reliability
- available in four measurement ranges
- maximum remains until rest
- optional accessories available (page 71)
- 3-year warranty; CE certified

#### electronic dynamometers

12-0340 50 lb./22.5 kg. 12-0341 100 lb./45 kg. 12-0342 250 lb./112.5 kg. 12-0343 500 lb./225 kg.

#### mechanical spring





- comes with accessory pack: 3 push pads (padded curved, padded straight, and 1cm² circular), 1 pull hook, and 1 snap-lock hook
- maximum remains until rest
- 1 year warranty; CE certified

12-0300 2.2 lb./1 kg. 12-0301 5.5 lb./2.5 kg. 12-0302 11 lb./5 kg. 12-0304 22 lb./10 kg. 12-0303 66 lb./30 kg.



# push-pull dynamometer accessories



### pads and hooks



### accessory pack push pads

12-0371 padded curve [a] 12-0370 padded straight [b] 12-0372 1cm<sup>2</sup> circular [c] pull hooks

12-0377 medium hook [d] 12-0379 oval snap [e]

### optional pull hook

12-0376 small hook [f] 12-0446 threaded oval [g]

7 piece set (1 each) 12-0470 7 piece set

### lift platforms



lift platforms - includes base, 5' chain, 1 snap oval, 1 threaded oval. Adjust chain length for test continuity.

### regular base

12-0406 (15x15") [a]


### large base

12-0407 (24x24") **[b]** 

### accessories

12-0443 chain 12-0445 snap oval [c] 12-0446 threaded oval [d] /ft

### wall anchor



### WalSlide™ anchor slides and locks to any position along 6' Walslide™ 10-5094 6' adjustable

### handles



12-0385 single grip [a] 12-0389 dual grip [b]

### pull cuffs













- measurement accessory is securely fastened to subject near muscle group being evaluated
- cable tensiometer or pushpull dynamometer is then connected to measurement accessory
- made with padded foam.

12-0420 ankle [a] 12-0422 thigh [b]

12-0424 wrist [c]

12-0428 foot stirrup [d] 12-0427 shldr harness [e]

12-0429 head harness [f]

### back-leg-chest dynamometers



- · measure strength of back, leg and chest muscles
- base provides secure footing
- chain length is adjusted to accommodate for height differences or to vary the point of force application
- comes fully assembled with 5 foot chain
- oversize base 24" x 24"; regular base 18" x 18"
- pointer remains at maximum until reset
- reads lbs. and kgs.

12-0400 660 lb. (300 kg.) adult with standard base 12-0403

12-0401 12-0402

660 lb. (300 kg.) adult with oversize base 330 lb. (150 kg.) adolescent with standard base 165 lb. (75 kg.) child with standard base

# Baseline® cable tensiometer




- · secure one cable to a stationary object; the other to the measurement accessory
- with 4 measurement accessories (ankle, thigh, wrist, head)
- with 2 five foot chains, 4 snap hooks
- · pointer remains at maximum until reset
- · reads lbs. and kgs.

12-0410 adult (660 lb./300 kg.) 12-0411 adolescent (330 lb./150 kg.) 12-0412 child (165 lb./75 kg.)



# DASE NE ALUA ION INSTRUMENTS



### lifting boxes

- test functional capacity evaluation of an individual's material handling activities within work hardening programs
- · with removable weight holders, hinged lid with stop mechanism

LxWxH (in.)		empty weight	weigl	nt capacity (lbs
55-1010	14 x 14 x 17	21½	100	regular
55-1011	8 x 8 x 12	9	25	small
55-1012	10½ x 10½ x 14½	14	50	medium
55-1013	13 x 13 x 161/4	19	80	large
55-1014	24 x 12 x 71/4	19	80	double
55-1015	4-piece SET: 2 [55	5-1010], 1 [55-	1011]	
	and 1 [55, 1012]			



# model 55-1021 and 55-1022 together

### heavy duty sleds

- for testing, evaluating and conditioning individuals in work hardening and industrial rehab programs
- · simulate everyday work motions
- control the resistance and strengthening treatment for shoulders, backs, hips, thighs and calves by adding weights
- anti-friction runners for a variety of surfaces
- sled specs: 18" wide x 36" long; handle adjusts from 32" to 44" high; foam covered grips; 4 removable weight holders; capacity of 250 lbs.
- remove weight holders and insert box to use with sled
- accessory box specs: 15" wide x 30" long x 6" high; inside dimensions 14" x 22½" removable rubber mat

55-1020 set: weight sled and accessory box

55-1021 weight sled

55-1022 accessory box for weight sled

# work hardening







### work device and push-pull weighted sled

- combination lifting box and push cart for work simulation, conditioning and evaluation
- strap, foam covered hand grips and hand slots for lifting
- cart promotes push/pull activities
- box is 13½ x 13½ x 12"; cart is 17 x 14 x 4¾" with 2" casters;
   T- handle adjusts 32-44" high; removable weight holder with weight keeper; 10 lb. empty, 80 lb. capacity
- companion sled can be used alone or with work device
- straight handle adjusts 36-48" high; interchangeable box; sled measures 16 x 13<sup>5</sup>/<sub>6</sub> x 6"; anti-friction runners

55-1030	cart with T-handle and accessory box
55-1032	cart with straight handle and accessory box
55-1031	cart with T-handle
55-1033	cart with straight handle
55-1034	sled with T-handle and accessory box
	sled with straight handle and accessory box
55-1035	sled with T-handle
55-1037	sled with straight handle
55-1038	sled and cart, box, both handles

# use Baseline push-pull dynamometer to generate quantitative measurements

see pages 70-71 for more information and accessories

### hydraulic analog (Dial) readout 12-0392 50 lb./22.5 kg. 12-0393 100 lb./45 kg. 12-0394 250 lb./115 kg. 12-0388 500 lb./225 kg.

### hydraulic digital (LCD) readout

12-0397 50 lb./22.5 kg. 12-0398 100 lb./45 kg. 12-0399 250 lb./115 kg. 12-0387 500 lb./225 kg.

### electronic digital (LCD) readout

12-0340 50 lb./22.5kg. 12-0341 100 lb./45kg. 12-0342 250 lb./112.5kg. 12-0343 500 lb./225kg.





# work hardening



# table top electrical assembly center

- table top center simulates multiple job tasks that involve assembly activities
- promotes upper extremities ROM as well as working in confined spaces
- basic electrical assembly supplies included - do not wire unit with live current

55-1054 assembly center





### work station

- wall mounted unit to improve manual dexterity
- height adjustable feature allows to recreate a variety of working postures
- widely spread pattern of bolt holes; includes a variety of fasteners
- specs: 48" long x 10" wide

55-1052 work station



### hand assembly device

- table top device is composed of 102 parts designed to include activities using bolts, washers, and nuts that require either hand tightening and/or a wrench
- specs: 24" x 8" x 12"
- hardware included

55-1053 assembly device



### shoveling work station

- for repetitive tasks. Multiple shelf heights allow for a varied range of arm, shoulder, leg and back motions which simulate on-the-job duties
- 3 angled shelves
- shovel included
- specs: 72" x 36" x 60"

55-1000 shovel station



# adjustable height shelf assembly and sorting bin

### sorting bin

- encourages cognitive, perceptual skills
- used to develop ROM
- 4 adjustable panels on top; 6 fixed-sized bins on bottom; 29½" x 13" x 16"

### adjustable height shelf and frame

- height adjustable work surface used to evaluate and develop ROM, strength and to learn proper lifting techniques
- 6" x 36" base; 12" x 34½" top shelf; 17¼" x 29¾" adjustable shelf; shelf height adjusts from 6" - 48" high in 6" increments

55-1040 shelf assembly with sorting bin

55-1042 sorting bin alone

55-1041 adjustable shelf assembly alone



### deluxe electrical assembly

- stand alone center simulates multiple job tasks and assembly activities that range from simple to complex
- encourages whole-body ROM and working in confined spaces
- basic electrical supplies included; do not wire unit with live current
- specs: 22½" wide x 10" deep x 72" high;
   28" wide x 20" deep base

55-1055 electrical assembly center





### lifting/bolt center

- upper extremity ROM assembly box and work capacity box all in one
- for vocational rehab and WCE
- removable weight holder, weight keeper, includes variety of fasteners; empty weight of 11½ lbs. with capacity of 80 lbs.
- specs: 24" x 131/2" x 12"

55-1050 lifting/bolt center



# strength and ROM evaluation

### MicroFET 2™ manual muscle tester



- ergonomic shape makes the MF2 easy to use
- versatile attachments ensure proper muscle isolation for increased testing validity
- recognized by medico-legal community for documentation accuracy and validity
- digital display located on gauge, for ease of use
- comprehensive software sold separately\*

12-0381 MF2 12-0381S w/ software

# MicroFET 3™ manual muscle tester / inclinometer





- perform two evaluation functions using a single device when diagnosing and treating muscular disorders
- the microFET 3 MT & ROM provides fast, accurate and reliable digital results for MT and ROM testing
- microFET 3 MT & ROM has the capability of using the latest wireless radio frequency (RF) system
- the microFET 3 MT & ROM with RF technology eliminates the need for cables, wires or cords, making it completely portable
- the single inclinometer eliminates the need for a secondary unit, freeing up one hand to stabilize or assist the patient through ROM and eliminates the need for an assistant

12-0382 MF3 12-0382S w/ software

### MicroFET 4™ grip/pinch tester



- combines grip and pinch testing into one device
- portable, flexible and costeffective option for hand rehab/upper extremity testing

### performable tests:

- maximum grip test
- 5 position grip test
- rapid exchange grip test
- hand grip fatigue test
- key pinch
- palmar pinch
- tip pinch

12-0298 MF4 12-0298S w/ software

### DynEx™ exerciser



maximum strength: ten consecutive maximum strength test results may be stored

rapid exchange test: provides a means to compare right/left hand exchange evenly spaced over a period of exchanges

**endurance test**: patient tries to maintain a given grip strength over time

In addition to the grip strength tests the **DynEx2** provides isometric grip strength therapy

12-0450 DynEx1 12-0451 DynEx2

### MicroFET 5™ single inclinometer



- complete with pre-programmed protocols for cervical, thoracic and lumbar range of motion tests
- uses simple keypad prompts to automatically calculate and display spinal segment ROM results
- can be combined with optional software for data collection and report generation

12-0299 MF5

### MicroFET 6™ dual inclinometer



The MF6 ROM is ideal for therapists, doctors and chiropractors because of its convenient, objective findings

- automatic shut-off feature
- accuracy within 1 degree
- less than one pound
- ergonomic design
- can be used as a single or dual site inclinometer
- easy-to-read LCD screen
- usage chart included

12-0289 MF6 12-0289S w/ software

### J-Tech Dualer IQ™ inclinometer



- stores data for 19 tests with up to six repetitions each
- · enhanced data recording
- when measuring large extremity ranges of motion, Dualer IQ generates test data without the alignment issues often associated with goniometers
- Dualer IQ can simulate many goniometer ROM measurement techniques

12-1059 dualer IQ



# strength and ROM evaluation



### J-Tech® commander

Each testing device automatically collects data, stores results, computes statistics and prompts for the next test. When combined with optional Easy Docs software, data can be transferred to PC word processors and spread-



Built-in protocols for 5-position, rapid exchange and sustained grip tests. Calculates CV, maximum force and fatigue for each repetition, plus test average and bilateral average.

> 12-0259 GripTrack



### **b** Pinch Track™ dynamometer

Rated to 50 lbs. with 0.2 lbs. resolution for accurate testing. Calculates CV, maximum force and fatigue for each repetition, plus left and right averages and bilateral deficit.

12-0211 PinchTrack



### © Power Track™ (manual muscle tester)

Accurately and objectively tests muscle strength and identifies areas of weakness. Stores/analyzes 40 tests with up to 4 repetitions per side. Records peak force and calculates CV and strength deficits.

12-0383 PowerTrack III



### d algometer (pain threshold tester)

Measures pressure thresholds. Protocols for evaluating trigger point tenderness and fibromyalgia. Comes with 2 12-0380 Lafayette MMT [h] tips. Calculates CV, test averages, and percent deficit.

> 12-1485 algometer



### e force gauge

Toggle from push to pull. Stores 40 tests with 1-4 reps. Records peak and average force and calculates CV. Comes with handles, snap hook and curved pad.

12-0288 force gauge



### f Tracker Easy Doc software

Software takes the work out of reporting information, questionnaires and objective musculoskeletal testing. Creates tables and graphs, formats text, interprets questionnaire scores, and performs calculations like CV, norms, averages and validity. With Narrative Writer you can select components to create reports that include charts and graphs.

Easy Doc software

12-0181 impairment module

12-0182 hand/upper extremity module



### commander downloader software and adapter

Download information into a computer text file for printing or inclusion in documentation software. (not shown)

12-0185 downloader software

> \*please view our website for technical specs and additional product information

### Lafayette® manual muscle tester



Select a wide range of options, such as, data storage, test times (1-10 seconds) and choice of high/low threshold settings. The large, LCD screen displays all information.

- dual measurement range, 0-300 lbs. or 0-50 lbs. for increased sensitivity
- selectable test time from 1-10 seconds
- · storage of up to 52 test results
- · records peak force, time to reach peak, and total time
- displays in pounds and kilograms



data transfer module allows wireless communication between manual muscle tester and computer

12-0380S software w/wireless module [i]

### Acumar dual inclinometer



- compact, handheld unit that features a large digital display for easy reading
- store measurements with the hold button, and review maximum, minimum and average values
- reduces examination time and enhances the objectivity of measurement and documentation
- with built-in wireless transmitter to optional computer interface

12-1061 [h,i,j] complete set

12-1063 **[h]** dual inclinometer

12-1064 ruler ſij

wireless interface 12-1065 **[j]** 


12-1062 [g] single inclinometer





# transparent plastic goniometers


### Baseline® plastic 360° ISOM (STFR)



- permits observation of joint's axis of motion, and range of motion
- 360° head has three scales calibrated to be used with the ISOM
- one arm has a linear scale in inches and centimeters
- scale reads in 1° increments
- · goniometer arm remains at measurement until rest

	1 each	25 each
12 inch ISOM	12-1000	12-1000-25
8 inch ISOM	12-1001	12-1001-25
6 inch ISOM	12-1002	12-1002-25

### Baseline® plastic pocket goniometers



- inexpensive clear plastic goniometers have 5-inch linear scale calibrated in inches and centimeters
- · goniometer fits in shirt pocket

25 each 1 each 12-1005 6" pocket 180° 12-1005-25 6" rulongmeter 360° 12-1006 12-1006-25


### plastic flexion-hyper extension



- clear plastic finger goniometer measures a range of 30° hyperextension through 120° flexion
- for ROM measurements of MCPs, PIPs and DIPS

1 each 25 each flexion/hyper-extension 12-1012 12-1012-25

### Baseline® HiRes™ plastic 360° ISOM (STFR)



(shown), 8 inch, and 6 inch

- white background behind the gradations, numbers, and text assures high contrast for added readability and resolution
- scale reads in 1° increments

	1 each	25 each
12 inch ISOM HiRes™	12-1000HR	12-1000HR-25
8 inch ISOM HiRes™	12-1001HR	12-1001HR-25
6 inch ISOM HiRes™	12-1002HR	12-1002HR-25

### Baseline® HiRes™ plastic pocket goniometers



white background behind the gradations, numbers and text assures high contrast for added readability and resolution

25 each 1 each 6" pocket 180° HiRes™ 12-1005HR 12-1005HR-25 12-1006HR-25 6" rulongmeter 360° HiRes™ 12-1006HR

### Baseline® 6-piece plastic goniometer set



- includes 1 each 12", 8" and 6" 360° goniometer, 1 each pocket 180°, rulongmeter 360° and flexion/hyper-extension gauge
- available in standard or HiRes™ (hyper- extension gauge is standard)

1 each 25 each standard set 12-1028 12-1028-25 HiRes™ set 12-1028HR 12-1028HR-25



# Baseline® absolute+axis™ goniometers DACLINE





### when guessing isn't good enough<sup>™</sup>

Many goniometric measurements require one arm of the goniometer to be positioned either horizontally or vertically. The other arm is then lined-up with the patient's body part.


Until now, the practitioner had to "assume" or "eyeball" the vertical or horizontal position. This interpolation is often incorrect and introduces measurement error. The Absolute+Axis (AA) goniometer removes the "assume" and "eyeball" by establishing a "true" vertical or horizontal.

The Absolute+Axis goniometer removes measurement variability and allows the practitioner to concentrate on correct placement of the other goniometer arm and the patient's body position.

Some measurements the AA goniometer is designed to improve are: shoulder and hip rotation (internal, external); cervical flexion, extension, rotation and side bending; hamstring flexibility; and active ROM.

### Baseline® Absolute+Axis™ goniometers (built-in level)

12-1025HR 360° Baseline® Hi-Res™ 12" goniometer (with built in Absolute+Axis™)



• integral absolute vertical and horizontal levels are permanently incorporated into both the AA Baseline (2-color) 12 inch and Baseline® Hi-Res™ (3-color) 12 inch goniometers

12" Absolute+Axis™ goniometer 12" Hi-Res™ Absolute+Axis™ goni. 12-1025HR

each 12-1025


25 each 12-1025-25 12-1025HR-25








### Absolute+Axis™ attachment



### convert your 12" goniometer to an Absolute+Axis™ goniometer



### Absolute+Axis™ attachment for 12" plastic goniometers

- for those who own a 12" 360° goniometer
- slides "snugly" onto one arm of the goniometer
- · does not interfere with the goniometer's function

25 each

Absolute+Axis™ piece 12-1026

12-1026-25

### Baseline® digital Absolute+Axis™ goniometer

- integral absolute vertical and horizontal levels are incorporated into the digital goniometer arm
- reads 0-180 degrees on LCD screen (separation angle between goniometer arms)
- powered by one 9V battery (included)
- powder-coated steel with inch/cm marks printed onto arms
- · freeze angle measurement
- horizontal and vertical levels to indicate level or plumb starting point

each

25 each

digital Absolute+Axis™ goniometer

12-1027

12-1027-25








# finger goniometers


### Baseline® 180° digit goniometer



• readings can be taken by the therapist using only one hand; this frees up the therapist's other hand to stabilize the patient during measurement or to record the measurement; the standard finger goniometer requires two hands to make a finger ROM measurement

stainless steel digit goniometer, each 12-1014 plastic (black) digit goniometer, each stainless steel digit goniometer, 25 each 12-1014-25 plastic (black) digit goniometer, 25 each 12-1016-25

### Baseline® finger (small joint) goniometer



 head has opposing 180° scales in 5° increments

• linear scale in inches and centimeters

• deluxe model has larger head

12-1010 standard finger, 6 inch 12-1011 deluxe finger, 6 inch 12-1015

12-1010-25 standard finger, 6 inch 12-1011-25 deluxe finger, 6 inch small finger, 31/2 inch small finger, 31/2 inch 12-1015-25

### finger motion gauge

- · clear plastic gauge simplifies and speeds measurement of composite finger flexion and finger-thumb opposition
- · scales measure up to 10 centimeters

12-1060 finger motion gauge



25 each

# stainless steel goniometers

DASELINE ALUA ION INSTRUMENTS


180° 8" goniometer

180° 14" goniometer

360° 14" goniometer

180° 6" pocket goniom.





 tension between the arms is controlled by the thumb knob  two opposing scales marked in 1° increments

extra-long arms for measuring all joints

 two opposing 180° scales and one 360° scale, all marked in 1° increments

 arm tension is controlled by thumb knob  Robinson pocket-sized goniometer has a 180° scale in 5° increments and a 5-inch (12½ cm) linear scale

12-1042 SS 6", 180° 12-1042-25 25 each

12-1040 SS 8", 180°

12-1041 SS 14", 180°

12-1050 SS 14", 360°

### 180° X-ray goniometer



- · document angles on X-ray negatives
- 180° head has two opposing scales marked in 2½° increments
- black aluminum with easy-to-read white radiopaque markings

12-1030 14", 180° 12-1031 8", 180° 12-1032 6", 180°


### large joint goniometer



- 18-inch stainless steel goniometer measures "large" joints (shoulders, hips, knees)
- two opposing 180° scales and a linear scale with inches and centimeters
- marked in 5° increments

12-1020 SS large joint


### 180° extendable goniometer



- easy-to-read, accurate instrument offers versatility for measuring large and small joints
- 4x magnifying lens to read the dual 180° high resolution scale and two stainless steel arms extending 9½" to 26"
- · smooth precision glide

12-1035 extendable 180°

### ISOM measurement chart



- 36"x30" wall chart illustrates measurement and recording rules from the American Academy of Orthopedic Surgeons
- goniometry, X-ray evaluation & anthropometric measurement

12-1095 ISOM chart



### 6-piece goniometer set

- 6 frequently requested Baseline<sup>®</sup> stainless steel and X-ray goniometers
- goniometers in a padded carrying case
  - 6" 180° Robinson pocket (12-1042)
  - 14" 180° stainless steel (12-1041)
  - 14" 360° stainless steel (12-1050)6" stainless steel finger (12-1010)
  - 8" 180° x-ray (12-1031)
  - 8" 180° stainless steel (12-1040)



12-1043 gonjometer set

### Baseline® wrist inclinometer



- measures range of motion of the wrist along all axes
- have patient grasp inclinometer handle; turn dial until scale reads 0; set first tab, take joint through range; set second tab, read range traveled directly from dial

12-0502 wrist inclinometer



### Baseline® bubble inclinometer



- standards codified in AMA Guide to the Evaluation of Permanent Impairment, third edition
- place inclinometer near joint to be measured; turn dial until scale reads 0; take joint through its range; read range traveled directly from dial
- some measurements require the simultaneous use of 2 inclinometers

12-1056

bubble inclinometer

12-1056-2

bubble inclinometer set (2 units)

12-1056-25

bubble inclinometer set (25 units)

### Baseline® digital inclinometer





- read ROM directly after the joint has been taken through its range
- place inclinometer near the joint to be measured; press "zero" button; move joint through its range; press the "hold" button

12-1057 digital inclinometer 12-1057-2 set (2 units)

### Baseline® gravity inclinometer



- place inclinometer near the joint to be measured; turn the dial to zero; take the joint through its range; read the total range directly from the dial
- · some measurements require the simultaneous use of 2 inclinometers

gravity inclinometer 12-1055 12-1055-2 set (2 units)

### inclinometers

### Baseline® AcuAngle inclinometer



- between feet
- pointer is dampened by fluid to assure accurate ROM measurements
- place inclinometer near joint to be measured; turn dial to zero; take joint through its range; read range-of-motion from dial
- adjustable feet along scale (side-to-side) adapt to body contours

12-1149 Baseline® AcuAngle

12-1149-2 set of 2

### Baseline® universal inclinometer



- measures ROM
- set moveable dial so that the weighted indicator points to zero, move joint through its range and read dial

12-1070 inclinometer w/ clip 12-1071 inclinometer w/ headband

### spinal assessment devices















- CROM Deluxe Cervical Range-of-Motion in-
- CROM Basic CROM main frame, rotation arm and video
- CROM 3
- **Deluxe Universal Inclinometer**
- BROM II Back Range-of-Motion instrument
- BROM Basic BROM rotation lateral unit plus two universal inclinometers
- PALM PALpation Meter (skeletal alignment and leg length discrepancy instrument)
- PAA Deluxe inclinometer

12-1151

a BROM basic

12-1155 12-1156

**b** CROM basic c CROM deluxe

12-1180

d PALM

12-1150 12-1097

12-1066 12-1157

**BROM II** 

PAA universal inclinometer

Deluxe universal inclinometer


CROM 3



# chiropractic measurement

### Baseline® large joint (arthrodial) protractor





- measure range-of-motion for all major articulations, cervical rotation, lateral flexion of the head and anterior-posterior cervical flexion
- bubble level assures measurement is made on horizontal plane
- two 180° opposing scales in 5° increments

12-1076 arthrodial protractor

### Baseline® posture evaluator



- provide posture evaluations that are more precise than those generated by the posture grid
- use the angle pointer to find and measure alignment

### Baseline® posture grid



- evaluate posture
- the plastic posture grid (25"x48") can be suspended from the ceiling or wall bracket
- set includes a suspension cord and a plumb bob to measure alignment

# BASELINE EVALUATION INSTRUMENTS

### Baseline® wall goniometer



- measure thoracic and lumbar ROM, lateral and forward flexion and extension of the spine
- height adjustable goniometer can be adjusted to align with the client's L4 and L5 axis points for measurements

12-1096 wall-mount goniometer

### Baseline® posture sets

- posture grid and posture evaluator can be used together or separately to measure and record posture
- three piece set including the posture evaluator, posture grid and arthrodial protractor is all you may need to perform complete posture evaluations and range-of-motion measurements

12-1077 grid and evaluator 12-1082 protractor, evaluator and grid

12-1078 posture evaluator

### 12-1080 posture grid

### Baseline® body level



- determine whether body parts are properly aligned
- place moveable arms on two sites
- read amount of tilt (in degrees) from the scale
- · Ideal for scoliosis screening

12-1090 body level



### Baseline® scoliosis meter



- measurements can detect scoliosis and abnormal anteropostero curve
- also measures the unstable lumbosacral, cervical and thoracic curves
- scales in centimeters, inches and degrees

scoliosis meter 12-1091

### **Scoliometer**



- measure the degree of rotation of a deformity of the back
- use to screen which persons to refer for further medical evaluation
- 3.5" x 7.5"
- with storage pouch

12-1099 scoliometer



### Baseline® Tactile™ Semmes-Weinstein™ type monofilaments



- · measure cutaneous sensory perception threshold
- each monofilament represents a unique amount of force
- fitted protective case enables safe storage and transport

sets normal 12-1660 2 piece: 2.83,3.61 protective 12-1668 2 piece: 4.56,5.07 12-1662 hand, set of 5 2.83,3.61,4.31,4.56,6.65 12-1664 foot, set of 6 2.83,3.61,4.31,4.56,5.07,6.65

set of 20 complete

12-1666

individual monofilaments		
	size	grams of force
12-1640	1.65	0.008
12-1641	2.36	0.02
12-1642	2.44	0.04
12-1643	2.83	0.07
12-1644	3.22	0.16
12-1645	3.61	0.4
12-1646	3.84	0.6
12-1647	4.08	1.0
12-1648	4.17	1.4
12-1649	4.31	2.0
12-1650	4.56	4.0
12-1651	4.74	6.0
12-1652	4.93	8.0
12-1653	5.07	10.0
12-1654	5.18	15.0
12-1655	5.46	26.0
12-1656	5.88	60.0
12-1657	6.10	100.0
12-1658	6.45	180.0
12-1659	6.65	300.0
cases		
12-1637	for 1 or	2 evaluators
12-1638	for 5 or	6 evaluators

for 20 evaluators

# monofilaments / aesthesiometers

### WEST™ sensory evaluators







Weinstein Enhanced Sensory Test (WEST) uses five monofilaments to measure skin sensitivity and peripheral nerve involvement. Monofilament tips are round and textured for reliable testing. Two or more filaments can be used for testing. Includes case (4½" by 6½"), test form and manual.

WEST-D provides an evaluation range from reduced tactile sensation on the face up to loss of protective sensation of the foot. Ideal instrument for testing tactile sensation at various locations.

WEST hand, foot and "D" monofilaments are individually calibrated.

12-1406 hand monofilaments 0.07, 0.2, 2, 4, 200

12-1407 foot monofilaments 0.5. 2. 10. 50. 200

12-1408 WEST-D monofilaments 0.02, 0.07, 0.5, 2, and 10

### Touch-Test™ sensory evaluators



- filament bends when peak force threshold is reached
- test sensory levels and obtain objective data on the status of diminishing or returning sensibility
- handle opens to a 90° angle
- finger indentations on handle provide a firm, comfortable hold
- color code shows pressure level
- each set comes in protective case with instructions included

### sets

12-1560 5 piece, hand 2.83,3.61,4.31,4.56,6.65

12-1562 6 piece, foot

2.83,3.61,4.31,4.56,5.07,6.65

12-1564 20 piece complete

accessories 12-1565 screening forms, hand 12-1566 screening forms, foot

\*accessories include screening forms and colored pencils.

### Baseline® folding monofilaments



- inexpensive monofilaments detect changes in neurological status
- · for screening peripheral nerve impairment
- used in diabetic foot screening (10 gm)

12-1395 1 gm 12-1396 10 gm 12-1397 75 gm 12-1398 3 piece set

### Baseline® retracting monofilaments



- damage resistant, filaments detect changes in neurological status
- for screening peripheral nerve impairment
- used in diabetic foot screening (10 gm)

12-1390 1 gm 12-1391 10 gm 12-1392 75 gm 12-1393 3 piece set

### Baseline® disposable monofilaments



- disposable, individual 10gm monofilaments
- choose from ADA LOPS or LEAP program

**ADA LOPS** 

12-1670 1 ea. 12-1670-20 20 ea. 12-1670-40 40 ea.

**LEAP** 

12-1671 1 ea. 12-1671-20 20 ea 12-1671-40 40 ea.



12-1639

# sensory evaluation

### Touch-Test™ 2-pt discriminator



- tests static, moving 2-pt discrimination
- consists of two sturdy, rotating, plastic disks that are joined together
- for testing following nerve repair, grafts and innervated tissue transfer for desensitization, or to determine the level of impairment
- includes instructions and carrying case
   12-1495 Touch Test® 2-point

### Dellon® 2-pt Disk-Criminator™



- for testing static and dynamic 1 and 2-pt discrimination
- includes 2 disks that quantify innervation density from 1 to 25 mm
- therapeutic aid for desensitization and home programs following nerve repair, graft and innervated tissue damage

12-1490 Dellon® 2-point

### Baseline® dolorimeter




- · evaluate pain threshold
- standardized (1.52cm²) flat circular probe is pushed against subject's skin until pain threshold is reached
- reads lbs. and kgs.

	_
12-1440	2.2 lb./ 1 kg.
12-1441	5.5 lb./ 2.5 kg.
12-1444	11 lb./ 5 kg.
12-1442	22 lb./ 10 kg.
12-1443	66 lb / 30 kg



### Baseline<sup>®</sup> Discrim-A-Gon™ 2-point discriminator

includes set of 2 octagons, each with 8 labeled 2-point distances




- 2 separate 2-point discrimination octagons (D1 and D2)
- each octagon measures a different range of 8 labeled, fixed 2-point intervals ranging from 1 mm to 25 mm
- easy-to-use, lightweight plastic wheel is the perfect sensory evaluation tool to test static and dynamic 1 and 2-point discrimination

12-1492 Discrim-A-Gon™ 12-1492-25 25 each

# DAJELINE EVALUATION INSTRUMENTS

### Baseline® 2-point discriminators



[b] plastic with third point



- evaluate cutaneous sensitivity and touch threshold
- third point permits alternating between single and double point stimulation without changing setting
- plastic points minimize influence of temperature
- measures up to 14 centimeters plastic

12-1480 [a] 2-point discriminator 12-1481 [b] with third point

stainless steel

12-1482 [c] 2-point discriminator

### pinwheels

### 12-1450 Wartenberg pinwheel



12-1451 Cleanwheel sterile neurological pinwheel



12-1455 Neuro-Aide neurological pinwheel



- consists of a rotating pinwheel attached to a handle
- used by orthopedists, neurologists, neurosurgeons, and other physicians to test a patient's ability to feel sharp contact with the skin
- test for nerve damage, or nerve dysfunction

12-1450 Wartenberg 12-1450-25 Wartenberg, 25 each 12-1451 Cleanwheel 12-1451-25 Cleanwheel 25 each

12-1455 Neuro-Aide 12-1455-25 Neuro-Aide 25 each



# sensory hammers and tuning forks

### Baseline® tuning forks



• use to evaluate hearing and vibratory sensation

### tuning fork sets

set of six with case 12-1460

(128,256,512,1024,2048,4096 cps)

sensory set (30 and 256 cps) 12-1461

Scrisory Set (60 and 200 ops)	, 12 1-	101
weighted tuning forks	each	25 each
30 cps weighted	12-1465	12-1465-25
64 cps weighted	12-1464	12-1464-25
128 cps weighted	12-1466	12-1466-25
256 cps weighted	12-1467	12-1467-25
unweighted tuning fork	S	
512 cps	12-1468	12-1468-25
1024 cps	12-1469	12-1469-25
2048 cps	12-1470	12-1470-25
4096 cps	12-1471	12-1471-25
student grade tuning fo	rks	
128 cps weighted	12-1475	12-1475-25
256 cps	12-1476	12-1476-25
512 cps	12-1477	12-1477-25
special tuning forks		
[a] Reider Seifer tuning fork	(64 and 128 cps)	12-1472
20 to 4096 cps, variable	frequency tuning	fork 12-1462
[b] 64 and 128 cps adjustab	le tuning fork	12-1517
The state of the s		



### Baseline® Taylor hammer





Use to elicit responses during examination. Can be combined with tuning fork.

	each	25 each
Taylor, red hammer	12-1500	12-1500-25
with 128 cps tuning fork	12-1501	12-1501-25
with 256 cps tuning fork	12-1502	12-1502-25
Taylor, yellow latex-free	12-1570	12-1570-25
Taylor, red latex-free	12-1571	12-1571-25
Taylor, green latex-free	12-1572	12-1572-25
. ,	12-1573	12-1573-25
Taylor, black latex-free	12-1574	12-1574-25

### Baseline® Buck neurological hammer





Two heads, point/brush [12-1510] for cutaneous sensation evaluation. Buck hammer [12-1516] has a pinwheel. Buck with tuning fork also has pinwheel.

	each	25 each
Buck hammer	12-1510	12-1510-25
w/ pinwheel, brush, needle	12-1516	12-1516-25
with 128 cps tuning fork	12-1511	12-1511-25
with 256 cps tuning fork	12-1512	12-1512-25
Buck, latex-free hammer	12-1514	12-1514-25

### Baseline® Babinski hammer



Adjustable head can be oriented in two directions. Use pointed bottom as a "gross" aesthesiometer.

	each	25 each
Babinski hammer	12-1520	12-1520-25
Babinski extra long body	12-1521	12-1521-25
Babinski, latex-free	12-1524	12-1524-25
extra long, latex-free	12-1526	12-1526-25



# blood pressure, oxygen level, pulse

### cuffs for sphygs. portable sphyg.

### wall mount sphyg.

### mobile floor sphyg.



• 2-tube cuff bladder system with bulb and valve for mercurial and aneroid units

12-2243 adult 12-2244 large adult 12-2242 child

12-2241 infant 12-2240 newborn



· high quality aneroid sphygmomanometer with adult inflation system, carrying bag

12-2250 adult w/case 12-2250-25 25 each



- 90° swivel mount and cuff
- · aneroid sphygmomanometer with adult inflation system

aneroid 12-2261 wall mount



- adjustable height caster base, double cuff container
- aneroid sphygmomanometer with adult inflation system

aneroid

12-2262 mobile floor

### electronic BP and pulse meters



semi-auto arm



auto arm



deluxe auto arm

all blood pressure monitors:

display reading of systolic and diastolic

• are latex free; batteries not included

standard auto arm cuff (4 AA batteries)

deluxe auto arm cuff (4 AA batteries)

· deluxe arm model has wide range cuff

semi-auto arm cuff (4 AA batteries)

· digital display; manual inflation

12-2273 manual inflation arm

12-2272 auto inflation arm

for arms 9" to 17'

comes with AC adapter jack

blood pressure and pulse; date and time • store 30 measurements for 2 users



auto wrist

### pulse oximeters

 pulse oximeters measure pulse (HR) and SpO<sub>2</sub> levels













### Nonin® oximeters

- fast, accurate SpO2 and pulse readings
- · adjusts from pediatric to adult finger sizes GO2 can be used outside the clinic, for
- high altitude and sports enthusiasts

### fingertip model

12-1924 G02 model [a] 12-1921 Onyx 9500, original 12-1920 Onyx II 9550, Next Gen [d]

### fingertip with wristwatch

- WristOx model has storage capacity and download capability
- wristwatch displays data in real-time 12-1910 WristOx 3100 [c]

### fingertip with handheld monitor

• handheld monitors have data options and alarms

12-1911 PalmSat 2500 [b] 12-1912 8500 series [e]

### pulse oximeter / BP combo

• portable unit monitors SpO<sub>2</sub>, HR, BP 12-1960 Avant 2120 [f]

### fingertip pulse oximeters



- measures pulse
- measures SpO<sub>2</sub>
- · compact and lightweight

### Choice® oximeter

12-1926 Choice® model oximeters

12-1922 pulse oximeter 12-1923 with clip

### stethoscopes







nurses

dual head

Sprague-Rappaport

nurses anodized aluminum chest piece, plastic ear tips, 22" tube


12-2200 nurses 12-2200-25 25 each

dual head turn to convert to non-child ring; external spring binaural, 22" tube

12-2210 dual-head

12-2210-25 dual-head (25 each) 12-2211 Littman type

Sprague-Rappaport rotating drum stem and external spring binaural: with 3 bell chest pieces and 2 diaphragms, 22" tube 12-2220 Sprague-Rappaport



detects irregular heartbeat

12-2274 deluxe wide auto arm

 irregular heart rate detection feature 12-2151 auto inflation wrist





### Baseline® skinfold caliper



- lightweight aluminum caliper
- "floating tips" to assure accurate and reproducible results
- 70mm scale is large and easy to read
- complete with booklet and tables
- protective carrying case

12-1110 skinfold caliper

### hand held BF analyzer



- measures body fat and water %
- 4%-45% body fat, 37.8%-66% body water with 0.1% resolution
- 8-person memory
- adult and child user modes
- lithium batteries included

12-1122 hand-held size, each

### BodyLogic® analyzer



- body fat measurement in just seconds
- no need to pinch skinfolds, just grip the lightweight analyzer with both hands
- stores 9 personal profiles
- five-button keypad

12-1121 BodyLogic 306

### Baseline® economy SF caliper



- plastic skinfold caliper is the inexpensive alternative for skinfold analysis
- spring-loaded arm assures reproducible results
- 50mm scale
- for screening large populations
- complete with booklet and tables

12-1112 economy skinfold caliper 12-1112-25 caliper (25 each)

### calibration block



- check calibration of any skinfold caliper
- 10mm-50mm

12-1116 calibration block

### digital fat tracker™



- caliper stores 3 measurements, calculates and displays body-fat percentage
- uses 2 AA batteries

12-1124 fat tracker 12-1129 deluxe fat tracker

### mini BF analyzer



- measures/monitors body fat %
- 5-50% body fat, 0.1% resolution
- 8-person memory
- 3V battery, included!
- features clock and alarm

12-1140

palm size, each

# skinfold calipers

### Baseline® slim guide SF



- for screening large populations
- 60mm scale
- spring-loaded arm
- with instructions and tables

12-1125 skinfold caliper 12-1125-25 caliper (25 ea.)

### Skyndex® skinfold caliper



 body fat % calculated on display; Jackson-Pollock is standard, Slaughter Lohmann for children, Durnin for research

### single formula

12-1100 Jackson-Pollock (J-P)

12-1101 Durnin

12-1103 Slaughter/Lohmann

### dual formula

12-1107 Durnin + J-P

12-1108 Durnin +S/Lohmann

12-1109 J-P+S/Lohmann

### credit-card size BF analyzer



- measures/monitors body fat %
- 5-50% body fat, 0.1% resolution
- 5-person memory
- 3V battery, included!
- features clock and alarm

12-1141

credit card size, each




# volume, circumference and breadth



### Baseline® volumetric edema gauges







### each set includes

- measurement cavity
- overflow beaker
- graduated cylinder
- designed to objectively measure volume of the body part by using the fluid displacement method
- body part is placed in a water-filled tank; the resultant overflow is collected and measured
- amount of edema can be measured and quantified

### upper extremity (LxWxH)

12-3500 hand set 5x3x9" 12-3501 arm set 7x7x30" 12-3504 forearm set 6x6x24"

### lower extremity (LxWxH)

12-3502 foot set 13x5x9" 12-3505 partial leg set 13x6x24" 12-3508 large partial leg set 13x8x16"

### Baseline® measuring tape



- reads in inches and centimeters and unfurls to 5 feet
- locks in place until reset button is pressed

12-1210 tape 12-1210-25 25 each

### Gulick measurement tape



- spring exerts constant force on woven fabric for accurate and repeatable results
- 5-foot tape in inches and cm
- push-button retractor

12-1201 Gulick 12-1201-25 25 each

### Martin breadth caliper



- for accurate breadth readings for chest, shoulder, pelvis, etc.
- 20-inch scale is marked in both inches and centimeters

12-1230 Martin caliper

### Collyer breadth caliper



- breadth caliper can be used to measure chest, shoulder, pelvis, etc.
- 15-inch scale is marked in both inches and centimeters

12-1231 Collyer caliper

### Baseline® circumference



- measure the circumference of a joint or digit
- measure by yourself: simply make loop, position loosely around waist, wrist or area being measured, and press button
- measuring tape will stay in place securely, while you have two free hands to make adjustments and record measurements



12-1205 circumference tape 12-1205-25 circumference tape, 25 each

# finger circumference gauge



- measure the circumference of a joint or digit
- wrap webbing around the finger and read circumference

12-1220 gauge, inches 12-1221 gauge, cm



# **Detecto® scales**

### (visit our website for more information)

### LCD floor scale



Space-saving physician's scale with digital display provides BMI (Body Mass Index) by conveniently selecting your height with easy-to-use scroll keys.

- 480 lb/220 kg capacity
- 0.2 lb/0.1 kg resolution

12-1348 floor scale



Scale with waist-high digital display provides BMI. Available with or without height rod.

- 500 lb/225 kg capacity
- 0.2 lb/0.1 kg resolution

12-1340 WH scale 12-1341 w/ height rod

### LCD waist-high scale LCD comfort-level scale



Scale with 48" high comfort-level digital display and built-in wheels. Provides BMI. Available with or

- 480 lb/220 kg capacity
- 0.2 lb/0.1 kg resolution

12-1342 comfort-level scale 12-1343 w/ mechanical HR 12-1344 w/ digital HR

# without height rod (HR).

500 lb/225 kg capacity

LCD eye-level

Scale with 59" high eye-level

digital display provides BMI.

Available with or without height

• 0.2 lb/0.1 kg resolution

rod (HR).

12-1345 eye-level scale 12-1346 w/ mechanical HR 12-1347 w/ digital HR

### remote indicator scale





Remote indicator features a large 1"/25mm display, units conversion, and tare. Lightweight, portable, perfect for mobile clinics, home care nurses.

- 400 lb/181 kg capacity
- .5 lb/.2 kg resolution

12-1349 floor scale

Dual reading beam. Available with mechanical height rod and wheels.

- 400 lb/175 kg capacity
- 4 oz / 100 g resolution

12-1350 eye-level scale 12-1351 w/ HR

12-1352 w/ HR and wheels

### stadiometer



Lightweight and portable wall mount height rod. Kit included for mounting on wall. Available with analog or digital readout.

12-1094 mechanical 12-1353 digital

### large dial display



Professional large dial scale with 61/2" dial in pounds (black) and kilograms (red), with nonslip surface, 330 lb/ 150 kg capacity with 1 lb/ 0.5 kg increments. Platform is 17x11".

12-1320 142KL scale

### Omron® body monitor scales

### HBF-500 body composition



- total body monitoring scale not only calculates visceral fat, but resting metabolism and skeletal muscle mass levels
- · displays body fat %, weight and BMI 12-1193 HBF-500

### HBF-400 fat loss monitor



· helps monitor overall weight while determining body fat %, a true indicator of healthy weight loss-the loss of fat, not muscle and body mass index (BMI)

12-1192 HBF-400



- portable large rubber deck safety plate platform with integral ramp
- · weigh patients in manual, power wheelchairs, dialysis patients, or walk-ons
- 1000 lb x .5 lb / 400 kg x .2 kg

12-1355 wheelchair/bariatric scale



### Baseline® body fat scales





- measures weight, body fat and water %
- max capacity 330 lbs. / 150 kgs
- 4%-45% body fat; 37.8%-66% body water
- division 0.2 lbs., 0.1% BF/water
- 8-person memory
- · adult and child modes

### tempered glass scale

12-1190 tempered glass\* \*requires 1 9V alkaline battery

### stainless steel scale

12-1191 steel\*
\*required 4 AA batteries

disposable additional

mouthpieces

### digital refractometer



- use to screen for drug compounds
- place a drop of fluid in the titanium well and press a button
- microprocessor delivers quick readout on one of five scales
- hundreds of different units of measure and scale combos
- · automatically calibrates to water

### Palm Abbe PA-202X

12-1902 Palm Abbe PA-202X

### Baseline® spirometer

- compact, lightweight spirometer for screening large populations
- subject blows into disposable mouthpiece
- indicator registers forced vital lung capacity up to 7 liters
- stays at maximum until reset
- comes with 50 plastic mouthpieces

### 12-1710 spirometer with 25 plastic mouthpieces

 plastic mouthpieces

 12-1711
 100 additional

 disposable mouthpieces

 12-1714
 250 additional

 12-1712
 500 additional

 12-1713
 1000 additional

### Baseline® sit-and-reach trunk flexibility test


- perform sit-and-reach test
- scales include: fitness, program, Cooper Institute/YMCA, AAHPERD
- scales printed in inches and centimeters
- sturdy construction with a built-in foot plate and a maximum stretch indicator that remains until reset
- powder-coated metal throughout

12-1085 standard

12-1086 deluxe



 deluxe model used for modified sit-and-reach test enabling the testing protocols to be used that vary the scale start point; this normalizes scores to account for variation in the arm and leg length of the person being tested





# heartrate watches, pedometers and timers

### Ekho® wrist watch heart rate monitor

















FiT-28 12-2044

standard watch functions - time of day, date, stopwatch, water resistant, interface with heart rate interactive equipment, light, user changeable batteries in transmitter and watch

standard heart rate functions - current heart rate, ECG accurate, transmitter and strap; additional functions below

E-10

standard basic HR and watch functions

12-2045 E-10

FiT-8, FiT-9

Target zone, out of zone alarm FiT 9 has thinner watch face.

12-2040 FiT-8 12-2041

FiT-28

Target zone, out of zone alarm, average HR, max HR, % of max HR, time in/above/below zone, calories burned, fat burned, countdown timer, 5-minute warm-up and recovery, key lock.

FiT-28 12-2044

E-15 -

Target zone, out of zone alarm, average HR, max HR, calories burned.

12-2046 E-15

FiT-18. FiT-19

Target zone, out of zone alarm, average HR, max HR, time in/above/below zone, calories burned. FiT 19 has thinner watch face.

12-2042 FiT-18 12-2043

# FiT-19

### transmitter / straps

standard with watches. 12-2050 TE-15 transmitter with strap

12-2052 elastic strap


### pedometer











12-1951 b distance traveled

12-1952 c steps taken and distance traveled

12-1953 d steps taken and distance traveled, credit card size

- large digital display
- hinged protective case
- · loss prevention strap
- · user changeable batteries

### Ekho® pedometers

12-1940 steps taken

12-1941 steps, distance, calories

12-1942 steps, distance, calories, activity time

### Polar® wrist watch heart rate monitor



FT1 12-2152





FT7F 21-2154F







RS300X 12-2157

RS400 12-2157

· come with strap, transmitter, and manual

- heart rate BPM; HR-based target zones with visual/audible alarms; manual target zone; plus items specified below
- watch features backlight, date/weekday, time, 30m water resistant

FT1 - shows heart rate; manual heart rate target zone; displays summary of workout; records average and maximum heart rate 12-2152 FT1 Black

FT4 - shows when you're improving fitness based on your heart rate; displays calories burned

12-2153M FT4 Male Black/Silver FT4 Female Bronze/Bronze

FT7 - EnergyPointer tells you if the main effect of your training is fitness improvement or fat burning; displays calories burned; compatible with FlowLink to load data to computer

12-2154M FT7 Male Black/Silver 12-2154F FT7 Female Black/Gold

FT60 - creates training program based on personal goals and sets new weekly training targets; gives feedback on the effect of training; measures aerobic fitness at rest and displays progress

12-2155M FT60 Male Black/White display

12-2155F FT60 Female Purple

RS300X - speed, distance, and heart rate; helps you train at the right intensity with personal sport zones; measures aerobic fitness at rest and displays progress; displays calories burned RS300X Running Watch

RS400 - for endurance athletes and runners who rely on systematic training; measures detailed performance information; with software RS400 Advanced Running Watch

12-2158 IrDA USB adapter links data to ProTrainer5 software

12-2159 FlowLink data transfer device

### stopwatch



- start-stop, timeout, split time, calendar modes
- registers 1/100th for first 30 seconds and 1 second thereafter up to 24 hours
- programmable alarm

12-2100 stopwatch

### timers



- rings after pre-set period elapses
- · long ring time

12-2004 60 minute

### metronomes



- 12-2031: speed adjusts 40-208 BPM; sweep LED; stand; requires 9V battery
- 12-2032: speed adjusts 40-250 BPM; includes watch battery; beep sound, swing out back; pocket size

12-2031 electric, analog dial input 12-2032 electric, digital input



# manipulation and dexterity tests

### Purdue<sup>®</sup> pegboard



- · used extensively to measure dexterity for two types of activity: gross movements of hands, fingers and arms; and fingertip dexterity
- complete with pins, collars, washers and examiner's manual with norms

12-3000 pegboard 12-3001 extra parts

12-3002 record blanks (25)

### Lafayette® pegboard



- measure a timed unilateral motor response
- subject places 25 pegs in round holes as quickly as possible across the rows
- scores should improve with successive trials

12-3005 Lafayette peg

O'Conner Dexterity Tests

· finger test requires hand placement of three pins per hole · has been successfully used as a predictor for skills involving

• tweezers test requires the placement of a single pin in a 1/16"

high score indicates manual aptitude for work involving precision

12-3030

12-3031

### grooved pegboard



- manipulative dexterity test consisting of 25 holes with randomly positioned slots
- pegs with a "key" along one side must be rotated to match hole before inserted
- requires more complex visual-motor coordination

12-3010 grooved

### 12-3011 extra pegs (32)

# Minnesota manual



- · manual dexterity test meas-
- · complete with carrying case,

12-3050 MN Dexterity 12-3051 extra pegs, each

# dexterity



- ures capability for simple eye-hand-finger movement
- record blanks and norms

12-3052 record blanks (50)

### extra pins (105) prehension bench

finger dexterity test

dexterity test

rapid manipulation of small objects

diameter hole with tweezers

board has 100 holes

12-3020

12-3031



- graded resistance and prehension for both single-hand and bilateral activities
- 36 springs deliver resistance up to 21 pounds
- different surface textures on the knob aid in decreasing hypersensitivity

10-1162 bench



### extra pins (105) work bench

dexterity test

tweezers dexterity test



- improve strength and endurance of shoulders, elbows, wrists and hands; fine motor coordination
- front allows bilateral work that duplicates a hidden, confined work area
- toolbox with 6" and 8" adjustable wrenches, 2 screwdrivers and 4 hex keys

10-1163 work bench

### hand-tool dexterity test



- measure proficiency in using simple mechanic's tools
- consists of tools and two uprights with bolts
- the object is to disassemble all bolts from one upright and reassemble them on the other upright
- includes norms 10-1164 hand-tool

### Roeder manipulation



- measures hand, arm and finger dexterity and speed
- · the test has two parts: sorting and assembling
- complete with instructional manual, norms and record blanks

12-3080 Roeder 12-3082 extra parts 12-3081 record blanks (50)

### complete Minnesota



- measures eye-hand coordination, speed and dexterity
- · test placing, turning, displacing, one-hand turning, and two hand turning and placing
- 2 folding boards, 60 blocks, carrying case and manual

12-3060 MN test 12-3051 extra pegs, each

### pipe-tree




- · exercise for bilateral coordination, problem solving
- includes 31 pieces of 1/2" PVC pipe
- cover incorporates a set-in PVC pipe and is used as base of construction; accommodates all the piping for convenient storage

10-0766 pipe tree

# 3B°

# 3B Scientific® anatomical models



### Stan the Classic Skeleton

This Classic skeleton model (we call him Stan) has been the standard of quality in hospitals, schools, universities, and laboratories for over 50 years. Ideal for active use in teaching and demonstrations, this is the most affordable full-size skeleton of this quality available anywhere! Don't settle for imitations which compromise quality in workmanship and materials. Stan's the man! Comes on waist mounted roller stand. 67" 17lbs.

12-4500 Stan the Classic Skeleton

### Sam The Super Skeleton Does it All!

Sam has all the standard benefits of Stan the Classic Skeleton. Plus a fully flexible vertebral column, hand-painted muscle origins and insertions, numbered bones, flexible joint ligaments, and an L3-L4 disc prolapse. Over 600 structures of medical/anatomical interest are identified. 67" 18lbs.

12-4503Stan the Super Skeleton12-4501Max the Muscle Skeleton12-4502Leo the Ligament Skeleton12-4504Phil the Physiological Skeleton





### Classic Skull, 3-part

Perfect for basic anatomical studies or an attractive medical gift. Natural cast, medical quality 3-part skull. 8x5x6" 2lbs.

12-4547 Classic 12-4550 Numbered 12-4549 Painted

### Skull, 22-part

The Skull Kit is a natural cast replica that can be disassembled into its 22 individual bones. The individual bones connect by means of inconspicuous, stable connectors attached at slightly simplified skull sutures. The bones are depicted in 9 didactic colors and each pair of bone plates has the same color. 8x5x6" 21lbs.

12-4553 Didactic Color 12-4554 Natural Bone



### Classic TMJ Skull with Masticatory Muscles, 2-part

Demonstrate the function of the masticatory muscles with jaw occlusion, the initial stage of jaw opening and the movements of the mandible to the side and front. Removable skullcap. 8x5x6" 2lbs.

12-4548 Natural Bone

### Skull on Cervical Spine, 4-part

Skull with movable jaw and seven cervical vertebrae flexibly mounted. Includes representations of the hindbrain, spinal cord, cervical nerves, vertebral arteries, basilar artery and rear cerebral arteries.

8x5x6" 3lbs.

12-4552 Didactic Color 12-4551 Natural Bone



### Shorty the Mini Skeleton

Made from durable, unbreakable plastic. Featuring a removable 3-part skull, removable arms and legs, and flexible spine with soft intervertebral discs. Muscled Shorty features painted muscle origins and insertions. 37" 41lbs.

12-4506 Shorty, on base 12-4505 Shorty, hanging 12-4508 Muscled Shorty, on base 12-4507 Muscled Shorty, hanging



### Loose Bones

Made from durable, unbreakable plastic.
All medical quality, natural casts.
Wire Mounted. Specify left (L) or right (R)
following product number.

12-4580R/L Hand Skeleton 12-4581R/L Hand with Ulna 12-4582R/L Arm Skeleton 12-4583R/L Arm with Scapula 12-4584R/L Foot Skeleton 12-4585R/L Foot with Ulna 12-4586R/L Leg Skeleton 12-4587R/L Leg with Hip Bone



# 3B Scientific® anatomical models

### Classic Flexible Spine

Our most popular spines are also our most affordable. Fully flexible and designed for hands-on demonstrations. Features full male pelvis and occipital plate, L3-L4 disc prolapsed, spinal nerve exits, cervical vertebral artery. Stand is not included, please see 12-4546. 29" 4lbs.

12-4529 Flexible Spine

12-4530 Flexible w/ Femurs

12-4531 Flexible w/ Femurs, Painted

12-4532 Spine with Female Pelvis

### **Didactic Flexible Spine**

The same features as our Classic Spine, additionally featuring Didactic painting to differentiate the five sections of the spinal column.

12-4536 Didactic Flexible Spine

12-4537 Didactic Spine w/ Femurs

### **BONElike™ Youth Vertebral Column**

A true-to-life model of a 5-year-old child's vertebral column. Flexible vertebral column includes occipital plate, pelvis and sacrococcyx. Within the spinal canal, the spinal cord with cauda equina and exiting nerve roots are represented in flexible material.

12-4545 Youth Vertebral Column

# Multifunctional Stand for Spines & Skeletons, 3-part

Nickel-plated steel. Can be placed on the floor or mounted on the wall

12-4546 Multifunctional Stand

### **Vertebral Degeneration Model**

Natural cast model illustrates degenerative changes to the vertebrae and intervertebral discs of the lumbar spine in five progressively serious degrees. From healthy disc and vertebrae all the way to very severe degeneration and bony constriction of the left intervertebral foramen with pressure on the spinal nerve root L5. 9" 1lb.


12-4542 Vertebral Degeneration



### Osteoporosis Model

Magnetic model compares osteoporotic and normal thoracic vertebrae using T11 and T12. Reproductions of sequential osteoporotic thoracic vertebrae with narrower intervertebral disc are located on the left of the stand. Two healthy vertebrae with discs are on the right side. A detailed illustration on the base depicts two 3D micro CT images obtained from bone biopsies. 6x4x4" 1lb.

12-4544 Osteoporosis Model





### **Cervical Section**

Consisting of occipital plate, the 7 vertebrae with discs, cervical nerves, vertebral arteries and spinal cord. 7" 1lb.

12-4539 cervical



### **Thoracic Section**

12 thoracic vertebrae with discs, thoracic nerves and spinal cord. Flexibly mounted on stand. 13" 1lb.

12-4540 thoracic



### **Lumbar Section**

5 lumbar vertebrae with discs, sacrum with flap, coccyx, spinal nerves and spinal cord. 13" 1lb.

12-4541 lumbar



Sports Shoulder

Comprises the upper half of the humerus, clavicle and shoulder blade. The muscles of the rotator cuff are displayed and the sites of origin (red) and insertion (blue) are highlighted in color. By removing the four individual muscles, all movements of the shoulder joint can be performed: abduction, adduction, inward rotation, outward rotation, anteversion, retroversion, elevation and circumduction. 7x7x9" 2lbs.

12-4525 Sports Shoulder

### **Hip Joint** with Removable Muscles, 6-part

Adult male hip joint with the individual muscles and color-coded muscle origins and insertions on the femur and the hip bone. The hip muscles are removable. 7x13x7" 2lbs.

12-4526 Hip Joint

### **Elbow Joint** with Removable Muscles, 8-part

Adult male elbow joint with individual muscles and color-coded muscular origins and insertions on the humerus, radius, and ulna. The muscles are removable. 9x16x10" 2lbs.


12-4528 Flhow Joint

### **Muscle Knee Joint** with Removable Muscles, 12-part

Muscle Knee shows different removable muscles and muscle portions of the knee area. Color-coded and raised areas indicate the muscle origin and insertion points on the femur, tibia, and fibula. 13x7x7" 2lbs.

12-4527 Knee Joint

### **Spinal Cord with Nerves**



the spinal nerve with its four branches, and the three layers of the spinal column. 10x7x5" 1lb.

12-4561 Spinal Cord with Nerves

# 3B Scientific® anatomical models



### Muscle Arm, 6-part

This model illustrates both the superficial and deeper muscles, five of which are removable. Tendons, blood vessels, nerves and bone components of the left arm and shoulder are shown in great detail. Regular is 3/4 lifesize, 24x7x7" 4lbs. Deluxe is lifesize 31x10x10" 9lbs.

12-4556 Regular Arm, 6-part 12-4558 Deluxe Arm, 6-part

### Muscle Leg, 9-part

The model illustrates both the superficial and deeper muscles, eight of which are removable. Tendons, vessels, nerves and bone components of the left leg and foot are shown in great detail. Regular is 3/4 lifesize 24x7x7" 4lbs. Deluxe is lifesize, 31x10x10" 9lbs.

12-4555 Regular Leg, 9-part 12-4557 Deluxe Leg, 7-part



Foot Skeleton with Ligaments and Muscles, 6-part

The model features not only the bones but also the muscles, tendons, ligaments, nerves, arteries, and veins. The frontal view features the extensor muscles of the lower leg. Tendons can be followed to their insertion points. On the dorsal portion of the model the gastrocnemius muscle is removable to reveal deeper anatomical elements. The sole of the foot is represented in three layers. 9x10x7" 2lbs.

12-4524 Foot Skeleton, 6-part

### **Hand Skeleton**

with Ligaments and Muscles, 4-part

Features the bones, muscles, tendons, ligaments, nerves, arteries, and veins of the hand and lower forearm. The dorsal side shows the extensor muscles as well as portions of the tendons at the wrist as they pass under the extensor retinaculum. The palmar face of the hand is represented in three layers. 13x5x5" 1lb.

12-4522 Hand Skeleton, 4-part



Hand & Foot Skeleton with Ligaments

All the anatomically important ligaments and tendons, large and small, are shown on these two detailed models. Hand includes the interosseous membrane between the radius and ulna, removable flexor retinaculum. Foot includes the Achilles and peroneus longus tendons of the ankle.

Hand Skeleton 12x6x4" 1lb. 12-4521 Hand

Foot Skeleton 9x7x12" 1lb. 12-4523 Foot

### **Joint Models**

High-quality, natural size joint models with ligaments show the anatomy and all possible physiological movements.

12-4513 Functional Shoulder

12-4509 Classic Shoulder

12-4519 Miniature Shoulder

12-4515 Functional knee

12-4511 Classic knee


12-4518 Miniature knee

12-4516 Functional elbow 12-4512 Classic elbow

12-4520 Miniature elbow


12-4514 Functional hip

12-4510 Classic hip 12-4517 Miniature hip





# 3B Scientific® anatomical models



### Classic Brain, 5-part

A midsagittally sectioned original cast of a real human brain. The left half has removable frontal and parietal lobes, temporal and occipital lobes, encephalic trunk and cerebellum. 5x6x7" 1lb.

12-4562 Classic Brain

12-4560 Neuro-Anatomical Brain

### Classic Eye, 6-part

This model, 3 times life-size, dissects into both halves of sclera with cornea and eye muscle attachments, both halves of choroid with iris and retina, lens and vitreous humor. 4x4x6" 1lb.

12-4563 Classic Eye, 6-part 12-4564 Deluxe Eye, 7-part



Our Classic Heart is highly detailed and at a price you will love. Just slightly smaller than life-size with exquisite detail throughout including ventricles, atria, valves, veins, and the aorta. The front heart wall is detachable to reveal the chambers and valves inside. 7x5x5" 1lb.

12-4567 Classic Heart 12-4568 Heart with Bypass

### Larynx, 2-part

Our functional model features all the structures of the larynx including epiglottis, vocal cords, arytenoid cartilage, hyoid bone.

12-4569 Functional Larynx 12-4571 Classic Larynx





### Classic Giant Ear, 3 times life-size, 4-part

This ear with representations of the outer, middle, and inner ear. Removable eardrum with hammer, anvil and stirrup. 2-part labyrinth with cochlea, and auditory/balance nerve. 13x6x7" 3lbs.

12-4565 Deluxe Ear, 6-part

# 12-4566 Giant Ear, 4-part


### 3B Scientific® anatomical charts



Human Brain 12-4600



Foot Deformities 12-4601



Acupuncture 12-4602



Ear Acupuncture 12-4603



Foot Reflex Zone



12-4575

12-4574

12-4572



Female Pelvis, 6-part

Details the topography of bones, liga-

ments, vessels, nerves, pelvic floor muscles and pelvic organs. The whole pelvic

floor is represented with partially remov-

able external anal sphincter, external

urethral sphincter, deep and superficial

transverse perineal and bulbospongio-

sus muscles. Rectum, uterus with fallop-

ian tubes and ovaries and vagina are

removable and can be disassembled. The right pelvic half demonstrates the divisions

and topographical anatomy of the com-

mon iliac artery, the external and internal

artery, the common iliac vein and the ex-

Female Pelvis, 6-part

Female Pelvis, 4-part

Female Pelvis, 3-part

ternal iliac vein. 7x11x7" 4lbs.

Arthritis 12-4605



Human Ear 12-4606



Human Eye 12-4607



12-4608



Anatomy 12-4609



Human Heart 12-4610





12-4611



Larynx 12-4612



Lymphatic System 12-4613



Musculature 12-4614



12-4615





Pelvis and Hin 12-4617



Rheumatism 12-4618



Shoulder and Elbow 12-4619



Human Skeleton 12-4620



Human Skull 12-4621



Spinal Column 12-4622



Sports Injury 12-4623



12-4624



laminated (L) 20.00 sticky back (S) 35.00

specify P (paper), S (sticky), L (laminated) when ordering







WaxWel™ unit can be ordered without paraffin and accessories

• ideal for large departments

11-1601 WaxWel™ unit only

### WaxWel™ paraffin refills (blocks)



- · convenient paraffin blocks for use in all paraffin units
- individually sealed in plastic bags, packaged six 1-pound bags to a box

scent	6 pound box	36 pound case
unscented	11-1720-6	11-1720-36
wintergreen	11-1722-6	11-1722-36
lavender	11-1714-6	11-1714-36
citrus	11-1716-6	11-1716-36
peach	11-1718-6	11-1718-36
rose blossom	11-1730-6	11-1730-36

# paraffin heat therapy

### WaxWel™ paraffin unit

- use to apply heat to hand or foot
- heavy duty for continuous use in high traffic clinics
- unit has oversized tank that can accommodate a man's size 12½ foot
- unit is lightweight and portable and can even be used in the home

# comes with washable terry cloth mitt and bootie, 65 hand/foot plastic liners and 6 pounds of wax

11-1600	unit with 6 pounds unscented wax
11-1603	unit with 6 pounds wintergreen wax
11-1604	unit with 6 pounds lavender wax
11-1605	unit with 6 pounds citrus wax
11-1606	unit with 6 pounds <b>peach</b> wax
11-1607	unit with 6 pounds rose wax
	11-1603 11-1604 11-1605 11-1606



### WaxWel™ paraffin refills (beads)



- when additional paraffin is needed, simply add proper amount of beads to fill tank
- beads are easy to pour (dispense) and quick to melt
- available unscented and in five favorite fragrances
- packaged in one-pound "zip-lock" plastic bags, 6 per box
- 3 pounds of beads in handy carrying bucket

scent	3 lb bucket	6 pound box	36 pound case
unscented	11-1750-3	11-1750-6	11-1750-36
wintergreen	11-1752-3	11-1752-6	11-1752-36
lavender	11-1744-3	11-1744-6	11-1744-36
citrus	11-1746-3	11-1744-6	11-1746-36
peach	11-1748-3	11-1748-6	11-1748-36
rose blossom	11-1760-3	11-1760-6	11-1760-36



# paraffin heat therapy

### Therabath®





with 6 lb. wintergreen wax

with 6 lb. unscented wax

- · hospital grade paraffin warmer with lid, grill, and guide
- inside unit dimensions: 6.75" W x 5" H x 12.9" L
- operating temperature of 130° F

11-1171 with 6 lb.unscented wax 11-1170 with 6 lb.wintergreen wax


### Dickson® paraffin baths





PB101





- Dickson PB-104: 14x7x9", 20lb. capacity, includes removable mobile stand. Includes timer-controlled sterilization and heated drain
- Dickson PB-105: 25x11.5x10", 50lb. capacity, includes removable mobile stand. Includes timer-controlled sterilization and heated drain

11-1150 PB-101 with 6 pounds paraffin

11-1151 PB-107 with 6 pounds paraffin

11-1152 PB-104 with stand with 20 pounds paraffin

11-1153 PB-105 with stand with 50 pounds paraffin

11-1161 mobile stand for PB-101 and PB-107

### Parabath<sup>®</sup>



- lightweight+portable stainless steel tank for clinic or home use
- for the treatment of patients suffering from pain associated with arthritis, strains, sprains, joint stiffness
- also used to increase range of motion
- deep oversized tank (6"W x 6"H x 12" L) allows easy immersion of hand and wrist, foot and ankle, or elbow
- stainless steel tank with insulated outer housing to retain heat
- operating temperature is 125°-134° F

11-1110 with 6 lb. unscented wax

11-1119 stand

### Whitehall® ParaTherapy



- Whitehall PT-06-s: hand paraffin bath 12"L x 6"W x 6"D, 6 lb capacity
- Whitehall PT-18-s: arm or hand paraffin bath 20"L x 11.5"W x 6"D, 18 lb capacity

11-1148 PT-6-s with 6 pounds paraffin 11-1149 PT-18-s with 18 pounds paraffin



# topical analgesics - hot and cold pain relief

### Point Relief® ColdSpot® pain relief gel and spray



- offers soothing relief
- all natural ingredients, including menthol and MSM
- available in gel, spray and roll-on applicators
- temporary relief of minor aches pains of muscles and joints associated with simple backache, arthritis, bruises, strains and/or sprains
- water-based formula minimizes grease and residue

water-based formula	a minimizes grease and residue		
each	box [qty]	case [qty]	
gel tube			
<b>4 oz</b> 11-0730-1	11-0730-12	11-0730-144	
gel dispenser bottle			
<b>16 oz</b> 11-0710-1		11-0710-24	
<b>32 oz</b> 11-0711-1		11-0711-16	
<b>1 gal</b> 11-0712-1		11-0712-4	
gel packet			
<b>5 gr</b> 11-0740-1	11-0740-100	11-0740-1000	
roll-on			
<b>3 oz</b> 11-0720-1	11-0720-12	11-0720-144	
spray			
<b>2 oz</b> 11-0700-1	11-0700-12	11-0700-144	
<b>4 oz</b> 11-0701-1	11-0701-12	11-0701-144	
<b>16 oz</b> 11-0702-1		11-0702-18	



ColdSpot™ Point Relief™ spray is also available in Be Better™ kits

Point Relief® HotSpot® pain relief and massage gel

- warming cinnamon scent
- all natural ingredients, including capsaicin
- available in gel application
- temporary relief of minor aches and pains of muscles and joints associated with simple backache, arthritis, bruises, strains and/or sprains
- warming gel offers soothing massage when rubbed in with hands
- water-based formula

gel	each	box [qty]	case [qty ]
4 oz	11-0780-1	11-0780-6	11-0780-24
8 oz	11-0781-1		11-0781-12
32 oz	11-0782-1		11-0782-8





# topical analgesics - hot and cold pain relief

### BioFreeze® pain relief

- contains an ILEX base, and a herbal extract for a quick onset, long lasting, pain relief
- BioFreeze® Cryotherapy® pain relief formula is greaseless, stainless, and has a vanishing scent
- lightly massage Bio-Freeze<sup>®</sup> into all sides of sore muscles
- use to relieve pain from sore muscles and muscle spasms, back,

shoulder and neck pain, arthritis, painful ankle, knee, hip & elbow joints, and muscle strains

use to relieve pain prior to massage therapy, soft tissue trigger point therapy, rehabilitation exercises, and pre and post workout stretch

	each	case	quantity ——
gel tube			
f 4 oz. tube	11-1031-1	11-1031-144	144 ea
4 oz. applicator	11-1030-1	11-1030-144	144 ea
gel dispenser bo	ttles		
<b>e</b> 16 oz.	11-1033-1	11-1033-24	24 ea
<b>a</b> 32 oz.	11-1034-1	11-1034-16	16 ea
c 1 gallon	11-1035-1	11-1035-4	4 ea
roll-on			
2 oz.	11-1039-1	11-1039-144	144 ea
<b>b</b> 3 oz.	11-1032-1	11-1032-144	144 ea
spray			
<b>g</b> 4 oz.	11-1037-1	11-1037-144	144 ea
<b>h</b> 16 oz.	11-1038-1	11-1038-18	18 ea
wipes			
24 / box	11-1029-24	11-1029-1152	48 boxes
packets			
<b>d</b> 1/4 oz. 100/box	11-1036-100	11-1036-1000	10 boxes


# Polar Frost® cold gel [A] [B] [C] [D] [E]

- cold gel is specially formulated to provide soothing relief for soft tissue stress injuries, muscle tension, inflammation and stiffness
- aloe vera has an anti-inflammatory effect and keeps the skin smooth and moist throughout treatment

each		case qu	ıantity		
		_	Cacii	case qu	acitity
	gel tu	be			
	[E]	5 oz.	11-0107-1	11-0107-12	12 ea
	gel di	spenser bottl	les		
	[B]	17 oz.	11-0106-1	11-0106-12	12 ea
	[A]	1 gallon	11-0103-1	11-0103-4	4 ea
	roll-o	n			
	[D]	2.5 oz.	11-0108-1	11-0108-24	24 ea
	spray				
		7.5 oz.	11-0109-1	11-0109-12	12 ea
	packe	ets			
	[C]	¼ oz. pack		11-0105-200	200 ea

### Prossage™ heat pain relief

- 100% natural therapeutic warming ointment for myofacial release, deep tissue work and trigger point therapy
- Prossage<sup>™</sup> heat developed to produce gentle penetrating heat to the skin, fascia, muscles & other soft tissue structures
- as tissues soften, more blood and oxygen penetrates congested areas reducing the ischemia responsible for many chronic pain syndromes, orthopedic conditions and common sports injuries



 Prossage™ heat is an area specific warming ointment, not an all over bodywork and massage oil

		each	— case quantity —
gel	bottles		
Α	3 oz.	11-0200-1	11-0200-24 24 ea
В	8 oz.	11-0201-1	11-0201-12 12 ea
С	32 oz.	11-0202-1	11-0202-8 8 ea

### Flexall® 454 pain relief





- help increase patient compliance and enhance ultrasound, cryotherapy, TENS, and massage therapies
- Flexall (7% menthol) and Maximum Strength Flexall (16% menthol) are aloe vera gel formulas providing fast, longlasting relief from arthritis, backache, bursitis, tendonitis, muscle strains and sprains, bruises, and cramps
- gels absorb quickly and completely, are greaseless, nonstaining, and gentle on the skin

Standard stre	ength each	 cas	e quantity ——
4 oz.	11-0220-1	11-0220-12	12 ea
16 oz.	11-0221-1	11-0221-6	6 ea
7 lb. w/ pump	11-0222-1		
packet 1/4 oz. 24/box	11-0223-24	11-0223-144	6 boxes
Max Strength gel bottle	each	 cas	e quantity ——
1½ oz.	11-0224-1	11-0224-24	24 ea
3 oz.	11-0225-1	11-0225-12	12 ea
16 oz.	11-0226-1	11-0226-6	6 ea
7 lb. w/ pump	11-0227-1		



# topical analgesics - instant hot and cold pain relief

### CryoDerm<sup>®</sup> cold and ThermoDerm<sup>®</sup> hot products

# CryoDerm® cold pain relie La Contains: lyptus, IL

pain relieving CRYOTHERAPY that Lasts Colder Longer™

contains: MSM, arnica, boswella, eucalyptus, ILEX, menthol, peppermint oil

gel each	box quantity	case quantity
E 4 oz 11-0404-1	11-0404-12 12	11-0404-144 12 boxes
<b>C</b> 16 oz 11-0405-1		11-0405-20 20 each
gel packet		
<b>G</b> ½ oz 11-0406-1	11-0406-108 108	11-0406-1080 10 boxes
spray		
<b>D</b> 4 oz 11-0401-1	11-0401-12 12	11-0401-144 12 boxes
<b>B</b> 16 oz 11-0402-1		11-0402-20 20 each
<b>A</b> 32 oz 11-0403-1		11-0403-12 12 each
roll-on lotion		
<b>F</b> 3 oz 11-0400-1	11-0400-12 12	11-0400-144 12 boxes

### ThermoDerm® hot



pain relieving warm therapy that Lasts Warmer Longer™

contains: MSM, arnica, boswella, menthol, green tea extract, aloe, yucca, tocopheryl acetate (vitamin E), ILEX

lotion each	box quantity	case quantity
I 4 oz 11-0244-1	11-0244-12 12	11-0244-144 12 boxes
<b>H</b> 16 oz 11-0245-1		11-0245-20 20 each
packet		
<b>K</b> 1/4 oz 11-0246-1	11-0246-108 108	11-0246-1080 10 boxes
roll-on lotion		
<b>J</b> 3 oz 11-0240-1	11-0240-12 12	11-0240-144 12 boxes

- · pain relieving therapy that lasts longer
- targets the pain spasm cycle, addressing the secondary effects of a subluxation, pain, ischemia, hypertonicity and inflammation
- use spray, gel or roll-on to relieve muscle spasm, reduce inflammation and decrease spindle cell response
- clear water-based formula never balls up or peels off skin and maintains glide
- · no artificial colors or isopropyl alcohol
- complements physical therapy, chiropractic and podiatric care, massage therapy, rehabilitation exercise, pre and post athletics

### 3 functional dispensers

### gel lotion

Massaging gel onto the patient enhances it's effects. Human Touch has been shown to speed up and enhance the healing process.

### roll-on

Many patients like the Roll-On application because they can apply pressure to painful trigger points and sore muscles.

### spray

Effective in acute conditions where manipulation or massage is contraindicated. Sprays at any angle, allowing patients to treat even the hardest to reach places without assistance

### topical skin refrigerants



Spray & Stretch



Instant Ice




Pain Ease Rx only

Topical skin refrigerant sprays which are used for temporary relief of pain associated with injections, minor surgical procedures, starting IV's, venipuncture, myofacial pain, restricted motion, muscle spasm, and minor sports injuries.

the state of the s		
	each	case quantity
Spray & Stretch cold spray, fine stream, 3.5 oz.	11-0800-1	11-0800-12 12 ea.
· · · · · · · · · · · · · · · · · · ·		
Instant Ice cold spray, mist, 3.5 oz.	11-0810-1	11-0810-12 12 ea.
Instant Ice cold spray, medium stream, 3.5 oz.	11-0811-1	11-0811-12 12 ea.
Pain Ease cold spray, mist 3.5 oz.	11-0830-1	11-0830-12 12 ea.
Pain Ease cold spray,	11-0831-1	11-0831-12 12 ea.

### Cat's Claw™ gel





 effective pain relief for arthritis, neck, shoulder and back pain, painful joints, sore muscles, muscle sprains & spasms

1186
Sale

	each	case	quantity ——
gel	Caon	Case	quantity
4 oz. tube	11-0951-1	11-0951-12	12 ea.
16 oz. pump	11-0952-1	11-0952-12	12 ea.
1 gallon bottle	11-0953-1	11-0953-4	4 ea.
packet			
5 gram packet	11-0954-100	11-0954-1000	1000 ea.
roll-on			
3 oz. roll-on	11-0950-1	11-0950-12	12 ea.



# heat therapy

### digital moist heat pad



- hand controller digitally adjusts temperature and treatment time
- · includes flannel cover

### digital pads

shoulder (pictured)

11-1135 king 26x14" 11-1136 medium 18x14" 11-1137 mini 19x7"

11-1138




analog moist heat pad

- electric pad provides moist heat without adding water (fomentation)
- hand controller adjusts temperature
- includes flannel cover

### analog pads

11-1120 king 13x27" 11-1121 medium 13x13" 11-1123 small 4x14"

# heavy duty UV IR and UV/IR lamps





### ultra violet (UV)

- 800 watt bulb with wire guard
- light is 64% UVA/ 36% UVB
- yoke construction for +/- 180° movement
- 10 minute timer
- #153 protective goggles
- 5 caster base 42-621/2" height

18-1220 [A] desk with analog timer 18-1225 [A] desk with digital timer 18-1280 [C] mobile with analog timer

18-1393 800 watt UV bulb

### IR 600 watt / UV 800 wat

- 600 watt IR / 800 watt UV
- UV: same specs as UV description
- IR: (2) 300 watt ceramic bulbs
- toggle from IR to UV

18-1320 [A] desk [C] mobile 18-1380

18-1393 800 watt UV bulb

18-1373 300 watt IR bulb 2ea

### infra-red (IR)

### economy

- 175 watt ruby bulb
- 19" flexible arm
- protective wire guard
- stationary base 10%"
- 5 caster base 13" dia.
- 35-58" height

18-1138 18-1168 18-1135

18-1181

18-1372

18-1370

18-1371

18-1165

### standard

- 250 watt ceramics bulb
- 12" flexible arm
- protective wire guard
- stationary base 10%"
- 5 caster base 13" dia.
- 35-58" height
- [B] economy stationary [E] economy mobile [B] standard stationary
- [E] standard mobile
- [D] professional mobile

175 watt ruby bulb 250 watt ceramic bulb 750 watt ceramic bulb

### professional

- 750 watt ceramic bulb
- protective wire guard
- 30 minute timer
- variable intensity
- 5 caster base 23" dia.
- 31-52" height

### economy heat pad



- 3 heat settings
- washable cover
- add water to the sponge for moist heat; remove sponge for dry heat

### dry heat pad

11-1130 dry, small 12x15" 11-1131 dry, large 12x24"

### moist or dry heat pad

11-1132 moist/dry, small 12x15" moist/dry, large 12x24" 11-1133

### hot water bottle



• hospital grade rubber hot water bottle with screw stopper

2 quart capacity

11-1140 each 11-1140-12 dozen

### UV / IR / vitamin D lamps

- digital timer
- 800 watt UV light is 70% UVA/30% UVB
- table-top has adjustable angle lamp
- tripod floorstand is height adjustable
- vitamin D flourescent lamp for individuals who don't receive enough sunlight or can't tolerate vitamin D supplements, 2 pairs of protective eyewear included





[A] table top [C] vitamin D

### UV lamp (800 watts)

18-1235 UV w/ digital timer, table top UV w/ digital timer, tripod 18-1234

### UV (800 watts) IR (600 watts)

18-1233 digital timer, table top 18-1232 digital timer, tripod floorstand

### vitamin D flourescent

18-1331 tabletop model

### accessories

18-1390 800 watt UV bulb (each) 18-1391 flourescent vitamin D bulbs (pair)

18-1395 protective eyewear, amber (each)



# **RELIEF PAK®** moist heat packs and covers



and cover

11-1313
halfsize pack
with
11-1363
halfsize cover

11-1364
pocketed standard
terry foam cover

ideal for at-home use;
heat pack in water on the
stove or in the microwave

shown above - heat in water-filled pan on stove

- simple and effective method of applying moist heat; provides 30 minutes of therapeutic moist heat
- soothe discomfort from stiff neck, backache, arthritis, sprains and bruises
- pack can be heated in water on stove or in microwave
- when treatment is complete pack can be stored in a plastic bag in freezer
- always apply moist heat pack with a Relief Pak® cover or a heavy towel
- terry cloth cover allows moist heat to pass without burning the skin; covers are filled with foam filling or all terry; covers are washable; hook and loop fasteners in terry cover keep moist heat packs in place
- new circular pack size is ideal for home use

Relief P	ak® moist heat	packs	Ro	elief Pak <sup>®</sup> cove	r
hot pack circular (10" diameter)	<b>each</b> 11-1205	<b>bulk (12 each)</b> 11-1205-12	foam filled terry circular (12" diameter)	<b>each</b> 11-1226	
standard (10x12") neck (24") oversize (15x24") half size (5x12") spine small (10x18") spine large (10x24") knee-shoulder (10x20") hand (7x13")	11-1310 11-1311 11-1312 11-1313 11-1314 11-1315 11-1316 11-1321	11-1310-12 11-1311-12 11-1312-12 11-1313-12 11-1314-12 11-1315-12 11-1316-12 11-1321-12	standard (20x24") neck (9x24") oversize (24x30") halfsize (12x15") standard with pocket oversize with pocket all terry microfiber circular (12" diameter)	11-1360 11-1361 11-1362 11-1363 11-1364 11-1365 each 11-1227	
	k <sup>®</sup> pack and co	over sets	standard (20x24") neck (9x24")	11-1366 11-1367	
hot pack AND foam-filled terry cover circular set	<b>each</b> 11-1308	<b>bulk (12 each)</b> 11-1308-12	oversize (24x30") halfsize (12x15")	11-1368 11-1369	
standard set neck set	11-1300 11-1301	11-1300-12 11-1301-12	hand (14x14") dual hand (23x14")	11-1376 11-1377	



bulk (12 each)

bulk (12 each)

11-1226-12 11-1360-12 11-1361-12 11-1362-12 11-1363-12 11-1364-12 11-1365-12

11-1227-12 11-1366-12 11-1367-12 11-1368-12 11-1369-12 11-1376-12 11-1377-12

# Hydrocollator® moist heat packs and covers







terry covers

- eight sizes designed to fit every body contour
- provides up to 30 minutes of deep soothing, moist heat
- brown tabs on most packs for easy rotation in busy clinics
- cost effective, provides hundreds of treatments from one pack

Hydrocollator <sup>®</sup> packs			
	each	bulk (12 each)	
e standard (10x12")	00-1006	00-1006-12	
<b>c</b> neck (24")	00-1002	00-1002-12	
h oversize (15x24")	00-1004	00-1004-12	
f half size (5x12")	00-1014	00-1014-12	
a spine small (10x18")	00-1008	00-1008-12	
<b>b</b> spine large (10x24")	00-1010	00-1010-12	
g knee-shoulder (10x20")	00-1012	00-1012-12	
<b>d</b> hand (6½x12")	00-1026	00-1026-12	

### Hydrocollator® pack and cover sets

_	•	
	each	bulk (12 each)
standard set	00-1064	00-1064-12
neck set	00-1062	00-1062-12

### Hydrocollator® covers foam-filled bulk (12 each) 00-1102-12 j standard (20x24") 00-1102 00-1120-12 k neck (9x24") 00-1120 00-1124-12 i oversize (24x30") 00-1124 standard with pocket 00-1108 00-1108-12 00-1110-12 oversize with pocket 00-1110 all terry bulk (12 each) 00-1118-12 i standard (20x24") 00-1118 j neck (9x24") 00-1104 00-1104-12 k oversize (24x30") 00-1126 00-1126-12 00-1131-12 hand (11x14") 00-1131 00-1132-12 dual hand (23x14") 00-1132

00-4016  $\,$  wall mounted towel rack 33"x 2" w/ six 13" hooks

### hot pack utility cart

# Hydrocollator® stationary heating units





- heating units are small enough to be kept on the counter-top
- water bath keeps packs at 168°F

<ul> <li>water bath keeps packs at 168°F</li> </ul>			
	units and	d packs	HxLxW
00-2102	E-1 with	4 standard	15" x12½" x8"
00-2102-2	E-1 with 3	3 standard, 1 neck	15" x12½" x8"
00-2102-3	E-1 with 2	2 standard, 2 neck	15" x12½" x8"
00-2802	E-2 with 2-OS, 3-std.,1-neck		16½" x14" x12"
11-1390 11-1391	cart, hold cart, hold		
		E-1 insulation wrap E-2 insulation wrap	
	11-1397	hot pack tongs	

### Hydrocollator® mobile heating units



- mobile units are large enough for a complete department
- no plumbing required; simply fill with water and plug in
- thermostatically controlled to maintain a temperature of 168°F

	units and packs	HxLxW
00-2302	SS-2 with 8 standard	33" x 20¾" x 15½"
00-2402	M-2 with 12 standard	33½" x 26½" x 15½"
00-2502	M-4 with 24 standard	33" x 203/4" x 20"
00-4010 00-4230	optional side rack extra shelf for side rack	



n h

> tŀ е

# RELIEF PAK® reusable cold packs and covers

• soothe discomfort from injury, sprains and strains, muscle soreness, bruises, swelling and tension headache









11-1000 standard

11-1002 oversize

11-1003 half size

11-1007 circular



11-1001 neck

- · re-usable cold pack is an effective cold therapy
- · place pack in a chilling unit or freezer
- pliable and soft even below freezing
- maintains its therapeutic level for up to thirty minutes

• urethane packs made of heavier material that can be chilled t	to
lower temperatures (do not chill below 0° F)	

- black urethane and black vinyl packs have 1 year warranty
- covers provide convenient cloth barrier between pack and skin

Relief Pak® black urethane cold packs

- maintaine its therapedate level for up to thirty minutes				each	bulk (12 each)
durable heavy-duty materials			circular 10" diameter	11-1257	11-1257-12
<ul> <li>apply with hand towel between skin or pack, or in cover</li> </ul>		standard (11x14")	11-1250	11-1250-12	
Relief F	Pak® blue vinyl co	ld packs	neck (6x23")	11-1251	11-1251-12
	each	bulk (12 each)	oversize (11x21")	11-1252	11-1252-12
circular 10" diameter	11-1007	11-1007-12	half size (7x11"	11-1253	11-1253-12
standard (11x14")	11-1000	11-1000-12	Relie	f Pak® cold pack	covers
neck (6x23")	11-1001	11-1001-12		each	bulk (12 each)
oversize (11x21")	11-1002	11-1002-12	circular 10" diameter	11-1009	11-1009-12
half size (7x11")	11-1003	11-1003-12	standard (13x16")	11-1010	11-1010-12
quartersize (5x7")	11-1004	11-1004-12	neck (9x25")	11-1011	11-1011-12
throat (3x11")	11-1005	11-1005-12	oversize (13x23")	11-1012	11-1012-12
eye (3x8")	11-1006	11-1006-12	half/quarter (9x13")	11-1013	11-1013-12

### Nylatex® wraps



### • for holding packs in place

	g paono piaco
00-1200	assorted 9 ea. (1 each excluding 00-1213, 00-1264)
00-1202 00-1264 00-1204 00-1206	2½x18" 3 ea. 2½x24" 3 ea. 2½x36" 3 ea. 2½x48" 3 ea.
00-1208 00-1210 00-1212 00-1213	4x18" 3 ea. 4x36" 3 ea. 4x48" 3 ea. 4x60" 3 ea.
00-1214 00-1216 00-1218	6x18" 3 ea. 6x36" 3 ea. 6x48" 3 ea.

### Relief Pak® chilling units





front loading chilling unit (left); top loading chilling unit (right)

front loading chilling unit: 4.6 cubic foot capacity. 2 full width door shelves and 2 full width wire shelves. Adjustable thermostat control. Defrost drain, and manual defrost. 22" deep x 27" wide x 39" high. Weighs 137 lbs.

11-0501

11-0506 front loading with 12 standard Relief Pak™ (11-1000)

top loading chilling unit: 5.0 cubic foot capacity with 1 removable basket. Adjustable thermostat control, defrost drain and manual defrost. 24" deep x 30" wide x 32" high. 110 lbs.

11-0500 top loading

top loading with 12 standard Relief Pak™ (11-1000) 11-0505



# cold packs and cold therapy

### Col-PaC® reusable cold packs







- re-usable cold pack for cold therapy place pack in a chilling unit or freezer
- remain pliable and soft even below freezing
- maintains therapeutic level for up to thirty minutes
- apply with hand towel between skin and pack
- · black urethane pack is more supple and maintains its flexibility at a lower temperature than the blue vinyl pack

blue vinyl cold packs			
	each	bulk (12 each)	
standard (11x14")	00-1500	00-1500-12	
neck (6x23")	00-1508	00-1508-12	
oversize (11x21")	00-1512	00-1512-12	
half size (7x11")	00-1506	00-1506-12	
quartersize (5x7")	00-1504	00-1504-12	
throat (3x11")	00-1502	00-1502-12	
eye (3x8")	00-1510	00-1510-12	

# black urethane cold packs

	each	bulk (12 each)
standard (11x14")	00-1552	00-1552-12
neck (6x23")	00-1554	00-1554-12
oversize (11x21")	00-1556	00-1556-12
half size (7x11")	00-1562	00-1562-12

### AirCast CryoCuff®









ankle



thigh

shoulder

• combines compression and Cryo/Cuff® (cuff and cooler) cold: compression minimizes

swelling, while cold minimizes pain • water cooler, pressurized ice water, and cuffs made for

 cooler holds enough ice and water for 8 hours of cryotherapy

specific body parts

· compression is gravity con-

cuffs are made of nylon and vinyl

 use with gravity feed (included) or continuous AutoChill® electric pump (sold separately)

11-1550 ankle

11-1552 medium foot 11-1556 medium knee 11-1557 large knee 11-1560 shoulder XL shoulder 11-1561 11-1563 XL thigh 11-1564 calf

11-1565 back/hip/rib 11-1566 elbow

11-1562 thigh 11-1567 hand/wrist

11-1588 cooler only 11-1589 AutoChill® pump

please inquire for additional cuffs and sizes; cuffs can also be purchased separately or with AutoChill® pump unit; units above come with cooler

### ColPaC® Freezer



- freezer with 12 standard size packs
- five cubic capacity
- stainless steel shelves are designed to quickly freeze ColPaC®s
- stainless steel trays and cup holders, perfect for ice massages
- measures 33 H x 21 L x 26" W

09-0910K ColPac freezer unit with 12 standard packs

### ColPaC® chilling units





- · mobile unit has large storage and chilling capacity
- · requires no plumbing
- plugs into any 110v outlet

### packs included HxLxW (in.) 00-3102 C-5 6 std. 6 half 33x27x16"

00-4010 optional side rack 00-4230 extra shelf for side rack

### specialty packs





• packs filled with non-toxic silica gel that delivers up to 30 minutes of soothing relief of acute pain, swelling and fever

poo-	-boo pac™	
	each	dozen
green	00-1533-1	00-1533-12
blue	00-1534-1	00-1534-12
SPO	RT-PAC™	
	each	10-pack
soccer	00-1517-1	00-1517-10



### Relief Pak® Torex™ hot and cold sleeve







- hot & cold therapy sleeves, gently roll on to cover limb or finger
- sleeve design is easy to use and stays in place hands free
- 360° coverage treats the entire circumference of the injured area, delivering hot or cold deep into soft tissues
- · Stays soft when chilled or heated to mold to the body

	each	12 each
Torex hot/cold sleeve, finger	11-1590	11-1590-12
Torex hot/cold sleeve, small	11-1591	11-1591-12
Torex hot/cold sleeve, medium	11-1592	11-1592-12
Torex hot/cold sleeve, large	11-1593	11-1593-12
Torex hot/cold sleeve, x-large	11-1594	11-1594-12


### Flexi-PAC™ reusable hot and cold



- Flexi-PAC<sup>™</sup> can be microwaved to provide soothing heat, or stored in freezer to provide penetrating cold therapy
- non-toxic flexible gel packs conform to any body contour at any temperature

	each	bulk
5x6"	00-4026-1	00-4026-48 48 ea
5x10"	00-4020-1	00-4020-24 24 ea
8x14"	00-4029-1	00-4029-12 12 ea.

### TheraTemp® microwavable moist heat therapy



- beads and cotton combine to provide soothing moist heat
- place pack in microwave to heat; no water needed!

(e)

 Thera-Temp® packs produce moist heat by drawing moisture from humidity in the air and retain it in the therapeutic beads

### rectangular pack with elastic velcro belt

11-1281	10 x 12"	(b)
11-1282	6½ x 12"	(d)
11-1283	12 x 14"	(a)
11 1201	61/ v 10"	(0)

9 x 18"

### specialized wraps

11-1285

11-1280 collar/ abdomen 11-1286 eye/ sinus 11-1287 wrist

### Relief Pak<sup>®</sup> cold n' hot™ compress





11-1297 circular (10" diameter)

- hot and cold compress has a comfortable soft fabric that conforms to the curves of your skin; fabric is gentle on the skin, frost-free, and washable
- heat pack in microwave or chill in freezer
- when heated, soothes discomfort from stiff neck, backache, arthritis, sprains and bruises
- when chilled, soothes discomfort from injury, sprains and strains, muscle soreness, bruises, swelling and tension headaches
- · non-toxic, non-crystallizing gel allows pack to contour
- use standard velcro hook to strap pack in place

circular (10" diameter	—— each —— 1) 11-1297	——— bulk ——— 11-1297-12 12 ea
small (3x5") medium (7x12") large (10x13") neck (7x22")	11-1291 11-1292 11-1293 11-1294	11-1291-12 12 ea 11-1292-12 12 ea 11-1293-6 6 ea 11-1294-12 12 ea
tri-sectional (8x16")	11-1298	11-1298-6 6 ea



# cold therapy products

# English-style ice cap





• English style ice bags are made of waterproof rubberized fabric with an aluminum screw cap

diameter	each	dozen
6"	11-1060	11-1060-12
9"	11-1061	11-1061-12
11"	11-1062	11-1062-12

### ICE-UP™ ice on the spot!™



- portable carry cooler keeps Ice-Up™ stick frozen up to 12 hours!
- cap design easy to open with wet hands
- patented ice grippers allow full range of controlled massage pressure
- · simple and effective: just fill it, freeze it and carry it for Ice on the spot™!

11-1088 Ice-Up 11-1088-12 (12 each)

### hot/cold packs economy



economic hot/cold pack alternative dozen case-qty.

4x6"	11-1066	11-1066-96
6x8"	11-1067	11-1067-24
6x10"	11-1068	11-1068-18
8x10"	11-1069	11-1069-12



# CryoStim™ ice massager



- application of ice massage therapy
- used over broad areas of the body or over a small spot to stimulate local "trigger" points for fast, effective pain relief
- apply CryoStim gel over treatment areas to enhance conductivity of probe and prevent skin irritation
- · includes two probes, two ounces of gel

11-1092 ice massager 11-1093 ael. 2 ounces

# instant cold packs



- in 3 seconds the pack becomes a 34°F by just squeezing the bag
- · cold is achieved by safe chemical reaction
- double lined bag; prices per dozen

11-1020 5x10" 11-1021 4x6"

# **CryoCup™**



- reusable, easy to use, cost-efficient icer
- plastic with cool retardant handle
- fill cup with water and freeze
- · ready to use for ice

11-1090 CryoCup (12 each) 11-1090-12

### Roller Ice® ice massager





- · targeted cold or hot therapy
- ball rolls in all directions
- · bullet provides additional weight
- neoprene "bulb" cover insulates liquid reservoir and protects hand from cold

11-1094 ball massager 11-1095 bullet 11-1096 ball with magnet 11-1097 reservoir cover

# Relief Pak® insulated ice bags

 refillable, posable with velcro band or string ties

 easy fill neck opening



large bag: 7"x13", small bag: 5"x13"

velcro belt small ice bag large ice bag	— <b>10-pack</b> — 11-1241-10 11-1240-10
tie strings small ice bag large ice bag	11-1243-10 11-1242-10

# shortwave diathermy

# Chattanooga® Intelect® SWD100 shortwave diathermy

- safe, effective deep tissue heating therapy of shortwave diathermy
- full-color touch-screen interface, onboard clinical resources, customizable treatment parameters and more

### ideal for:

• pain relief

00-1600

- · reduction of muscle spasm
- decreasing joint stiffness

Intelect SWD100 shortwave diathermy

- contractures
- increased blood flow
- · chronic inflammatory conditions
- bursitis

# Mettler Auto\*Therm® 390 diathermy

Auto\*Therm® 390 features continuous and pulsed shortwave diathermy. Moves easily between treatment rooms when soft-rubber capacitive pads are being used. Can be mounted to optional cart to provide a solid platform for inductive coil drum electrode.



Input: 100-240 VAC, 50-60Hz, Frequency: 27.12 MH, Continuous mode: 100 W Average Power, Pulsed mode: 200 W Peak Power, Pulse frequency: 10, 20, 50, 100, 400 Hz, Pulse duration 60, 100, 200, 300, 400  $\,\mu$ s, Treatment time: 1-30 minutes, dimensions: 390:15 pounds, 7 in. (H) x 14 in (W) x 16.5 in (D) Cart: 45 pounds, 301/4 in (H) x 19 in (W) x 19 in (D)

Auto\*Therm® 390 packages: 13-3062 ME390: 2 4.7x7.1" soft rubber plate applicators with flexible cables, 6 felts spacer, 2 cloth covers, and 2-48" elastic velcro straps

ME391: standard ME390 PLUS instrument cart, 13-3063 four-joint folding arm, 1-14cm diameter coil applicator, and 1 special field coil applicator

### cart

13-3088 ME97 mobile cabinet with plastic doors and shelves 13-3089 ME96 cart for use with arm and coil applicator

### plate applicator options

13-3070 12x18cm rubber plate applicator with cables 13-3071 14x20cm felt spacer for soft rubber plate applicators 13-3072 18x26cm cloth cover for soft rubber plate applicators

### coil applicator options

13-3073 14 cm diameter coil applicator without cable 13-3074 special coaxial cable for (13-3073) coil applicator 13-3075 four-joint arm

13-3077 strap for securing 14cm coil applicator (13-3073) to patient

13-3076 lanyard for patient safety switch

# Mettler Auto\*Therm® 395 diathermy

The Mettler® Auto\*Therm® 395 features continuous and pulsed shortwave diathermy. Automatic tuning unit has eight different treatment applicators: three inductive coils, three spaced condensers and two soft, flexible plate electrodes. Multi-jointed arms lock in place for ease of use and safety. Membrane panel and intensity control knob make treatment set-ups quick and easy. Weight: 99 pounds. Dimensions: 33.5 in. (H) x 15 in. (W) x 15.4 in. (ID).

# 395 package

13-3061 Auto\*Therm 395 shortwave diathermy (includes 2 arms, 1-14cm coil with cable, and 2-13cm 200W condensers with cables)

13cm condenser, 200W

condenser cable

999.99

# 395 applicators / accessories

13-3252	soft plate applicator set #1 includes 2-8x14cm soft-rubber applicators w/ flexible cables, 2 6 felt spacers, 2 cloth covers; 2 perforated rubber bands with	
13-3253	soft plate applicator set #2 includes 2-12x18cm soft-rubber applicators w/ flexible cables, 6 felt spacers, 2 cloth covers; 2 perforated rubber bands with	
13-3254 13-3255 13-3256	8cm coil, w/o cable 30W 14cm coil, w/o cable 90W 15x35cm tripole coil w/ cable, 200W	
13-3260	inductor coil cable	- J.3
13-3257	4.2cm condenser, 21W	
13-3258	8.5cm condenser, 80W	







13-3259

13-3261

# Fluidotherapy®

- increase patient's circulation and ROM with the soothing dry heat of the Fluidotherapy® unit...without having to deal with the messy cleanup and wetness associated with other therapies
- programmable; standard with automatic timer, wake-up preheater, and locking casters
- adjustable controls include: temperature (maximum 125°), airflow (velocity of particulate medium), pulse mode (pressure bursts)



single FT110

- single extremity unit for hand, wrist, elbow, ankle or foot
- includes 20 pounds of Cellex medium

FLU110D

single FT110



double ULTRA 115

- treats two extremities
- removable from stand for easier lower extremity applications
- includes 40 pounds of Cellex medium

FLU115D double ULTRA 115



cellex medium

- 10 pounds of medium available in two 5 pound containers
- can be used with both Fluidotherapy<sup>®</sup> and Thermotherapy® systems

MFD0001

Cellex-10 lbs.

# thermo-therapy

- granulated medium is kept in suspension by heated air
- thermostatically controlled heater can quickly heat medium to to 125° F; temperature indicated on the unit
- built in electric timer can control treatment time
- variable stream blowers 1720RPM) control heated streams
- · all units come standard with granulated medium and instructions



# single hand TT101

- 10-12 lb capacity
- 15"H x 19"D x 11"W (shown with optional stand)

11-1900 11-1920

single TT101 [a] TT101 stand

# • 12 lb capacity

15"H x 22"D x 11"W (shown with optional stand)

single hand TT101L

11-1901

single TT101L [b] 11-1920 TT101L stand

# dual hand TT202

- 40-50 lb capacity
- 39"H x 29"D x 20"W

11-1911 double TT202 [c]



# self-contained iontophoresis systems

# **IontoPatch**





IontoPatch SP™

- self-contained wireless technology eliminates need for external batteries or wires
- allows patient to return to daily activities while receiving time-released iontophoresis
- · made of flexible, breathable material
- · technician applies medication directly to the electrode
- shuts off automatically when prescribed dose has been administered
- · STAT offers the same benefits as IontoPatch with shorter wear time

	<u>IontoPatch 80</u> ™	<u>IontoPatch SP</u> ™	<u>IontoPatchSTAT</u> ™
active area	15.5 cm <sup>2</sup>	8.0 cm <sup>2</sup>	12.0 cm <sup>2</sup>
dosage	80 mA-min	40 mA-min	80 mA-min
fill volume	1.2-1.6 cc	1.2 - 1.6 cc	1.2 - 1.4 cc
wear time	24 hours	24 hours	4 hours

13-5220	IontoPatch 80™	6 treatment kits per box (patch, saline
13-5221	IontoPatch SP™	vial, alcohol prep)
13-5223	IontoPatch Stat™	viai, alconoi prepj

# **lomed Companion 80**



- reserve battery capacity: compensates for patients with higher skin resistance
- hypoallergenic adhesive: adheres well even in the shower
  - high tech design conforms well to treatment areas and fits comfortably beneath clothing

Companion oc	Com Com	<u>panion</u>	80
--------------	------------	---------------	----

active area dosage fill volume wear time 6.45 cm<sup>2</sup> 80 mA-min 1.1cc 24 hours


6 treatment kits per box (patch, saline vial, alcohol prep) 50-0004-4 Companion 80

# single use iontophoresis electrodes and controllers

# **IOGEL** electrodes

- clear adhesive allows precise placement on target site
- silver-silver chloride conducts more efficiently than carbon and stabilizes pH without using chemic





- GelSponge<sup>®</sup> ensures complete contact between the electrode and the skin and provides a conductive medium similar to ultrasound gel
- EZ Fill Window™ protects adhesive during hydration and handling and provides no-stick handling

12 treatment kits per box\*

50-0002-1 small 1.5cc 50-0002-2 medium 2.5cc 50-0002-3 large 3.5cc

# OptimA™ electrodes

- pH defining agent adjusts pH to optimize transport when using a carbon conductive element
- pH buffering agent maintains ideal pH
- viscosity enhancing element creates uniform current distribution for greater patient comfort
- adhesive foam stretches and adheres well when under stress

12 treatment kits per box\*

50-0005-0 small 1.5cc 50-0005-1 medium 2.0cc 50-0005-2 large 3.0cc

# **TransQ E electrodes**

- 2 shapes: rectangular and circular for various treatment sites
- adhesive foam stretches and adheres well under stress

12 treatment kits per box\*

50-0000-9 TransQE small 1.5 - 2.0cc 50-0000-8 TransQE medium 2.5 - 3.0cc

# TransQ Flex electrode

- clover shape: surrounds highly contoured treatment area
- extra adhesive: adds more surface area to increase adhesion when under stress



12 treatment kits per box\*

50-0003-3 TransQ Flex 2.5cc

# Ionto+ Hi-Performance electrodes

- new EZ-fill window
- greater adhesive area eliminates pad movement
- silver-silver chloride conductor (superior to carbon) for higher conduction (more drug delivery)
- buffering agents for additional pH stabilization

12 treatment kits per box\*

13-5250 wings 13-5251 sm. square 13-5252 med. square 13-5253 large square



# Chattanooga lonto controller

- dual channel: treat 2 sites or deliver 2 medications simultaneously
- constant current provides precise delivery regardless of skin impedance
- current can be set in 0.1 mA increments between 0.5 mA and 4.0 mA



- automatic 30 second current ramp up and down during power on/off for patient comfort
- beeping sounds if electrode faults for open current or high impedance, low battery, treatment completion or power left on
- previous dosage can be set as default to maintain frequently used treatment parameters

00-1335 Ionto controller 01-3214 lead wire replacements

# id<sup>3</sup> digital controller

- single channel
- 40 mA adjustable;
   1 minute settings
- 1,2,3,4 mA adjustable settings
- three 9 volt batteries included
- automatic complete circuit setting which measures the amount of current going through the dispersive pad
- synthetic buffering provides consistent pH balancing to prevent the effects of hydrolysis

13-5230 id<sup>3</sup> controller



<sup>\*</sup> treatment kit includes ionto electrode, dispersive electrode and alcohol prep

# electrodes and gel

# economy electrodes

- self-adhesive
- · economic and functional
- made with high quality Tyco gel
- made with white foam
- conveniently packed in bags of 4, 10 bags per case (40/case)

13-1287 2" round 13-1288 3" round 13-1291 2" square 2" x 3.5" rectangle 13-1293

# economy conductive spray

- conductive electrolyte spray
- · for electrotheray stimulation devices to increase conductivity

13-1294 2 oz. spray 13-1295 8 oz. spray



# **Dura-Stick® electrodes**

- electrodes and packaging with patient usage tracking
- · re-sealable packaging with fields for patient name, first-use, and number of uses
- electrodes sold in 40/case (10 packages of 4)

# Dura-Stick® premium

- stainless steel mesh for enhanced con-
- blue gel series for special needs skin

04-2171 2" round 04-2174 2" square 2" x 3.5" rectangle 04-2175 04-2172 1.5" x 2.5" oval 04-2173 3" x 5" oval 04-2176 2" square blue

1.5" x 3.5" rect. blue 04-2177

# Dura-Stick® plus

- blue foam backing
- · tabs for ease of removal

04-2182 2" round 2" square 04-2183 04-2181 2" x 3.5" rectangle 04-2159 1.5" x 2.5" oval 04-2160 2" x 4" oval

# Dura-Stick® standard

- · white cloth backing
- carbon film conductive material

04-2178 2" square 04-2179 2" x 3.5" rectangle



premium

TENS/NMES/

FES electrodes



TENS/NMES/

FES electrodes

# Dura-Stick® single use

- · white cloth backing
- · carbon film conductive material
- single use disposable
- ideal for TENS applications

04-2180 2.25"x2.5" disposable

# Intelect® ultrasound gel



- · ultrasound and electronic coupling agent
- · hypo-allergenic, water-soluble, greaseless and non-staining gel packed in 5 liter container
- with reusable 8 ounce bottle

00-4266-1 5 liter (1.3 gal) 00-4266-4 5 liter 4/case

# Conductor™ Gel Highly conductive gel

00-4248-1 8½ oz 00-4248-12 81/2 oz 12/case 00-4248-24 81/2 oz 24/case 00-4238-1 5 liter (1.3 gal) 00-4238-4 5 liter 4/case





- · for massage or as an ultrasound coupling agent
- menthol, hexachlorophene, or polysorbates
- non-allergenic, water soluble, greaseless
- with aloe vera

00-4262-1 gallon 00-4262-4 gallon 4/case 00-4260-1 8½ oz 00-4260-12 8½ oz 12/case 00-4260-72 81/2 oz 72/case

# Myossage® lotion



- for soothing massage
- non-alcoholic, non-greasy, non-staining
- water-based

00-4210-1 gallon 00-4210-4 gallon 4/case 00-4208-1 8½ oz 00-4208-12 81/2 oz 12/case 00-4208-72 81/2 oz 72/case



# Chattanooga® Intelect® **LEGEND XT**™ Therapy Systems

The first electronic "platform" capable of providing multiple modalities in one therapy

station. You can start simple (stim and ultrasound) and add modalities (additional stim, laser, battery module and patient management software) when needed.

- · stim unit offers 6 clinical waveforms including: interferential, Russian, bi-phasic, high volt, microcurrent, and premodulated
- full functioning ultrasound delivers 1 and 3.3 MHz through a 5cm<sup>2</sup> applicator. Multiple waveforms available in combination mode
- built-in patient data card device can record patient modality and pain profile information
- Intelect Legend XT Therapy System features user-defined protocols to customize any waveform to your specific need.

Each unit ships with 5 Patient Data Cards, 4 each 3" inch diameter carbon electrodes (2 red, 2 black), 4 each DuraStick II 23/4 inch diameter electrodes and 1 21/2 x 24 inch Nylatex 

• optional mobile cart has 6 storage bins to minimize cord and cable hassle

00-2780 mobile cart



units (stim left, combo right) shown above with optional cart.

Sum umi	5	
00-2763	2-channel stim	<b>~</b>
00-2793K	2-chan. stim w/ cart	<b>~</b>
00-2786K	4-channel stim	<b>~</b>
00-2797K	4-chan. stim w/ cart	<b>2</b>



# combo units (stim plus ultrasound)

 standard with 5cm<sup>2</sup> (1 and 3.33MHz) applicator, 9 oz. bottle conductor gel

00-2760K	2-channel combo	2
00-2791K	2-chan. combo w/ cart	2
00-2788K	4-channel combo	<b>**</b>
00-2795K	4-chan. combo w/ cart	<b>~</b>

# modular components

• stim/ultrasound combo device is enhanced with the addition of one or more of the following therapy modules: 2-channel stim module, laser therapy module



2-channel stim module

• upgrades base unit from 2 to 4 channels of stimulation

00-2781 additional stim



# laser therapy module

- independent control over parameters
- includes protective eyewear (2x)

00-2758 laser module





carry bag

- ergonomically designed bag for easy transport of Legend XT model
- holds unit and accessories securely 02-7133 Legend XT carry bag



# battery module

• automatically charges when the base unit is plugged in

00-2767 battery module

# lead wires for Legend XT

02-7312	lead wires for channel 1
02-7313	lead wires for channel 2
02-7314	lead wires for channel 3
02-7315	lead wires for channel 4

# additional applicators for Legend XT and Transport

02-7381	1cm <sup>2</sup>
02-7382	2cm <sup>2</sup>
02-7383	5cm <sup>2</sup>
02-7384	10cm <sup>2</sup>



# patient interrupt switch

02-7469



 patient data cards allow modality and pain profile information to be recorded from the unit and



reviewed on the PC, or filed with the patient's treatment records for use in future treatment sessions

02-7465 25 data cards

# Chattanooga® Intelect® TRANSPORT™ electrotherapy

# Transport™ 2-channel **stimulation**

- maximum portability with easy to use operation
- 2-channel output with four standard waveforms: interferential, premodulated, HV, and Russian
- logical and simple user interface with 10 programmable memory positions
- built-in positional stand allows desktop or wall-mount use
- standard with 4 ¾ inch diameter Dura Stick II electrodes. 4 3" diameter carbon electrodes (2 red, 2 black), lead wires and 1 2½" x 24" Nylatex wraps

00-2783K2-channel stimulation00-2894Kstim with cart00-2783Bstim with bag and battery



Transport™ stim unit

# Transport™ 2-channel stimulation / ultrasound combination

- use in combination, electrotherapy and ultrasound work seamlessly to heat and stimulate muscles for more rapid healing and better patient outcome
- features Transport 2-channel stimulation and ultrasound

00-2738K Intelect Transport combo: 2-channel stim and ultrasound w/ 5cm² head 00-2738KC Intelect Transport combo with mobile cart 00-2738KB Intelect Transport combo with with bag and battery pack



Transport™ combo unit

# Transport™ ultrasound

- transportable ultrasound delivers 1 and 3.3 MHz frequency in continuous (100%) or pulsed (10%, 20%, 50%) format
- includes 10 programmable memory positions
- powered by battery or 110 volt AC
- output is a clean, culminated beam with a BNR that is of the lowest available (max 5:1)
- built-in positional stand facilitates desktop or wall-mount use
- standard with 1-9 oz. bottle of ultrasound gel and a 5cm<sup>2</sup> ultrasound applicator

00-2782K ultrasound w/ 5cm² app. 00-2893K ultrasound with cart

**~** 

Transport™ ultrasound unit

### ■MAP (minimum advertised price applies to this product)

# Transport battery pack



 self contained nickel metal hydride battery stores conveniently in the bottom of Intelect Transport™ units

02-7478 battery 00-2911K battery & bag

# Transport carry bag



- ergonomically designed bag for easy transport of *Trans*port<sup>™</sup> models
- holds unit securely in place, and fits accessories

02-7467 carrying bag

# additional soundheads



 additional applicators for Legend XT™ and Transport

02-7381 1cm<sup>2</sup> 02-7382 2cm<sup>2</sup> 02-7383 5cm<sup>2</sup> 02-7384 10cm<sup>2</sup>

# mobile cart





 optional mobile cart and adapter for Intelect<sup>™</sup> Transport<sup>™</sup> unit

Transport™ cart
00-2891K cart+adapter
Transport™ adapter for XT™
00-2884 adapter



# Chattanooga® Intelect® ultrasound

# Chattanooga® Intelect® ultrasound / stim combination



# Intelect® Ultrasound

- dual frequency clinical ultrasound system
- 1 and 3.3 MHz frequencies, 10%, 20%, 50% and continuous duty cycle selections with head warming; bright back-lit LCD display
- includes soundhead and a 9 oz. bottle of Conductor<sup>™</sup> gel

INT001K-5 ultrasound w/ 5cm² head INT001K-2 ultrasound w/ 2cm² head INT001K-10 ultrasound w/ 10cm² head

additional soundheads (for use with Intelect® ultrasound and combo units)

07-8048 2cm² sound head only 07-8047 5cm² sound head only 07-8046 10cm² sound head only



# Intelect® 4-channel Combo

- dual frequency four-channel can be used independently or in combination with Pre-modulated, Interferential, Russian and High Volt waveforms
- 1 and 3.3 MHz frequencies, 10%, 20%, 50% and continuous duty cycle selections with head warming; brightlit back LCD display
- includes soundhead and a 9 oz. bottle of Conductor™ gel

00-7550K-5 w/ 5cm² head 00-7550K-2 w/ 2cm² head 00-7550K-10 w/ 10cm² head



# Intelect® 2-channel Combo

- dual frequency 2-channel can be used independently or in combination with Pre-modulated, Interferential, Russian and High Volt waveforms
- 1 and 3.3 MHz frequencies, 10%, 20%, 50% and continuous duty cycle selections with head warming; bright-lit back LCD display
- includes soundhead and a 9 oz. bottle of Conductor™ gel

INTOCBK-5 w/ 5cm² head INTOCBK-2 w/ 2cm² head INTOCBK-10 w/ 10cm² head

# Chattanooga® Intelect® stim



### Intelect® 2-channel and 4-channel stimulator

- two or four channel electrotherapy system
- four waveforms: Interferential, Premodulated, High Volt and Russian; HV available on all channels
- bright back-lit LCD display
- includes Nylatex® wraps, both red and black 3" round electrodes and lead wires

INT002 00-7560 INT02, 2 channel stimulator 4 channel stimulator

# Chattanooga® D-HVP





# Intelect® 1-channel digital portable D-HVP

- one-touch digital programming access, custom touch-pad interface, pause feature and 29 custom or standard presets
- portable stimulator for multiple patient conditions which include maintaining or increasing ROM, relaxing muscles in spasm, preventing or retarding muscle disuse atrophy, muscle re-education and increasing local blood circulation

00-7570 D-HVP 1-channel

# MB-TD MB-TD

# Chattanooga® storage

# stainless steel mobile stands

- 2 shelf design 191/2" high x 183/4" wide x 163/4" deep
- 3 shelf design 29½" high x 18¾" wide x 16¾" deep

00-4002 2 shelf MB

00-4004 3 shelf MB-T

00-4018 3 shelf w/drawer MB-TD



# Amrex® ultrasound / stim combos

# US/50 with LV



- ultrasound features 20 watts of power with pulsed or continuous output modes
- single channel LV with surge, tetanize and pulsation

13-3144 w/10cm head 13-3156 w/5cm head 13-3157 w/both heads

# US/54 with LV



- utrasound features 20 watts of power with pulsed or continuous output modes
- dual channel LV with surge, tetanize and pulsation

13-3146 w/10cm head 13-3158 w/5cm head 13-3159 w/both heads

# US/752 with HV



- · ultrasound has 20 watts of power with pulsed or continuous output
- dual channel HV with var. freq.1-160 pps, reciprocal & ramped surge (2-90 seconds), patient interrupt switch and adjustable "on/off"

13-3153 w/10cm head 13-3154 w/5cm head 13-3155 w/both heads

# US/HVG50 w/ LV&HV



- ultrasound has 20 watts of power with pulsed or continuous output
- dual channel HV with reciprocal and ramped surge
- LV galvanic with variable surge and reverse polarity

13-3145 w/10cm head 13-3160 w/5cm head 13-3161 w/both heads

# Amrex<sup>®</sup> stim units

# **MS/322 AC LV**



- single channel
- LV AC current with surge, tetanize and pulse modes

13-3120 1 MS322

# LVG/325A DC LV



- single channel
- LV DC galvanic current
- iontophoresis/phonophoresis
- adj. surge and alternating polarity surge rates
- · normal or reverse polarity

13-3121 LVG325A

# MS/324 LV



- dual channel
- LV AC current with surge, tetanize, reciprocal and pulse modes

13-3123 MS324A

• (AB) has battery for use when AC not used

13-3124 MS324AB

# **MS/324C LV**



- dual channel
- LV AC current with surge, tetanize, reciprocal and pulse modes
- 110/220V 50/60 Hz power system

13-3125 MS324C

Amrex<sup>®</sup> ultrasound units

# **HV/752 DC HV**

# MS/401B AC LV



- dual channel
- with independent meters
- HV DC current with surge, tetanize, reciprocal and pulse modes
- reciprocal and surge modes ramped
- bi-polar application

13-3133 HV/752





- quad channel
- bi/mono phasic stimulation
- · LV AC current with surge, and reciprocal modes

13-3126 MS401B

# U/20 ultrasound unit



- · continuous output mode
- quick disconnect head
- 12 watt 1.0 MHz

13-3103 w/10cm head 13-3106 w/5cm head 13-3107 w/both heads

# U/50 ultrasound unit'



- continuous and pulsed output modes
- can be combined with Stim unit\*
- quick disconnect head
- 20 watt 1.0 MHz

13-3100 w/10cm head 13-3101 w/5cm head 13-3105 w/both heads

\*use with any Amrex low-volt or high-volt stim for combination therapy



# Sonicator® Plus 992 Combo



# Sonicator® Plus 994 Combo



# ultrasound + stim

# Sonicator® Plus 930 Combo



each Sonicator Plus comes standard with a 5cm<sup>2</sup> dual frequency 1 and 3 MHz. applicator

- microprocessor-controlled 2-ch. combo
- six treatment waveforms: Interferential, Premodulated, Medium Frequency, BiPhasic, HV & Microcurrent.
- 3 different protocols can be set up simultaneously using different timers

13-3001 Sonicator 992

- microprocessor-controlled 4-ch. combo
- six treatment waveforms: Interferential, Premodulated, Medium Frequency, BiPhasic, HV & Microcurrent
- 5 different protocols can be set up simultaneously using different timers

13-3000 Sonicator 994

- therapeutic ultrasound and three forms of waveforms: Interferential, Premodulated and Medium Frequency. Membrane panel
- select treatment parameters and adjust intensity by twisting dial

13-3014 Sonicator 930

# Sonicator® 940 Plus combination therapy

- each Sonicator Plus comes standard with a 5 cm<sup>2</sup> applicator at 1 and 3 MHz
- a 10 cm² / 1MHz applicator or a 1 cm² / 3.3 MHz applicator may be ordered to expand the therapeutic ultrasound capabilities
- microprocessorcontrolled Sonicator Plus 940 is a 4 channel combination unit





# • combines 1 and 3 MHz ultrasound with 8 treatment waveforms: interferential, premodulated, medium frequency (Russian), EMS, HV, TENS, microcurrent and direct current

- all waveforms are available on all 4 channels.
- a dual-frequency 5.5 cm², 1/3 MHz applicator comes standard with the unit
- membrane panel and back-lit LCD touch panel for quick setup of treatment parameters
- 61 preset programs and space for 80 user-defined programs.
- an optional 0.9 cm² dual-frequency applicator is available

13-3019 940 Plus Ultrasound

# Sonicator® 740, 740x ultrasound



- 2-freq. (1, 3.3 MHz) delivered in continuous or 10, 20 or 50% pulsed modes
- applicator is contoured for easy application, has a LED in the handle that lights up when adequate coupling achieved
- 740 has 5 cm<sup>2</sup> applicator, 740x has 3 applicators (1, 5, 10 cm<sup>2</sup>)

13-3300 Sonicator 740 13-3301 Sonicator 740x 13-3320 battery (optional)

# Mettler® soundhead applicators



 soundheads fit specified Mettler ultrasound unit(s), can be interchanged without calibration

area/frequency			fits Sonicator® model					
		cm <sup>2</sup> /MHz	740	930	940	992	994	
13-3009	ME7310	10/1				х	Х	
13-3012	ME7513	5/1&3.3		Х		Х	Х	
13-3010	ME7331	1/3.3 p*				Х	Х	
13-3310	ME7431	1/3.3 p*	Х					
13-3311	ME7413	5/1&3.3	Х					
13-3312	ME7410	10/1	Х					
13-3313	ME9402	10/1&3.3			Х			
13-3314	ME9401	5/1&3.33			Х			
	P* pencil s	tyle applicator.						



# neuromuscular stimulators



# Sys\*Stim® 208 stimulator



- one-channel stimulator with Surge, Pulse, Tetanize modes
- · easy-to-use knobs that control stimulation intensity and frequency, a digital timer and a remote patient cut-off switch

13-3015 Sys\*Stim 208

# Sys\*Stim® 208A stimulator



- two-channel stimulator with Surge, Pulse, Tetanize Reciprocation modes
- · easy-to-use knobs controls stimulation intensity and frequency, digital timer and remote patient cut-off switch

13-3016 Sys\*Stim 208

# Sys\*Stim® 294 stimulator



- four-channel neuromuscular stimulator with six stimulation waveforms: interferential, pre-modulated, medium, biphasic, HV and microcurrent
- 4 different treatment protocols may be run simultaneously

13-3020 Sys\*Stim 294

# Sys\*Stim® 206 stimulator



- one-channel stimulator with 4 wave-
- narrow pulse is suited for TENS while wide and AC pulses can be used for muscle stimulation
- · direct current can be used to stimulate denervated muscles or to perform iontophoresis

13-3021 Sys\*Stim 206

# Sys\*Stim® 226 stimulator



- Interferential, Pre-modulation, Medium Frequency, and Biphasic with adjustable phase duration and frequency waveforms
- · surge, reciprocation, surge with reciprocation and vector rotation modes
- two timer mode allows for 2 simultaneous separate treatment protocols

13-3032 Sys\*Stim 226

# Trio\*Stim® 215 stimulator



- · battery operated neuromuscular electrical stimulator with three modes: TENS, muscle stimulation and microcurrent
- 2-channels allows user to provide a device for pain control with TENS or Microcurrent stimulation, or use for muscle rehabilitation, re-education, spasm reduction

13-1315 Trio\*Stim 215

Trio\*Stim AC adapter 13-1316

# plastic therapy cart

plastic 3 shelf cart has casters and easy to grab handles on each side

34"x24"x18"





# metal therapy cart

 metal 3 shelf cart has casters and easy to grab handles on each side

34"x24"x18"

13-3087 ME74 13-3086 ME75 with 3 hospital outlets

rugged



- padded tote complete with shoulder strap
- fits unit and accessories



13-3064 for Sonicator 740 (shown) 13-3091 for Sys-Stim/Sonicator

13-3065 for 992,994,294

# Z-Stim IF150 interferential Z-Stim IF250 interferential

# Quattro IF interferential

# **Quattro EMS stimulator**



- dual channel
- 6 pre-modulated frequency ranges and 6 continuous frequency selections

Amrex® MC-1

13-3180 IF150



- dual channel
- Quad-Polar, Bi-Polar and Russian modalities
- · frequency, time and milliamp digital displays

13-3181 IF250

# Intelect® IFC interferential



- single channel
- Spectrum MC1 micro-current stimulator provides easy to use treatment protocol
- intensity (20-1000 mA), frequency (0.5 - 620MHz) and polarity are variable
- used with pads or probe (both included)
- uses AC power or battery pack

galvanic 1-channel

portable

battery operated dual chan-

 for 2 or 4 electrode use; adj. pulse amplitude, interference frequency and difference frequency parameters; 4 stim modes: continuous, 1/1, 8/8, and 10/10 ramped; 120V AC adapter

07-7719 IFC

nel IFC stim

13-3171 MC1

# digital HV pulsed




output mode: synchronous or alternation

alternation mode-LED alternation between green and red color

13-1304 galvanic 1-chan

power: 9V battery or AC adapter



galvanic stim, 1-chan.

- ous, alternating
- · pulse: variable rate, fixed duration
- system includes: stimulator, lead wires, 4 self stick electrodes, 9V battery and AC adapter

13-1305 digital galvanic



- 4-channel, 4-mode (constant, adjustable, auto sweep (15 seconds), frequency shift (30% below to 60% above set frequency) interferential unit
- symmetrical balanced sine wave w/ adj. pulse amplitude (0-70mA), rate (4-160bps), duration (to 125µs); 4000 MHz; modulating frequency from 4004 to 4160 MHz
- AC powered

13-3360 Quattro IF

# clinical EMS unit features Asymmetric Bi-Phasic square pulse wave with

4-channel, 3-mode (constant,

synchronous)

- adjustable pulse amplitude (0-100mA) and rate (1-160MHz) with constant width (300μs); adjustable ramp, and on/off cycle time
- with 9 volt adapter

alternation.

13-3350 Quattro FMS

# Intelect® NMES



 battery operated portable dual channel NMES stimulation unit

• includes adjustable pulse amplitude, pulse frequency and contraction/relation time parameters

07-7715 NMES

# Intelect® digital NMES



- battery operated 2-chan digital NMES stim
- includes adi. pulse amplitude, width,and frequency and contraction/relaxation time; 3

stimulation modes: synchronous, constant and alternate; east to read LCD display and convenient digital controls

07-7717 digital NMES

# IF 2-channel w/ AC adapter



- 4 Modes: continuous, 1/1 abrupt, 8/8 abrupt, 10/10 ramped
- · frequency shift: from 25% below set frequency to 55% above
- carrier frequency: 4000Hz
- modulating frequency: 1-150Hz
- wave form: symmetrical biphasic square
- stimulator, lead wires, 4 self stick disposable electrodes, 9V battery, AC adapter

13-1303 interferential

# digital interferential 2-channel



- 9 Modes
- shift: 30% below to 60% above
- carrier frequency: 4000Hz
- •modulating freq.: 4004-4160Hz
- symmetrical

balanced sine wave

- adjustable pulse rate
- timer and compliance meter
- stimulator, lead wires, 4 self stick electrodes, 9V battery

13-1302 digital IF



# portable TENS, IFC, NMES, IF, GALVANIC and HV

# **EMS 2-channel**



- variable pulse: rate 5, 30, 100(Hz), width 250(µS)
- variable ramp time: 1,3 or 5 S
- variable contraction/relaxation time
- wave: 3-types
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery, instructions

13-1306 EMS

# **EMS 2-channel**



- 3 Modes
- adjustable pulse rate and width
- variable ramp time: 1, 3 or 5 S
- variable contraction/relaxation time; timer and compliance meter
- wave: 3-types
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery, instructions

13-1309 digital EMS

# Intelect® economy TENS



- dual channel modified square waveform with adjustable amplitude (0-80mA), pulse frequency (2-150 Hz) and pulse width (60-25 uS) and timer
- includes 9-volt battery, Dura-Stick II self adhesive electrodes, carrying case

07-7710 TENS economy

### accessories

07-7610 battery recharger 07-7609 rechargeable NiCd 9V battery

# Intelect® standard TENS



- dual channel with adjustable parameters including frequency, pulse duration and 3 modulation modes and timer
- includes 9-volt battery, Dura-Stick II self adhesive electrodes, carrying case

07-7600 TENS standard

### accessories

07-7610 battery recharger 07-7609 Rechargeable NiCd 9V battery

# Mettler® TENS Sys\*Stim



- two independently controlled channels deliver burst or modulated Symmetric Bi-Phasic square pulse wave output
- intensity (1-80mA, 0-100 volt), duration (30-260µs) and rate (2-150MHz) of the output is controlled
- 210T has timer
- TENS 212 has digital display of treatment parameters

13-1317 Sys\*Stim 210 13-1318 Sys\*Stim 210T 13-1314 Sys\*Stim 212

# Intelect® digital TENS



- battery operated portable dual channel digital TENS stimulation unit
- adjustable frequency and pulse duration and timer; 5 modes: burst, normal, modulated rate, and width, strength duration and B-Pulse; and easy to read LCD display and convenient digital controls

07-7712 digital TENS

# NeuroTech® Aviva TENS XP



- clinically proven, safe and effective means of providing a drug-free method of pain relief
- portable so patients can use it at home or on-the-go

13-1366 AvivaTENS XP 13-1368 2 electrode pkg 13-1369 3 electrode pkg

# NeuroTech® Aviva Stim XP



- effective muscle stimulation device that re-educates and strengthens atrophied, weakened or immobilized muscles
- portable so patients can use it at home or on-the-go

13-1367 AvivaStim XP 13-1368 2 electrode pkg 13-1369 3 electrode pkg



# portable TENS/EMS (Transcutaneous Electrical Nerve Stimulation)

# **TENS 2-channel**



- · adjustable pulse rate and width
- includes carry case, stimulator, lead wires,4 self stick electrodes, 9V battery

1 Mode (constant), no timer (A) 13-1313 TENS - 1 mode

3 Modes (Burst, Modulation, Constant), timer (B)

13-1323 TENS - 3 modes

# digital EMS/TENS 4-channel



- 3 TENS modes, 3 EMS modes:
- variable pulse rate and width
- variable ramp time
- variable contraction/relaxation time
- wave: Symmetrical Bi-phasic square pulse
- 2 timers and compliance meter
- comes with 8 electrodes

13-1349 digital EMS/TENS

# **TENS** with heat

- freq: 1-180 Hz; pulse width: 100-300; intensity: 100 amp; wave: alternating
- treatment modes: shoulder, waist, elbow, foot, pound, push, press, knead, auto
- heating: 96-108° F
- timer: 10 / 20 / 30 minutes; unit will shut off if skin temperature reads 114° F
- adaptor: DC 9V/1.2A; battery: 1.5 V x 3 (AAA) (load 500)

13-1365 TENS with heat

# digital TENS 2-channel



- 5 Modes
- · adjustable pulse rate and width
- wave: Asymmetrical Bi-Phasic Square Pulse
- timer (5~60 mins or continuous) and compliance meter
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery

13-1308 TENS - 5 modes

# microcurrent TENS 2-channel




- dual channel microcurrent / micro TENS
- 3 Frequencies: .3Hz, 8Hz, 80Hz
- · changes polarity each second
- wave: Modified square Biphasic
- timer: constant, 20, 40 min.
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery, instructions

13-1322 Microcurrent TENS

# **EV-803P digital TENS**

- 9 preset programs
- channel:dual, isolated between. Pulse amplitude: 0-80mA peak to peak. Pulse

0



- mode: B,N,M and P (from 1-9)
- voltage: 0-40V peak; timer: 30, 60 min. or continuous; compliance meter

13-1361 EV-803P

# **EMS/TENS 2-channel**



- 3 TENS modes, 2 EMS modes:
- variable pulse rate and width
- variable ramp time
- variable contraction/relaxation time
- wave: Asymmetrical Bi-phasic square pulse
- timer and compliance meter
- includes carry case, stimulator, lead wires, 4 self stick electrodes, 9V battery

13-1348 EMS/TENS



# wireless TENS/EMS



- single channel output, 3V CR2025 battery (included), output 0 ~ 90 mA
- density: change of 16 steps
- timer: 30 min. backward to auto turn off
- wave form: dual phase square wave
- LCD display w/ timer/mode/density
- mode: part A P1~P5 are 5 TENS programs; part B P6~P0 are 5 EMS fitness program; contraction (on time) 2s. or 4s.; relaxation (off time) 2s. or 4s

13-1360 mini unit

# BakWel TENS w/ wrap





- temporary relief of pain associated with sore/aching muscles in the lower back due to strain from exercise or normal activities
- includes battery, belt, electrodes and cables

13-1362 Bakwel TENS



# Chattanooga® Vectra® GENiSYS Laser



- · easy-to-read graphical LCD display
- Continuous and pulsed treatment modes
- · real time feedback of dosage delivered
- useful selection of clinical indications & user defined protocols
- choice of 8 interchangeable laser diode and cluster applicators w/ power ranging from 100-1440mW
- electronic signature recognition allows each laser applicator to be used in a "plug-and-use" format
- laser therapy can be controlled from the laser applicator
- lightweight design, with battery-powered option and custom designed mobile carrying bag, make transport easy
- Unit designed for table top, wall mount or mobile use
- laser protective eyewear included (2x)

### Vectra® GeniSys Laser

00-2784 GeniSys Transport Laser (w/out applicator) 00-2896 GeniSys Laser w/ cart and adapter (w/out applicator)

## Single Laser Diode Applicator

02-7840 850 nm Laser Diode 100 mW

02-7841 850 nm Laser Diode 200 mW 02-7805 820 nm Laser Diode 300 mW

9 Diode Cluster Applicators (Laser and LED) 540 mW

02-7811 (5) 850 nm Laser Diodes @100 mW each (4) 670 nm LED Diodes @10 mW each

9 Diode Cluster Applicators (Laser and LED) 1040 mW

02-7812 (5) 850 nm Laser Diodes @200 mW each

(4) 670 nm LED Diodes @200 mW each

13 Diode Cluster Applicators (Laser, SLD, LED) 415 mW

02-7814 (3) 850 nm Laser Diodes @100 mW each

(3) 950 nm SLD Diodes @15 mW each

(7) 670 nm LED Diodes @10 mW each

# laser and applicator sold separately





# 13 Diode Cluster Applicators (Laser, SLD, LED) 715 mW

02-7816 (3) 850 nm Laser Diodes @200 mW each

(3) 950 nm SLD Diodes @15 mW each

(7) 670 nm LED Diodes @10 mW each

### 33 Diode Cluster Applicators (Laser, SLD, LED) 1440 mW

02-7808 (5) 850 nm Laser Diodes @200 mW each

(8) 880 nm SLD Diodes @25 mW each

(8) 950 nm SLD Diodes @15 mW each

(12) 670 nm LED Diodes @10 mW each

### accessories

02-7525 Laser Protective Eyewear

02-7478 Transportable Battery Pack

02-7133 Vectra® GeniSys Therapy System Carry Bag

02-7467 Vectra® GeniSys Transportable Carry Bag

00-2780 Vectra® GeniSys Cart 00-2883 Cart Adapter (Blue)

00-2758 Laser module

# METTLER ≤LECTRONICS® Laser Sys Stim 540









left to right: Sys\*Stim 540; laser applicator over laser output sensor showing visible proof of laser activity; laser applicator showing one laser diode at 785 nm surround by 3 LEDs at 470nm; cluster applicator showing 12 SLDs at 950 nm with seven LEDs at 660 nm

- 785 nm laser probe emits visible blue LEDs that illuminate the treatment field
- allows the clinician to set treatment time, continuous/pulse modulation modes and laser output energy
- treatment parameters are displayed on the LED displays and controlled via the membrane keypad
- store applicators in integrated holders when not in use
- output controls on both the laser and cluster applicators

Weight: 5.5 pounds, (2.5 kg)

Dimensions: 4.7" H x 9.4" D x 9.4" W 12cm H x 24cm D x 24cm W

# Laser Sys\*Stim 540

13-3200 Sys\*Stim 540L w/ 785 nm applicator

13-3201 Sys\*Stim 540C w/ 950 nm SLD,660 nm LED cluster

13-3202 Sys\*Stim 540LC w/ 785 nm applicator, 950 nm SLD 660 LED cluster

# accessories

13-3210 Sys\*Stim 540 785 nm applicator

13-3211 Sys\*Stim 540 660 nm SLD, 950 nm LED

cluster applicator

13-3212 Sys\*Stim 540 protective UVex glasses (1 pair)



# cervical and lumbar traction devices

# overdoor cervical traction



- mounts "over" door
- includes universal velcro head halter, 20
   lb. water bag, hardware
- machine washable

50-1000 50-1000-25 overdoor cervical bulk set (25 each)

# NECKPRO

- uses a compression spring coupled with a ratcheting device that delivers a precise amount of cervical traction
- each click audibly signals an increase in the tension being applied
- NeckPro<sup>™</sup> overdoor cervical traction system includes NeckPro<sup>™</sup> compression spring and ratchet assembly, head halter and overdoor bracket (fits doors measuring up to 1½" (3.8cm) in thickness
- commercial size bracket available for doors measuring up to 13/4" (4.4cm) in thickness

00-4395 NeckPro™

# optional accessories

04-3907 thick door-mount (13/4") 04-3911 replacement halter

# TracCollar™ pneumatic cervical traction



- portable and lightweight traction collar gently stretches and relaxes muscles slightly for pain relief
- developed by orthopedic engineers for comfortable, effective, daily treatment
- for professional and home use
- does not require assembly or batteries
- fastens easily with buckles
- includes carrying case
- collar is constructed from soft, plastic foam; internal inflatable bladders contain natural rubber latex
- 97% cotton and 3% Lycra® fabric covering can be easily removed and cleaned with warm water and mild soap
- not for recent sprains/strains, injuries or acute conditions

50-1095 TracCollar S/M, 14-16" red 50-1096 TracCollar L, 16-18" blue

# Posture-Pro™ cervical traction









- up/down motion promotes disc lubrication/joint nutrition
- portable, easy to operate

50-1040 Posture Pro 1000 cervical 50-1041 Posture Pro 2000 lumbar 50-1045 Posture Pro 4100 combo

50-1046 Posture Pro 6100 electronic combo

# supine C-Trax® cervical traction






- used to administer static or intermittent cervical traction while patient is in supine position
- mounts to any door frame
- includes head halter, 20 lb. water bag and mounting hardware

50-1001 supine C-Trax® cervical w/ regular halter 50-1002 supine C-Trax® cervical w/ TMJ halter



# cervical and lumbar traction devices











# Saunders® Cervical HomeTrac® deluxe

- · directs traction force towards occiput
- applies up to 50 pounds of traction
- user has total control of force at all times
- provides a comfortable and cost-effective option to continuing clinical traction treatments
- includes user's guide and deluxe carrying case

00-7044 Saunders® Cervical HomeTrac®
50-1007 optional extension foot for HomeTrac®
50-1008 replacement pump/gauge assembly
50-1080 instructional video

### ComforTrac® cervical traction

- · adjustable neck cushions slide in and out for custom fit
- quick release allows user to immediately release all pressure; increasing safety and convenience
- allows three angles of traction at 10, 15 and 20 degrees with no additional parts, attachments or accessories required
- width of neck wedges are independently adjustable
- unit comes standard with a wheeled trolley bag carrying case

50-1070 ComforTrac® Cervical Traction

# Pronex® pneumatic cervical traction

- even distraction in the anterior and posterior cervical discs
- patient controls the amount of traction applied
- squeezing an inflator bulb increases the pressure, a release knob gently reduces the pressure

 50-1050
 regular
 (14-16" neck)

 50-1051
 large
 (16-18" neck)

 50-1052
 wide
 (18" and up neck)

# Saunders® lumbar HomeTrac® deluxe

- actively moving treatment surface, and pump/gauge
- effectively applies up to 200 pounds of traction
- friction-free surface ensures smooth motion

50-1005 lumbar HomeTrac®
50-1009 replacement pump/gauge
50-1081 instructional video

☎ MAP (minimum advertised price applies)

# ComforTrac® lumbar traction

- unit comes standard with a trolley bag carrying case
- automatic carriage return increases ease of use, and aids in intermittent traction applications
- innovative hinged design requires no assembly or set-up
- quick release device on pump allows user to immediately release all pressure; increasing safety and convenience
- hand pump displays in both pounds and kilograms

50-1075 ComforTrac® lumbar



# Chattanooga® TX® traction tables, units and accessories

# TX TTFT-200 table

# TX TTET-200 table

# TX ADP-400 tables





4 section ADP-400 with TX traction unit (00-6870)


- high-low
- hand control
- optional foot control
- optional factory installed casters



# TX unit on mobile TX stand

TX (00-6870)

with mobile traction stand (00-3600)



- easy to operate
- · digital touch screen interface
- intermittent, static and cyclic traction programs
- user-defined hold, rest, treatment times
- real-time digital treatment tracking
- built-in safety features
- patient interrupt switch
- · cervical maximum poundage warning · audible signal at end of treatment and when patient interrupt switch is
- activated
- traction tension parameters 0 200 lbs.
- progressive and regressive settings
- store up to 10 user protocols
- store up to 14 patient treatment sessions on patient data card
- patient pain mapping

### TX traction unit

00-4759 TX unit

### TX traction stand

00-3600 mobile traction stand

### TX traction tables

00-6840 TTFT fixed height 00-6870 TTET height adj. AD-P300 ADP300 height adj. AD-P300C ADP300 w/ front casters AD-P400 ADP400 height adj. AD-P400C ADP400 w/ front casters

### TX integrated package

00-6878K TTFT package\* 00-6876K TTET package\* 00-6874K ADP400 package\* 00-6874KC ADP400 w/ front casters\*

\* package contains table with TX traction unit, flexion stool (TXS-1), deluxe universal belt system and Saunders cervical system

# ADP400 options

08-6049 retractable front casters 06-5988 hi-low foot switch 04-8082 patient interrupt switch





TX can pivot

# Saunders® cervical system



- pulls at base of the occiput for more effective cervical traction
- performs horizontal traction
- clevis included, choose from list at right

# Saunders cervical unit

00-7040K traction device\*

# clevis choices\*

06-2449 TX/Triton mobile stand (prior to 2006) 04-8038 TX/Triton (after 2006)

06-2596 Triton (prior to 2006)

06-2587 Tru-Trac 02-0575 Midland

additional clevis' available at \$ each



all clevis'

# Chattanooga® TX® traction accessories

# pivotal therapy system

components may be used individually to address specific problems or integrated to form a complete system to treat the entire spine

### individual pieces

**Occipivot** (a) treat tension and stress syndromes, headaches, TMJ, cervical syndrome and tension at base of skull. Promotes the release of the intrinsic/extrinsic muscles of the cervical spine.

**Cervical pivot point** (b) useful as a substitute for gentle cervical traction. It can also be used as a self-mobilization and strengthening tool. Promotes self-release for sternocleidomastoid, upper trapezius and 1st ribs.

Short Thoracic Pivot (c) use in combination with cervical (b) and lumbar (e) pivots, for craniosacral self treatment, as a sacral wedge for mobilization, as a rib release tool, and for sternal and anterior rib mobilization.

Standard Thoracic Pivot (d) treat lower thoracic/upper lumbar dysfunction use anterior release techniques or specific rib-release techniques. Ideal posture re-training tool for patients with forward head/kyphosis slumping.

**Lumbar Pivot** (e) Promote correct lordosis for static positioning or pivot in figure-eight along with bi-directional pelvic tilts. Add a sandbag to lower the abdomen for additional extension, self mobilization, and strengthening exercises for the abdomen.

00-4502 occi (a) 00-4504 cervical (b) 00-4506 short thoracic (c) 00-4508 standard thoracic (d) 00-4510 standard lumbar (e)



# a – occipivot

b - cervical pivot

c - short thoracic

pivot



d – standard thoracic pivot



e – lumbar pivot

# pivotal therapy set (a,b,c,d,e)

- for combined treatment of spine, surrounding muscles
- includes one each: occipivot

   (a), cervical pivot (b,d), short
   and standard thoracic (c)
   pivot and standard lumbar
   (e) pivot

00-4500 5-piece set



# Cambion® foot care



 shock dampening, cushioning and pressure distributing visco-elastic foot care products



- open cell structure adds breathability and reduces heat
- antibacterial treatment reduces mold and mildew
- Poron® foam top responds to changes in your foot while walking
- smooth, non-piling, easy-to-clean

Heel Pad [1] 01-3105 A 01-3106 B 01-3107 C 01-3108 D	Posted Heel Cush. [3] 01-3126 A 01-3127 B 01-3128 C 01-3129 D
Heel Spur Pad [2]	Full insoles [4]
ricci opui i uu [L]	
01-3109 A	01-3100 A
01-3109 A	01-3100 A
01-3109 A 01-3110 B	01-3100 A 01-3101 B

# face cradle pillow



- provides maximum patient comfort when lying in the prone position
- elevated base and nose cutout allow patient to breathe freely

00-4206 pron-pillo®

# Hi-Lo stool



 adjustable height (13-20") feature allows therapist to place lumbar spine into neutral or flexed position for patient comfort during lumbar traction

00-7050 stool


# knee bolster



- enhance patients positioning during traction therapy
   knee holster measures 10 y
- knee bolster measures 10 x 10 x 12"

04-2901 knee, 2 each

# TXA-1 traction package



traction package (TXA-1) has adjustable cervical halter, thoracic restraint, pelvic belt, 2 restraining straps, 2 cervical pillows, and spreader bar

**Traction accessory package** 00-1440 TXA-1 package

TXA-1 package components 00-1405 adj. cerv. halter [a] 00-1301 TX pillow (ea) [d] 00-1341 pillow case (ea) [d] 06-0030 17" spreader bar [e] 00-1433 HD pelvic pkg. [b]

HD pelvic package components 00-1432 2 thoracic straps [c]



# G5® percussor massage machines

# Thumper Pro



- variable speed dial controls the anatomically designed massaging spheres pivot up and down one quarter inch up to 35 times per second
- · optional foot stand is available

14-1080

14-1083 foot stand 14-1084 carry bag

# Thumper mini-Pro



- powerful 1/4" percussion stroke, 2 massage spheres to "HUG"
- elongated and angled handle to reach down your back
- 3 power settings

14-1081 mini-pro 14-1085 carry bag

### Thumper Versa-Pro



- five-in-one lower body massager for feet, calves, hamstrings, quadriceps, lower back
- 7 auto-program LED indicators
- 6 power level LED indicators

14-1082 versa-pro



- height adjustable trolley from 34-
- penetrating massage using frequencies 15 to 55 CPS
- ideal for sports medicine and physical therapy clinics
- includes PRO-PKG applicators (14-1512 below)

14-1500 Plus

masseur®

length: 21 in.

width: 21 in

height: 45 in.

comes with

trollev

speed ranges 15 to 60 CPS

ideal for high volume PT clinics

complete with caster stand and

tray. Can be used with pole, table-top or wall-mounted (mounting hardware extra)

includes PRO-PKG applicators

weight: 25 lbs.

Gemini™ lenath: 18 in. width: 22.5 in. height: 43 in. weight: 48 lbs. comes with

ideal for physical, massage and respiratory therapy

trolley

- speed ranges 15 to 60 CPS
- · offers a macro and micro head
- custom cart with built-in bins and PRO-PKG and MCT-PKG applicators (14-1512 and 14-1511



speed ranges 15 to 60 CPS

trolley

- stroking percussion apparatus designed as an adjunct for massage
- complete with caster stand and tray. Can be used with pole, table-top or wall-mounted (mounting hardware extra)
- includes PRO-PKG applicators (14-1512 below)

### 14-1504 workout mas 14-1505 GK-3 applicator and applicator cover packs



(14-1512 below)

MCT-PKG



PRO-PKG



PRO-PORT-PKG

- micro right angle adapter, micro soft sponge rubber, micro fine sponge with back, 10 micro fine sponge replacements, 10 micro course sponge replacements, micro blunt tipped applicator, micro skin and scalp applicator, micro round firm applicator, 50 extra small disposable applicator covers 14-1511 micro pack
- large sponge applicator 3.5" diameter, multiple prong applicator, fourball applicator, pointed-tip applicator, large round firm rubber applicator 5" diameter, curbed sponge applicator, 50 medium and large disposable applicator covers

14-1512 pro pack

four and two-ball applicator, pointed-tip applicator, large round firm rubber 5" diameter, curved sponge applicator, 25 medium and large disposable applicator covers

14-1513

pro port pack

# Therassist® length: 21 in. width: 21 in. height: 42 in. weight: 25 lbs comes with trolley

- speed ranges 15 to 60 CPS
- ideal for stress reduction, myofacial release, sports rehabilitation, postural drainage
- complete with caster stand and tray. Can be used table-top or wall-mounted (mounting hardware extra)
- includes PRO-PKG applicators (14-1512 below)

14-1503 Therassist

# ProPower®



- · portable hand-held massager
- penetrating massage using frequencies 15 to 55 CPS
- use before or after activity; deal for self massage at home
- includes PRO-PORT PKG applicators (14-1513 below)

14-1506 ProPower

# **GBM**®



- · portable hand-held massager
- penetrating massage using frequencies 15 to 55 CPS
- ideal for light duty clinics and
- includes PRO-PORT PKG applicators (14-1513 shown left)

14-1501 GBM

individual applicators and applicator covers available



# Chattanooga® compression therapy

# PresSsion® multi-3 sequential





PresSsion® intermittent

compression sleeves

- intermittent compression sleeves can be used with 3-chamber sequential unit to provide intermittent compression
- must use with included adapter

# arm

04-3013 half 17" 04-3014 full 29" full wide 29" 04-3015

leg

foot/ankle 10" 04-3016 04-3017 half 20" 04-3018 three quarters 26"

04-3019 short 29" 04-3020 med 32" 04-3021 long 35" 04-3023 long wide 35"

# • use PresSsion® intermittent sleeves with the sequential pump adapter for intermittent compression

and 3 chamber sequential compression

• use the 3 chamber sleeves for sequential compression • adj. pressure (30-100mmHG) and fixed

cycle (inflation 35 sec., deflation 10 sec.)

• two treatment outputs

### pump only

3 chamber sequential 00-4331

# PresSsion® multi-6 sequential



- 6-chamber sequential compression
- the sleeve inflates sequentially from the distal to the proximal chambers
- distal chamber receives the full selected pressure; pressure to the proximal chamber is at a reduced level
- adjustable pressure (30-100mmHG) and cycle (inflation 5-120 sec., deflation 5-60 sec.). 15-120 minute timer
- two treatment outputs

00-4333 6 chamber sequential

# PresSsion® sequential compression sleeves




- sequentially inflates from distal to proximal chambers
- use on 3 or 6 chamber compression units

arm	type	size	chambers
04-3034	half	17"	3
04-3093	full	30"	3*
04-3035	full	30"	6
leg			
04-3036	half	20"	3
04-3194	3/4	26"	3*
04-3149	3/4	26"	6
04-3092	full	32"	3*
04-3032	full	32"	6

\* for use with multi-3 only

# Nylatex® wraps





- heavy-duty elastic wrap functions very well as a compression wrap
- · velcro closure assures that the wrap remains tight
- Nylatex® wraps can also hold and position cold and hot packs

00-1200 assorted 9 ea. (1 each excluding 00-1213, 00-1264)

00-1202	21/2x18" 3 ea
00-1264	21/2x24"3 ea.
00-1204	21/2x36" 3 ea
00-1206	21/2x48" 3 ea
00-1208	4x18" 3 ea.
00-1210	4x36"3 ea.
00-1212	4x48" 3 ea.
00-1213	4x60" 3 ea.
00-1214	6x18"3 ea.
00-1216	6x36" 3 ea.
00-1218	6x48" 3 ea.
00-1219	6x60" 3 ea.



# Kinetek<sup>™</sup> continuous passive motion (CPM)

# Kinetek® Spectra™ knee CPM machine



- long-life DC motor, anatomically correct alignment to increase patient compliance and decrease joint stress, wide carriage that accommodates entire spectrum of users, ROM from hyperextension (-10°) through full extension (120°), patient lockout which prevents program changes, waterproof hand control, pause settings for up to 15 minutes of rest ad flexion or extension, and a 24 month limited warranty
- weighs only 24 lbs

13-1185 Spectra knee CPM

13-1186 Replaceable Patient Pad Kit 13-1187 Spectra CPM knee trolley 13-1188 Spectra CPM bed lift

# Kinetek® Performa™ precision knee CPM machine




- UL approved for hospital use, this CPM machine features a leg carriage designed for automatic motion
- adult and pediatric capability
- features waterproof computer, quiet operation, programmable muscle stim capability, carrying case for easy transportation, built-in therapy timer, attached sliding hip bar, manual set-up mode, pad system and straps

13-1170 Performa knee CPM

13-1171 Replacement Patient Pad Kit

# Kinetek® Maestra® hand CPM machine



- cut set-up time by incorporating a low-profile rigid splint
- lightweight and easy to operate; patient can easily and quickly receive comprehensive motion therapy for variety of indications
- Clinical benefits: anatomically shaped drive bar enhances abduction and adduction of the web spaces as joints are being mobilized; up to 15 minutes of pause at extension and flexion allows for a controlled stretch and rest period; composite flexion, intrinsic plus and minus motion included; convenient digital display of angle, ROM settings, speed and force
- Indications: indications include the treatment and prevention of intra-articular adhesions, extra-articular contractures, and excessive post-operative swelling

13-1165 Kinetek Maestra hand CPM

# Kinetek® Centura™ shoulder CPM machine



- provides post-operative passive range of motion
- helps prevent joint stiffness, soft tissue contractures, atrophy
- patient set-up is anatomical, comfortable, and quick with a colorcoded set-up system
- programmable hand control allows pre-programming sixteen different protocols
- · change for right and left is easy
- breaks down for convenient transportation; also includes a carrying case for detachable parts
- includes pad kit, carrying case and reusable chair cover

13-1150 Centura shoulder CPM
13-1151 replacement patient pad kit
13-1152 Centura reusable chair cover
13-1153 Centura disposable chair cover



# Chattanooga® knee continuous passive motion (CPM)





# **Opti-Flex 3 knee CPM**

 designed to provide anatomical motion for all types of CPM knee patients: pediatric, athletic and bariatric

### Features:

- Progressive ROM automatically increasing daily ROM
- Oscillation Zone<sup>™</sup> allows more time spent in "active" ROM
- Fast Back™ accelerates through nonworking ROM while knee spends more time in "active" ROM
- Comfort Zone™ addresses patient pain threshold by providing therapy at a reduced flexion angle limit

## **Artromot K1 knee CPM**

weighs just 23 lbs, making it the lightest knee CPM on the market

- only 38 inches (97 cm) long
- wide low-profile, closed-base frame for stability
- faster motor than the K3
- accommodates patients 4 ft to 6.6 ft (1.2 m 2 m)
- reversible footplate for short leg length
- on-board power supply, ROM -10° to 120°

# knee CPM units

00-2090K 00-2026

Opti-Flex K3 knee CPM Artromot K1 knee CPM

patient kits

02-0533

Opti-Flex patient kit, tan

01-1432 01-1438 Breg patient Kit Kinetek patient kit

straps

04-0128

Opti-Flex 34" strap kit

# **Opti-Flex knee CPM trolley**

- OptiFlex® HI-LO knee CPM trolley
- OptiFlex® HI-LO is a multifunctioned CPM trolley design
- ease of storage and transportation
- allows for proper positioning of CPM
- adjustable height from 18.25" to 28"

# **positioning trolley** 00-4520 hi-lo knee trolley







# Artromot K3 and K4 knee CPM

- German designed and engineered to provide effective post-operative ROM for the knee and hip joint with minimal adjustments for ease of operation
- K4 features a true anatomical hip axis providing precise knee and hip motion
- lifting, transporting and storage are easy with the "grab easy" base design

# K3 and K4 features

- K4 features true anatomical hip axis providing precise knee and hip motion
- handheld programming of device parameters, i.e. extension, flexion, pause, force and speed
- four hidden programming features: warm-up, patient runtime, device runtime, full speed and motion
- wide, deep leg orthosis to accommodate larger patients
- reliable solid-state digital electronics
- wide, stable skid resistant base
- removable power cord

### knee CPM units

00-2043 Artromot K3 knee CPM 00-2044 Artromot K4 knee CPM

# patient kits

02-0710 02-0711 Patient Kit K3 fleece Patient Kit K4 fleece



# Chattanooga® shoulder, elbow, ankle continuous passive motion (CPM)



# **Artromot S3 shoulder CPM**

- provides anatomical motion for the shoulder, providing enhanced patient comfort, compliance and more beneficial outcomes
- addresses full range of motion capabilities

### **Features**

- hand control for storing treatment data on memory chipcard
- fully synchronized motors allow treatment parameters to be customized to meet patient's specific therapy protocols
- ergonomic healthy armrest and fully adjustable chair
- patented universal left / right shoulder design
- folds easily for transport / shipping
- lightweight design 55 lbs (25 kg)

# **Opti-Flex S shoulder CPM**

- easy-to-use hand control to store treatment data on memory card
- fully synchronized motors allow treatment parameters to be customized to meet patient's specific therapy protocol
- · ergonomic rest for healthy arm and fully adjustable chair

# shoulder CPM units

00-2029 OptiFlex S

### patient kits and accessories

Patient Kit Optiflex S 02-0660 02-0719 Patient Kit S2/S3 ARM

02-0715 Patient Kit S2

opti-card data storage

02-0661 5-card pack



card reader













# Artromot E2 elbow

- provides anatomical motion for the elbow, and to enhance patient comfort and compliance
- Artromot® elbow addresses full range of motion capabilities pronation, supination, flexion and extension



00-2038 E2 Elbow CPM

# patient kits

02-0723 Patient Kit

opti-card data storage

02-0661 5-card pack





# Artromot SP2 2M ankle

- provides anatomical motion for the ankle, providing enhanced patient comfort and compliance
- addresses full ankle range of motion capabilities dorsi flexion/plantar flexion and inversion/eversion
- easy adjustable motors allow treatment parameters to be customized to meet patient's specific therapy protocols
- adjustable base allows for use in bed or chair

# OptiFlex ankle CPM

- easy to adjust motors so treatment parameters can be customized to meet patient's specific therapy protocols
- · easy-to-use hand control can store data on memory card
- · adjustable base for use on bed or chair
- patented universal left/right design

# ankle CPM units

00-2027 OptiFlex ankle CPM

### patient kits

02-0716

patient kit



# biofeedback

# biofeedback Stabilizer™





- Stabilizer™ pressure biofeedback provides valuable information to ensure quality and precision in exercise performance and muscle testing
- provides biofeedback to determine if the patient is able to selectively isolate and maintain contractions of the cervical or lumbopelvic core stabilization muscles
- measuring range is 0-200 mmHg analog pressure with an accuracy of +/- 3mmHg pressure
- use with exercise to protect and stabilize joints
- prevent and treat low back/neck pain
- improve core stabilization exercises
- enhance home exercise program
- · complete with display gauge, pressure cell,instruction booklet

00-9296 Biofeedback Stabilizer™, complete

# GSR (galvanic skin response)





- GSR module is self contained
- readings can be taking using the fingertip pads on the module, or the two electrodes can be placed on any body part
- GSR module gives an audible feedback (earphone comes standard and can be plugged in to the unit) or a visual feedback (optional meter can be plugged into the unit)
- optional temperature module (right) comes with a temperature probe that attaches directly to the GSR for either audio or visual meter feedback (temperature meter comes with unit)

15-4266 GSR monitor

15-4267 visual feedback monitor 15-4268 Temp 2X (right)

# surface EMG biofeedback system

- includes highly sensitive Myoscan active sensor; electrodes sold separately
- single channel surface EMG system allows freedom of movement for ambulatory exercises and provides a bright light bar and proportional tone biofeedback
- three sensitivity ranges for muscle strengthening, rehabilitation or relaxation therapy, and provides stable readings in less than a second
- accurate and versatile biofeedback for muscle re-education
- · two modes of biofeedback LED and 5 way tone
- 5 way tone feedback includes proportional, above and below threshold and 4-second delayed alarm
- 13-segment LED bargraph offers precise visual feedback
- wide and narrow band pass filters provide precise monitoring of normal and fatiguing muscles and eliminate heart-beat artifact
- 0-2 Volt DC output allows connection to third party equipment; permits monitoring
- low current consumption provides long battery life

### **EMG**

13-1442 EMG unit, with Myoscan active sensor

### accessories

450mm shielded, low noise extender cables for use with 12-0051 12-0053 and 12-0054

### sEMG electrodes


[A] triode electrode with triangular configuration; 2cm spacing 12-0052 Ag-AgCl electrode backed with nickel plated brass snaps; 100 ea.

12-0053 [B] strip electrodes; use as strip or separate; 100 strips

12-0054 [C] pre-gelled electrode for placement on sensitive, dry skin; 100 ea.







12-0052





# WHITEHALL® stainless steel whirlpools

- whirlpools have a single ½ horsepower motor to drive the agitation, spot spray and emptying operations
- dial varies the degree of agitation by controlling the amount of air introduced into the water
- the agitator can be raised and lowered on its shaft to change the directional flow; can be removed from the shaft if service is needed
- available with a polished satin finish or an electrostatically applied color
- · comes with thermometer
- tank can be emptied through the gravity drain or by using whirlpool turbine assembly and drainage hose









mobile lo-boy (66x24x18")

mobile hi-boy (48x24x28")







stationary hi-boy (48x24x28")

mobile extremity unit (25x13x15")

podiatry (22x13x12")

- mobile models mounted on 4 heavy-duty five inch casters
- stationary models designed for permanent installation with direct plumbing connections to the water supply and drain
- lo-boy long and shallow units are designed for easy access
- lower height for pediatric and geriatric patients requiring hydrotherapy treatment
- can provide treatment of various extremity regions along hip and back areas

		LxWxD**	height**	cap.**	
42-1214	mobile slant back SB-100-S	76x40x21	24	100	
42-1215	mobile L-75-M	52x24x18	24	75	
42-1216	mobile L-90-M	60x24x18	24	90	
42-1217	mobile L-105-M	66x24x18	24	105	
42-1211	stationary L-90-S	60x24x18	24	90	
42-1212	stationary L-105-S	66x24x18	24	105	

- hi-boy complete coverage to the back, hip and shoulder areas
- if full immersion is not desired, an optional adjustable suspension seat can be utilized
- when coupled with the optional tank top seat, unit can be used to treat lower extremities

		LxWxD**	height**	cap.**
42-1225	mobile H-60-M	36x20x28	34	60
42-1228	mobile H-105-M	48x24x28	34	105
42-1223	stationary H-105-S	48x24x28	34	105

• extremity - compact design saves space; mobile units(\*) work well with Hydrolift (right)

		LxWxD**	height**	cap.**	
42-1245*	mobile podiatry P-10-M	22x13x12	18	10 *	
42-1255*	mobile foot/ankle/arm E-15-MU	25x13x15	28	15 *	
42-1256*	mobile foot/ankle/calf/arm E-22-MU	28x15x18	34	22 *	
42-1258	mobile foot/ankle/calf/knee/arm E-45-M	32x15x25	31	45	
42-1251	stationary foot/ankle/calf/arm E-22-SP	28x15x18	34	22	
42-1253	stationary foot/ankle/calf/knee/arm E-45-S	32x15x25	31	45	

\*these units work alone or with Hydrolift™ (right) \*\*dimensions in inches, capacity in pounds



### Whitehall Hydrolift™

- Hydrolift™ (attached to whirlpool above) is essentially two whirlpools in one, allowing you to treat lower and upper extremities with the same unit
- raises the <u>mobile extremity</u> and <u>podiatry</u> whirlpool to any height up to 3 feet off the floor (select model to raise 10, 15 or 22 gallon tank)
- turn handle raises the whirlpool quickly and easily, even when completely filled

# Hydrolift unit

42-1060 for 10-gal. unit 42-1061 for 15-gal. unit 42-1062 for 22-gal. unit

# Hydrolift with whirlpool

42-1070 w/10-gal. unit 42-1071 w/15-gal. unit 42-1072 w/22-gal. unit



# whirlpool accessories for physical therapy

- lower extremities can be treated easily and comfortably with the adjustable whirlpool chairs
- chairs are designed to permit a close fit to the end of stationary or mobile whirlpool baths
- height adjustable to fit over the rim of the whirlpool, where they can be locked into place
- chrome plated steel frame
- Padded black vinyl upholstered seat and back





### hi-boy chairs

- adjustable whirlpool chair with casters has four 3" swivel casters, two with locks; a safety belt; height adjusts from 33-46½"; 14x19"non-slip rubber platform; 18x20" seat, 10x16" back
- model without casters has the same features except rubber tips replace the casters, and the height adjusts from 29½-43½"

42-1052 hi-boy whirlpool chair with casters 42-1051 hi-boy whirlpool chair without casters

# lo-boy chairs

- lo-boy whirlpool chair comes with or without casters. Height adjusts from 22½-28½" on model with casters, and 19½-25½" on model with rubber tips
- seat measures 18x20" and back measures 10x16"

42-1053 lo-boy with casters 42-1054 lo-boy without casters

# adjustable suspension seat

• for lower extremity treatment


42-1410 20" hi-boy 42-1411 24" hi-boy 42-1412 24" lo-boy

# reclining seat

 floor seat with reclining back for optimum patient comfort

42-1420 low boy

suspension seat



reclining seat

# whirlpool table with side bench seatrim height seat for lower extremity

treatments; side benches for hand/arm treatments from either side of tank

42-1405 20" seat hi-boy 42-1406 24" seat hi-boy 42-1407 24" seat lo-boy

# tank top seat

- attaches securely onto tank leaving area underneath unobstructed
- 12" wide covered seat extends around front and sides of whirlpool

42-1400 20" seat hi-boy 42-1401 24" seat hi-boy



whirlpool table w/ side benches



tank-top seat

# other accessories

42-1430 head rest 42-1431 arm support (ea.) 42-1432 hand grips (pair)

42-1450 washout hose

# thermostatic mixing valves

42-1440 15 GOPURA ½" tube 42-1441 25 GPM ¾" tube

# Lil' Champ™ whirlpool



- for table top, floor, or in stand
- variable speed turbine
- 23.25L x 13.75W x 19.25H
- insulated lid

42-1080 AW-606



# portable bathtub whirlpool

- electric turbine ejector-aerator, mounted on a "U" shaped stainless steel base, is designed to fit inside a standard sized bathtub
- 1/3 Hp motor, moves 60 gal/min.
- portable

42-1260





# **DIPSTERS** patientwear

- limited-use/disposable garment is ideal for hydrotherapy and other situations requiring limited body coverage
- Dipsters® have proven effective in over 5,000 hospitals, clinics, facilities, training rooms... and even hotel pools since 1974
- limited-use/disposable patient drapes are tough enough to reuse, inexpensive enough to use only once
  - can be washed and gas autoclaved
  - made from DuPont's Tyvek® (polyethylene)
  - comfortable
  - · they dry immediately

# men's drawstring waistband boxer shorts

• boxer shorts with adjustable drawstring waistband and full supporter

men		dozen
20-1010 20-1011 20-1012	small medium large	
20-1013 20-1014 20-1015	x-large xx-large xxx-large	



# men's elastic waist and boxer shorts

• boxer shorts with elastic waistband and full supporter

men		dozen
20-1000	small	
20-1001	medium	
20-1002	large	
20-1003	x-large	
20-1004	xx-large	
20-1005	xxx-large	
boys		dozen
20-1020	small	
20-1021	medium	
20-1022	large	
	-	



# bib top and bikini



- drawstring ties top around neck and back · top is held in place with
  - bikini has elastic waistband and legs

wome	n	dozen
20-1030	small	
20-1031	medium	
20-1032	large	
20-1033	x-large	
20-1034	xx-large	
20-1035	xxx-large	

20-1065 xxx-lg girls dozen 20-1070 small 20-1071 medium 20-1072 large

stretch middle and legs

tie-strina

women

20-1060 small

20-1062 large

20-1064 xx-lg

20-1063 x-large

20-1061 medium

elastic trim and neck

dozen

girls dozen 20-1050 small 20-1051 medium 20-1052 large

# halter top and bikini



- · top has drawstring that ties around back and
- bikini has elastic waist and legs

women		dozen
20-1080	small	
20-1081	medium	
20-1082	large	
20-1083	x-large	
20-1084	xx-large	
20-1085	xxx-large	

# coverall and shorts



- · coverall top fits over head
- open sides closed with tie-string around waist
- · shorts have drawstring waist women

dozen

20-1090	small	
20-1091	medium	
20-1092	large	
20-1093	x-large	
20-1094	xx-large	
20-1095	xxx-large	

bib top and shorts

- bib top with boxer shorts
- drawstring top ties around neck and at back
- · Shorts have draw-string waist

women		dozen
20-1040	small	
20-1041	medium	
20-1042	large	
20-1043	x-large	
20-1044	xx-large	
20-1045	xxx-large	



# aquatic therapy products



- hand bars and ankle cuffs are made of buoyant and durable soft and non-abrasive closed-cell foam
- material weighs little on land, but when in water, these low-impact exercisers add resistance to your arm and leg movements
- for aquatic exercise routines and rehabilitation
- made specially so colors do not fade or leach

# aquatic ankle cuffs



- superior comfort and durability
- resistance for leg workouts
- adjustable

20-4040B ankle cuff blue, pair 20-4040R ankle cuff red, pair

# aquatic pull buoy



- made of smooth, contoured foam
- provides excellent lower body buoyancy
- does not cause chaffing and irritation

20-4050B/R pull buoy small 20-4051B/R pull buoy large

# aquatic exercise kit



- standard kit includes: jogger belt, hand bars
- deluxe kit includes: jogger belt, hand bars, ankle cuffs

20-4210B/R standard kit, small 20-4211B/R standard kit, medium 20-4212B/R standard kit, large 20-4200B/R deluxe kit, small

20-4201B/R deluxe kit, medium

20-4202B/R deluxe kit, large

Fabrication Enterprises

# aquatic hand bar



- two 5 1/2"x 2 1/2" floats on 12" bar
- for upper-body workouts
- non-slip hand grips

20-4021B hand bar blue, pair 20-4021R hand bar red, pair

# aquatic swim belt



- smooth, comfortable fit
- adjustable strap with dual locking buckle

20-4001B/R swim belt 2 floats 20-4002B/R swim belt 3 floats 20-4003B/R swim belt 4 floats

# kick boards



- embossed non-slip surface and contoured handgrips. "Classic" will help you maintain proper body positioning
- "Champion" promotes correct body positioning and form during training. Has several hand positioning locations, side handgrips

20-4100B/R classic junior 20-4101B/R classic adult 20-4110B/R champion

# aquatic swim bar



- improves balance and stabilization while increasing muscle tone and flexibility
- two 6" floats on 30" cushioned hand bar

20-4030B swim bar blue, each 20-4030R swim bar red, each

# aquatic jogger belt



- keeps you buoyant & properly positioned
- adjustable

20-4010B/R small (60-160 lbs) 20-4011B/R medium (160-220 lbs) 20-4012B/R large (220+ lbs)

# noodles



- buoyant enough to hold 200 lb. person
- superior comfort and durability
- noodle measures 3" diameter x 54" long
- noodle bends and contorts with ease

20-4220 exercise noodle - 1 20-4220-32 exercise noodle - 32 20-4220-96 exercise noodle - 96

# splint pan SP-1501



- dual thermostatic controls (175°F high and 160°F low)
- · integral digital thermometer
- insulated lid

max sheet 24-4090 SP-1501 24"x18"

# splint pan SP-1502



- dual thermostatic controls (190°F high and 150°F low)
- integral digital thermometer
- insulated lid

max sheet 24-4091 SP-1502 19"x11"

• for table top, floor, or in stand

splinting products

Lil' Champ whirlpool

- · variable speed turbine
- 23"L x 14"W x 19"H

42-1080 AW-606


# Orfit® splinting material

- Ortho/Ease/Drape provide superior finish/feel splinting material sold individually or in case of 4
- does not contain rubber; not sensitive to UV light comes in 3 sizes with or without macro perforation

# Orfit® Ortho® splinting material



- maximum resistance to stretch and makes excellent rigid and large splints
- adheres to itself when hot; good for uncooperative patients; fingerprint resistant
- limited memory will hold its shape well during remolding
- molded directly on the patient's skin without protection
- hot water bath, kitchen oven, microwave or hair dryer can be used





# Orfit® Ease® splinting material

- compares to Rolyan Ezeform, Tailor-Splint and Kay-Splint
- medium resistance to stretch and has draping and contouring features
- adheres to itself when dry heated; great for the fabrication of large or small splints
- controlled stretch even with large splints; maximum rigidity (flexural modulus 850mPa)
- not tacky but temporarily self-bonding, giving the splint maker a helping hand
- permanent bond by dry heat and hard pressure; not coated; no need for solvent

<b>Ease</b> ®	18 x 24 x 1/16"		18	x 24 x 1/8"	24 x 36 x 1/8"	
	each	case	each	case	each	case
non perf	24-4120	24-4120-4	24-4121	24-4121-4	24-4122	24-4122-4
macro perf	24-4130	24-4130-4	24-4131	24-4131-4	24-4132	24-4132-4



# Orfit® Drape® splinting material

- compares to Polyform and Polyflex II
- coated material to prevent accidental bonding
- water based coating allows temporary bonding (when warm and toweled)
- excellent drape for intimate fit with minimal handling (gravity molding)
- good rigidity even when stretched out. Excellent for bases for dynamic hand splints
- no recoil when reheated; splints hold their shape when reheated so that no fit is lost
- self-sealing edges hold the seal even when reheated
- permanent bonding possible by dry heat after physical removal of the coating

<b>Drape</b> ®	18 x 24 x 1/16"		18	x 24 x 1/8"	24 x 36 x 1/8"	
	each	case	each	case	each	case
non perf	24-4140	24-4140-4	24-4141	24-4141-4	24-4142	24-4142-4
macro perf	24-4150	24-4150-4	24-4151	24-4151-4	24-4152	24-4152-4



# taping products

# water resistant Kinesio® tape

 gentle and effective approach to re-education of the neuromuscular system while also improving blood/lymph circulation, relieving pain and providing comfort and stability for your patient care needs

		each		i .		case —	
color	size	sku	price		sku (case)	qty.	price
beige	1"x5.5yds	24-4900			24-4900-6	6 rolls	
beige	2"x5.5yds	24-4910			24-4910-6	6 rolls	
blue	2"x5.5yds	24-4912			24-4912-6	6 rolls	
red	2"x5.5yds	24-4914			24-4914-6	6 rolls	
black	2"x5.5yds	24-4916			24-4916-6	6 rolls	
beige	3"x5.5yds	24-4920			24-4920-4	4 rolls	
beige	2"x34yds	24-4942			1		

### **Kinesio Taping Materials**

24-4960	Kinesio tape, pro scissors
24-4961	Kinesio tape, illustrated taping book
24-4962	Kinesio tape, perfect manual
24-4963	Kinesio tape, clinical therapeutic applications
24-4964	Kinesio tape, clinical taping DVD
24-4965	Kinesio tape, for lymphedema and chronic swelling
24-4966	Kinesio tape, in pediatrics (fundamentals and whole body)




# KT® treatment tape

 breathable, water-resistant, elastic tape that has multiple uses for clinicians as well as many other medical and health care professionals

	İ	each		Г		case —	
color	size	sku	price		sku (case)	qty.	price
beige classic	2"x20'	25-3410			25-3410-8	8 rolls	
blue classic	2"x20' 2"x20'	25-3412 25-3414			25-3412-8 25-3414-8	8 rolls 8 rolls	
pink classic black classic	2 x20 2"x20'	25-3414 25-3416			25-3414-8 25-3416-8	8 rolls	
black Mega Roll		25-3416MR			25-54 10-0	0 10115	
beige pre-cut	2"x20'	25-3420			25-3420-8	8 rolls	
blue pre-cut	2"x20'	25-3422			25-3422-8	8 rolls	
pink pre-cut	2"x20'	25-3424			25-3424-8	8 rolls	
black pre-cut	2"x20'	25-3426	I	ı	25-3426-8	8 rolls	I

# spider tech pre-cut kinesiology tape

- tape-based therapeutic support developed to facilitate a functional approach to pain modulation and myofacial dysfunction
- "ready to apply" applications meet the needs of today's rehabilitative and performance enhancing professionals by providing a coupled, synergistic influence on pain reduction, injury prevention and performance enhancement called 3P





# AcuForce® 2.5



- uses it's own weight (2.5 lbs.) in applying pressure
- saves thumbs, wrists and elbows from fatigue and injury
- performs muscle stripping, cross-fiber friction and trigger point therapy
- · works well with lotions

14-1393 AcuForce 2.5

# AcuForce® Massage Star



- 3 applications, the wedge, on the top of the tool, to perform muscle-stripping and cross-fiber friction
- wide and narrow points to perform trigger point therapy and reflexology
- · XL for larger hands

14-1380 massage star 14-1381 massage star XL

# AcuForce® 7.0



- use it's own weight (7 lbs.) in applying pressure
- saves thumbs, wrists and elbows from fatigue and injury
- performs muscle stripping, cross-fiber friction and trigger point therapy
- works well with lotions

14-1390 AcuForce 7.0 14-1391 carry case

# massagers

# AcuForce® 3.0i



- highly effective cyro-therapy
- 3 pound unit forms a large ice surface under the tool to perform long-lasting therapy
- weight of the tool removes the need to apply pressure during treatment

14-1392 AcuForce 3.0i

# Point Relief® mini massager



- compact, portable, small enough to fit in purse or pocket
- battery operated (uses 1 "C" battery, included)

	each	25 each
massager	14-1050	14-1050-25
massager with heat	14-1051	14-1051-25

# facilitation massager



- battery operated vibrator with one speed (105cps)
- two "C" cells fit into plastic handle

14-1040 each 14-1040-25 25 each

# CAT® adjusting tool



- chiropractic thrust adjustment device for spine and extremities
- adjusts from 20-45 APR (average force range measured in pounds of force)

14-1400 CAT

# rotary massager

- large pad hand-held massager
- large handles support rotary (orbital) massager when used as a floor unit to massage feet

# massagers

single speed
dual speed
variable speed

14-1004 variable speed with access post14-1003 2-speed with foot attachment

# accessories

14-1008 foot/leg orbital massager holder 14-1009 sheepskin cover for 14-1008 14-1005 accessory kit for large pad 14-1006 carry bag

14-1006 carry bag 14-1007 wall hanger

# dual head massager



 twin head design allows simultaneous massage to both sides of the neck, shoulder blade and arm

14-1020 dual speed

# single head massager



- two-speed action with a large or small head
- soothing heat available on large head model

14-1030 small head (1) 14-1031 large head (2) 14-1032 large w/heat (3)




# manual massagers - deep-muscle therapy tools

big bend

14-1337

# Backnobber®



- hand crafted from a steel bar, bent in the shape of an "S" with wooden balls, one relatively small and one somewhat larger, permanently affixed to each end
- hook device over shoulder or under arm and use its leverage to apply direct, deep, steady pressure to trigger points

14-1335 regular: 19" length, curves 9" across 14-1336 long: 29" length, for people over 6 foot tall 14-1337 big bend: 26" length, for people over 200 lbs.

# Backnobber® II



- reinforced polymeric resin
- hooking device over the shoulder or under the arm
- 35-page users guide pictured above

14-1330 Backnobber II

## tips for locating trigger points, deep muscle pressure techniques, muscle stretching, strengthening exercises and frequently asked questions

Jacknobber® II

instruction card

14-1370 instruction card

# Index Knobber®



- 7/8" diameter wooden ball
- assist the "hands-on" therapist engaged in deep muscle therapy to avoid overuse injuries to the hands, wrists and forearms

14-1311 Index Knobber

# Index Knobber® II



- for clinic or home use
- durable polymer allows user to apply deep, sustained and precise pressure to trigger points, and sore, stiff muscles

14-1310 Index Knobber II

# Jacknobber<sup>®</sup>



- trigger point therapy tool as well as a general purpose, stress relieving massager
- apply to muscles around the shoulder blades for comforting backrub

14-1301 Jacknobber



- four knobs used to apply pressure to trigger points
- small knobs provide deeper, more intense pressure
- large knobs provide broader, gentler pressure

14-1300 Jacknobber II

# Palmassager<sup>®</sup>



- hold firmly in the palm of the hand and use 1, 2 or 3 of its smooth, round knobs to massage tense, tired or sore muscles of the neck, back, shoulders, arms or legs
- ergonomically placed finger holds for optimum comfort and control

14-1320 Palmassager

# Knobble®



- "cap" is held in the palm, and the smoothly rounded stem is used to apply pressure to trigger points in soft tissue
- reduces wear-and-tear to therapists hands

14-1340 Knobble (wood) 14-1341 Knobble II (plastic)

# Vibrassager®



- rolling massager equipped with 2 knobbed wheels that vibrate as they are rolled across the back for relaxing massage
- rounded handle end used as trigger point tool for those especially sore muscles

14-1360 Vibrassager

# Orbit Massager®



- give stress releasing massage to back, neck, shoulders and other muscles
- use on yourself or a partner
- freely rotating ball kneads the muscles to loosen and relax soft tissue and encourage circulation

14-1350 Orbit massager



# massage balls and rollers





# Cando® massage balls

- massage balls and rollers are used to help achieve muscle and tension release for all parts of the body
- "bumpy" textured surface also promotes sensory stimulation and is perfect for manual dexterity exercises
- colors may vary

	eacn	aozei
7cm 2.8in.	30-1995	30-1995-12
8cm 3.2in.	30-1996	30-1996-12
9cm 3.6in.	30-1997	30-1997-12
10cm 4.0in.	30-1998	30-1998-12
15cm 6.0in.	30-1999	30-1999-12

# Cando® massage rollers (pair)

· effective for massaging feet and large areas

pair

12 pair

6.5x16cm (2.6x6.3in.) 30-1994 30-1994-12

# massage lotion warmers

# Chattanooga® bottle warmers



Chattanooga mini Chattanooga TM1

- TM1 electric keeps massage lotion and ultrasound gel at a comfortable temperature
- easily attaches to any treatment table or wall. Mini-Lotion, (non-electric) warmer hangs on the front of the Hydrocollator E-1 or E-2 heating units

00-4012 TM1, for 3 8oz,btls, electric 02-2100 mini, for 3 8oz. bottles non-electric (conductive)

# NAImco® gel/lotion warmers



NAImco 3 bottle NAImco 5 bottle

- gel/lotion warmer are ergonomically designed heating units for massage gels and lotions
- comfortably warm up to three or five 250 ml bottles
- convenient, tabletop design is available
- designed for continuous 110V operation with NO water needed

50-6068 3 bottle warmer 50-6069 5 bottle warmer

# Ideal® lotion warmers

Electric gel warmers maintain desired temperature for the ultimate in patient comfort. Variable temperature range of 90° to 120°. Includes dispenser bottles. Comes with 6' grounded power cord. Wall mountable.

holds 1 8oz. bottles 50-6055 50-6056 holds 2 8oz. bottles 50-6067 holds 3 8oz. bottles 50-6057 holds 4 8oz. bottles 50-6061 holds 2 12 or 16oz. bottles



Ideal® gel warmer

# massage lotions



# Myossage® lotion

- professional-grade, soothing lotion is specifically designed for the clinician
- water-based, alcohol free formula is a great medium for giving massages

00-4210-1 gallon 00-4210-4 gallon 4/case 00-4208-1 8½ oz 8<sup>1</sup>/<sub>2</sub> oz 12/case 00-4208-12 00-4208-72 81/2 oz 72/case

# Ultra Myossage® with Aloe Vera

- contains aloe vera; does not contain menthol, hexachlorophene, or polysor-
- use as a massage lotion or as an ultrasound coupling agent
- non-allergenic, water soluble, no grease

gallon 00-4262-1 00-4262-4 gallon 4/case 00-4260-1 8½ oz 00-4260-12 81/2 oz 12/case 00-4260-72 81/2 oz 72/case



# portable massage tables

# economy table



- starter package for new graduates or those on a tight budget
- has many of the same great features as more expensive tables and includes piano hinges and interlocking legs
- includes carry case and stationary face rest

height range: 24" - 32" top: 30" x 73" weight: 34 lbs. capacity: 450 lbs.

15-3731 economy massage table

# massage chair

- arm height and angle adjustments
- adjustable chest pad; removable sternum pad; adjustable headrest; seat height and angle adjustments
- folds and totes easily
- chair, carrying case, removable sternum pillow, eyeglass / jewelry pouch

weight: 18 lbs. capacity: 300 lbs.

# table with adjustable back



- multi-purpose portable table that allows you to expand your menu of choices to include all forms of traditional massage, reflexology, facials and waxing
- tilt back adjusts to 10 different levels

height range: 24" - 32" top: 30" x 73" weight: 34 lbs. capacity: 450 lbs.

15-3733 table w/ adjustable back



- smaller piece of equipment for ease and mobility
- standard with carrying case

weight: 17 lbs. capacity: 300 lbs.

15-3734 desktop support

# deluxe table



- rounded corners, a stronger, flatter working surface, and Reiki end panels on each end of the table increasing your mobility all around the table
- with the Shiatsu release cables, this table can easily lay flat on the floor

height range: 23" - 34" top: 30" x 73" weight: 35 lbs. capacity: 450 lbs.

15-3732 deluxe massage table

# adjustable face cradle



- 9" on center
- tilts and elevates above surface of table
- adjustable face rest also elevates to support head in face up positioning
- lightweight and folds flat for storage

15-3735 base and cushion 15-3736 base only (no cushion)

# SaniCover



- crescent-shaped face rest covers blend the feel of cotton flannel with the convenience of a fitted disposable overlay
- soft, comfortable fabric mesh provides effective layer of sanitary protection

15-3737-50 bag of 50 15-3737-500 case of 500

15-3730 massage chair





# treatment tables

# Plynth™ treatment tables

- contemporary styled Plynth® has a solid welded steel frame for superior stability
- · powder coated paint not only makes the table more attractive, but adds to the durability, making the steel surface scratch resistant and impervious to rust
- high density, seamless cushions are comfortable and strong, eliminating snags, splits and frays
- Naugahyde® upholstery completes the stylish effect, available in your choice of colors (p. 143)








00-6480 1 section 28x78x30" 00-6484 2 section 28x78x30" 00-6485 3 section 28x78x30"

# electric hi-lo tables





footswitch control

power assisted caster system

400 lb. lifting weight capacity

15-1500 1 piece 27x76" top, 17-36" Hi-Lo no casters

15-1502 2 piece 27x76" top, 18-33" Hi-Lo head 15x27", 45°up / 90°down; center 61x27"

15-1506 3 piece 27x76" top, 18-37" Hi-Lo head 15x27", 45°up / 90°down w/ nose hole center 22x27", 30°up; foot 39x27", 73° up

15-1515 5 piece 27x76" top, 24x36" Hi-Lo head 24x36", 45°up / 90°down; center 22x27", 30°up; foot 39x27", 73° up armrest 13x6", 4" casters



15-1515

# electric pedestal tables

15-1536

- single pedestal design
- radiused, seamless corners on upholstered components
- footswitch control
- power assisted caster system
- 400 lb. lifting weight capacity
- one step caster system
- 15-1536 1 piece 27x76" top, 24-36" Hi-Lo no casters
- 2 piece 27x76" top, 24-36" Hi-Lo 15-1537 head 15x27", 45°up / 90°down w/ nose hole; center 61x27'
- 15-1539 3 piece 27x76" top, 24-36" Hi-Lo head 15x27", 45°up / 90°down w/ nose hole center 22x27", 30°up; foot 39x27", 70° up
- 5 piece 27x76" top, 24x36" Hi-Lo head 15x27". 45°up / 90°down w/ nose hole: center 22x27", 30°up; foot 39x27", 70° up armrest 15x6", 4" casters





- undercarriage access allows for sitting or standing upright during treatment
- height adjusts from 18" to 38" to accommodate wheelchairs and gurneys

ADP300



#### treatment tables

#### ErgoBasic™ Hi-Lo

## shown full loaded with options tilt headpiece

- accelerator pelvic drop
- paper holder and cutter
- seamless cushion
- height adjustable 18-24"

00-9040 soft foam top 00-9050 firm foam top

options

08-9295 pelvic seat 08-9288 EZ tilt headrest 08-9493

forward & toggle drop with tilt and elevation

#### TRE-CH3™ Hi-Lo table



- electric elevation from 20-38"
- tilting headpiece
- paper holder and cutter
- extending ankle rest
- 300 lb. capacity

TRE-CH3 (120V) 00-6630 06-2358 foot switch 06-2357 hand switch

#### ErgoWave™ IST



- 3 independent rollers
- height adjustable
- ultra-quiet operation
- 350 lb. capacity

00-9080 Ergowave 00-9080N w/ nose hole

options

08-7896 bolster leg pillow 08-8152 head pillow 06-7149 nose hole

#### upholstery colors:

note: color swatches may vary from actual upholstery color

black mocha burgundy

**I-Skin colors** Midnight Pacific Blue Terra Cotta Chocolate Sand Maroon

forest green slate midnight gray

dove gray imperial blue stone blue

#### Adapta® Mesa treatment tables



- ergonomically designed so no handles or knobs get between you and your patient
- scissor lift brings patients to your level using any of 6 easy-to-reach pedal con-
- gas assist moveable sections
- EZ lift retractable casters
- 450 lb. weight capacity

2-section Hi-Lo 18-40" top 28x78" 00-3319K Adapta I-skin

## Adapta® Summit Hi-Lo treatment tables



3 section

- PostureFlex "peaks" the center of the table to increase lumbar flexion/extension, postural elevation and positional traction
- ergonomically designed so no handles or knobs get between you and your patient
- scissor lift brings patients to your level using any of 6 easy-to-reach pedal controls
- gas assist moveable sections
- EZ lift retractable casters
- 450 lb. weight capacity

3-section Hi Lo 18-40" top 28x78"

00-3320K Adapta I-Skin

00-3323K Adapta I-Skin w/PostureFlex


7-section Hi Lo 18-40" top 28x78"

00-3321K Adapta I-Skin

00-3322K Adapta I-Skin w/PostureFlex



## ood color options mat color options hardwood treatment tables



#### plain treatment table

#### with H-brace

#### with shelf



standard treatment table

15-1010	24x72x30"
15-1011	24x78x30"
15-1012	30x72x30"
15-1013	30x78x30"



• H-brace construction for stability

15-1000	24x72x30"
15-1001	24x78x30"
15-1002	30x72x30"
15-1003	30x78x30"



- modified H-brace for stability
- full storage shelf is easily accessible

24x72x30"
24x78x30"
30x72x30"
30x78x30"

#### with adjustable back



- modified H-brace construction
- adjustable back raises to 85°
- storage shelf is easily accessible

15-1030 30x78x30"

#### enclosed treatment table



- table with enclosed cabinet, sliding doors and drawer
- · choose upholstered or raised rim
- 30 x 78 x 30"

15-1083 upholstered table 15-1082 raised rim treatment table 15-2123 mat for raised rim table

#### enclosed exam table



- exam table with enclosed cabinet, drawer and sliding doors
- space saving table has fold down 3-position head and foot rests
- 24" wide, 30" high, 52" folded, 72" extended; stirrups, 7x24" armboard

15-1073 enclosed exam table

#### Hi-Lo treatment tables

#### manual



• manual hydraulic tables: 30x78"; height adjusts 25-33"

#### plain top with shelf 15-1200 with drawer

15-1200 with drawer 15-1202 without drawer

adjustable back with shelf 2 section top (adjustable backrest) 15-1201 with drawer 15-1203 without drawer

#### raised rim top

15-1210 with drawer 15-1211 w/out drawer

#### mat for raised rim top

15-2123 30x78x2"

# 15-1220



- electro-mechanical drive offers smooth and quiet operation
- 78x30x2" top, adjusts from 27-39" high; foot control; 110V 60Hz
   6.5 amp 720 watt; thermal overload protection; UL hospital grade plug; power cord 120" long

#### plain top

15-1220 plain top 15-1221 plain w/shelves

#### 2 section adjustable top

2 section top (adjustable backrest) 15-1230 2 section top 15-1231 2 sec w/shelves



#### mat tables

#### fixed height mat tables and platforms

- 2" polyurethane mat is upholstered with vinyl coated nylon for patient comfort
- mat surface is 21" high to accommodate wheelchair transfer
- top has raised-rim to hold mat in place
- platform surface 19" high without the mat
- optional 2" mat sold separately

#### upholstered tables

#### raised-rim platforms



#### mats for raised rim platforms

wood finish is natural

mat colors

yellow

white

black

natural

navy

red

green

camel

- 2" polyurethane mat is upholstered with vinyl coated nylon for patient comfort
- fits snugly on raised rim platform



3x6 feet
3x7 feet
4x7 feet
5x7 feet
6x8 feet

#### height adjustable Hi-Lo mat platforms and tables

- 2" polyurethane mat is upholstered with vinyl coated nylon for patient comfort
- top has raised-rim to hold mat in place
- optional 2" mat sold separately

#### manual tables

#### manual platforms



15-2016	3x7 feet	15-2046	3x7 feet
15-2017	4x7 feet	15-2047	4x7 feet
15-2018	5x7 feet	15-2048	5x7 feet
15-2019	6x8 feet	15-2049	6x8 feet

#### wall mounted tables


up

## economy electric tables

#### economy electric platforms







15-2026	3x7 feet	15-2056	3x7 feet
15-2027	4x7 feet	15-2057	4x7 feet
15-2028	5x7 feet	15-2058	5x7 feet
15-2029	6x8 feet	15-2059	6x8 feet



#### deluxe electric tables

#### deluxe electric platforms



- heavy-duty, long lasting motor •
- foot pedal activator is fully electric

15-2082	4x7 feet	15-2072	4x7 feet
15-2083	5x7 feet	15-2073	5x7 feet
15-2084	6x8 feet	15-2074	6x8 feet

- space saving mat table
- when not in use, table can be folded and "locked" into place against the wall
- extends only 8" into the room leaving valuable floor space for other uses
- 2 " mat surface is filled with high-density polyurethane and is upholstered with heavy-duty vinyl coated nylon

15-2060	3x6 feet, 21" high
15-2062	4x7 feet, 21" high
15-2063	5x7 feet 21" high



## treatment table positioning

Skillbuilders® premium rolls and wedges have integral waterproof polyurethane coating that is easily washed for years of use.


#### Skillbuilders® wedges



- · positioning shapes for exercise therapy
- use under the back, arm or leg of a patient for comfortable and functional positioning
- use to provide the necessary leverage to promote an exercise regime; i.e. under the knee to enable weighted ankle lifts

width x length x height		incline	
30-1010	20x22x4"	small	11°
30-1011	20x22x6"	small	16°
30-1012	20x22x8"	small	20°
30-1013	20x22x10"	small	25°
30-1014	24x26x6"	medium	
30-1015	24x26x8"	medium	13°
30-1016	24x26x10"	medium	า 17°
30-1017	24x26x12"	medium	1 22°

#### Skillbuilders® rolls



used to foster exercise therapy and general positioning

length x diameter		
30-1000	24x4"	small
30-1001	24x6"	small
30-1002	24x8"	small
30-1003	36x8"	medium
30-1004	36x10"	medium
30-1005	36x12"	medium
30-1006	48x12"	large
30-1007	48x14"	large
30-1008	48x16"	large


#### Skillbuilders® roll rack





#### Skillbuilders® support wedge





wood

quad boards


- used on a treatment or mat table or on the floor
- supports resistance weight exercises
- folds flat when not in use

#### steel

- each section can be adjusted independently to match exact femoral and tibial lengths
- · chrome plated steel, folds flat

10-1140 wood board 10-1141 wood board, padded 10-1142 steel board






- supports knees and lower legs when user is in supine position
- wedge is useful for stretching and lower extremity exercises

		WxLxH
0-1250	wedge	22x13x17
0-1251	crescent	22x13x9"

3

#### Skillbuilders® half rolls



- combine stability of a wedge with the positioning versatility of a roll
- use velcro hook on flat bottom to anchor roll for additional security
- stable base for knee flexion, extension

wid	th x length	elevation
30-1230	4x18"	2"
30-1231	4x18"	3"
30-1232	6x24"	41/2
30-1233	8x30"	6"

#### face cradle pillow





- · comfort when patient is in prone position
- elevated base and nose cutout allow patient to breathe freely

00-4206 Pron-Pillo®

#### Chattanooga® Bolsters



- washable, heavy-duty integral foam bolsters are ideal positioning products for almost any use
- strong enough to support exercise modalities, comfortable to use for traction and versatile enough for many other uses

04-2901 knee crescent 10x10x12" (2ea.)



## pillows, supports and wedges

#### foam positioning and support rolls



- small profile makes it easy to achieve proper neck or back support without extra bulk
- · roll fits inside pillowcase to restore cervical curve
- precision-cut foam in a removable, washable cotton/poly cover

50-1210 roll, 19"L x 31/2"D 50-1210-25 roll, 19"L x 31/2"D 25 ea.

additional roll cover for 19"L x 31/2"D roll 50-1211

50-1215 roll, 19"L x 5"D 50-1215-25 roll, 19"L x 5"D 25 ea.

lumbar back roll

additional roll cover for 19"L x 5"D roll 50-1216

#### cervical neck pillow



- back support
- · adds comfort to seat
- attaches and stays in place with elastic strap
- standard with washable navy poly/cotton cover

full round

50-1217 11x5'

half round (D-shaped) 11x3" 50-1218

- · rounded cervical pillow
- polyester filled
- non-removable poly/cotton
- removable cover available

50-1200 7x17" 50-1200-25 7x17", 25 ea.

adt'l. cover 50-1201 50-1201-25 adt'l. cover, 25

#### bed wedges

- covered with heavy-duty vinyl-coated nylon and filled with high compression urethane foam
- shapes are ideal for institutional and home use for positioning patients in bed and in a
- · Ideal for leg elevation



	(LxWxH)		(LxWxH)
31-2000	22"x20"x4"	31-2005	28"x24"x10"
31-2001	22"x20"x6"	31-2006	28"x24"x12"
31-2002	22"x20"x8"	31-2050	20"x30"x8"
31-2003	28"x24"x6"	31-2051	30"x30"x16"
31-2004	28"x24"x8"	31-2052	40"x30"x16"

#### lumbar foam positioning and support



- · small profile makes it easy to achieve proper neck or back support without extra bulk
- · precision-cut foam in a removable, washable cotton/poly cover
- color choice: white, blue, red and black (blue is default)

50-1212 foam w/ cover, 14"x13"x3" 50-1212-25 foam w/ cover, 14"x13"x3", 25 ea. 50-1213 foam w/ cover, 18"x13"x3" 50-1213-25 foam w/ cover, 18"x13"x3", 25 ea. 50-1214 memory foam w/ cover, 13"x14"x3' memory foam w/ cover, 13"x14"x3", 25 ea. 50-1214-25

#### sleeping pillows



· pillows made with high-resiliency fiber and blended cotton cover

#### pillows

standard (24x16") 00-4280 00-4281 std. gentle (24x16") 00-4282 mid-size std. (22x15") 00-4283 petite (19x12")

made with 100% polyester

00-4284 economy (24x16")

#### cervical support



- · cradles and supports head and neck in any sleeping position
- combines comfort and the therapeutic benefits of hot/cold therapy (17x11½")
- ice pack not included

00-4272 Hot/Cold pillow

#### traveller pillow



- compact size to fit in suitcase
- use while sitting (lobes down, neck in center) or lying down (lobes up or down)
- · optional pillow case

00-4285 18"x9"

00-4295 blue pillow case



## treatment room stools and lamps


#### mobile revolving stool



- height adjustable (18"-24") "screw-type" revolving stool with round footrest
- 2" hooded carpet casters
- 14" diameter, 4" thick padded seat
- 250 pound weight capacity

16-1102 stool, black 16-1100 specify color







- 3" padded seat
- 2" hooded carpet casters.
- 18" fixed height
- 250 pound weight capacity
- available with and without carpet caster

#### mobile

16-1602 mobile, gray 16-1600 specify color

#### stationary

16-1603 stationary, gray 16-1601 specify color



- therapy stool with round, 3" upholstered swivel seat allows easy access for the patient and to equipment
- gas spring for easy raising and lowering of stool (20"-24")

#### without back

07-7063	black
07-7060	blue
07-7061	burgundy
07-7062	grav

#### with back

black
blue
burgundy
gray

#### exam lamp, goose neck style





- non-skid ribbed rubber platform
- reinforced rubber tips
- steel welded silver vein construction
- 9" height x 11" depth x 14" width
- 300 pound weight limit
- 34" hand rail optional

16-1701 footstool 16-1701-2 footstool 2 each 16-1710 with hand rail




- 16" flexible goose neck arm adjust 360°
- goose neck style allows the lamp to be adjusted in any direction
- adjustable height (48" to 72")
- UL Listed, for use with incandescent light bulbs up to 60 watts

#### stationary weighted steel base

17-1101 3-prong plug

#### caster base

17-1121 3-prong plug



- square armrests make sitting down and getting up safe and easy
- soft, comfortable and water-resistant vinyl seat cover
- · adjustable foot rest
- · attractive solid maple construction
- 300 pound capacity
- 26" floor to seat,24" width, 41" height, 23.5" depth

#### hip high chair

16-1800 maple high chair



#### mirrors

#### Ultra-Safe™ glassless mirrors

- distortion-free image
- shatterproof "glassless" mirrors cannot crack, shatter or chip but can be punctured
- lightweight mirrors can be securely mounted to a wall or ceiling surface with just a small piece of self-adhesive velcro

#### caster base



 features metal stabilization legs with 2" carpet casters

19-1014 24x72" 19-1015 24x96" 19-1016 36x72" 19-1017 36x96" 19-1018 48x72" 19-1019 48x96" 19-1044 - 60x72" 19-1045 60x96"

#### portable stationary



· features metal stabilization legs with floor levelers

19-1002	2	16x48"	
19-1004	1	24x72"	
19-1005	5	24x96"	
19-1006	3	36x72"	
19-1003	3	36x72"	horiz
19-1007	7 🖵	36x96"	
19-1008	3🖵	48x72"	
19-1009		48x96"	
19-1042	2	60x72"	
19-1043	3₽	60x96"	

#### wall mount



• comes with self-adhesive mounting pads 16 1 2 1 2 1

19-1032	10140
19-1033	24x72"
19-1038	24x96"
19-1034	36x72"
19-1037🛶	36x96"
19-1035🗪	48x72"
19-1035 <b>□</b> 19-1036 <b>□</b>	
.0 .000-	48x96"
19-1036	48x96" 60x72"

10.1032

#### portable triple-panel



- free-standing triple panel needs no additional support
- each panel is attached with heavy-duty hinges

19-1051 16x48" \* 19-1052 24x72" \*

#### portable double-panel



19-1072 24x72" \*

#### glassless speech therapy mirrors





- 12x16" Ultra-Safe™ glassless mirrors are designed for use by speech therapists and their patients
- single panel portable mirror has a protective cover that opens to form a desk stand
- triple panel mirror allows the patient to view each sound formation from any angle

19-1060 portable 12x16" 19-1061 triple panel 12x16" \*

\* symbol - size of each individual panel

Truck Shipments - these mirrors ship by truck and may take 3 weeks to arrive

#### plate glass mirrors







- distortion-free 1/4" plate glass mirrors are mounted in hardwood frames
- each panel of the adult mirror measures 28x75" glass size, each panel of the pediatric mirror measures 22x601/2" caster base triple wall mount panel

#### caster base

• mounted on hardwood legs with 2" carpet casters

vertical

19-1101 28x75" 19-1111 26x60" horizontal

19-1102 28x75" 19-1112 26x60"

 mirrors may be hung vertically or horizon-

19-1100 28x75" 19-1110 22x58" • free-standing triple panel need no additional support

• each panel is attached with hinges.

19-1103 28x75"\*

19-1113 22x60" \*



## tilt tables - upholstered

- tilt tables are built to meet the demands of continuous use in busy therapy departments
- top is 28" wide x 78" long with 1" vinyl covered padding
- metal walk off foot plate, 14"x171/2" with safety tread
- three wide safety straps
- special 4" easy roll, shock absorbing, semi pneumatic casters, with individual step on brakes
- welded steel construction
- black vinyl is standard on 15-3041, 15-3042, 15-3043
- space blue vinyl is standard on 15-3040

choose upholstery color from the following color list:



#### economy manual

#### professional manual

#### professional electric

## professional Hi-Lo electric



- affordable manual unit designed for home use
- adjusts from horizontal to vertical
- horizontal height is 32" off ground
- black metal frame, blue vinyl is standard

15-3040 economy



- tilts from full vertical to 9 degrees below the horizontal
- horizontal height is 31" off ground
- adjustment handle is retractable
- chrome plated frame, black vinyl is standard

15-3041 manual



- powered by quiet, electric motor that assures smooth, vibration-free adjustment
- adjusts from vertical to 9 degrees below the horizontal
- pneumatic hand-held control provides safety and control on either side of the table and comes with 7 ft. cord
- horizontal height is 31" off ground
- UL hospital grade plug, 10 ft. long power cord
- chrome plated frame, black vinyl is standard

15-3042 electric



- horizontal height adjusts from 26" to 38"
- powered by two heavy duty electric motors
- single pneumatic hand control provides safety and control for both vertical and horizontal positions and comes with an 7 ft. cord
- adjusts from vertical to 9 degrees below the horizontal
- UL hospital grade plug, 10 ft. long power cord
- chrome plated frame, black vinyl is standard

15-3043 electric hi-lo




#### work tables

## SKILLBUILDERS

#### work and O.T. tables

#### circular Hi-Lo work table

#### rectangular Hi-Lo work table



- foot operated hydraulic pump raises 11/4" butcher-block top from 28"-36"
- large area underneath table provides enough room for four wheelchairs

#### circular

15-3300	36" diameter, 23" base	
15-3304	42" diameter, 23" base	
15-3301	48" diameter, 27" base	
15-3302	60" diameter, 27" base	

- foot operated hydraulic pump raises 1¼" butcher-block top from 28"-36"
- H-frame base provides stability
- worktable is wheelchair accessible

#### rectangular

15-3310	30"x48" H-frame
15-3311	30"x60" H-frame
15-3313	36"x48" H-frame
15-3314	36"x60" H-frame
15-3316	48"x60"H-frame





 heavy duty work tables; O.T style tables have apron cutouts for wheelchair

O.T.	work	WxL (in.)
15-3260	15-3290	30x30
15-3261	15-3291	30x48
15-3262	15-3292	30x60
15-3263	15-3293	30x72
15-3264	15-3294	30x96
15-3281	15-3271	36x48
15-3282	15-3272	36x60
15-3283	15-3273	36x72
15-3284	15-3274	36x96





#### economy Hi-Lo work tables

- accommodates wheelchairs on all sides
- manual hydraulic tables adjust 28-351/2"
- built in floor levelers
- black enamel steel constructed frame

model 15-3250 work table: english oak laminate 15-3250 - 48"W x 60"L x  $1\frac{1}{4}$ " thick

model 15-3251 work table: butcher block laminate 15-3251  $48\text{"W} \times 96\text{"L} \times 1\%\text{"}$  thick



#### pedestal Hi-Lo work tables

- adjustable work table powered by electro-mechanical drive
- 21/4" maple butcher block
- · heavy duty black and chrome steel base
- · foot control is standard, hand control is optional
- operating range: 25-42"
- 110V 60Hz 4.9 amp 520 watt
- thermal overload protection
- · hospital grade plug
- power cord 120" long

15-3255 rectangular table 48x60x2<sup>1</sup>/<sub>4</sub>" 15-3256 circular table 60" diameter



## ambulation training

#### parallel bars

- space saving and adjustable; folding parallel bars with 42" base require only 14" of floor space when not in use
- . wood base unit comes with one stabilizing bar to lock parallel bars in open position; steel base unit comes with two stabilizing bars
- to store the space-saving bars, simply unlock and slide until handrails are nearly touching; add 3' to length when folded

#### floor base



#### wood base:

7'W x 42"L x 11/4" thick height adjusts 22-36" width adjusts 16-24" handrail size is 11/2" 15-4000 7-foot 15-4001 10-foot

#### steel base:





#### wall mount

- · adult size wall mount bars solve space problems by folding against a wall
- single plunger/lock keeps bars steady when in use
- · wood base plate included
- wood base measures 5" wide x 42" long x 3/4" thick; height adjusts from 28-42"; length 7'; handrail size 1½"; width between handrails 221/2"

15-4010 7' wall mount

## training stairs

steel base shown folded



- · straight staircase with platform enclosed on three sides
- overall length 55"; overall height 54"; height to platform is 24"; width is 30"; 4 steps; riser size is 6x10"; platform size is 30x24"; comes complete with bus step

15-4200 training stairs

#### training stairs



- · straight staircase for departments with low ceilings
- 2 and 3 step sides, bus step
- overall length 55"; overall height 40"; height to platform 12" width 30"; 2-step: 6x10" riser, 3-step: 4x10" riser; 30x24" platform size


15-4202 training stairs

#### training stairs



- · straight staircase
- overall length 96"; overall height 593/4"; height to platform 24"; width 36"; 8 steps; riser size 6x10"; platform size

15-4203 training stairs




#### wheelchair training aids



15-4220

- wheelchair training stools for teaching to get in and out of wheelchairs without assistance
- secure together in any sequence
- top size 20x103/4": height of stools 20,16,12,8,4 inches

15-4221 training stool set



15-4220 climbing stool set



- · curb and ramp training set
- · 3 curbs with handles can be nested and stored under ramp platform
- 44"x84" platform; 3/4"-8" incline; curb sizes 40x40x6", 36x36x4", and 32x32x2"; with anti-slip tread

15-4210 curb/ramp training set

#### convertible stairs





(above) stairs are easily converted

- 2 sections which lock together and can be converted from a straight line to right angle
- anti-slip treads
- handrails 18½, 27, 35" high
- overall height 54"; height to platform 24". Long section 5'10"

15-4204 24x24" platform; short section 4'6"

**15-4201**: 30x30" platform; short section 5'

15-4205 36x36" platform; short section 5'6'





## ambulation training


## professional Hi-Lo electric parallel bars platform mounted



- smooth, quiet and fast height adjustment with convenient hand-held control unit
- 28" width permits wheelchair access
- height automatically adjusts from 29-44"; width manually adjusts from 15-28"; stainless steel handrail size 1½"; Platform 10' or 15' x 38"
- 110V 60Hz 2.5 amp control box is double insulated; UL hospital grade plug

15-4140 10-foot long 15-4142 15-foot long

## adjustable height and width floor/platform mounted bars



 8" width adjustment (17½" to 25½" on platform mounted bars, height adjusts 31" to 41")

#### floor mounted bars

15-4090 7' long 15-4091 10' long 15-4092 12' long

#### platform mounted bars

15-4080 7' long, 9'x49" platform 15-4081 10' long, 12'x49" platform 15-4082 12' long, 14'x49" platform

#### adjustable height floor/ platform mounted bars



height adjusts from 31-41"; width between handrails is 22½" on platform mounted bars

#### floor mounted bars

15-4150 7' long 15-4151 10' long 15-4152 12' long

#### platform mounted bars

15-4040 7' long, 9'x47" platform 15-4041 10' long, 12'x47" platform 15-4042 12' long, 14'x47" platform

#### platform mounted parallel bars



- height and width adjustable height adjusts 26-44"; width adjusts 19-26"
- width permits wheelchair access
- platform 10 foot or 12 foot x 42"

15-4050 10' handrail; ht, width adj. 15-4051 12' handrail; ht, width adj.

Accessories for all platform mounted units:

15-4057 10' child railing

15-4058 12' child railing

15-4037 10' divider board

15-4038 12' divider board

15-4039 chrome platform hole cap (2)

#### platform mounted parallel bars



- height adjustable height adjusts 26-44", handrail width 25½"
- platform 10'/12' x 30"
  15-4030 10' handrail; ht adjustable
  15-4031 12' handrail; ht adjustable





#### floor mount parallel bars



- height and width adjustable
- economic parallel bars are ideal where space allows a permanent installation
- width and height adjustments accommodate school aged through adult sizes
- features: 1½" diameter steel, chrome plated, uninterrupted handrail

15-4060 10' length 15-4061 12' length 15-4062 14' length 15-4063 16'length 15-4064 18' length 15-4065 20' length 15-4066 22' length 15-4067 24' length

folds easily

Nimbo™ Roller

31-3650

31-3651

31-3652

31-3653

31-3730

31-3731 31-3732

31-3733

31-3655

31-3656

31-3657 31-3658

accessories 31-3654

#### pediatric walkers and rollers


#### Nimbo™ walker

#### Glider™ walker

#### Kaye™ posture control walker







- model; 45° handle on pediatric model • 5" flat-free tires allow both indoor and outdoor use, and tool free adjustments
- wide base provides stability and balance
- weight limit 175 lbs.

	height	red	blue
pediatric	12-21"	31-3010R	31-3010B
junior	22-33"	31-3012R	31-3012B
standard	29-37"	31-3014R	31-3014B

accessories	
platform attachment -ped	31-3020
platform attachment- jr	31-3021
platform attachment- std	31-3022
uni-directional front wheel- ped	31-3023
uni-directional front wheel- jr	31-3024
uni-directional front wheel- std	31-3025
45° handle, medium- junior	31-3027
45° handle, large- std	31-3028



- · 4-wheeled walkers, front swivel wheels
- · folds for easy storage
- pictured with seat
- · height measured to top of handle (floor to mid-buttocks)

walker	neight	
tyke	14.5-18.5"	31-3670
junior	16.5-22"	31-3671
youth	19-25"	31-3672
adolescent	23-30.5"	31-3673
young adult	28.5-36"	31-3674
walker with	seat	
tyke	14.5-18.5"	31-3690
junior	16.5-22"	31-3691
youth	19-25"	31-3692
adolescent	23-30.5"	31-3693
young adult	28.5-36"	31-3694

#### PVC platform walker

young adult seat harness

tyke, goldenrod yellow

youth, cornflower blue

tyke, goldenrod yellow junior, lime green

youth, cornflower blue

pelvic stability attachment

young adult, blue ice

forearm platforms

mounting bracket tyke/youth seat harness

young adult, blue ice

junior, lime green

optional seat attachment only



- · adjustable height walker has five 4" casters for extra stability (2 are locking)
- walker includes adjustable height padded support (table) platform

31-3580 pediatric (22"-30") 31-3581 child (30"-42")

#### Cricket™ walker



- · training walker stimulates walking move-
- for children (2-8 years) with motor and balance dysfunction

31-3560 small 31-3561 large

#### mobile child walker



- padded body ring opens in rear for easy entry and exit
- used walking or sitting; removable saddle seat adjusts 10-15" from floor

31-3570 mobile walker



#### pediatric gait trainers

#### Comet<sup>™</sup> anterior gait trainer



- anterior (push) trainer is height adjustable and folds for convenient storage and transport
- speed and direction be controlled with the variable resistance tabs, swivel wheel locking brackets and the one-directional rear wheel
- standard features include: ankle prompts, seat harness, pelvic stabilizer and a posterior positioning bar for adjustable depth.

31-3660 tyke (ht 17-20", w 22") 31-3661 ped. (ht 22-28", w 23½") 31-3662 adult (ht 29-36", w 26.½")

optional accessories

31-3610, adj. forearm platforms 31-3618 chest harness

#### Star™ posterior gait trainer



- posterior (pull) trainer is height adjustable and folds for convenient storage and transport
- speed and direction can be controlled with the variable resistance tabs, swivel wheel locking brackets and the one-directional rear wheel
- standard features include: ankle prompts, seat harness, vertical handgrips, pelvic stabilizer and a posterior positioning bar for adjustable depth

31-3680 tyke (ht 17-20", w 22") 31-3681 ped. (ht 22-28", w 23½") 31-3682 adult (ht 29-36", w 26½") optional accessories

31-3610 Adj. forearm platforms 31-3618 chest harness

31-3010 chest harness

#### Luminator™ gait trainer



- seat harness, pelvic stabilizer, positioning bar, vertical handgrips, variable resistance tabs, one-directional rear wheels, swivel wheel locking bracket
- · additional accessories available


#### **Anterior Gait Trainer**

31-3000	tyke	17-20"H
31-3001	junior	22-28"H
31-3002	youth	29-36"H

#### **Posterior Gait Trainer**

31-3003	tyke	17-20"H
31-3004	junior	22-28"H
31-3005	youth	29-36"H

#### forearm crutches



- vinyl coated, ergonomically contoured arm cuffs molded for comfort and stability; leg and forearm sections adjust independently for optimal sizing; vinyl hand grips are comfortable and durable; extra large tips for added stability
- sold in pairs

model			height	weight cap. lbs.	pair
43-2061	D		211/2-301/2	300	pan
43-2065	С		21½-30½	300	
43-2066	В	blue	14-21	250	
43-2067	Α	red	14-21	250	

#### pediatric folding rollers

#### anterior (push) safety rollers





- adjustability and modular components provide substantial built-in growth and allow for many years of use
- optional accessories: basket, adj. forearm platform, IV pole, oxygen tank holder, variable resistance tabs for rear wheels

	size	W (in)	H (in)	D (in)
31-3600	tyke	18	17-22	25
31-3601	ped.	20	22-28	25
31-3602	ped wB	221/2	24-30	25
31_3603	ned vi	221/2-271/2	23-20	25

#### posterior (push) safety rollers



- supports user from behind; requires the child to adopt a more upright posture and promotes trunk extension
- accessories (call for prices) include a basket, vertical handgrips, pelvic stabilizer, adjustable forearm platforms, posterior positioning bar, sling seat, ankle prompts, chest support, IV pole, oxygen tank holder and variable resistance tabs for the rear wheels

	size	W (in)	H (in)	D (in)
31-3630	tyke	21	17-20	22
31-3631	ped	221/2	22-28	25
31-3632	adult	27	29-36	27

curved

## canes and crutches

#### cane clamp





· attaches anywhere with adhesive backed hook and loop

cane clip

- accommodates canes up to 1" (2.5cm) in diameter
- accommodates canes up to 1" (2.5cm) in diameter

43-1160 cane clamp

quad

43-1161 cane clip, 2 ea.



- attaches with plastic screw to surfaces up to 11/8" (2.86cm)





- small plastic fitting mounts on cane to hold it on a table edge when not in use
- best when used in pairs to contact both bottom and top of table
- small rubber band (included) placed around holder end maintains a secure grip on cane shaft

43-1162

cane holder, 2 ea.

#### each

• canes are made of sturdy 1" diameter metal tubing

curved handle silver 43-2000 43-2000-6 offset handle silver vein 43-2010 43-2010-6 43-2012 43-2012-6 black 43-2014 43-2014-6 bronze

handle is height adjustable 29" to 38" (folding cane 33" to 37")

height adjustable canes

folding

offset

43-2020 folding black

		quad canes	
small base	chrome	43-2030	43-2030-4
large base	chrome	43-2040	43-2040-4

#### height adjustable crutches



- aluminum height adjustable crutches come in 3 sizes: adult 5'0" to 6'2"; tall adult 5'10" to 6'6"; youth 4'2" to 5'2"
- leg, forearm sections of forearm crutch adjust independently

#### underarm crutch case (8 pairs) 43-2050 43-2050-8 43-2052 43-2052-8

43-2054-8

pair

43-2060 adult youth 43-2061 43-2062 tall adult

43-2054

adult

youth tall adult

#### cane and crutch accessories

cane tip (pair)	43-2070
crutch tip, forearm (pair)	43-2071
crutch tip, underarm (pair)	43-2072
crutch handgrip, underarm open (pair)	43-2075
crutch handgrip, underarm closed (pair)	43-2076

#### cane hand loop



- no more dropping canes and walking sticks
- one end wraps around the cane/walking stick shaft while the larger loop, on the opposite end, wraps around the wrist
- made of strong webbing

43-1163

cane hand loop



## walkers and rollators

#### height adjustable folding walkers



- press the buttons to fold for easy transport or storage
- available with rubber tips or front wheels and rear glider tips
- height adjustable: adult 32-39"; junior 25-32"; oversize 32-39"

	each	case - qty.
adult, no wheels	43-2100	43-2100-4
junior, no wheels	43-2102	43-2102-4
oversize, no wheels	43-2104	43-2104-2
adult, 3 wheels	43-2110	43-2110-4
junior, 3 wheels	43-2112	43-2112-4
adult, 5 wheels	43-2120	43-2120-4
junior, 5 wheels	43-2122	43-2122-4



- 6" casters for indoor/outdoor use
- fold for easy transport or storage
- 4-wheel model has a seat for resting
- available with push or loop braking system

4-wheel rollator, blue 4-wheel rollator, red	<b>with loop brake</b> 43-2150 43-2152	with push brake 43-2140 43-2142
3-wheel rollator, blue 3-wheel rollator, red	43-2160 43-2164	

## optional basket 4-wheel rollator 43-2170 3-wheel rollator 43-2171



## platform attachment

#### walker platform attachment

- provides the ability to bear weight on the arm and relieve pressure for those with hand injuries or arthritis
- lightweight aluminum. No tools needed to attach to a walker
- may be mounted on inside or outside of walker side frame
- same unit can be used on either left or right side of walker

43-2130 platform attachment 2 ea.





- easy-to-install
- pre-cut to fit most walkers
- protects floor from scuff marks

43-2136 walker balls 2 ea.

#### adjustable glides



- easy-to-install
- for use with all 1" folding walkers
- durable composite construction ensures long life

43-2135 glides 2 ea.



#### walker rest seat



- comfortable, fold-down padded seat
- two inch foam padded vinyl cover cleans easily
- two hook and loop fasteners straps secure the seat to the uprights or horizontal rails

43-2180 rest seat

#### cane holder

- double-clip holder allows user to keep cane handy
- snap one side of the clip onto tubular structure such as a table leg, walker, commode, wheelchair, or lawn chair
- opposite side of clip will hold any cane
- soft plastic gripping rollers will not damage canes
- holder dims: 3x2x1¾" (7.6x5.1x4.4 cm)
   43-1164 cane holder



#### reacher clip

- keeps reachers handy and firmly attached to walker or wheelchair frame
- larger side of double clip snaps onto wheelchair, walker frame or bed rail, smaller side onto the reacher

43-2181 reacher clip

#### transport chair



- 8" composite wheels in front and rear
- · light, strong and durable
- rear wheel locks
- swing-away footrests
- back folds for storage and transport
- · seat belt for added safety
- · carry pocket on back of rest
- weighs 19 lbs: capacity 300 lbs

blue frame, black upholstery\*
17" seat 43-2200
19" seat 43-2210

\* other colors and styles available; please inquire

#### Silver Sport 2 - dual axle wheelchair

- powder coated steel frame is lightweight, strong and durable
- padded armrests, plastic footplates
- vinyl upholstery is durable, and easy to clean
- urethane tires on composite wheel; push to lock brakes
- 8" front casters, adjustable in 3 positions
- carry pocket on backrest provides additional convenience

43-2260

- swing-away footrests or swing-away elevating legrests with adjustable length riggings
- fixed or detachable desk arms
- · capacity: 250 lbs

18" seat





detachable arm
swing-away swing-away+elevating
footrest footrest
16" seat 43-2230 43-2231

#### oxygen tank holder for wheelchairs

 attach to handles and frame; tank position does not interfere with user's head or back

• heavy duty, waterproof nylon



Mini tank holder

 2 pockets hold 2 M6 tanks or supplies

43-2280 mini holder

#### leg lift

 helps lift leg or foot onto wheelchair footrest

43-2295 leg lift

 webbing with loop at each end for hand and foot

Standard tank holder

 fits size "D" and "E" oxygen tanks

43-2281 std. tank holder

#### laptop wheelchair desk

43-2261



- ideal for laptop computer users
- storage area below lid

43-2290 24x22"

#### cup holder with snack tray

- double-well cup holder
- attaches at any angle to tubular frame from 1/2 to 1½" diameter
- deeper well has pop-out bottom to accept glasses and cups

43-2285 cup holder 8"x6"x2", 31/8" well



## wheelchairs

#### bean bag lap desk



- plastic top with removable "styrofoam" bean bag base
- bean bag provides form fitting and comfortable support for the surface on your lap

43-2291 lap desk, 12x16"

## lightweight, molded rimmed wheelchair tray



- ideal for eating, writing, working and playing games
- smooth flat top with raised "spill" edge all around
- body/rim cut-out is secured to wheelchair arms with straps
- ultra-rigid plastic tray maintains straightness, doesn't sag at corners

43-2292 tray 24x20x1"

## wheelchair drink holder, clamp on



- slotted on both sides to hold a mug or cup with handles
- fits bottles and cans
- clamp has easy-to-tighten lever

43-2286 drink holder



## wheelchair cushions and accessories

#### wheelchair accessories

#### transfer boards



 tapered hardwood boards for those with lower extremity disabilities

#### wood

50-3000 8x24" no hole 50-3001 8x30" no hole

50-3004 8x24" 2 hole 50-3005 8x30" 2 hole

#### plastic

50-3007 8x24" no hole

#### wheelchair trays



- raised rim, cut-out design, attaches with velcro straps
- Skillbuilders<sup>®</sup> tray has non-slip
   fits into all standard wheelfinish

#### wood finish

50-1300 economy 24x20" 50-1301 standard 24x20"

#### clear acrylic

50-1302 w/ rim, straps

#### Skillbuilders® arm elevator



- coated arm elevator can be kept sanitary
- chair arm support trays

#### Skillbuilders® knee spreader



- coated knee spreader can be kept sanitary
- comfortably maintains 4" of knee separation to prevent hip abduction
- · impervious to liquid
- hook and loop straps hold wedge in place

30-1411 knee spreader

#### Amara wheelchair cushions with removable anti-microbial black cover

Amara 300

#### Amara 100

## general use:

- for low-moderate risk of skin breakdown
- minimal postural support

50-1800 16x16x3" 50-1801 18x18x3" 50-1802 20x18x3"

#### Amara 250



- for moderate/high risk skin breakdown; w/ foam insert
- minimal postural support

50-1820 16x16x3" 50-1821 18x18x3" 50-1822 20x18x3'

#### Amara 400



- · for moderate-hi risk of skin breakdown; w/ gel insert
- significant postural support

50-1840 16x16x3" 50-1841 18x18x3" 50-1842 20x18x3"

30-1410 arm elevator

#### Amara solid sea inserts



- · for moderate-hi risk of skin breakdown, w/ foam insert
- moderate postural support

50-1810 16x16x3" 50-1811 18x18x3" 50-1812 20x18x3"

#### · for moderate-hi risk of skin break down; w/ gel insert

· significant postural support

50-1830 16x16x3" 50-1831 18x18x3" 50-1832 20x18x3"



· wood insert helps eliminate sagging and hammocking of wheelchair sling seats

50-1850 16x16" 50-1851 18x18" 50-1852 20x20"

#### high seat cushion



- use to elevate seating level of normal chair to reduce stress on knees and legs
- · vinyl covered, wood base
- over the shoulder carry strap for portability
- available in black, burgundy,

specify color when ordering

50-1880 16x16x4"

#### economy foam cushions

 high density foam for seating comfort and a removable and washable cover



#### foam cushions

50-1325 16x18x2" 50-1326 16x18x3" 50-1327 16x18x4"

## economy gel cushions

- for decubitus ulcer prevention, control
- gel for displacement of pressure
- · removable cover

#### gel cushions

50-1360 16x18x2" 50-1361 16x20x2"



#### economy gel/foam cushions

- · for decubitus ulcer prevention, control
- · gel for displacement of pressure
- · foam for added comfort
- removable cover

#### gel/foam cushions 50-1362 16x18x2'

50-1363 16x20x2'



## bathroom products - toilet

#### height adjustable commodes







- commodes come with commode bucket, cover and splash shield
- made with steel or aluminum tubing (specified below)
- drop-arm commode has easy-release for safe lateral transfer to and from the commode
- commodes have 300 pound capacity, extra-wide bariatric commode can support up to 650 pounds

#### height adjustable (18-22"), folding [a]

43-2330 steel, plastic armrests

43-2330-4 steel. 4 ea

43-2332 steel extra wide (650 lb cap) 43-2332-2 steel extra wide, 2 ea

aluminum, padded armrests 43-2336-4 aluminum, 4 ea

43-2336

#### drop arm, height adjustable (18-24"), [b] with wheels

43-2350 steel, padded arm 43-2350-2 steel, padded arm, 2 ea


#### drop arm, height adjustable (19-23") [c]

43-2340 KD, steel, standard seat 43-2340-2 KD, steel, standard seat, 2ea

43-2342 KD, steel,padded seat

#### Columbia<sup>®</sup> toilet support





optional footrest/armrest sold separately

#### low back (safety belt & reducer ring)

45-2220 small

45-2221 medium

large

#### high back (H-brace & reducer ring)

45-2230 small

45-2231 medium

45-2232 large

#### accessories

45-2223 footrest/armrest 45-2260 reducer ring

#### Columbia® positioning commode



 freestanding commode chair includes bowl, footrest, armrest and wrap-around height-adjustable back support

#### low back (safety belt)

45-2240 small

45-2241 medium

45-2242 large

#### high back (H-harness)

45-2250 small

45-2251 medium 45-2252 large

accessory

45-2257 caster set (not included)

#### TumbleForms2® pediatric mobile potty seat



- · use as freestanding commode or position over toilet
- · positioning features of the Carrie seat plus utilitarian features of commode chair
- includes seat, adjustable neck support, safety harness, adjustable footrest and tray, mobile frame, splash guard and commode bucket

30-3380 medium (elementary)

30-3381 large (junior)

x-large (small-adult) 30-3382



## bathroom products - toilet seats

#### brackets



#### Slip-In bracket

- slides under toilet seat hinge
- no tools required

43-2580 Slip-In



#### bolt-down bracket

- · use toilet seat hinge bolts
- unscrew toilet seat hinge nuts to slide in bracket

43-2581 Bolt-Down

#### economy elevated seat



- elevates seat by 4"
- wider than typical donut, slightly contoured
- · hand grips for portability

43-2520 economy 43-2521 Slip-In 43-2522 Bolt-Down

#### standard elevated



- elevates seat by 5"
- openings in front and back for left side sloped easy access
- supports up to 300 lbs.
- easy to install without tools

43-2530 standard

#### Arthro® Tall-Ette seats





- elevates seat by 4"
- slope on one side for those requiring extended leg positioning

#### right side sloped

43-2540R standard Slip-In 43-2541R elongated

**Bolt-Down** 

43-2540L standard Slip-In 43-2541L elongated

Bolt-Down

contoured elevated toilet seats

#### · provides for comfort and support

· available in three heights

standard toilets: all heights come with a Slip-In bracket for fast and safe installation

elongated toilets: all sizes come with the Bolt-Down bracket for a more permanent attachment

standard seat

2" Slip-In

4" Slip-In

6" Slip-In

ALC:			2
		-8	


standard 6 inch

-	100	
Sections		

	elongated seat
43-2510	2" Bolt Down
43-2511	4" Bolt Down
43-2512	6" Bolt Down

#### hinged elevated toilet seats

- seat elevates by 4"
- two part seat allows for normal raising and lowering of lid
- uses existing toilet seat and lid
- does not appreciably alter the appearance of toilet
- hinged seat allows for easy cleaning





43-2570 standard 43-2571 elongated

#### elevated seats with arms



43-2500

43-2501

43-2502

• elevates seat by 4"

- for use with regular toilet seat and lid
- width between arms is 19" arm height is
- · 300 lb. capacity
- hardware included

43-2550 standard 43-2551 elongated



- elevates seat by 5"
- openings in front and back for easy access
- width between arms is 19", arm height is 41/2"
- 300 lb. capacity
- no tools required

43-2555 Lock-On™ seat w/arms



#### extra wide elevated seats with arms

- elevates by 4"
- · arm rests and Bolt Down bracket
- 221/2" seat width, 18" width between arms
- · for standard and elongated seats

43-2560 with arms



- elevates by 4"
- arm rests and Bolt Down bracket
- 221/2" seat width, 18" width between arms
- for standard and elongated seats
- seat with aluminum legs for stability
- · seat slides up legs to adjust height

43-2565 with arms and legs

## bathroom products - bathing



- · seat and back angle adjust independently
- fits all standard tubs
- includes two adjustable positioning belts and headrest for safety, comfort, support
- back support has 10 angle positions...including trendelenberg; seat has 5 tilt positions
- · quick drying sea green mesh fabric covers aluminum frame
- seat folds for easy storage, and is adjustable in depth (13-18")
- length: 39.5-44.5", width: 17", height from floor: 6", weight capacity: 100 pounds
- optional height adjustable elevation base is easily attached to the Dolphin
- adjusts in 1" increments (11-18") for easy transfer; can be used in bathtub or shower

45-2290 bathing chair 45-2291 elevation base

# Otter® bathing chair

- · used by children with physical, neurological, mobility/balance disability or spinal cord injury
- upholstered seat has a folding plastic frame with removable machine-washable mesh fabric
- includes two positioning belts and adjustable head cushions
- · seat and back have five angle adjustments from 0-90°

small to 36" tall 45-2190 45-2191 medium 33-50" large 46-68" 45-2192 45-2195 tub stand (optional)

#### Columbia® reclining bath chair



- tilts from 30° to 70°. Gives full trunk and head support
- PVC frame, mesh fabric and rubber feet
- · head pad, leg extension set that raises chair and additional 9" or 15" and caster set are sold separately

#### bath chair (tilt-in space)

45-2200 small to 65 lbs. 45-2201 medium to 130 lbs. large to 150 lbs. 45-2202


#### bath chair (articulating)

45-1840 small 45-1841 medium 45-1842 large

#### accessories (not included)

45-2205 leg extension 9" leg extension 15" 45-2208 45-2206 head pad 45-2207 caster set

#### wrap-around support



· wrap around support gives children freedom while seated upright in bathtub, wading pool, or on the floor

#### low back (safety belt)

45-2210 small 45-2211 medium 45-2212 large

#### high back (H-harness)

45-1850 small 45-1851 medium 45-1852 large

#### bath bench



- height adjustable (15-20") bath bench has drainage holes in seat and back.
- bench, back provide comfort, strength
- lightweight aluminum frame is durable and corrosion proof
- choose fully assembled or knock down

without back assembled, folding 43-2400

with back, removable arms assembled, folding 43-2410

without back, knock/down

43-2402

43-2402-4 4 each

with back, knock/down

43-2412

43-2412-4 4 each



## bathroom products - bathing

#### adjustable shower seats





- adjustable legs 16-21" with suction feet
- removeable back and arms available
- built-in hand held shower spray holder
- · durable and rust-proof
- 300 lb. capacity

45-2300 no arms, no back 45-2301 with arms only 45-2302 with back only 45-2303 with arms and back


#### adjustable transfer bench



45-2340

- provides easier access for getting in and out of shower
- adjustable legs 16-21" with suction feet
- removable back and arms
- built-in hand held shower holder
- durable, rust proof
- 350 lb. capacity

#### corner shower seat



- safe seating for people who have difficulty standing
- fits into shower stall corner
- built-in soap holder
- 300 lb capacity

45-2310 corner seat

#### shower stools



- for people who have difficulty standing
- adjustable legs 16" to 21" with non-skid feet
- all plastic; won't rust
- rotating model swivels to give caregiver easy access to patient
- 300 lb. capacity

45-2320 round, 15½" seat dia. 45-2330 rotating, 16" seat dia.


#### grab bars

#### **Groovy Grip™**



- plastic, wall-mounted bars with grooved surface for secure non-slip grasp
- mounted brackets conceal screws
- distance from wall: 2", 11/2" bar diameter
- hardware and instructions included

45-2350 16" length 45-2351 18" length 45-2352 24" length







- plastic, wall-mounted bars with grooved surface for secure non-slip grasp
- wall clearance is 3", ¾" rod diameter
- hardware and instructions included

45-2360 single level, 10" length 45-2361 double level, 17" length 45-2362 triple level, 25" length

#### sponges



transfer bench



#### swiveling back scrubbers [a]:

- · rotating head scrubber helps to wash back and lower body
- for people with reduced range of motion

45-2370 334" round sponge, straight handle 45-2371 334" round sponge, curved handle 45-2375 replacement sponge, bag of 2

#### straight arm sponges [b]:

- extends reach to wash lower back and extremities
- sturdy 22" long plastic handles

45-2380 3" thick, round sponge, 5" diameter 45-2381 3" thick, rectangular sponge, 5"x3"



#### lotion applicator



[c]



[d]

#### lotion applicator [c]:

- apply lotions and medications to hard to reach places
- 12" angled handle is ribbed w/ thumb rest
- replacement sponges available

45-2390 lotion applicator 45-2395 replacement sponge

#### roll on applicator with folding handle [d]

- comes with two pivoting roller heads that massage while applying lotion
- inverted roller hugs body's curves; ideal for arms, legs and neck
- round roller is ideal for the back
- handle folds for storage and travel: 14½" long, 8" folded 45-2400 applicator w/ folding handle

#### sof-seam™ knee support

- sof-seam™ neoprene sleeve in 3mm and 4mm; available in a variety of styles to provide compression, comfort and warmth
- provide stability and still let the patient maintain daily activities while providing support and easing the discomfort
- royal blue nylon inner lining that reduces skin irritation and black nylon exterior for easy cleaning; fits both right and left knee



[A] \*
S: 13-14"
M: 14-15"
L: 15-16"
XL: 16-18"
XXL: 20-22"

\* measurement taken around patella



[B/C] \*
S: 12-14"
M: 14-16"
L: 16-18"
XL: 18-20"
XXL: 20-22"



## [A] basic knee support with open patella

 3mm neoprene support for sprains, strains, reduction of edema and arthritis of knee joint

24-2600	small
24-2601	medium
24-2602	large
24-2603	x-large
24-2604	xx-large

## [B] 13" sport knee support with anterior closure and functional hinge

 4mm neoprene support for preventative rehabilitation, long term management, ACL, PCL or combined instabilities

24-2640	small
24-2641	medium
24-2642	large
24-2643	x-large
24-2644	xx-large
24 2044	AX large

## [C] 17" sport knee support with anterior closure and functional hinge

 4mm neoprene support for preventative rehabilitation, long term management, ACL, PCL or combined instabilities

24-2634	small
24-2635	medium
24-2636	large
24-2637	x-large
24.2638	vy Jargo

#### PneuGel® wrist, arm and elbow support other supports available. please inquire.

- PneuGel bladder is made from an inflatable material that encapsulates a freezable gel for cold and compression therapy
- outer pneumatic bladder compartment can be inflated to support a comfortable and secure fit and to provide measured com-

24-4541 tennis elbow



24-4542 elbow wrap

pression to reduce swelling and edema

• fabric layer next to skin reduces the potential for skin irritation and condensation build-up during cold therapy; one size fits all.



specify left or right (L/R) 24-4540L/R wrist wrap



24-4543 ankle support

#### Sroufe® support, splint, and spica



#### 7" elastic wrist

 cotton elastic, removable and malleable palm stay w/ hook and loop closures

wrist circ.	specify left or right (L/R)
5-6.5"	24-4550L/R
6-7"	24-4551L/R
6.5-8"	24-4552L/R
7-9"	24-4553L/R
8.5-10"	24-4554L/R



8" soft wrist splint

 low cut in the palmer crease area w/ wide thumb hole to reduce irritation

wrist circ.	specify left or right (L/R)
5-6.5"	24-4570L/R
6-7"	24-4571L/R
6.5-8"	24-4572L/R
7-9"	24-4573L/R
8.5-10"	24-4574L/R



## 8" soft wrist splint with abducted thumb

 nylon outer fabric, soft moisture-wicking lining, breathable foam padding

wrist circ. specify left or right (L/R) 5-6.5" 24-4575L/R 6-7" 24-4576L/R 6.5-8" 24-4577L/R 7-9" 24-4578L/R

24-4579L/R

8.5-10"



#### universal thumb spica

- soft, supportive materials deliver gentle protective positioning of MP and CMC joints
- universal sizing fits both left and right hand

24-4545 universal



#### Aircast® Air-Stirrup® braces - the "standards of care" for the treatment of ankle injuries

#### Air Stirrup® ankle brace



- for functional management of ankle sprains
- anatomic shell for left or right leg, Duplex™ aircells

left right
10.5" large 24-2650L 24-2650R
9" medium 24-2651L 24-2651R
8¾" small 24-2652L 24-2652R
6" pediatric 24-2656L 24-2656R

#### AirHeel™



- applies intermittent pulsating compression to tendon and fascia to help minimize swelling, promote circulation
- for plantar fasciitis, achilles tendonitis, and heel pain
- fits left and right foot24-2701 small24-2702 medium24-2703 large

#### sport Stirrup®



- based on the Air-Stirrup ankle brace, but with a narrower design for less bulk
- one size fits all left right
  one size 24-2653L 24-2653R

#### AirSport® ankle brace



 combines Air-Stirrup features with "step-in" design and automatic heel

 Left
 right

 x-small
 24-2710L
 24-2710R

 small
 24-2711L
 24-2711R

 medium
 24-2712L
 24-2712R

 large
 24-2713L
 24-2713R

 x-large
 24-2714L
 24-2714R

#### ankle support sizing

\* fits both right and left ankle small: M 5.5-7, W 5.5-8.5 medium: M 7.5-11, W 9-12.5 large: M 11.5-13, W 13-14.5

#### leg brace



 management of stress and stable fractures of lower leg

left right 13" small 24-2668L 24-2668R 15.5" med. 24-2666L 24-2666R

 includes anterior panel for added tibial protection

left right 13" small 24-2669L 24-2669R 15.5" med. 24-2667L 24-2667R

#### pneumatic tennis-elbow



- concentrates compression on extensor muscle, not around the arm, providing more support, less constriction
- available in beige and black
- one size fits all

50-5500 one size beige 50-5510 one size black



#### Ankle Sprain Care Kit™ by AirCast® with Cando® band



- video and booklet instruct patient regarding injury care
- practitioners time and cost is reduced
- kit includes: Air-Stirrup ankle brace, ankle wrap, cold pack, Cando® exercise band (4foot strip, green, medium), Instruction video and booklet


 left
 right

 small
 24-2721L
 24-2721R

 medium
 24-2722L
 24-2722R

 large
 24-2723L
 24-2723R

#### drop foot brace



le I (3-5)\* 24-1600L

small (3-5)\* 24-1600L medium (6-8)\* 24-1601L large (9-11)\* 24-1602L right

24-1600R 24-1601R 24-1602R

#### Support™ AFO orthosis



- low arch and open heel fit any shoe
- plastic material, design give AFO support, flexibility
- provides static dorsiflexion assistance and lateral stability for entire ankle/foot area
- use a scissor to trim insole for exact shoe fit; may be modified with a heat gun

24-1606R 24-1607R 24-1606L 24-1607L

small / medium large / x-large small / medium large / x-large right (3-7)\* right (8-12)\* left (3-7)\* left (8-12)\*


\*men's shoe size




## SKILLBUILDERS vestibular training and play therapy swing seats for

Skillbuilders® swing seats are lightweight, durable, and fun!

#### swing seat with pommel




#### swing seat without pommel



seat small	L x W x H 24" x 8" x 11"	weight cap.
medium	30" x 12" x 14"	150 lb
large	35" x 15" x 16"	200 lb

#### complete seat with pommel



- plastic seat is available with and without a built-in abductor (pommel)
- safety harness and crotch strap standard for proper positioning
- · adjustable suspension mechanism to fit any frame up to ten feet
- optional headrest, legrest, seat cover and swing tie-down straps

	with rope	with chain
small	30-1630	30-1590
medium	30-1631	30-1591
large	30-1632	30-1592

	with rope	with chain
small	30-1634	30-1593
medium	30-1635	30-1594
large	30-1636	30-1595

- complete full support swing seat comes with pommel, head rest and leg rest
- · small is for child, medium for teenager, and large for adult

n

#### bouncer accessory



- bounce adapter adds a new-dimension of vestibular challenge to swinging
- · child can bounce while swinging back-andforth, or without swinging

#### HD (125 lbs. each)

heavy-duty, 1 each 30-1692 30-1692-2 heavy-duty, pair

#### rotator accessories



 rotator adapter adds a new dimension and vestibular challenge to swinging

30-1695 single attachment (not shown) 30-1696 double attachment 30-1697 quadruple attachment 30-1698 swivel rotator (not shown)

#### swing seat accessories





headrest 30-1633 small/medium 30-1629 large leg rest 30-1643 small/medium 30-1644 large covers for seat with pommel 30-1645 small 30-1646 medium 30-1647

replacement seat harnesses

large

30-1630H small 30-1631H medium 30-1632H large

#### tie-down retention straps


• anchor swing to ground to limit motion

30-1637 straps (2 each)



## swing seats for vestibular training and play therapy SKILLBUILDERS

swing while experiencing the comfort of lying down!



- offers the special needs child the chance to experience the exhilaration and fun of swinging in the park or at home
- fully adjustable safety harness secures shoulders and waist
- · contoured design provides head, body and leg support
- for children up to 150 lbs

#### seat

30-1684 full-body swing with rope 30-1686 full-body swing with 8 foot chain 30-1687 full-body swing with 10 foot chain



Vestibulator® II

- · all benefits of ceiling suspension system
- 360° vertical rotation system for simultaneous vertical/horizontal stimulation
- complete contains all accessories


30-3601 complete unit

30-3602 frame with rope and ascender

#### accessories

30-3603 platform swing 30-3604 flexidisc 30-3605 roll swing 30-3606 prone net swing 30-3607 net swing with seat 30-3610

30-3609 rope and ascender



- [A] galvanized steel pipe construction, 77"H x 83"W x 70"D
- [B] zinc plated harware and T-clamps, 96"H x 96"W x 96"D
- [C] high carbon steel frame with galvanized safety "J" hooks; comes with blow molded bench seat with ropes and ground anchors; ndoor leg conversion kit sold separately 85"H x 78"W x 68"D

#### frame only

30-1671 [A] portable indoor frame

30-1672 [B] indoor / outdoor frame

30-1670 [C] outdoor frame w/ board swing 30-1673 indoor conversion kit for 30-1670

#### Flaghouse seat

[D] comes with adjustable harness / seatbelt strap

30-1650 toddler 24x8x11" 30-1651 youth 30x12x14" 150 30-1652 large 35x15x16" 200

#### bench seat

30-1580 blow molded bench seat

#### strap seat

30-1585 strap seat

#### Skillbuilders® protective helmets





- protects the special needs person from the dangers of physical activity, potential strike hazards (floor, wall, furniture, etc.), selfabuse, seizures and spasticity
- washable soft-top helmet protects front, back and side of the head
- filled with shock absorbent foam material
- adjustable velcro chin-strap keeps helmet in place

32-2200	x-small	(17½-18½")	
00 0004		,	
32-2201	small	(19-20")	
00.0000		` ,	
32-2202	medium	(20½-21½")	
22 2222	1	(22.221)	
32-2203	large	(22-23")	
32-2204	v lorgo	(23½-24½")	
32-2204	x-large	(23/2-24/2)	

size measured around temple



medium

## adaptive positioning wedges

- · adaptive positioning provides stability and support for individuals lacking some degree of body control
- prerequisite for achieving greater gross and fine motor skills
- wedges perform equally well for prone and supine activities
- choose the length and height of the wedge by considering what weight should be borne on the forearms and what extension the head, neck and trunk should elicit



- · versatile shape for gross motor activities like tumbling, rolling or
- ideal positioning shape for physically involved children
- suitable as physical therapy positioning wedges
- used with optional abductor wedge (30-1018)
- abductor wedge can be used with any Skillbuilders® or Tumbleforms® wedge or other positioning products
- hook on the abductor wedge can mate with any velcro loop material

	$W \times L \times H$		incline
30-1010	20x22x4"	small	11°
30-1011	20x22x6"	small	16°
30-1012	20x22x8"	small	20°
30-1013	20x22x10"	small	25°
30-1014	24x26x6"	medium	9°
30-1015	24x26x8"	medium	13°
30-1016	24x26x10"	medium	17°
30-1017	24x26x12"	medium	22°

#### 5x5x4" abductor

30-1018

#### Skillbuilders® strap wedge



- 5" wide velcro positioning straps provide support when positioning the child
- without the removable straps, the deluxe wedge functions just like the standard Skillbuilders® wedge
- · used with optional leg abductor wedge (30-1018)

30-1020 30-1021 30-1022 30-1023	W x L x H 20x22x4" 20x22x6" 20x22x8" 20x22x10"	small small small small	incline 11° 16° 20° 25°
30-1024	24x26x6"	medium	
30-1025	24x26x8"	medium	
30-1026	24x26x10"	medium	
30-1027	24x26x12"	medium	

#### Tumbleforms® wedges

abductor wedge



- · for positioning and gross motor activities
- optional abductor comes with velcro

	$W \times L \times H$		incline
30-3010	20x22x4"	small	11°
30-3011	20x22x6"	small	16°
30-3012	20x22x8"	small	20°
30-3013	20x22x10"	small	25°
30-3014	24x26x6"	medium	9°
30-3015	24x26x8"	medium	13°
30-3016	24x26x10"	medium	17°
30-3017	24x26x12"	medium	22°
30-3018	abductor		

#### Valu-Form® wedges



 Valu-Form® wedges are a cost effective line of positioning products

	$H \times W \times L$		incline
30-3180	4x20x22"	small	11°
30-3181	6x20x22"	small	16°
30-3182	8x20x22"	small	20°
30-3183	10x20x22"	small	25°
30-3184	6x24x28"	medium	9°
30-3185	8x24x28"	medium	13°
30-3186	10x24x28"	medium	17°
30-3187	12x24x28"	medium	22°

#### Tumbleforms® strap wedges



 removable positioning straps provide flexibility and support

	$W \times L \times H$		incline
30-3020	20x22x4"	small	11°
30-3021	20x22x6"	small	16°
30-3025	24x26x8"	medium	13°
30-3026	24x26x10"	medium	17°
30-3027	24x26x12"	medium	22°
30-3018	abductor		



## adaptive positioning wedges

#### Skillbuilders® bi-height wedge



- position at a higher or lower height (use on side 1 or side 2)
- · high molded sides give lateral support

hi-height

- positioning strap provides support
- Accommodates abductor wedge (30-1018)

	V V A L A	ויסוו-ווע	giit
		(side 1/	side 2)
30-1240	18x22 x	4"	6"
30-1241	18x22 x	8"	10"
30-1242	22x26 x	6"	8"
30-1243	22x26 x	10"	12"

## Skillbuilders® adjustable wedge system



- users can freely move arms and hands while lying on this adjustable wedge
- both the front and rear legs can be adjusted to the proper height with single button locks. 22" wide x 36" long
- adjustable: front 2-15", rear adjust 2-7"
- removable abductor adjusts 4-8"
- side rail for carrying, securing positioning straps
- wedge folds flat
- vinyl covered foam padded top

31-1050 adjustable wedge 22x36"

## Skillbuilders® therapist kit; wedge, roll and half roll



- therapy kits offer Skillbuilders positioning shapes including wedges, rolls and half rolls in one convenient package
- wedge (LxWxH") W
  roll (length x diameter) R
  half roll (width x length + elevation) HR

kit	contents	
30-0990	30-1011 <b>W</b> 30-1002 <b>R</b> 30-1232 <b>HR</b>	
30-0991	30-1017 <b>W</b> 30-1004 <b>R</b> 30-1233 <b>HR</b>	
30-0992	30-1011, 30-1016 <b>W</b> 30-1002, 30-1008 <b>R</b> 30-1233 <b>HR</b>	

#### Tumbleforms® Bi-Form® wedge




- position at a higher or lower height
- high molded sides give lateral support

WxL x bi-height

- positioning strap provides support
- accommodates abductor

		(side 1	/ side
30-3240	18x22x	4"	6"
30-3241	18x22x	8"	10"
30-3242	22x26x	6"	8"
30-3243	22x26x	10"	12"

#### Tumbleforms® Thera-wedge®



- large 2 sided wedge: one smooth, one contoured with high molded sides for passive support ("nestles" patient into proper alignment)
- smaller wedge "interlocks" with larger wedge
- straps provide support
- large 23x38x10"/16"; small 23x38x12" and abductor

30-3210 dual wedges

#### Tumbleforms® slotted wedge



- wedges offer a variety of positioning options for children and adults
- wedges help alignment of the hip, preventing rotation and providing adduction of lower extremities
- promotes comfort and optimal positioning during treatment sessions
- turn the wedge over, and the leg troughs offer support to lower extremities

30-3250 S 24x14x9" 30-3251 M25x14x12" 30-3252 L 26x16x13" 30-3253 XL 29x19x16"




## adaptive positioning rolls

#### Skillbuilders® rolls





- · adaptive shape for rolling and rocking movements, coordination and vestibular training and positioning
- all except the 4-inch diameter roll have a solid core to prevent "bottoming-out"
- store rolls on rack when not in use (6 51/2" long hooks)



48x16"


roll rack

24x4"

small

small

large



Skillbuilders® half rolls

- · half rolls offer stability of a wedge with positioning versatility of a roll
- loop can be added to flat bottom to anchor half roll for additional security

	width x length	+	elevation
30-1230	4x18"		2"
30-1231	4x18"		3"
30-1232	6x24"		4 1/2"
30-1233	8x30"		6"

# storage rack

#### **Tumble Forms® rolls** VALU-FORM® rolls

30-1000

30-1008

30-1009



- · adaptive shapes for positioning, neurodevelopment training, and rolling and rocking movements
- hang on roll rack with webbing strap

diameter x length							
30-3000	4x24"	small					
30-3001	6x24"	small					
30-3002	8x24"	small					
30-3003	8x36"	medium					
30-3004	10x36"	medium					
30-3005	12x36"	medium					
30-3006	12x48"	large					
30-3007	14x48"	large					
30-3008	16x48"	large					
30-3009	roll rack						



- Valu-Form® rolls and wedges from TumbleForms® are a cost effective line of positioning products
- versatile and economical
- · rolls do not have webbing strap
- seamless burgundy colored coating is urine-resistant and easily cleaned

#### diameter x length 30-3190 4x24" 30-3191 6x24" small 30-3192 8x24" small 30-3193 8x36" medium 30-3194 10x36" medium 30-3195 12x36" medium 30-3196 12x48" large 30-3197 14x48" large 30-3198 16x48" large

#### Tumble Forms® raised rolls



- raised rolls offer stability of a wedge with positioning versatility of rolls
- velcro on flat bottom can be used to anchor roll for additional security

	width x length	x elevation
30-3230	4x18"	2"
30-3231	4x18"	3"
30-3232	6x24"	41/2"
30-3233	8x30"	6"



# Tumble Forms2

#### positioning

#### sidelyer positioning system



- includes padded base with backrest and positioning straps
- optional two-piece block set can be used for leg abduction or as a headrest

30-3111 universal (42x23")

30-3112 2-block set

#### **Grasshopper®**



- achieves almost any positioning alternative by using any combination of the 17 standard components: wedges, adductors, rectangles, trapezoid, quarter roll, log, and straps
- mobile padded base allows mobility without re-positioning
- · includes accessory bag

30-3101 grasshopper (56x23")

#### Tadpole®



- creatively positions pediatric patient using any combination of standard components: log, half roll with semi-circular supports, contoured wedge and mirror/tray
- fits into carry-bag with shoulder straps

30-3105 complete (34x19")

30-3106 half roll 30-3107 tray with mirror 30-3108 wedge

#### Versa-Form® positioning pillows



- provides customized, semi-permanent (regular: 1 month, plus: 3 months) positioning
- shape/reshape in seconds by opening air valve and molding to desired shape and then evacuating air with vacuum pump
- Versa-Form® pillows with terry cloth covers pictured

	flat size	Versa-Form	Versa-Form Plus	covers
	16x20	30-2024	30-2004	30-2034
Account	22x12	30-2020	30-2000	30-2030
	22x26	30-2021	30-2001	30-2031
The state of the s	22x34	30-2022	30-2002	30-2032
	32x40	30-2023	30-2003	30-2033





vacuum pump

Single handed operation allows the user to form the pillow.

30-2040 pump



## Skillbuilders® modular seating system

- special needs children must be fed, hugged, transported and positioned; modular seating system makes all this possible in a dignified and humane manner
- built-in abductor, contoured interior for posturally correct seating and a safety harness

sizing chart on page 175

#### Skillbuilders® feeder seat





- basic positioning; seat has built-in abductor, posturally correct contoured interior and a safety harness
- · convert to floor sitter by adding wedge
- · add base to make mobile

#### feeder seat

30-1070 small 30-1071 medium 30-1072 large 30-1074 x-large

#### Skillbuilders® floor sitter



- · feeder seat and floor sitter wedge
- seat is secured to wedge with velcro
- adjust from upright to reclining position

#### floor sitter (seat and wedge)

 30-1080
 small

 30-1081
 medium

 30-1082
 large

 30-1083
 x-large

#### floor sitter wedge only 30-1085 for s, m, l 30-1086 x-large

#### head/neck support



- wrap-around hinged collar supports jaw line and occipital region for improved neck flexion, head control and midline positioning
- 1-piece neck support maintains head in mid-line and neutral position
- attaches to seat with velcro

30-1390 wrap collar

30-1395 collar, sm/med. 30-1396 collar, large/x-l

#### Skillbuilders® 3-piece mobile floor sitter



- with feeder seat, floor sitter wedge and flat mobile base
- seat, wedge and base are secured with velcro
- individual may be moved without re-positioning

#### 3-piece mobile floor sitter

30-1090 small

30-1091 medium

30-1092 large

30-1093 x-large

#### 3-piece mobile base only

30-1095 for small, medium, large

30-1096 x-large

#### Skillbuilders<sup>®</sup> 2-piece mobile floor seat



 2-piece sitter comes with feeder seat and mobile base

#### 2-piece mobile floor sitter

30-1170 small, w/ wood base 30-1171 medium, w/ wood base 30-1172 large, w/ wood base 30-1173 x-large, w/ wood base 30-1174 x-large, w/ metal base

#### 2-piece mobile base only

2-piece mobile base of 30-1088 small, wood 30-1089 medium, wood 30-1097 large, wood 30-1094 x-large, metal

## Skillbuilders® feeder seat cover



- soft, washable cover provides comfort and absorbs moisture
- easy to put on and remove

 30-1075
 small

 30-1076
 medium

 30-1077
 large

 30-1078
 x-large



# Tumble Forms2

## Feeder® Seat

#### 2-piece mobile floor sitter



- feeder seat and mobile base
- move individual without repositioning
- medium and large come with a wood box base, XL with a tubular steel base.

30-3171 medium 30-3172 large 30-3174 x-large mobile base **only** 30-3089 medium

large

x-large

30-3097

30-3094



Feeder seat (in red) plus wedge (in blue) make up the floor sitter

**color choice**: choose blue, red, teal or purple simply by putting a B, R, T or P after the sku code!

 large pneumatic tires allow travel over rough

 rear wheels fold for easy storage and transportability

#### feeder seat

- provides basic positioning
- seat has built-in abductor, posturally correct contoured interior and a safety harness

#### Sizing (height/hip)

small: 36/7", medium 48/10", large 60/13" x-large 72/16"

30-3070 small 30-3071 medium 30-3072 large 30-3074 x-large

#### floor sitter

- feeder seat and floor sitter wedge
- seat is secured to wedge with velcro allowing adjustment from upright to reclining position

30-3080 small 30-3081 medium 30-3082 large 30-3083 x-large

floor sitter wedge only

30-3085 s,m,l 30-3086 x-large

#### Rover® stroller with feeder seat



Rover® frame with feeder seat

30-3521 medium 30-3522 large

30-3523 x-large

Rover® frame only 30-3504 small/medium 30-3500 large

30-3505 x-large

Tray for Rover®

30-3506 medium frame

#### adjustable tray



- use with feeder seat, floor sitter, classroom seat and Carrie<sup>®</sup> seat
- adjustable tray height, angle
- non-slip large work surface (23x23") ideal for all activities
- large rubber feet keep tray from slipping

30-3079 adj. tray

#### feeder seat covers



 soft, washable cover provides comfort, absorbs moisture

30-3075 small 30-3076 medium 30-3077 large 30-3078 x-large

for feeder seat sizing chart, see page 175



# Tumble Forms2

#### Carrie® seat

- built-in abductor and contoured interior provide lateral trunk support and shoulder protraction
- adjustable neck stabilizer maintains head in mid-line and neutral position
- · harness securely positions child
- **sizing** (height/hip) small 36"/7", medium 48"/10", large 60"/13", x-large 72"/16"

 30-3326
 small (preschool)

 30-3327
 medium (elementary)

 30-3328
 large (junior)

 30-3329
 x-large (small adult)

#### Carrie® seat w/ footrest and tray

- seat, adjustable footrest, adjustable tray
- tray provides a surface for feeding and other activities
- footrest adjusts to maintain hip, knee and ankle flexion and to position child for optimum comfort

 30-3310
 small (preschool)

 30-3311
 medium (elementary)

 30-3312
 large (junior)

 30-3313
 x-large (small adult)

- Carrie<sup>®</sup> seat may be used as car seat, it meets F.M.V.S.S.213
- secure seat using standard vehicle seat belt

medium seat approved for use on airplanes





- Carrie seat, adjustable footrest and tray, and box-type mobile base
- base raises seat to preferred height. Individual may be moved without repositioning

Carrie® seat w/ mobile base, footrest, tray

30-3330 small (preschool) 30-3331 medium (elementary)

base only 30-3320 small/medium

#### seat covers & footrests



- soft, washable cover provides comfort, absorbs moisture seat cover
- 30-3305 small 30-3306 medium 30-3307 large 30-3308 x-large foot rest 30-3290 small 30-3291 medium 30-3292 large 30-3293 x-large trav 30-3280 small 30-3281 medium 30-3282 large

x-large

#### Carrie® collar



- wrap-around hinged collar supports jaw line and occipital region for improved neck flexion, head control and midline position
- neck support maintains head in mid-line and neutral position
- attach to seat with velcro

30-3390 collar

#### Rover® stroller with Carrie® seat







folds for storage

- Carrie® Rover® stroller combines the Carrie seat with a stroller base for maximum mobility indoors or out
- stroller has front swivel casters and large rear wheels

#### Rover with Carrie seat

30-3510 small (preschool) 30-3511 medium (elementary) 30-3512 large (junior) 30-3513 x-large (small adult)

#### Rover Carrie seat frame only

30-3501 s/m (pre-school, elem.) 30-3502 large (junior) 30-3503 x-large (small adult)



30-3283

## modular seating guide

for seating, refer to pages: Skillbuilders® feeder seat: 172 Tumble Forms 2® feeder seat: 173 Tumble Forms 2<sup>®</sup> Carrie<sup>®</sup> seat: 174





			nall school		dium entary		rge nior		rge adult
_	t (inches) t (pounds)		to 36 )-40		-48 -60		-60 -100		-72 120
<u>J</u>	inside height (inches)	FS CS	18 22	FS CS	23 26	FS CS	29 30	FS CS	35 33
	inside depth (inches)	FS CS	8	FS CS	9½ 10	FS 1 CS	31/4 15	FS CS	161/2 15½
	inside width inches	FS CS	7¾ 8	FS CS	9 20	FS CS	12½ 13	FS CS	15 14
	outside width inches	FS CS	10½ 12	FS CS	14 14½	FS CS	17 16¼	FS CS	21 18½
	outside height inches	FS CS	21 25	FS CS	27 30	FS CS	34 35	FS CS	41 39

#### Tumble Forms® corner chair



- large base prevents tipping; built-in abductor helps leg extension and reduces extensor thrust
- adjustable harness support and hip strap securely position child
- removable tray (28"x28") can lie flat for feeding or tilt to 30°
- sizing: up to 60"

30-3050 corner chair

#### Tumble Forms® square chair





- use as chair (top) or table (bottom)
- position neurologically impaired children (up to 60")
- shapes modify the seating position
- 8-piece shape set contains 2 large rectangles, 2 small rectangles, 2 triangles, 1 ring and 1 disc

30-3062

seat with shape set

30-3061

shape set



#### Outstander™ multi-positioning stander



- tray can be moved for each position
- adjustable headrest provided for use in supine position
- pads can be removed for very small or short children and then added to the frame as they grow
- angle adjustment varies amount of weight carried by lower body
- included accessories: foot restraint and heel cups, knee blocks, curved rigid hip supports, lateral supports, trunk support, headrest, tray

31-3251 small stander 16-35" 100 lbs 31-3250 medium stander 23-43" 150 lbs 31-3255 replacement height bar\*

\*the stander can grow as the child grows -- call with needed dimensions.

## Tumble Forms® Tristander® multi-positioning stander



- · for prone, supine and vertical standing
- pneumatic tilt adjusts stander from 15° to 90°
- each cushioning pad adjusts to securely accommodate a wide range of positioning needs
- included accessories: frame, shoes with bar, hip, knee and trunk supports, headrest and 3 inch double locking casters
- interchangeable modules between 45 Tristander (child) and 58 Tristander (adolescent)

child (32-45")		adolescent (43-58")
31-3201	unit only	31-3224
31-3200	unit with tray	31-3221
31-3230	tray only	31-3210

## Marvel™ mobile free stander



- Marvel free stander provides vertical alignment and upright support
- easy to assemble, adjust and reposition
- four points of support: foot, knee, pelvic and trunk level
- included accessories: foot restraints and heel cups, knee blocks, posterior pelvic stabilizer, trunk support, tray, wheels
- small 18-29" 75 lbs medium 25-44" 125 lbs

31-3461 small 31-3460 medium

#### TUGS™ standers





Horizontal

## TUGS mobile supine stander

- board tilts from vertical to horizontal and can be secured at any in between angle
- supports child at head, trunk, pelvis, knees and feet
- headrest and armrest fold down, 30-50", 100 lbs.
- · adjustable tray is included

31-3420 supine

## TUGS mobile prone stander

- position with adjustable main board, footboards, support blocks, foot plates, 25-48", 100lbs.
- height/angle adjustable tray is included

31-3430 stander

#### **HUGS™** standers



## **HUGS** stationary vertical stander

- durable stander has 4 adjustable support points: chest, pelvis, knee, foot
- tray is included
- little 32-42", 50 lbs. big 43-53", 65 lbs.

31-3440 little 31-3441 big

#### prone mobile stander



- benefits of standing with convenient mobility
- padded supports provide flexibility and positioning
- foot plate accommodates leg length discrepancies
- child 30-37", up to 50 lbs.

31-3450 mobile stander



## positioning system

# **SKILLBUILDERS**®



## Skillbuilders® MSS Tilt and Recline positioning system







- designed for children with moderate to severe physical involvement
- basic seating module can be used as chair sitter, floor sitter, mobile floor sitter or mobile high chair
- tilts to 45° and reclines 80° to 140°
- basic seating module includes:
  - back: adjustable height and angle small: 10-12½", medium: 12½-16"
  - seat: adjustable depth: 5-14"
     width: small: 9", medium: 10½"
  - footrest: adjustable height and depth
  - headrest: adjustable height, depth and angle; rotates and can be placed off-center
  - abductor: adjustable height and depth; position between thighs or knees
  - lateral supports: adjustable height and width
  - armrest: adjustable height and angle
  - trunk support: zippered H-harness
  - cushion covers: removable and washable

#### chair sitter (a)

- remove the seating module from the mobile floor sitter or mobile high chair and secure it to the chair sitter attachment bolster kit
- use the straps to secure the chair sitter to a standard chair
- 31-3705 chair sitter attachment

#### mobile floor sitter (c)

- mobile base allows the caregiver to move the child from place to place without the hassle of repositioning
- the low mobility base can be purchased separately
- 16" floor to seat height

31-3710 small mobile floor sitter 31-3711 med. mobile floor sitter

31-3706 add'l. low mobile base

#### floor sitter (b)

- to use as a floor sitter simply remove the mobile base from the mobile floor sitter or mobile high chair
- 13.5" floor to seat height
- cannot be ordered separately (order as mobile floor sitter or mobile high chair and remove base)

#### mobile high chair (d)

- high base will bring the child to table height to facilitate eating or other activities
- mobile base allows the caregiver to move the child from place to place without the hassle of repositioning
- add'l high mobility base purchased separately; 22.5" floor to seat height

31-3720 small mobile hi chair 31-3721 med, mobile hi chair

31-3707 add'l. high mobile base

#### 





- Traveler stroller accommodates the MSS tilt and recline seating module
- lightweight stroller folds easily with the seating module in place; has three tilt positions, oversized rear wheels and front swivel wheels

31-3490

traveller stroller base only (left)

accessories

31-3497

canopy for traveller stroller (middle)

31-3498

basket for traveller stroller





## chairs for school and early intervention facilities

#### Skillbuilders® first class school chair



- · use in classroom, clinic or home
- · for physically challenged children with mild to moderate involvement; multiple adjustments and accessories meet the needs of each child as they grow and develop
- seat is height and depth adjustable with a 15-degree anterior or posterior tilt; armrests are height and width adjustable
- optional support kit contains trunk H-harness for anterior chest support, adjustable depth abductor and adjustable height/width lateral support; pelvic strap helps maintain proper hip position

school chair		mobile school chair
31-3800 small 31-3802 large		31-3810 small 31-3812 large
Accessories Tray Mobility base Extension legs Support kit (see above)	small 31-3832 31-3830 31-3838 31-3834	31-3840 3 31-3848
Footrest Headrest Hip guides Anti-tip supports	31-3835 31-3833 31-3836 31-3837	31-3843 31-3846

#### roll chair



- padded roll helps to maintain proper hip abduction
- · front and rear legs may be adjusted to various angles
- seat belt included
- chestnut vinyl is standard

31-1120 small (to 50" tall) 31-1121 medium (to 68" tall)

#### accessories

31-1122 adj. clear tray small 31-1123 adj. clear tray med. 31-1124 adj. lateral supports



- adjustable upholstered seat, sling-style adj. depth back and adj. footboard
- seat belt included
- · chestnut vinyl is standard
- optional wood grained tray has raised rim: 24"x24", 7"x10" cutout

31-1130 adj. chair accessories 31-1131 tray

#### Kettler® Herlag® high chairs



specify: "N" Natural, "M" Mahogany, "B" Blue, "W" White or "R" Red

- Tipp Topp high chairs have adjustable seat and footrest, 3 point safety harness, and quick release removable feeding tray
- made of solid European Beechwood
- recommended for ages 8 months and up

Tipp Topp (to 150 lbs) specify N,M,B,R 31-3880

- Junior high chairs have adjustable seat and footrest, removable safety bar, 3 point safety harness, and resin floor protectors
- made from solid hardwood
- recommended for ages 10 months and up.

junior (to 125 lbs) specify M,N

- Mini high chairs are 3-position height adjustable
- · made of solid hardwood

mini (to 100 lbs) specify M,N,W

#### multi-use chairs













D





- seating and positioning needs of any neurologically delayed or orthopedically involved child can be met
- back and footplate are height adjustable; seat depth adjustable
- seat belt included; chestnut vinyl is standard child chair

31-1140

31-1141 adol	escent chair	
accessories	child	adolescent
(a) adj. clear acrylic tray	31-1122	31-1123
(b) adj. lateral supports	31-1124	31-1124
(c) knee abductor	31-1145	31-1145
(d) adj. pelvic support	31-1146	31-1147
(e) adj. foot bracket	31-1148	31-1149



# seating system

Skillbuilders® products for balance and stabilization training activities

#### Seat2Go positioning seat





with: abductor and headrest accessory.

- lightweight padded positioning seat is designed for children that need support while sitting on the floor, at the table or at play
- contoured seat prevents child from sliding forward and improves seating posture (95° seat to back angle)
- double strap trunk support keeps child securely positioned
- in addition, the Seat2Go can be used as a stroller insert

31-3850 small seat (50 lbs. capacity) 31-3851 medium seat (75 lbs. capacity)

#### **Options**

31-3860 headrest (adjustable height) 31-3861 abductor (adjustable height)

#### Columbia<sup>®</sup> TheraPedic™ IPS™ car seat

- safely transport special children
- shell side and positioning pads provide laterial support

#### standard features:

- high profile sides for optimum position and support
- vacuform ABS plastic shell and chromed steel side tubes
- durable, machine-washable, soft plush cover
- latex and silicone free
- 5-point harness system
- meets FMVS213, crash test and fire code 302
- · certified for use in aircraft
- available colors: blue, rose





31-1010 IPS™ 2000 child seat: 20-102 lbs, up 60" 31-1016 IPS™ 2500 small adult seat: 40-130 lbs, 54-66"








- lightweight padded back support promotes improved posture and trunk control
- · securely straps to chairs, wheelchairs and strollers
- · contoured back maximizes support and comfort
- two-strap system with buckle closure secures back support to an existing back surface
- straps can be stored in back pocket when not in use
- removable lateral supports are individually height and width adjustable
- zippered vinyl trunk harness for anterior chest support

31-3870 14" wide x 16" high 31-3860 headrest

#### Columbia® TheraPedic™ IPS™ sroller transport base

 lightweight folding stroller for use with Thera-Pedic™ IPS™ car seat



 perfect for shopping at the mall, going to the playground, or traveling

#### standard features:

- five angles of tilt: 20, 25, 30, 35, or 40 degrees
- height adjustable footrest
- large rear wheels for a smooth ride on any surface
- foldable frame for convenient storage and transport
- seat locks to safely secure car seat to mobility base
- lightweight, yet strong, durable frame

31-1011 IPS 2100 stroller for child seat 31-1017 IPS 2520 stroller for s-adult seat 31-1013 umbrella stroller for child seat

#### seat with transport base

31-1012 child seat + IPS™ 2100 base 31-1018 s-adult seat +IPS™ 2520 base 31-1014 child seat + umbrella stroller



31-1012

31-1011





# mobility products

• children who can't pivot, crawl, creep or walk often benefit from equipment that helps them move and explore their environment

#### CrawlAbout® crawlers







- height adjustable aluminum frame
- foam padded suspended harness
- four easy glide pivoting casters

32-2050 small (purple), 22-25"H, 50 lb. capacity 32-2051 medium (blue), 26-29"H, 75 lb. capacity

- height adjustable child crawler promotes crawling movements when patients need added body support
- used to treat children who have difficulty tolerating movement

32-2015 crawler 32-2016 positioning belt

#### Whirl-a-Wheel™ extension scooter



- mobility unit is low to the ground
- children with limited upper body strength can get in and out without assistance

32-2011 child, 15" seat width 32-2013 optional back support



# Tumble Forms<sup>®</sup> scooter board

- round board has padded nonslip top with 2 hand loops for security and to help prevent "pinched" fingers
- 2" casters, 100 pound weight capacity

30-3153 24" diameter

#### Tumble Forms® Jettmobile® extension scooter



- abductor, chest positioner, 2 half-moon shapes and positioning strap
- wheels position child
- · abductor keeps legs from scissoring
- push or pull with arms or legs.

30-3140 child 28" complete 30-3141 teen 38" complete 30-3142 child 28" w/o accessories 30-3143 teen 38" w/o accessories



#### Tumble Forms® Ready Racer®

- pre-wheelchair mobility aid for children 18 months to 5 years old (to 55 lbs.)
- improves perceptual-motor skills, strengthens arms and upper trunk
- optional handle adjusts (24-36")

30-3160 racer

30-3161 adjustable handle



# sensory and vestibular ball pits and group activity units

#### free standing ball tower



- steel reinforced, fully padded structure with 1000 balls
- structure even includes canopy!
- nylon netting keeps toddler and balls inside pool
- snap together unit easy to put together
- dimensions: 4' hexagon (54" corner to corner) x 15" deep x 72" high. Weighs 64 pounds

#### indoor/outdoor ball tower

32-2490 pit structure, 1000 balls

#### panel-sided ball pit



- sensory environment pit gives child visual and tactile feedback, provides the "buoyancy" of water without the mess
- expandable, portable and safe
- side panels (48" L x26" H) solid with foam padding, reinforced vinyl cover
- · design allows for shape flexibility

	pane	els size	balls
32-2408	4	4'x4'	1000
32-2400	4	4'x4'	2500
32-2401 🖵	<b></b> 5	6'x6½'	3500
32-2402 🖵	<b>6</b>	5'x7'	5500
32-2403 🖵	8	10'x10'	9000
32-2432 ac	ditior	nal panel	4'x2'x3'

products ship by truck

#### Tortoise® shell

- adds dimension to individual or group vestibular therapy programs
- use lid as recessed base to prevent rocking
- deluxe model includes net with balls.
- hang from Vestibulator<sup>®</sup> using optional suspension kit





30-3640 standard (shell with lid) 30-3641 deluxe with net, balls 30-3643 deluxe + suspension kit

30-3644 suspension kit 30-3642 net with balls

#### additional balls for sensory pit

- brightly colored balls support the child and encourage gross motor activity
- 500 balls can be stored (and washed) in the nylon mesh ball bag

32-2416 ball bag (each)

 	, p		
3" (73mm	) 500/case	250/case	12/case
red	32-2410R-500	32-2410R-250	32-2410R-12
yellow	32-2410Y-500	32-2410Y-250	32-2410Y-12
blue	32-2410B-500	32-2410B-250	32-2410B-12
green	32-2410G-500	32-2410G-250	32-2410G-12
purple	32-2410P-500	32-2410P-250	32-2410P-12
orange	32-2410O-500	32-2410O-250	32-24100-12
multi	32-2410M-500	32-2410M-250	32-2410M-12

#### stair steps



 develop climbing skills and other gross motor functions

	steps	WxDxH
31-2300	2	22x24x14
31-2301	3	22x32x21

#### barrel roll



 multi-colored barrel roll is ideal for development of balance, body awareness and coordination, as well as for a games

31-2600 38"x32" L x d

#### convertible crawl box



- 2-in-1 product converts from a 2 foot by 8 foot folding mat to a 2-foot cube
- use to develop gross motor skills


31-2670 crawl box



# pediatric vestibular physical therapy and play

#### Tumble Forms?

# Turtle™ system



- vestibular stimulation, sensory integration and mobility
- use dome curved-side down as rocker board, or flat-side down as vestibular trainer
- · add scooter to dome to create vestibular mobility trainer
- use with adj. leg as T-stool
- deluxe w/ net, rope, ascender

30-3221 standard 30-3220 deluxe scooter/stool 30-3222

# Turt-L-Shell





- inflatable half ball
- · use with either the flat or curved side down to perform vestibular exercises
- adjust difficulty level by adding or removing air
- 24"(60cm) diameter

30-1940 Turt-L-Shell

#### Tumble Forms?



- board has padded non-slip top with recessed sides to help prevent pinched fingers
- use board for vestibular activities involving orthopedically and neurologically involved children
- 18x24" top, 6" height

30-3152 board

#### Tumble Forms?



- 6-foot reversible beam is covered with TumbleForm® nonslip surface
- position to have 3" or 6" width
- top is 7" from floor. Weight capacity is 200 lbs.

32-2111 balance beam

#### Tumble Forms?

# barrel roll



- · for vestibular training, position child on top of barrel and "roll"
- use barrel as a tunnel for play therapy or to practice crawling
- · soft outer shell for comfort with solid core support to prevent
- roll is 2 feet long with 22" outside/11" inside diameter

30-3040 barrel roll

#### **Tumble Forms**



- for vestibular, balance training and exercise activities
- soft outer shell for comfort with solid core support to prevent "bottoming out"

#### diameter

30-3030	11"
30-3031	16"
30-3032	22"

#### Cando<sup>®</sup> wedge



- inflatable wedge requires active participation to assume correct sitting posture
- · combines dynamic seating and static positioning
- inflate or deflate wedge to vary degree of difficulty

30-1925 large (15x15") 30-1926 small (10x10")

#### Cando® jump ball



- jump ball is a fun way to proproprioceptive vide vestibular stimulation
- heavy-duty grip handles
- holds up to 300 lbs.

30-1825 yellow 16" 18" 30-1826 red 30-1827 green 30-1828 blue



#### junior mat

- lightweight and portable
- · comes with carry bag; has die-cut handles
- "air channel" construction allows mat to conform to body position

36x17x1.2" 38-0551 38-0551-10 case of 10 mats





# pediatric vestibular physical therapy and play

#### plastic ball stool/trainer



- versatile exercise ball chair holds a 50cm ball (included!) securely in place
- use as a chair, or use as an inflatable exercise dome
- casters have locking mechanism to make chair stationary
- · lock casters when using during vestibular exercises . . . you do not want exercise ball chair to roll away
- includes ball

30-1796 mobile exercise ball stool

#### plastic ball chair



- · mobile ergonomic seating system combines the benefits of ball therapy with the convenience of a chair base
- ball held securely in place by the base
- features removable back, without arms
- child and adult sizes available
- casters glide on carpet or hard surfaces
- includes ball colors may vary
- light assembly required

#### plastic chair without arms

adult size w/ 55 cm ball 30-1792 30-1795 child size w/ 35 cm ball

#### optional locking casters

30-1797 2 each



#### Cando® inflatable exercise balls - regular






- used to improve balance, coordination, flexibility, strength and even just for fun
- · used by therapists as an aid for vestibular movement and equilibrium therapy
- · non-slip surface is ribbed for extra security
- under inflate to give a soft, mushy feel or fully inflated for a firm, bouncy feel
- support 300 pounds (136kg)

#### nolyhanned

polybaggeu			
30-1800	30cm	11.8in. b	lue
30-1801	45cm	17.7in. y	ellow
30-1802	55cm	21.6in. d	range
30-1803	65cm	25.6in. g	jreen
30-1804	75cm	29.5in. r	ed
30-1805 30-1806	85cm 95cm	33.5in. b 37.4in. r	

#### Cando® inflatable saddle rolls





- · easier to control because motion is limited to forward and backward direction
- · aids in developing muscles, coordination and balance
- the "saddle" offers extra support for the patient
- supports up to 300 pounds

	cm-HxL	in-HxL	
30-1725	40x90	16x35	yellow
30-1726	50x100	20x39	orange
30-1727	60x110	24x43	green
30-1728	70x120	28x47	red
30-1729	80x130	32x51	blue



# **Oycem**Non-Slip Solutions

#### netting on roll



- · soft mesh non-slip material
- improves grip and provides padding
- mesh allows air to circulate, avoiding clamminess and maintaining comfort
- works well placed on chairs or wheelchairs to help prevent cushions or people from sliding off
- can be easily cut to size
- Ivory only

50-1550 24" x 61/2 foot roll

#### circular table mats/pads



#### rectangular table mats/pads



- mats for healthcare, home and leisure provide a secure surface that anchors items such as cups, plates, tools, note pads and telephones to a table
- mats are great for trays and any dry surface which needs extra grip
- during therapy the extra stability can facilitate a variety of coordination and spatial awareness exercises
- Y-yellow, R-red, B-blue, G-forest green, MG-metallic green, BLK-black

50-1595	5½" diameter	50-1590	10" x 71/4"
50-1596	7½" diameter	50-1591	10" x 14"
50-1597	81/2" diameter *	50-1592	15" x 18"
50-1598	10" diameter *		

\*8½" and 10" diameter circular pads available in yellow, red and blue only

#### floor mat



- used in hospitals, residential homes, and dry surface where there might be a risk of slipping and where there is a need for increased stability
- gray only

50-1641 18" x 24"

#### walking mat



- assists in feet stabilization during a number of movements such as transferring from bed to chair
- gray only

50-1640 24" x 12'

#### animal-shaped play mat





- provide secure surface that anchors items
- animals help retain attention
- Y-yellow, R-red, G-forest green, B-blue

50-1620 hippo 50-1621 tortoise

#### animal seat mat



- adds fun and retains child's attention during therapy
- added stability can be used to facilitate coordination and spatial awareness exercises
- blue only

50-1630 hippo

#### jar and bottle opener



- jar opener fits in the palm of the hand; effective on screwtop jars including child resistant and tamper-proof lids
- bottle opener fits bottles, caps
- B-blue, Y-yel

50-1650 jar opener 50-1651 bottle opener

#### tray liner



- place on any tray to provide a non-slip surface
- prevents cups and plates from sliding when the tray is tilted
- can be cut to fit smaller tray
- F-floral, C-clear, GH-gingham

50-1610 15x10½"

#### cup/can/glass holder



- fits snugly onto bottom of a mug or tumbler
- non-slip base helps to prevent movement and spillage
- can be used with hot drinks
- blue only

50-1652 holder

#### cleaning wipes



- cleaning wipes are an alternative to soapy water
- wipes leave your product instantly dry and ready for use
- contains 30% alcohol and 70% de-ionized water

50-1660 wipe 10-pack



- non-slip material solves many stabilization and gripping problems encountered by people with limited dexterity, restricted mobility and reduced wrist or upper limb strength
- grips on both sides to prevent movement
- ideal for providing support, stability and confidence during everyday activities and exercise

color - Specify a color by adding the color abbreviation from the list at the end of each style, for example: 50-1500R is a red original 8"x10yd roll. IF NO COLOR IS SPECIFIED, BLUE WILL BE DELIVERED.

#### original Dycem® on rolls



- effective anti-slip material
- original material is 0.4mm with a thick non-woven center
- can be cut to any size and is malleable so it can be fitted around various objects
- Y-yellow, R-red, B-blue, BLK-black, Wwhite, G-forest green

50-1500	8" x 10 yard roll
50-1501	8" x 6½ foot roll
50-1502	8" x 3¼ foot roll
50-1508	16" x 16 yard roll
50-1505	16" x 10 yard roll
50-1506	16" x 6½ foot roll
50-1507	16" x 3¼ foot roll

#### grip-it Dycem® on rolls



- 0.6mm thick polyester-film base material
- more rigid than the original and should be used where a stiffer surface is required, such as under writing paper
- F-floral, GH-gingham, C-clear

50-1510	8" x 10 yard roll
50-1511	8" x 6½ foot roll
50-1512	8" x 3¼ foot roll
50-1515	16" x 10 yard roll
50-1516	16" x 6½ foot roll
50-1517	16" x 3¼ foot roll

#### comfort Dycem® on rolls



- 1.1mm thick with center carrier
- soft foam texture has been designed to give extra padding during a variety of everyday activities and can help make therapy more relaxed and comfortable
- blue only; shown with optional dispenser

50-1521 8" x 6½ foot roll 50-1522 8" x 3¼ foot roll

#### roll rack dispenser

• fits up to 16"x 16 yards of any Dycem®

50-1509 roll rack dispenser

#### self-adhesive discs



- put on any piece of equipment to provide a secure base
- prevents movement especially on angled surfaces
- simply peel off the backing from the disc and press on to object to be secured
- black only

50-1580 5/8" diameter, 12 each 50-1581 7/8" diameter, 12 each

#### self-adhesive strips



 original: Dycem® strips are made with an adhesive backing to permanently secure the strip; B-blue, Y-yellow

50-1560 16" x 1-1/8" strip (3 ea.)

 comfort: Dycem® self-adhesive strips are available with foam for enhanced cushioning, protection, and vibration absorption; blue only

50-1570 16" x 21/4" comfort strips (3 each)

#### self-adhesive material panels



original: self-adhesive material permanently secures the Dycem® piece in place; B-blue, W-white, R-red, Y-yellow

50-1530 16" x 1 yard original SA roll

 comfort: self-adhesive material with foam for enhanced cushioning, protection, and vibration absorption; blue only

50-1540 8" x 1 yard comfort SA roll



# hook and loop fastener

hook or loop rolls; pressure sensitive adhesive back

# hook or loop rolls; standard non-adhesive webbing back





• standard woven-back is sewn or glued to object

aggages	0

• self-adhesive back is peeled off and pressed against an object

	10 yard			25 yard		10 yard			25 yard
hook	loop	width	hook	Іоор	hook	loop	width	hook	loop
24-7066	24-7067	1/2"	24-7064	24-7065	24-7056	24-7057	1/2"	24-7054	24-7055
24-7002	24-7003	1"	24-7000	24-7001	24-7012	24-7013	1"	24-7010	24-7011
24-7006	24-7007	11/2"	24-7004	24-7005	24-7026	24-7027	11/2"	24-7024	24-7025
24-7022	24-7023	2"	24-7020	24-7021	24-7032	24-7033	2"	24-7030	24-7031

specify color desired by applying the following after SKU (standard and adhesive back)
W = white, BLK = black, T = tan; white is default color

#### elastic (stretchable) loop roll



- secure an object with this stretchable elastic loop
- loop can be applied to standard and adhesive back hook

24-7017 1" x 10 yd 24-7015 1" x 25 yd 24-7018 2" x 10 yd 24-7016 2" x 25 yd

specify W = white, BLK= black

#### back-to-back hook and loop roll



- · hook and loop are bonded back-to-back
- create circular closures for positioning, restraining and strapping

10 yd length	width	25 yd length
24-7070	1/2"	24-7074
24-7071	1"	24-7075
24-7072	11/2"	24-7076
24-7073	2"	24-7077


#### double sided loop roll



- soft, durable, non-stretchable strapping
- loop fabric is laminated to both sides of padding and can mate with any hook

24-7255	1" x 10 yd
24-7256	1" x 25 yd
24-7257	2" x 10 yd
24-7258	2" x 25 yd

# hook circles and tabs (self-adhesive)



- self-adhesive tabs
- use to attach splint straps, or other applications that require small pieces of hook

24-7245 rectangular 2"x1" (175/pack) 24-7246 circular 7/8" dia. (850/pack) 24-7247 circular 1-3/8" dia. (600/pack) 24-7248 circular 1-7/8" dia. (450/pack)

# loop strap with self adhesive hook



regular loop with 2 inches of self adhesive hook on one end

24-7265 1" x 12", 10 each
24-7266 1" x 18", 10 each
24-7267 1" x 24", 10 each
24-7268 2" x 12", 10 each
24-7269 2" x 18", 10 each
24-7270 2" x 24", 10 each

#### D-ring hook and loop strap



- D-ring hook and loop strap
- with or without self adhesive material

10 each	regular	adhesive back
1" x 12"	24-7100	24-7106
1" x 18"	24-7101	24-7107
1" x 24"	24-7102	24-7108
2" x 12"	24-7103	24-7109
2" x 18"	24-7104	24-7110
2" x 24"	24-7105	24-7111



# dressing aids

# sock and stocking aid with ridges





- · makes putting on socks easier
- three "fingers" to accommodate heel and foot contours
- loop strap can be used with one hand
- length is adjustable up to 30"

86-0001 sock and stocking aid

#### slip on dressing aid





- a notch helps remove shoes and socks
- handles have hook-and-loop closure at one end to hold onto clothing
- squeeze top of handles to open sock for entry of foot and leg

86-0020 slip on dressing aid, 19'

#### double button loop

- two loops fit a large variety of button sizes
- handle provides several grasping positions
- loops fold when not in use
- handle 5" long; button loops 4½" long 86-0300 double button loop

#### elastic shoe laces

- for those who require better support
- can be tied in usual manner
- size: 3/16" x 24"



86-1125 elastic shoe laces, 2 pr., black 86-1126 elastic shoe laces, 2 pr., white 86-1127 elastic shoe laces, 2 pr., brown

# deluxe flexible sock and stocking aid





- three "finger" design holds stockings
- inside nylon material reduces friction
- outside terrycloth holds stocking
- strong molded insert; hand washable
- two long loop handles for people with limited hand function

86-0002 deluxe sock and stocking aid







- combination hook/pusher assists in putting on shirt and jacket, pulling up slacks or removing socks
- hook on the opposite end aids in pulling zippers and shoelace hoops

86-0030 dressing stick, 24" long

#### button loop/zipper pull combo



 built-up foam rubber handle with buttoning loop at one end, and zipper pull hook at other end

86-0310 button loop / zipper pull

#### rotating button loop/zipper pull

- complete buttoning process with one hand
- built-up foam rubber handle with buttoning loop at one end, and zipper pull hook at other end
- button loop rotates for flexibility

86-0311 rotating button loop / zipper pull

# compression stocking aid heel guides





- rigid, extra wide design
- extra long handles, 17" with 3 grip levels
- · comes apart for traveling
- for compression hose up to 40 mmHg
- heel guide keeps foot in perfect position

86-0010 compression sock aid

#### hip kit



- ideal for anyone recovering from hip, knee or back surgery
- includes 4 popular assistive devices for independent dressing and bathing
- contains a round sponge, reacher, "T" handle shoehorn, sock and stocking aid 86-0090 hip kit

#### plastic "T" handle shoehorn

 unique handle allows wrist to remain in neutral position while putting on shoes

86-0350 "T" handle, 20"

#### metal shoehorns



- extra long shoe and sock removers
- notch to remove shoes and socks

86-0360 24.5" shoehorn 86-0361 27" shoehorn

#### plastic shoehorn

• flexible shoehorn

86-0370 plastic coated handle



# household gadgets

#### key holders and turners

#### Hole-in-One key holder



- can be used with fingers inserted through hollow center or around outside
- has loop for wrist strap, key chain, key ring or to hang on hook.

60-0100 key holder, 5" long

#### **Key Turner for one key**


- · neon yellow for visibility
- improves grip and turning leverage when using a key
- key is secured to turner with binding post
- handle has 2 finger holes for control 60-0101 key turner



#### Key Holder for two or four keys

- rigid handle provides extra leverage to
- keys fold back into plastic holder when

for two keys, 3½" x 1"



- turn key in lock
- not in use

60-0103 for four keys, 3½" x 1"

#### electric light switches

#### **Big Lamp Switch**

- · easy to see, grip and turn
- 2"dia. three-spoked knob creates more leverage to turn lamp switches on and off
- unscrew and remove the original knob and replace with Big Lamp Switch

60-1100 Big Lamp Switch 60-1100-3 bag of 3



#### wall switch extension handle

- provides 12" extension to any switch
- puts wall switches within easy reach of children, wheelchair users or people with limited reach
- hardware is included 60-1101 bag of 2

#### plug puller

- · attaches permanently to plugs
- secure with supplied plastic ties
- two large finger holes make it easy to hold and manipulate

60-1102 plug puller



#### door grippers

- door knob extender
- 5" extension handle provides extra leverage
- fits over standard door knob. Allows complete access to keyholes
- cord or string can be inserted though slot for people with limited reach

60-1110 doorknob extender, 2 ea.

#### door knob / faucet gripper

- inner ribs grip knob
- · outer surfaces allow hand to turn knob
- installed by slipping device over knob 60-1111 door knob/faucet gripper, 2 ea.



#### door knob gripper

- stretch over door knob to produce a soft ribbed surface
- reduces strength required to open a door
- can also be used on outdoor water facet handles and kitchen and bathroom taps
- latex free

60-1112 door knob gripper

#### T-handle knob turner

- can be used on any non-round or ribbed knob or appliance handle
- grip is secured when pressed onto knob
- measures 21/4" in diameter but can grip larger irregular-shaped knobs
- T-shaped handle measures 51/4" long 60-1113 knob turner



#### jar / bottle openers

#### easy twist jar opener

- · select the appropriate opening size, grip firmly and twist
- · soft rubber gasket creates a firm grip easy twist jar opener



#### Penguin-Pop jar opener

- penguin shaped jar and bottle openers fit jar lids and bottle caps small and large
- place around the jar or bottle top, hold firm and turn counter clockwise.

60-0001 penguin-pop jar opener

#### 4-way opener

- · soft, comfortable non-slip grips
- · opens bottle caps
- · pops jar lids; opens soda cans 60-0002 4-way opener





#### Hot Hand® hand protector / jar opener

- pockets fit thumb and fingers
- works as a gripper to unscrew tight jars
- secure grip on hot, cold, or slippery items 60-0020 hand protector / jar opener



#### reachers

#### Matey® reachers



- · provide a tight grip
- can pick up small flat items such as dimes or bulky items such as newspapers, clothing or cans
- magnet on the tip helps pick up small steel parts
- short post at the end can be used as a dressing aid
- handle is ergonomically designed to fit the hand
- maximum grip opening is 21/2"
- handle has a hole for hanging
- · made of aluminum and plastic

#### A • standard length reaching aid

standard model, 26- $\frac{1}{2}$ " long, weighs 6 oz.

50-1120 standard length reaching aid

#### b • extended length reaching aid

use when extra reach is needed, 32-1/2" long, weighs 8 oz.

50-1121 extended length reaching aid


#### C • pocket-fold reaching aid

folds to 12-½" long. Can be carried in a pocket or bag Extended length 22" long, weighs 6 oz.

50-1122

Pocket-Fold reaching aid

#### E-Z reachers



- pick up a variety of objects ranging in size and weight
- full hand trigger allows for more strength to be used than a conventional finger trigger

50-1130	40" E-Z reacher
50-1131	30" E-Z reacher
50-1134	20" E-Z reacher
50-1132	30" E-Z reacher, locking
50-1133	30" E-Z reacher, folding
50-1139	replacement suction cup, each (includes suction cup and acorn nut)

#### handi-reachers



- lightweight reacher has an ergonomic trigger that is operated with whole hand
- · constructed with wired-driven jaws
- frame resists twisting and bending
- includes magnet for picking up metal objects and a hook to aid with dressing

50-1140 standard reacher, 24" 50-1141 long reacher, 30" 50-1142 x-long reacher, 35"

#### Raptor® reachers



- a low cost, durable, lightweight, reacher
- The jaw opens a full 2½" and provides a tight grip on even the smallest items
- A hook-like extension on the front can be used as a dressing aid
- 24" long. Weighs only 6 oz.

50-1150 Raptor reacher

#### reaching aids





- aids patient with limited reach and hand strength
- pistol grip has magnetic tip to retrieve objects
- drag lug helps position products for easy pickup
- open-jaw (closes when squeezed)

#### metal frame

50-1100 26½", 6oz. 50-1101 32", 8oz. 50-1102 folding 22", 6oz.

wood frame with pistol grip 50-1104 32", 8oz.



#### Inner-Lip™ plate



- deep inner lip keeps food from sliding off the plate
- bring utensil to the edge of the plate and push the food onto the utensil
- 9" (23 cm) diameter
- dishwasher, autoclave and microwave safe, except non-skid version

62-0100 red, plastic
62-0110 light blue, plastic
62-0111 sandstone, plastic
62-0111-12 sandstone, plastic, 12 ea.

62-0120 non-skid plastic, sandstone

#### partitioned scoop dish



- three 1-½" deep compartments keep food separated
- two smaller sections hold 6.75 oz. The larger one has a 12.6 oz. capacity
- high compartment walls may be used to push food onto forks and spoons
- stackable with clear plastic covers for easier transporting and storing of food
- diameter measures 81/4" (20 cm)
- dishwasher and microwave safe

62-0130 partitioned dish w/lid, sandstone 62-0131 partitioned dish w/lid, red

# plates, bowls, etc.

#### non-skid scoop dish



- non-skid base
- combination plate/bowl has a contoured lip that can be used to help push food onto cutlery
- measures 6 x 7½" x 1¾" deep on the high side
- not recommended for microwave and dishwasher

62-0140 non-skid dish, sandstone 62-0141 non-skid dish, red

#### scoop bowl with suction cup base



- high rim and reverse curve on one side aids in scooping food onto a utensil without spilling
- bottom of dish has a rubber suction base to prevent skidding
- hemispherical, 4½" (11.4 cm) diameter
- dishwasher safe

62-0150 bowl with suction cup base to prevent skidding

#### scoop plates



- high rim and reverse curve on one side aids in scooping food onto a utensil without spilling
- bottom has a rubber suction base to prevent skidding
- 6-3/4" diameter
- dishwasher safe62-0160 plate w/suction cup



- flat bottom and high rim with a reverse curve on one side aids in scooping food onto a utensil without spilling
- 6-3/4" diameter
- dishwasher safe 62-0161 plate

#### stainless steel plate guard



- for people who push food off the plate
- · clip onto any circular dinner plate
- plastic bumpers fit plate diameters 9 -11"
- dishwasher and autoclave safe

62-0170 translucent plastic 62-0171 white plastic



- prevents food from being pushed off plate
- three spring action clips secure to plate
- fits from 9-11" diameter plates

62-0172 plate guard

#### bib clip and cord



- make "do-it-yourself" bibs out of napkins, towels, plastic film, or other materials.
- consists of braided nylon cord with a spring tension pinch clip at each end
- long version has a sliding plastic bead to make it length adjustable

62-0200 short bib clip, pair 62-0201 long bib clip



# eating utensils

#### built-up handle stainless steel utensils

- comfort handle easy to hold
- hole at one end can be used to attach a wrist strap
- top rack dishwasher safe

#### straight utensils



61-0000 knife 61-0001 fork 61-0002 teaspoon 61-0003 soup spoon

61-0009 complete set

#### **UBend-It™**



 a twist in the shaft allows bending either side at any angle

61-0014 fork 61-0015 tablespoon

# weighted cutlery



- each utensil weighs 7.3 oz.
- comfortable grip
- soup spoon has recessed edge and deep bowl to prevent spilling
- straight stainless steel

61-0020 knife 61-0021 fork

61-0022 teaspoon 61-0023 soup spoo

#### teaspoon uter



- white food grade polypropylene cutting board
- stainless steel spikes hold food in place for cutting
- 1/2" high corner guards
- 4 suction cups anchor board

61-0200 board, 12x12"





- each utensil weighs 8 oz.
- use to reduce hand tremors

61-0030 straight knife

61-0031 straight fork

61-0032 straight teaspoon 61-0033 straight soup spoon

specify left (L) or right (R) with sku

61-0036L/R fork

61-0037L/R teaspoon

61-0038L/R soup spoon

#### utensil holders



61-0100 hand clip



- webbing with utensil pocket
- adjustable velcro closure

61-0110 pediatric, beige 61-0111 sm., med.,beige

61-0112 large, beige

#### comfort grip

- lightweight, 3 oz.
- soft vinyl coating with finger indentations for secure grip
- knife can also be used as a rocker knife, and is suitable for one handed use



straight

61-0040 knife 61-0041 fork 61-0042 teas

61-0042 teaspoon 61-0043 soup spoon 61-0044 small teaspoon

left / right handed



specify left(L) or right(R) with sku 61-0046L/R fork 61-0047L/R teaspoon 61-0048L/R soup spoon

#### featherlite



- very light weight 1.7 oz.
- handle is 1" in diameter and slightly tapered for easy grip
- knife can be used regularly or as a rocker

#### straight

61-0050 straight knife 61-0051 straight fork 61-0052 straight teaspoon 61-0053 straight soup spoon

#### specify left (L) or right (R) with sku

61-0055L/R knife

61-0056L/R fork

61-0057L/R teaspoon

61-0058L/R soup spoon

#### soft handle



- large, foam grip handles
- foam grip for easy cleaning
- lightweight

#### straight

61-0060 knife

61-0061 fork

61-0062 teaspoon

61-0063 soup spoon

#### specify left(L) or right(R) with sku

61-0066L/R fork

61-0067L/R teaspoon

61-0068L/R soup spoon

#### knife / fork combo



 rocker knife and fork with comfort grip; cut food with knife then use fork end for eating

61-0070 knife /fork

#### roller knife / fork



 roller knife / fork with comfort arip

61-0071 knife /fork

#### rocker knife



- cut meat with rocker motion using only one hand
- curved blade, solid handle

61-0072 rocker knife

#### meat cutter knife



 cut meat with rocking motion using only one hand

61-0073 meat cutter knife

#### Kennedy J- cup



- lightweight spillproof drinking cup
- ideal for people of all ages
- ergonomic "J" handle
- · won't spill even when turned upside down
- useable with any style straw
- dishwasher safe

60-1000 7 oz. cup 60-1001 replacement lid





- · designed with a nose cutout to encourage correct head position
- allows user to drink with little or no shoulder flexion, wrist extension or head/neck movement even when wearing a cervical collar
- dishwasher safe plastic

60-1050	[c] 8oz., clear
60-1051	[b] 8oz., light blue
60-1052	[a] 8oz., sandston


# [g]

60-1040 [d] 4 oz. [e] 8 oz. 60-1041 60-1042 [f] 12 oz.

[g] 2-handled, 8 oz. 60-1043

#### spillproof cups and lids





- cups and mugs have up to 8 ounce, capacity
- · useable with or without antisplash lid
- · dishwasher safe plastic

60-1010	tumbler
60-1011	1 handle
60-1012	2 handle



#### spout lid

- · fits comfortably in mouth
- flow of liquid can be regulated by placing a finger over hole
- dishwasher safe plastic

60-1030	spout
60-1031	anti-splash
60-1032	vacuum

#### Thumbs-Up™ cup



- design aids in raising cup without wrist motion
- double-wall design protects hands from heat and cold while retaining liquid temper-
- holds 8 ounces, thumbs-on grip reduces spilling due to trembling
- top rack dishwasher and microwave safe

60-1060 cup and lid 60-1061 replacement lid

#### Hand-to-Hand™ mug



- · contoured grip; easy to hold with one or two hands
- · double-wall protects hands from heat and cold while maintaining inside tempera-
- holds 10 ounces
- top rack dishwasher and microwave (cup only) safe

60-1070 cup and lid 60-1071 without lid replacement lid 60-1072

#### mugs and lids



- insulated cup and mug keep food/drink warm or cold
- weighted base helps reduce hand tremors
- mug color is granite, cup is blue

#### insulated

60-1080 insulated mug, no-spill lid, 12 oz. 60-1085 insulated cup, no-spill lid, 8 oz. 60-1086 insulated cup, tube lid, 8 oz.

#### weighted



weighted base cup

60-1200 weighted mug, no-spill lid, 12 oz. 60-1205 weighted cup, no-spill lid, 8 oz. 60-1206 weighted cup, tube lid, 8 oz.

#### lids

60-1090 no-spill lid for cup/mug, pkg/3 60-1091 tube lid for cup/mug, pkg/3

#### No-tip cup keeper



- · prevents cups from tipping
- fits cups and glasses with 21/4" base
- 5" diameter at base

60-1250 no-tip cup keeper

#### vacuum feeding cup



- · permits dribble-free drinking without raising head
- · mouthpiece large enough to use with straw
- graduated at 4, 6, and 8 oz.
- dishwasher and autoclave safe


60-1210 vacuum cup 60-1210-10 box of 10




#### bedroom

#### bed assists









- easy to grasp handles can be used for repositioning and turning in bed
- · bars slide between mattress and box spring
- handles are 20" wide and height adjustable from 16-19"
- side pockets provide convenient storage areas

single\* bar extends underneath mattress; includes strap that attaches to bed frame for additional stability; fits twin to king size beds 86-0100 single **double\*** bar installs on both sides of bed; fits twin to queen size beds 86-0101 double

\* not for use on mechanical beds

#### blanket lift bar



- protects sensitive feet from weight of a blanket
- creates free space of about 30x16x12"
- one side slips under mattress; the other side supports blanket
- plastic coated steel, 16x12x16"

86-0110 blanket lift bar

#### bed tray



- has 16" wide leg room with flat top surface of 13x10"
- large side pockets
- side handles makes carrying easy

86-0120

bed tray, 22x12x7"

#### bed rope ladder



- provides assistance in sitting up from a lying position
- has 5 rungs at 5" intervals
- looped at end for easy attachment
- total length is 57"

86-0130 rope ladder

#### shampoo rinse basin



- large, rigid plastic rinse basin with recessed neck rest
- side drain hose with a shut-off clamp
- open for continuous flushing; close to retain water

86-0200 rinse basin, 21x18x6"

#### head wash system





- double walled basin has recessed neck rest and a plug that allows used water to drain into reservoir below basin
- after washing, secure drain plug then carry the basin away to be emptied
- 2 gallon capacity

86-0230 head wash system

#### shampoo basin



- inflatable tray provides support for head, neck and shoulders during shampooing
- 8" double tube design prevents spilling
- drain hose conveniently drains water from basin

86-0210 shampoo basin

#### inflatable crescent basin



- at shampoo time remove pillow and replace with the crescent-shaped basin; no repositioning is necessary
- large drain tube for emptying
- inflatable plastic basin

86-0220 crescent basin



#### medication

#### jumbo pill box

#### pill crusher

#### tablet crusher





- each compartment holds 65 ASA size pills
- Braille on lids

85-0100 jumbo pill box



- · small enough to carry in purse
- crusher mechanism doubles as a small water cup
- dishwasher safe

85-0110 pill crusher



- crush any tablet
- mix with liquid to swallow
- crushing cavity holds a 1" diameter plastic cup (12 supplied)

85-0120 tablet crusher 85-0121 cups, 100 each

#### Autodrop® eyedrop guide



- · correctly positions bottle over eye
- prevents blinking by keeping lower lid open
- unique pinhole directs eyesight away from drops

85-0130 Autodrop®



#### tube squeezer

- · accepts any metal or plastic tube 2"wide
- twist key to squeeze out contents

85-0150 tube squeezer 85-0150-2 bag of 2

#### Autosqueeze™ bottle squeeze



- provides leverage to make squeezing bottles easier
- can be used with Autodrop®

85-0140 Autosqueeze™

#### Easy Out squeezer




- for one-handed use
- squeeze contents complete out from bottom up

85-0151 bag of 2

#### hygienic

#### Bottom Buddy™ toilet tissue aid





- designed for self-cleaning when reaching is difficult
- soft, flexible head grips any tissue or
- soiled paper is released into the bowl with a push of a button

85-0200 Bottom Buddy™ aid

#### Self Wipe® toilet aid



- rounded design allows for more accurate placement and pressure
- tissue is rolled then placed around the angled clamp
- after use the tissue is released by pressing button on handle
- rotating handle allows correct amount of tissue to be dispensed
- 10" handle

85-0201 Self Wipe® aid



### transfer

#### padded transfer belts



side release buckle

- these padded transfer belts assist patient transfers
- allows for the comfort and dignity of patient
- safety for both patient and caregiver
- pads are 6" wide with four handle loops spaced for safe grip points
- sizes color-coded for quick identification in a rehab department
- available with plastic side release SR, or the extra safe auto buckle



auto buckle

side release buckle

50-5120S	small, blue	24"-30"
50-5120M	medium, green	28"-48"
50-5120L	large, black	40"-64"


#### auto buckle

50-5121S	small, blue	24"-30"
50-5121M	medium, green	28"-48"
50-5121L	large, black	40"-64"

#### gait belts

• made of heavy duty 2" wide cotton webbing to provide secure hold for therapists and caregivers while working with patients

#### metal buckles



- nickle chrome plated steel buckle
- available in six colors and eight sizes

50-5130-32	32" belt
50-5130-36	36" belt
50-5130-40	40" belt
50-5130-44	44" belt
50-5130-48	48" belt
50-5130-50	50" belt
50-5130-60	60" belt
50-5130-72	72" belt

#### quick release plastic buckles



• press top and bottom buttons to release

50-5131-32 50-5131-36 50-5131-40 50-5131-44	32" QR belt 36" QR belt 40" QR belt 44" QR belt
50-5131-48	48" QR belt
50-5131-50	50" QR belt
50-5131-60	60" QR belt
50-5131-72	72" QR belt

#### safety quick release buckles



- · cover plate prevents accidental opening
- acetal plastic buckle

50-5132-36	36" SQR belt
50-5132-40	40" SQR belt
50-5132-48	48" SQR belt
50-5132-50	50" SQR belt
50-5132-60	60" SQR belt
50-5132-72	72" SQR belt

#### standing / sitting turntable

- $\bullet\,$  rotating pivot disk makes standing or seated transfers easy and safe
- user can sit on turntable and turn on bed or chair
- useful for getting in and out of car seats
- non-skid disk at bottom provides a secure controlled motion
- can be used in conjunction with a gait belt
- 15" diameter, 280lb. weight capacity

50-5150 turntable







#### sensors

# A 59-0100 B 59-0101 C 59-0102


#### seat belt sensors

- alarm sounds when patient opens buckle
- click-release (automobile type) seatbelt buckle
- Early Warning belt has extra 12" of Velcro to remove before removing belt
- · alarm easily attaches to wheelchair back
- an alternative to chair sensor pads

59-0100 [A] buckle release seat belt

59-0101 **[B]** EZ release seat belt 59-0102 **[C]** Early warning EZ release belt

59-0109 12" belt extender



#### pad sensors

- thinner and more sensitive making them better for lighter patients
- less compression foam for faster recovery when pressure is released
- no air holes make pad liquid resistant
- · attaches with sticky strips on pad
- 6 month warranty

59-0110 chair pad 10x15" 59-0110-10 chair pad 10x15" 10 ea 59-0111 bed pad 11x30" 59-0111-10 bed pad 11x30" 10 ea

# cable

proximity sensor

patient safety devices

#### proximity pad sensors

- works with patient on chair or bed
- · provides earlier fall warning
- works on all air support surfaces
- no weight requirements; one size fits all, cut to desired length
- undetectable to patients, thin in design

59-0120 Mylar proximity sensor 59-0120-10 proximity sensor, 10 ea 59-0121 universal connector 59-0122\* cable (specify suffix)\*

\*specify cable when ordering (030, 030R, 040, 050R, 060R, 070, 110R, 120)

#### alarms

#### TR2 alarm



- tamper resistant
- alarm always active (no on/off switch)
- connect to sensor pads or nurse call
- select 85 or 100 db volume alarm tone or silent
- ABS mounting holster with mounting clip
- · battery check system
- 4 AA batteries with low battery indicator 59-0200 TR2 alarm 59-0200-10 TR2 alarm 10 ea

#### IQ Duo plus alarm



- functions as pull cord and sensor alarm
- silent setting that can be used in conjunction with the nurse call connector
- select 85 or 100 db volume/alarm tone
- tamper resistant: no on/off switch; hidden reset button
- mounting holster with strap
- 3-AAA batteries; low battery indicator

59-0210 IQ Duo Plus alarm59-0210-10 IQ Duo Plus alarm 10 ea

#### IQ mini pull cord alarm



- magnetic pull cord releases in any direction
- ideal for single patient use
- easily attaches to wheelchair or bedside
- non replaceable battery included

59-0220 mini pull cord alarm 59-0220-10 mini pull cord alarm 10 ea

#### basic magnetic pull cord alarm



- magnetic cord, pin or connect to nurse call
- select 85 or 100 db volume alarm tone
- comes with mounting clip
- 2-AAA batteries with low battery indicator 59-0230 basic pull cord alarm

59-0230-10 basic pull cord alarm, 10 ea

#### basic pull pin alarm



- use magnetic cord, pin or connect to nurse call system
- alarm sounds locally and at station when using nurse call feature
- select 85 or 100 db volume alarm tone
- comes with mounting clip
- 9 volt battery (not included) with low battery indicator

59-0240 basic pull pin alarm 59-0240-10 basic pull pin alarm 10 ea

#### voice alarm



- plays pre-recorded message of up to 20 seconds when alarm is triggered
- select 85 or 100 db vol. or silent alarm tone
- use magnetic pull cord or connect to sensor pads, seatbelts or nurse call system
- comes with mounting holster
- 4-AA batteries with low battery indicator

59-0250 voice alarm 59-0250-10 voice alarm 10 ea



			sku-sor	ted pag	ge inde	X		
item page	item page	item page	item page	item page	item page	item page	item page	item page
00-1002 103	00-2783 113	00-7570 114	04-2182 111	10-0215 26	10-0601 32	10-0833 19	10-0967 22	10-1183 59
00-1004 103 00-1006 103	00-2784 121 00-2786 112	00-9040 143 00-9050 143	04-2183 111 04-2901 125	10-0216 26 10-0217 26	10-0602 32 10-0603 32	10-0834 19 10-0835 19	10-0968 22 10-0969 22	10-1184 43 10-1185 43
00-1008 103	00-2788 112	00-9080 143	04-3013 127	10-0218 26	10-0604 32	10-0836 19	10-0970 22	10-1191 58
00-1010 103 00-1012 103	00-2791 112 00-2793 112	00-9296 131 01-1432 129	04-3014 127 04-3015 127	10-0219 26 10-0250 26	10-0606 32 10-0616 32	10-0837 19 10-0838 19	10-0994 22 10-0995 22	10-1192 58 10-1193 58
00-1012 103	00-2795 112	01-1432 129	04-3016 127	10-0250 26	10-0640 32	10-0839 19	10-1000 16	10-1193 58
00-1026 103	00-2797 112	01-3100 125	04-3017 127	10-0252 26	10-0660 32	10-0840 24	10-1001 16	10-1230 24
00-1062 103 00-1064 103	00-2802 103 00-2883 121	01-3101 125 01-3102 125	04-3018 127 04-3019 127	10-0253 26 10-0254 26	10-0661 32 10-0662 32	10-0841 24 10-0842 24	10-1002 16 10-1003 16	10-1236 16 10-1237 16
00-1102 103	00-2884 113	01-3103 125	04-3020 127	10-0255 26	10-0663 32	10-0843 24	10-1004 16	10-1238 16
00-1104 103	00-2891 113	01-3104 125	04-3021 127	10-0256 26	10-0664 32 10-0667 32	10-0844 24	10-1006 16	10-1239 16
00-1108 103 00-1110 103	00-2893 113 00-2894 113	01-3105 125 01-3106 125	04-3023 127 04-3032 127	10-0257 26 10-0258 26	10-0667 32 10-0668 32	10-0845 24 10-0850 24	10-1007 16 10-1008 16	10-1240 16 10-1270 24
00-1118 103	00-2896 121	01-3107 125	04-3034 127	10-0300 27	10-0669 32	10-0851 24	10-1009 16	10-1271 24
00-1120 103 00-1124 103	00-2911 113 00-3102 105	01-3108 125 01-3109 125	04-3035 127 04-3036 127	10-0301 27 10-0302 27	10-0673 32 10-0674 32	10-0852 24 10-0853 24	10-1010 16 10-1016 16	10-1272 24 10-1273 24
00-1124 103	00-3319 143	01-3110 125	04-3092 127	10-0302 27	10-0691 32	10-0854 24	10-1017 16	10-1273 24
00-1131 103	00-3320 143	01-3111 125	04-3093 127	10-0304 27	10-0710 33	10-0855 24	10-1018 16	10-1275 24
00-1132 103 00-1200 104	00-3321 143 00-3322 143	01-3112 125 01-3126 125	04-3149 127 04-3194 127	10-0305 27 10-0350 29	10-0711 33 10-0712 33	10-0856 24 10-0857 24	10-1019 16 10-1036 16	10-1276 24 10-1290 21
00-1202 104	00-3323 143	01-3127 125	04-3907 122	10-0351 29	10-0713 33	10-0859 24	10-1037 16	10-1291 21
00-1204 104	00-3600 124	01-3128 125	04-3911 122	10-0352 29	10-0714 34	10-0860 24	10-1041 16	10-1292 21
00-1206 104 00-1208 104	00-4002 114 00-4004 114	01-3129 125 01-3214 110	04-8038 124 04-8082 124	10-0353 29 10-0354 29	10-0715 34 10-0716 34	10-0861 24 10-0862 24	10-1042 16 10-1074 17	10-1293 21 10-1310 17
00-1210 104	00-4010 103	01-8010 35	06-0030 125	10-0355 29	10-0720 33	10-0863 24	10-1075 17	10-1311 17
00-1212 104	00-4012 140	01-8025 34 01-8030 34	06-2357 143	10-0356 29 10-0357 29	10-0722 33 10-0724 33	10-0864 24 10-0865 24	10-1076 17	10-1312 17
00-1213 104 00-1214 104	00-4016 103 00-4018 114	01-8030 34 02-0533 129	06-2358 143 06-2449 124	10-0357 29 10-0358 29	10-0724 33 10-0725 33	10-0865 24 10-0866 24	10-1077 17 10-1078 17	10-1313 17 10-1314 17
00-1216 104	00-4020 106	02-0575 124	06-2587 124	10-0359 29	10-0733 18	10-0870 24	10-1080 16	10-1315 17
00-1218 104 00-1219 127	00-4026 106 00-4029 106	02-0660 130 02-0661 130	06-2596 124 06-5988 124	10-0360 29 10-0361 29	10-0734 18 10-0735 18	10-0871 24 10-0872 24	10-1081 16 10-1089 17	10-1316 17 10-1320 17
00-1219 127	00-4029 100	02-0001 130	06-7149 143	10-0362 29	10-0736 18	10-0872 24	10-1009 17	10-1320 17
00-1301 125	00-4208 111	02-0711 129	07-7060 148	10-0363 29	10-0737 18	10-0874 24	10-1091 17	10-1322 17
00-1335 110 00-1341 125	00-4210 111 00-4230 103	02-0715 130 02-0716 130	07-7061 148 07-7062 148	10-0383 29 10-0384 29	10-0738 18 10-0739 18	10-0875 24 10-0876 24	10-1092 17 10-1093 17	10-1323 17 10-1324 17
00-1341 125	00-4238 111	02-0719 130	07-7063 148	10-0440 29	10-0740 18	10-0880 24	10-1099 42	10-1325 17
00-1432 125	00-4248 111	02-0723 130	07-7064 148	10-0441 29	10-0741 18	10-0881 24	10-1100 42	10-1326 17
00-1433 125 00-1440 125	00-4260 111 00-4262 111	02-2100 140 02-7133 112	07-7065 148 07-7066 148	10-0442 29 10-0443 29	10-0742 18 10-0743 18	10-0882 24 10-0883 24	10-1101 42 10-1102 42	10-1327 17 10-1328 17
00-1500 105	00-4266 111	02-7312 112	07-7067 148	10-0444 29	10-0744 18	10-0884 24	10-1103 42	10-1330 21
00-1502 105	00-4272 147	02-7313 112	07-7600 119	10-0445 29	10-0745 18	10-0885 24	10-1104 43	10-1331 21
00-1504 105 00-1506 105	00-4280 147 00-4281 147	02-7314 112 02-7315 112	07-7609 119 07-7610 119	10-0470 24 10-0471 24	10-0746 18 10-0747 18	10-0886 24 10-0900 22	10-1110 32 10-1119 32	10-1332 21 10-1333 21
00-1508 105	00-4282 147	02-7381 112	07-7710 119	10-0472 24	10-0749 18	10-0901 22	10-1120 58	10-1334 21
00-1510 105 00-1512 105	00-4283 147 00-4284 147	02-7382 112 02-7383 112	07-7712 119 07-7715 118	10-0473 24 10-0540 19	10-0750 18 10-0756 20	10-0902 22 10-0903 22	10-1128 43 10-1129 43	10-1341 21 10-1342 21
00-1512 105	00-4285 147	02-7384 112	07-7717 118	10-0541 19	10-0757 20	10-0905 22	10-1129 43	10-1342 21
00-1533 105	00-4293 146	02-7465 112	07-7719 118	10-0542 19	10-0766 91	10-0906 22	10-1131 58	10-1344 21
00-1534 105 00-1552 105	00-4295 147 00-4331 127	02-7467 113 02-7469 112	07-8046 114 07-8047 114	10-0543 19 10-0544 19	10-0775 20 10-0776 21	10-0907 22 10-0908 22	10-1132 43 10-1133 43	10-1350 25 10-1351 25
00-1554 105	00-4333 127	02-7478 113	07-8048 114	10-0545 19	10-0777 21	10-0910 22	10-1140 57	10-1352 25
00-1556 105 00-1562 105	00-4395 122	02-7525 121	08-6049 124	10-0546 19	10-0778 21 10-0779 21	10-0911 22	10-1141 57	10-1353 25
00-1562 105 00-1600 108	00-4500 125 00-4502 125	02-7805 121 02-7808 121	08-7896 143 08-8152 143	10-0550 29 10-0551 29	10-0779 21 10-0784 20	10-0912 22 10-0913 22	10-1142 57 10-1150 58	10-1370 59 10-1371 59
00-2026 129	00-4504 125	02-7811 121	08-9288 143	10-0552 29	10-0785 20	10-0918 22	10-1152 60	10-1372 59
00-2027 130 00-2029 130	00-4506 125 00-4508 125	02-7812 121 02-7814 121	08-9295 143 08-9493 143	10-0553 29 10-0554 29	10-0786 20 10-0787 20	10-0919 22 10-0920 22	10-1153 60 10-1154 60	10-1373 59 10-1374 59
00-2029 130	00-4510 125	02-7816 121	09-0910 105	10-0555 29	10-0790 21	10-0920 22	10-1155 60	10-1374 59
00-2043 129	00-4520 129	02-7840 121	10-0150 26	10-0556 29	10-0791 21	10-0923 22	10-1156 60	10-1380 59
00-2044 129 00-2090 129	00-4759 124 00-6480 142	02-7841 121 04-0128 129	10-0200 26 10-0201 26	10-0557 29 10-0558 29	10-0792 21 10-0793 21	10-0924 22 10-0925 22	10-1157 60 10-1158 60	10-1385 59 10-1460 22
00-2090 129	00-6484 142	04-2159 111	10-0202 26	10-0559 29	10-0794 21	10-0926 22	10-1159 60	10-1461 22
00-2102 103	00-6485 142	04-2160 111	10-0203 26	10-0560 29	10-0796 21	10-0940 22	10-1160 58	10-1462 22
00-2302 103 00-2402 103	00-6630 143 00-6840 124	04-2171 111 04-2172 111	10-0204 26 10-0205 26	10-0561 29 10-0563 29	10-0797 21 10-0798 21	10-0941 22 10-0942 22	10-1161 60 10-1162 91	10-1463 22 10-1464 22
00-2502 103	00-6870 124	04-2173 111	10-0206 26	10-0564 29	10-0800 19	10-0943 22	10-1163 91	10-1466 22
00-2738 113	00-6874 124	04-2174 111	10-0207 26	10-0565 29	10-0801 21	10-0944 22	10-1164 91	10-1467 22
00-2758 112 00-2760 112	00-6876 124 00-6878 124	04-2175 111 04-2176 111	10-0208 26 10-0209 26	10-0566 29 10-0575 29	10-0812 21 10-0813 21	10-0948 22 10-0956 22	10-1167 58 10-1168 60	10-1468 22 10-1469 22
00-2763 112	00-7040 124	04-2177 111	10-0210 26	10-0577 29	10-0815 21	10-0957 22	10-1170 59	10-1470 22
00-2767 112	00-7044 123	04-2178 111 04-2179 111	10-0211 26 10-0212 26	10-0578 29 10-0579 29	10-0816 21	10-0958 22	10-1173 59	10-1471 22
00-2780 112 00-2781 112	00-7050 125 00-7550 114	04-2179 111 04-2180 111	10-0212 26	10-0579 29 10-0581 29	10-0819 24 10-0821 21	10-0959 22 10-0960 20	10-1174 59 10-1179 59	10-1473 22 10-1480 22
	00-7560 114			10-0600 32		10-0962 20	10-1182 59	10-1481 22



						<b>S</b>	sku-s	sor	ted p	oag	ge in	de	X					
10-1484   22   01-1612   28   01-1690   91   01-2291   20   02-2643   22   01-2640   23   01-3180   30   01-50965   14   01-5267   7   10-1494   21   01-1613   28   01-1627   28   01-1627   29   01-1	item	page	item	page	item p									age	item p	age	item	page
10-1486   22   01-6167   28   01-820   19   10-2292   20   02-664   22   02-871   35   03-3181   30   10-5096   14   10-5267   7   10-1468   22   01-871   35   03-3181   30   10-5096   14   10-5267   7   10-1468   22   01-871   35   03-3181   30   10-5096   14   10-5267   7   10-1468   22   01-871   35   03-3181   30   10-5096   14   10-5267   7   10-1469   21   01-1462   23   01-825   19   02-220   16   02-568   22   02-287   36   03-3191   30   05-5100   14   10-5268   7   10-1462   24   01-825   19   02-220   16   02-568   22   02-287   36   03-3191   30   05-5100   15   10-5268   7   10-4691   23   03-485   30   03-5100   15   03-5268   7   10-4691   23   03-485   30   03-5100   15   03-5268   7   10-4691   23   03-485   30   03-5100   15   03-5268   7   10-4691   23   03-485   30   03-5100   15   03-5268   7   10-4691   23   03-485   30   03-5100   15   03-5268   7   10-4691   23   03-485   30   03-485   30   03-5100   15   03-5268   7   10-4691   23   03-485   30   03-485   30   03-5100   15   03-5268   7   10-5268																		
10-1469   21   0-1616   22   0-1621   9   10-2293   20   0-2646   22   0-2877   35   10-3182   30   10-5997   14   10-5020   7   10-1499   21   10-1619   28   10-1623   19   10-2295   20   10-2695   22   10-2673   35   10-3183   30   10-5998   14   10-2632   7   10-1491   21   10-1619   28   10-1623   19   10-2295   20   10-2695   22   10-2673   35   10-3183   30   10-5998   14   10-2633   7   10-1494   21   10-1619   28   10-1825   19   10-2201   16   10-2685   22   10-2293   36   10-3193   30   10-5103   16   10-2635   7   10-1494   21   10-1622   28   10-1855   19   10-2201   16   10-2654   22   10-2293   36   10-3193   30   10-5103   16   10-2635   7   10-1495   21   10-1622   28   10-1855   19   10-2201   16   10-2654   22   10-2293   36   10-3193   30   10-5103   16   10-2655   7   10-1495   21   10-1622   28   10-1855   19   10-2203   16   10-2658   22   10-2293   36   10-3193   30   10-5103   16   10-2657   7   10-1495   21   10-1622   28   10-1855   19   10-2203   16   10-2658   22   10-2293   36   10-3193   30   10-5103   16   10-2626   7   10-1622   28   10-1623   28   10					l .													
10-1494   21   10-1616   28   10-1823   19   10-2286   20   10-2650   22   10-2873   35   10-3184   30   10-5096   14   10-26283   7   10-1495   21   10-1616   28   10-1824   39   10-2300   16   10-26285   22   10-2827   38   10-3181   30   10-5101   16   10-26286   7   10-1495   21   10-1616   28   10-1825   39   10-3200   16   10-26285   22   10-2827   38   10-3181   30   10-5101   16   10-26286   7   10-1495   21   10-1610   25   10-1622   28   10-1825   39   10-2620   24   10-1610   25   10-1622   28   10-1825   39   10-2620   24   10-1611   25   10-1622   28   10-1825   39   10-2308   16   10-2308   16   10-2685   22   10-2825   36   10-3183   30   10-5101   15   10-2620   24   10-1611   25   10-1622   28   10-1825   39   10-2308   16   10-2308   16   10-2680   22   10-2825   36   10-3183   30   10-5101   15   10-2620   14   10-1611   25   10-1622   28   10-1825   39   10-2308   16   10-2308   16   10-2680   22   10-2845   16   10-2308   16   10-2620   24   10-1611   25   10-1622   28   10-1825   39   10-2308   16   10-2680   22   10-2845   16   10-2308   16   10-2680   22   10-2845   16   10-2308   16   10-2680   22   10-2845   16   10-2308   16   10-2680   22   10-2845   16   10-2308   16   10-2680   12   10													10-3182				10-5261	7
10-1492   21   10-1617   28   10-1627   19   10-2296   20   10-2651   22   10-2650   23   10-2650   30   10-5100   15   10-2264   7   10-1618   28   10-1620   31   10-1620   31   10-1620   32   10-1620   31   10-1620   32   10-16																		
10-1499   21   10-1619   28   10-1825   19   10-2300   16   10-2505   27   10-2505   27   10-1406   21   10-1620   28   10-1825   19   10-2302   16   10-2505   22   10-2502   36   10-1813   30   10-5101   15   10-2505   77   10-1610   28   10-1620   28   10-1																		
10-1469   21   10-1620   28   10-1801   10   10-2302   16   10-2503   16   10-2505   22   10-2505   26   10-2505   27   10-2505   27   10-1513   25   10-1602   28   10-1805   28   10-1																		
10-1410   25   10-1622   28   10-1825   29   10-2303   16   10-2305   16   10-2655   22   10-2262   36   10-318   30   10-5106   15   10-2220   4   10-1617   25   10-1623   28   10-1825   19   10-2305   16   10-2660   22   10-2640   15   10-3186   30   10-5106   15   10-2220   4   10-1617   25   10-1623   28   10-1825   19   10-2305   16   10-2660   22   10-2640   15   10-3186   30   10-5106   15   10-2220   4   10-1617   25   10-1618   26   10-1620   27   10-1618   26   10-1618																		
10-151  25   10-1622   28   10-1825   19   10-2905   16   10-2660   22   10-2940   16   10-3195   30   10-5106   15   10-5284   4																		
10-1512   25   10-1623   28   10-1826   19   10-2206   16   10-2661   22   10-2941   15   10-3197   30   10-5106   15   10-5281   7   10-1513   25   10-1630   28   10-1806   19   10-2308   16   10-2662   22   10-2441   15   10-3198   30   10-5101   15   10-5291   7   10-1513   25   10-1640   28   10-1806   19   10-2309   16   10-2662   22   10-2441   15   10-3198   30   10-5101   15   10-5292   7   10-1610   28   10-1610																		
10-1514   25   10-1630   28   10-1806   19   10-2309   16   10-2663   22   10-2443   15   10-3199   30   10-1510   15   10-5293   7						19												
10-154  25   10-164  28   10-1860   29   10-2809   27   10-2664   22   10-2944   15   10-399   30   10-5120   15   10-5293   77																		
10-1516   25   10-1641   28   10-189   21   10-250   27   10-266   22   10-294   15   10-3201   10   10-512   15   10-5295   7   10-1519   25   10-1643   28   10-1892   21   10-2502   27   10-267   23   10-289   28   10-3202   10   10-5152   10   10-5295   7   10-1519   25   10-1643   28   10-1892   21   10-2502   27   10-267   23   10-289   28   10-3202   10   10-5152   10   10-5296   7   10-1612   25   10-1650   28   10-1893   21   10-2503   27   10-267   23   10-289   28   10-3203   10   10-5162   11   10-5301   12   10-5162   10   10-5162   11   10-5301   12   10-1612   12   10-					l .													
10-15   25   10-1643   28   10-1893   21   10-2502   27   10-2670   23   10-2980   28   10-3204   10   10-5152   10   10-52697   7   10-1521   25   10-1650   28   10-1894   21   10-2504   27   10-2670   23   10-2982   28   10-3205   10   10-5162   11   10-5301   10   10-1522   25   10-1652   28   10-1895   28   10-189	10-1515	25	10-1641	28	10-1890	21	10-2500	27	10-2664	22	10-2944	15	10-3201		10-5121		10-5294	7
10-15 2  25   10-1645   28   10-1849   21   10-2504   27   10-2676   23   10-2864   28   10-3205   10   10-5160   11   10-5300   12   10-1522   25   10-1651   28   10-1896   21   10-2506   27   10-2686   23   10-2896   28   10-3205   10   10-5161   11   10-5300   12   10-1523   25   10-1652   28   10-1896   21   10-2506   27   10-2686   23   10-2896   28   10-3205   10   10-5161   11   10-5300   12   10-1523   25   10-1652   28   10-1896   21   10-2506   27   10-2686   23   10-2896   28   10-3225   10   10-5161   11   10-5320   6   10-1523   25   10-1653   28   10-1904   40   10-2506   27   10-2686   23   10-2896   28   10-3225   10   10-5161   11   10-5320   6   10-1523   25   10-1706   42   10-1932   17   10-2516   27   10-2716   23   10-3065   30   10-3234   10   10-5168   11   10-5335   62   10-1706   42   10-1934   17   10-2516   27   10-2716   23   10-3065   30   10-3234   10   10-5168   11   10-5336   12   10-1531   25   10-1706   42   10-1941   17   10-2516   27   10-2716   23   10-3096   30   10-3246   40   10-5168   11   10-5336   12   10-1534   25   10-1706   42   10-1944   17   10-2516   27   10-2716   23   10-3096   30   10-3607   36   10-5196   5   10-5366   13   10-1536   25   10-1707   42   10-1941   17   10-2516   27   10-2712   23   10-3096   30   10-3607   36   10-5196   5   10-5366   13   10-1536   25   10-1707   42   10-1941   17   10-2516   27   10-2712   23   10-3096   30   10-3607   36   10-5196   5   10-5366   13   10-1536   25   10-1710   42   10-1941   17   10-2516   27   10-2712   23   10-3096   30   10-3607   36   10-5196   5   10-5366   13   10-1546   25   10-1710   42   10-1941   17   10-2516   27   10-2712   23   10-3096   30   10-3607   36   10-5196   5   10-5366   13   10-1546   25   10-1710   42   10-1941   17   10-2516   27   10-2712   23   10-3096   30   10-3607   36   10-5196   5   10-5366   13   10-1546   25   10-1710   42   10-1941   17   10-2516   27   10-2712   23   10-3096   30   10-3607   36   10-5196   5   10-5366   13   10-1546   25   10-1716   42   10-1946   17																		
10-1522   25   10-1650   28   10-1894   21   10-2506   27   10-2606   27   10-2606   28   10-2905   28   10-2305   10   10-5162   11   10-5300   12   10-1523   25   10-1652   28   10-1904   01   10-2507   27   10-2608   28   10-2905   28   10-2221   10   10-5163   11   10-5300   6   10-1525   25   10-1654   28   10-1904   01   10-2507   27   10-2608   27   10-2608   28   10-2905   28   10-2221   10   10-5163   11   10-5320   6   10-1525   25   10-1654   28   10-1904   01   10-2507   27   10-2608   27   10-260					l .													
10-1523   25   0-1652   28   0-1994   01   0-2506   27   10-2881   23   0-2985   28   0-3222   10   10-5165   11   10-5320   6   10-1526   25   10-1666   28   0-1994   01   0-2507   27   10-2808   23   10-2806   28   0-3222   10   0-5166   11   0-5333   6   10-1526   25   0-1666   04   0-1931   17   10-2508   27   10-2710   23   10-3006   30   10-3231   10   10-5166   11   10-5333   6   10-1508   25   10-1700   42   10-1932   17   10-2510   27   10-2710   23   10-3006   30   10-3232   10   10-5166   11   10-5333   12   10-1533   25   10-1700   42   10-1933   17   10-2512   27   10-2713   23   10-3006   30   10-3232   10   10-5166   11   10-5333   12   10-1533   25   10-1700   42   10-1934   17   10-2512   27   10-2713   23   10-3006   30   10-3233   10   10-5169   11   10-5340   12   10-1533   25   10-1706   42   10-1943   17   10-2514   27   10-2715   23   10-3006   30   10-3234   10   10-5195   5   10-5346   13   10-1533   25   10-1706   42   10-1943   17   10-2514   27   10-2715   23   10-3006   30   10-3203   30   10-3610   36   10-5195   5   10-5346   13   10-1533   25   10-1706   42   10-1944   17   10-2516   27   10-2716   23   10-3003   30   10-3600   36   10-5195   5   10-5346   13   10-1543   25   10-1716   42   10-1944   17   10-2516   27   10-2721   23   10-3003   30   10-3600   36   10-5195   5   10-5366   13   10-1542   25   10-1716   42   10-1945   17   10-2516   27   10-2721   23   10-3003   30   10-3600   36   10-5195   5   10-5366   13   10-1542   25   10-1716   42   10-1945   17   10-2516   27   10-2721   23   10-3003   30   10-3600   36   10-5195   5   10-5366   13   10-1542   25   10-1716   42   10-1945   17   10-2516   27   10-2722   23   10-3100   31   10-3600   36   10-5195   5   10-5366   13   10-1542   25   10-1716   42   10-1945   17   10-2516   27   10-2723   23   10-3100   31   10-3600   36   10-5201   4   10-5376   13   10-1542   25   10-1716   42   10-1945   17   10-2516   27   10-2726   23   10-3100   31   10-3600   36   10-5201   4   10-5376   13   10-1400   10-1402   25   10-																		
10-1542   25   0-1656   28   0-1904   40   0-2507   27   0-2682   23   0-2966   28   0-3222   10   0-5166   11   0-5325   6   0-1526   25   0-1668   20   0-1920   13   0-2508   27   0-2680   23   0-3020   0-3230   10   0-5166   11   0-5336   12   0-1566   11   0-5336   12   0-1566   11   0-5336   12   0-1568   11   0-1568   11   0-5336   12   0-1568   11   0-1568   11   0-5336   12   0-1568   11   0-1568   11   0-5336   12   0-1568   11   0-156																		
10-1526   25   10-1664   28   10-1920   18   10-2508   27   10-2709   23   10-3080   30   10-3230   10   10-5166   11   10-5330   12   10-1530   25   10-1700   42   10-1932   17   10-2510   27   10-2710   23   10-3086   30   10-3231   10   10-5168   11   10-5334   12   10-1532   25   10-1702   42   10-1931   17   10-2511   27   10-2712   23   10-3086   30   10-3231   10   10-5168   11   10-5336   12   10-1532   25   10-1702   42   10-1934   17   10-2512   27   10-2713   23   10-3086   30   10-3233   10   10-5189   11   10-5346   12   10-1533   25   10-1706   42   10-1941   17   10-2512   27   10-2713   23   10-3086   30   10-3233   10   10-5189   15   10-5346   13   10-1533   12   10-1533   25   10-1706   42   10-1942   17   10-2514   27   10-2715   23   10-3090   30   10-3080   36   10-5193   5   10-5356   13   10-1535   25   10-1707   42   10-1941   17   10-2516   27   10-2723   31   10-3090   30   10-3080   36   10-5193   5   10-5356   13   10-1535   25   10-1710   42   10-1941   17   10-2516   27   10-2723   23   10-3090   30   10-3080   36   10-5193   5   10-5356   13   10-1540   25   10-1717   42   10-1942   17   10-2516   27   10-2723   23   10-3090   30   10-3080   36   10-5193   5   10-5356   13   10-1540   25   10-1717   42   10-1942   17   10-2516   27   10-2723   23   10-3090   30   10-3607   36   10-5195   5   10-5356   13   10-1541   25   10-1717   42   10-1942   17   10-2518   27   10-2723   23   10-3090   30   10-3607   36   10-5195   5   10-5356   13   10-1543   25   10-1717   42   10-1942   17   10-2518   27   10-2723   23   10-3100   30   10-3607   36   10-5195   5   10-5356   13   10-1543   25   10-1717   42   10-1940   17   10-2518   27   10-2723   23   10-3100   31   10-3637   36   10-5200   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386   38   10-5408   4   10-5386																		
10-1551   25   10-1710   42   10-1932   17   10-2510   27   10-2711   23   10-3066   30   10-3232   10   10-5168   11   10-5336   12   10-1532   25   10-1702   42   10-1934   17   10-2511   27   10-2713   23   10-3066   30   10-3234   10   10-5161   5   10-5346   13   10-1532   25   10-1705   42   10-1942   17   10-2514   27   10-2715   23   10-3060   30   10-3605   36   10-5191   5   10-5346   13   10-1533   25   10-1706   42   10-1942   17   10-2514   27   10-2715   23   10-3090   30   10-3605   36   10-5193   5   10-5356   13   10-1535   25   10-1707   42   10-1942   17   10-2516   27   10-2720   23   10-3090   30   10-3605   36   10-5193   5   10-5366   13   10-1535   25   10-1710   42   10-1944   17   10-2516   27   10-2720   23   10-3090   30   10-3617   36   10-5195   5   10-5365   13   10-1546   25   10-1711   42   10-1951   17   10-2518   27   10-2723   23   10-3090   30   10-3617   36   10-5195   5   10-5365   13   10-1543   25   10-1716   42   10-1952   17   10-2518   27   10-2723   23   10-3091   30   10-3605   36   10-5200   4   10-5375   13   10-1543   25   10-1716   42   10-1954   17   10-2518   27   10-2725   23   10-3100   31   10-3605   36   10-5200   4   10-5385   36   10-1544   25   10-1726   42   10-1954   17   10-2516   27   10-2762   23   10-3101   31   10-3663   36   10-5200   4   10-5382   8   10-1548   25   10-1723   42   10-1961   25   10-2555   27   10-2765   23   10-3131   31   10-3645   36   10-5200   4   10-5382   8   10-1548   25   10-1723   42   10-1961   25   10-2556   27   10-2760   23   10-3131   31   10-3645   36   10-5200   4   10-5480   12   10-1546   25   10-1725   42   10-1970   18   10-2556   27   10-2760   23   10-3131   31   10-3645   36   10-5200   4   10-5480   12   10-1546   25   10-1725   42   10-1970   18   10-2556   27   10-2760   23   10-3131   31   10-3645   36   10-5200   4   10-5480   12   10-1546   25   10-1725   42   10-1970   18   10-2556   27   10-2760   23   10-3136   31   10-3665   36   10-5200   4   10-5480   12   10-1546   25   10-1525   10-1525																		
10-1532   25   10-1701   42   10-1933   17   10-2511   27   10-2712   23   10-3066   30   10-3233   10   10-5169   11   10-53446   13   10-1533   25   10-1706   42   10-1941   17   10-2513   27   10-2714   23   10-3090   30   10-3600   36   10-5192   5   10-53446   13   10-1535   25   10-1707   42   10-1943   17   10-2515   27   10-2716   23   10-3091   30   10-3605   36   10-5193   5   10-5356   13   10-1536   25   10-1707   42   10-1943   17   10-2516   27   10-2721   23   10-3092   30   10-3610   36   10-5193   5   10-5366   13   10-1536   25   10-1710   42   10-1945   17   10-2516   27   10-2721   23   10-3093   30   10-3610   36   10-5194   5   10-5366   13   10-1540   25   10-1716   42   10-1951   17   10-2516   27   10-2721   23   10-3094   30   10-3620   36   10-5198   5   10-5306   13   10-1541   25   10-1716   42   10-1953   17   10-2516   27   10-2721   23   10-3094   30   10-3620   36   10-5198   5   10-5306   13   10-1541   25   10-1716   42   10-1953   17   10-2516   27   10-2722   23   10-3110   31   10-3620   36   10-5201   4   10-5380   88   10-1544   25   10-1716   42   10-1954   25   10-2551   27   10-2761   23   10-3110   31   10-3630   36   10-5201   4   10-5382   88   10-1544   25   10-1722   42   10-1960   25   10-2551   27   10-2761   23   10-3112   31   10-3640   36   10-5203   4   10-5382   88   10-1545   25   10-1724   42   10-1962   25   10-2553   27   10-2762   23   10-3112   31   10-3640   36   10-5203   4   10-5384   88   10-1545   25   10-1726   42   10-1976   18   10-2555   27   10-2762   23   10-3130   31   10-3640   36   10-5203   4   10-5384   88   10-1545   25   10-1726   42   10-1976   18   10-2555   27   10-2762   23   10-3130   31   10-3640   36   10-5205   4   10-5384   88   10-1555   16   10-1726   42   10-1977   18   10-2555   27   10-2762   23   10-3130   31   10-3640   36   10-5205   4   10-5361   88   10-1555   16   10-1726   42   10-1978   18   10-2555   27   10-2762   23   10-3140   31   10-3660   36   10-5206   4   10-5561   88   10-1555   89   10-1555   16   1																		
10-1532   25   10-1702   42   10-1944   17   10-2513   27   10-2713   23   10-3066   30   10-3234   10   10-5191   5   10-53447   12   10-1534   25   10-1706   42   10-1942   17   10-2516   27   10-2715   23   10-3091   30   10-3605   36   10-5193   5   10-5356   13   10-1535   25   10-1710   42   10-1944   17   10-2516   27   10-2721   23   10-3092   30   10-3605   36   10-5195   5   10-5356   13   10-1536   25   10-1710   42   10-1944   17   10-2516   27   10-2721   23   10-3093   30   10-3610   36   10-5195   5   10-5365   13   10-1540   25   10-1715   42   10-1952   17   10-2516   27   10-2721   23   10-3093   30   10-3610   36   10-5195   5   10-5365   13   10-1542   25   10-1716   42   10-1952   17   10-2516   27   10-2723   23   10-3100   31   10-3620   36   10-5196   5   10-5365   13   10-1542   25   10-1716   42   10-1954   17   10-2516   27   10-2723   23   10-3110   31   10-3625   36   10-5200   4   10-5380   8   10-5434   25   10-1717   42   10-1954   17   10-2550   27   10-2762   23   10-3110   31   10-3627   36   10-5200   4   10-5380   8   10-5434   25   10-1723   42   10-1961   25   10-2556   27   10-2761   23   10-3131   31   10-3647   36   10-5202   4   10-5382   8   10-1544   25   10-1723   42   10-1961   25   10-2556   27   10-2761   23   10-3131   31   10-3641   36   10-5204   4   10-5384   8   10-1545   25   10-1723   42   10-1961   25   10-2556   27   10-2761   23   10-3131   31   10-3641   36   10-5204   4   10-5480   10-1552   16   10-1726   42   10-1967   48   10-2556   27   10-2762   23   10-3131   31   10-3640   36   10-5204   4   10-5480   10-1552   16   10-1726   42   10-1967   48   10-2556   27   10-2765   23   10-3135   31   10-3640   36   10-5206   4   10-5480   40-5480																		
10-1533   25   10-1705   42   10-1941   17   10-2516   27   10-2716   23   10-3090   30   10-3600   36   10-5192   5   10-5356   13   10-1535   25   10-1707   42   10-1943   17   10-2516   27   10-2716   23   10-3092   30   10-3610   36   10-5194   5   10-5365   13   10-1536   25   10-1710   42   10-1943   17   10-2516   27   10-2721   23   10-3093   30   10-3610   36   10-5194   5   10-5366   13   10-1540   25   10-1711   42   10-1951   17   10-2517   27   10-2721   23   10-3093   30   10-3617   36   10-5196   5   10-5366   13   10-1541   25   10-1716   42   10-1953   17   10-2517   27   10-2721   23   10-3094   30   10-3620   36   10-5196   5   10-5366   13   10-1541   25   10-1716   42   10-1953   17   10-2518   27   10-2721   23   10-3110   31   10-3620   36   10-5201   4   10-5376   13   10-1541   25   10-1717   42   10-1953   17   10-2519   27   10-2725   23   10-3110   31   10-3620   36   10-5201   4   10-5380   88   10-5201   4   10-5380   88   10-5201   4   10-5380   88   10-1543   25   10-1717   42   10-1960   25   10-2551   27   10-2760   23   10-3120   31   10-3640   36   10-5203   4   10-5382   88   10-1543   25   10-1722   42   10-1960   25   10-2551   27   10-2760   23   10-3121   31   10-3640   36   10-5203   4   10-5384   88   10-1543   25   10-1724   42   10-1962   25   10-2553   27   10-2762   23   10-3132   31   10-3640   36   10-5205   4   10-5450   12   10-1551   16   10-1725   42   10-1976   18   10-2555   27   10-2762   23   10-3132   31   10-3640   36   10-5205   4   10-5450   12   10-1551   16   10-1725   42   10-1976   18   10-2555   27   10-2762   23   10-3132   31   10-3640   36   10-5205   4   10-5450   12   10-1551   16   10-1725   42   10-1976   18   10-2555   27   10-2762   23   10-3143   31   10-3640   36   10-5205   4   10-5510   88   10-1551   18   10-1720   41   10-1978   18   10-2555   27   10-2762   23   10-3143   31   10-3660   36   10-5205   4   10-5510   88   10-1553   18   10-1720   41   10-1978   18   10-2555   27   10-2762   23   10-3143   31   10-3660   36   10-52																		
10-1585   25   10-1770   42   10-1943   17   10-2516   27   10-2720   23   10-3093   30   10-3610   36   10-5195   5   10-5306   13	10-1533	25	10-1705	42		17	10-2513	27	10-2714	23	10-3090		10-3600		10-5192	5	10-5347	12
10-1580   25   10-1710   42   10-1984   17   10-2516   27   10-2721   23   10-3094   30   10-3617   36   10-5198   5   10-5306   13					l .													
10-1540   25   10-1711   42   10-1951   17   10-2517   27   10-2723   23   10-3094   30   10-3602   36   10-5198   5   10-5370   13   10-1541   25   10-1714   42   10-1954   17   10-2519   27   10-2723   23   10-3110   31   10-3602   36   10-5201   4   10-5380   8   10-1542   25   10-1717   42   10-1954   17   10-2519   27   10-2725   23   10-3112   31   10-3603   36   10-5201   4   10-5380   8   10-1544   25   10-1712   42   10-1960   25   10-2552   27   10-2761   23   10-3103   31   10-3604   36   10-5203   4   10-5384   8   10-1544   25   10-1723   42   10-1962   25   10-2552   27   10-2761   23   10-3131   31   10-3641   36   10-5203   4   10-5384   8   10-1546   25   10-1725   42   10-1962   25   10-2553   27   10-2761   23   10-3131   31   10-3641   36   10-5203   4   10-5384   8   10-1565   16   10-1725   42   10-1972   18   10-2555   27   10-2762   23   10-3132   31   10-3645   36   10-5206   4   10-5460   12   10-1552   16   10-1727   42   10-1974   18   10-2555   27   10-2764   23   10-3135   31   10-3645   36   10-5206   4   10-5481   15   10-1553   16   10-1727   42   10-1972   18   10-2555   27   10-2764   23   10-3143   31   10-3660   36   10-5206   4   10-5481   15   10-1553   16   10-1727   42   10-1972   18   10-2556   27   10-2764   23   10-3140   31   10-3660   36   10-5206   4   10-5481   15   10-1553   16   10-1731   41   10-1974   18   10-2556   27   10-2776   23   10-3140   31   10-3660   36   10-5206   4   10-5511   8   10-1570   15   10-1733   41   10-1976   18   10-2565   27   10-2770   23   10-3140   31   10-3660   36   10-5211   4   10-5513   8   10-1570   15   10-1733   41   10-1976   18   10-2600   22   10-2772   23   10-3140   31   10-3660   36   10-5211   4   10-5513   8   10-1570   15   10-1733   41   10-1987   18   10-2600   22   10-2773   23   10-3140   31   10-3660   36   10-5215   4   10-5513   8   10-1570   15   10-1733   41   10-1987   18   10-2600   22   10-2773   23   10-3140   31   10-3660   36   10-5215   4   10-5513   8   10-1570   15   10-1733   41   10-1987   18																		
10-1542   25   10-1716   42   10-1953   17   10-2519   27   10-2725   23   10-3112   31   10-3663   36   10-5201   4   10-5380   8   10-1545   25   10-1712   42   10-1960   25   10-2550   27   10-2765   23   10-3112   31   10-3663   36   10-5203   4   10-5384   8   10-1545   25   10-1723   42   10-1960   25   10-2552   27   10-2761   23   10-3112   31   10-3661   36   10-5203   4   10-5384   8   10-1545   25   10-1723   42   10-1960   25   10-2553   27   10-2761   23   10-3113   31   10-3661   36   10-5205   4   10-5480   15   10-1551   16   10-1725   42   10-1971   18   10-2555   27   10-2764   23   10-3135   31   10-3664   36   10-5205   4   10-5481   15   10-1553   16   10-1727   42   10-1972   18   10-2555   27   10-2764   23   10-3135   31   10-3665   36   10-5206   4   10-5481   15   10-1553   16   10-1727   42   10-1972   18   10-2555   27   10-2764   23   10-3143   31   10-3665   36   10-5206   4   10-5481   15   10-1553   16   10-1733   41   10-1974   18   10-2555   27   10-2770   23   10-3143   31   10-3665   36   10-5208   4   10-5511   8   10-1572   15   10-1733   41   10-1976   18   10-2601   22   10-2772   23   10-3142   31   10-3665   36   10-5211   4   10-5513   8   10-1572   15   10-1733   41   10-1976   18   10-2603   22   10-2773   23   10-3143   31   10-3667   36   10-5213   4   10-5513   8   10-1572   15   10-1733   41   10-1980   18   10-2603   22   10-2776   23   10-3143   31   10-3667   36   10-5213   4   10-5515   8   10-1572   15   10-1733   41   10-1981   18   10-2603   22   10-2776   23   10-3143   31   10-3667   36   10-5213   4   10-5515   8   10-1572   15   10-1735   41   10-1981   18   10-2603   22   10-2776   23   10-3145   31   10-3667   36   10-5213   4   10-5515   8   10-1573   15   10-1740   41   10-1981   18   10-2603   22   10-2776   23   10-3146   31   10-3667   36   10-5215   4   10-5516   8   10-1572   15   10-1735   41   10-1981   18   10-2605   22   10-2786   23   10-3145   31   10-3665   36   10-5215   4   10-5516   8   10-1576   15   10-1745   41   10-1981   18		25		42	l .	17	10-2517	27	10-2722	23	10-3094	30		36		5		13
10-1543   25   10-1717   42   10-1964   17   10-2550   27   10-2762   23   10-3120   31   10-3637   36   10-5202   4   10-5382   8   10-1545   25   10-1723   42   10-1961   25   10-2551   27   10-2761   23   10-3121   31   10-3640   36   10-5204   4   10-5460   12   10-1561   16   10-1724   42   10-1961   25   10-2553   27   10-2763   23   10-3131   31   10-3645   36   10-5206   4   10-5480   15   10-1551   16   10-1726   42   10-1970   18   10-2555   27   10-2763   23   10-3132   31   10-3645   36   10-5206   4   10-5481   15   10-1555   16   10-1726   42   10-1971   18   10-2555   27   10-2763   23   10-3136   31   10-3645   36   10-5206   4   10-5481   15   10-1553   16   10-1727   42   10-1972   18   10-2556   27   10-2763   23   10-3136   31   10-3650   36   10-5206   4   10-5481   15   10-1553   16   10-1730   41   10-1973   18   10-2556   27   10-2776   23   10-3141   31   10-3660   36   10-5206   4   10-5411   8   10-1555   16   10-1731   41   10-1974   18   10-2558   27   10-2770   23   10-3141   31   10-3660   36   10-5210   4   10-5513   8   10-1571   5   10-1733   41   10-1976   18   10-2601   22   10-2773   23   10-3143   31   10-3660   36   10-5214   4   10-5516   8   10-1571   5   10-1733   41   10-1968   18   10-2601   22   10-2773   23   10-3145   31   10-3690   36   10-5214   4   10-5516   8   10-1571   5   10-1734   41   10-1983   18   10-2603   22   10-2776   23   10-3145   31   10-3690   36   10-5214   4   10-5516   8   10-1573   15   10-1741   41   10-1983   18   10-2603   22   10-2778   23   10-3145   31   10-3690   36   10-5214   4   10-5516   8   10-1573   15   10-1741   40   10-1985   18   10-2603   22   10-2778   23   10-3145   31   10-3690   36   10-5214   4   10-5516   8   10-1573   15   10-1741   40   10-1985   18   10-2603   22   10-2786   23   10-3145   31   10-3690   36   10-5214   4   10-5516   8   10-1573   15   10-1741   40   10-1985   18   10-2605   22   10-2786   23   10-3145   31   10-3690   36   10-5216   4   10-5516   8   10-1576   15   10-1761   40   10-1985   18   10																		
10-1544   25   10-1722   42   10-1960   25   10-2551   27   10-2760   23   10-3121   31   10-3640   36   10-5203   4   10-5304   8   10-5364   1																		
10-15646   25   10-1724   42   10-1962   25   10-2553   27   10-2762   23   10-3132   31   10-3645   36   10-5205   4   10-5481   15   10-1552   16   10-1726   42   10-1971   18   10-2555   27   10-2765   23   10-3136   31   10-3650   36   10-5206   4   10-5481   15   10-1553   16   10-1726   42   10-1971   18   10-2555   27   10-2765   23   10-3136   31   10-3650   36   10-5207   4   10-5511   8   10-1555   16   10-1730   41   10-1973   18   10-2555   27   10-2766   23   10-3140   31   10-3665   36   10-5207   4   10-5511   8   10-1555   16   10-1730   41   10-1973   18   10-2557   27   10-2776   23   10-3141   31   10-3665   36   10-5210   4   10-5511   8   10-1575   10-1732   41   10-1975   18   10-2600   22   10-2772   23   10-3143   31   10-3665   36   10-5211   4   10-5514   8   10-1571   15   10-1733   41   10-1976   18   10-2600   22   10-2773   23   10-3143   31   10-3676   36   10-5213   4   10-5514   8   10-1571   15   10-1735   41   10-1980   18   10-2600   22   10-2773   23   10-3144   31   10-3676   36   10-5213   4   10-5516   8   10-1572   15   10-1735   41   10-1980   18   10-2600   22   10-2775   23   10-3146   31   10-3685   36   10-5213   4   10-5516   8   10-1573   15   10-1741   41   10-1981   18   10-2604   22   10-2776   23   10-3146   31   10-3685   36   10-5215   4   10-5516   8   10-1576   15   10-1741   41   10-1983   18   10-2604   22   10-2780   23   10-3146   31   10-3685   36   10-5215   4   10-5516   8   10-1576   15   10-1741   41   10-1983   18   10-2604   22   10-2780   23   10-3146   31   10-3685   36   10-5215   4   10-5516   8   10-1576   15   10-1741   41   10-1983   18   10-2604   22   10-2780   23   10-3146   31   10-3685   36   10-5215   4   10-5516   8   10-1576   15   10-1741   41   10-1983   18   10-2604   22   10-2780   23   10-3146   31   10-3685   36   10-5215   4   10-5516   8   10-1576   15   10-1743   40   10-1983   18   10-2604   22   10-2780   23   10-3146   31   10-3685   36   10-5215   4   10-5516   8   10-1576   15   10-1745   41   10-1983   18   10-260						25												
10-1551   16																		
10-1552   16																		
10-15554   16   10-1730   41   10-1974   18   10-2557   27   10-2770   23   10-3141   31   10-3660   36   10-5210   4   10-5511   8   10-1555   16   10-1731   41   10-1975   18   10-2500   22   10-2771   23   10-3143   31   10-3665   36   10-5211   4   10-5514   8   10-1570   15   10-1732   41   10-1976   18   10-2600   22   10-2772   23   10-3143   31   10-3670   36   10-5213   4   10-5514   8   10-1571   15   10-1733   41   10-1976   18   10-2600   22   10-2773   23   10-3145   31   10-3677   36   10-5213   4   10-5515   8   10-1572   15   10-1735   41   10-1980   18   10-2602   22   10-2774   23   10-3145   31   10-3677   36   10-5215   4   10-5516   8   10-1573   15   10-1740   41   10-1981   18   10-2603   22   10-2776   23   10-3145   31   10-3665   36   10-5215   4   10-5517   8   10-1577   15   10-1741   41   10-1982   18   10-2604   22   10-2781   23   10-3145   31   10-3665   36   10-5215   4   10-5517   8   10-1578   15   10-1742   41   10-1982   18   10-2605   22   10-2781   23   10-3145   31   10-3665   36   10-5215   4   10-5519   8   10-1579   15   10-1742   41   10-1984   18   10-2606   22   10-2781   23   10-3145   31   10-3665   36   10-5215   4   10-5519   8   10-1580   15   10-1742   41   10-1985   18   10-2606   22   10-2781   23   10-3148   31   10-3667   36   10-5215   4   10-5519   8   10-1581   15   10-1745   42   10-1985   18   10-2610   22   10-2781   23   10-3150   30   10-3720   36   10-5218   4   10-5520   8   10-1581   15   10-1751   41   10-1990   18   10-2612   22   10-2785   23   10-3153   30   10-3725   36   10-5221   5   10-5523   8   10-1584   15   10-1751   41   10-1992   18   10-2612   22   10-2785   23   10-3154   30   10-3737   36   10-5221   5   10-5525   8   10-1585   10-1															10-5207			
10-1555   16   10-1731   41   10-1974   18   10-2558   27   10-2771   23   10-3142   31   10-3665   36   10-5211   4   10-5514   8   10-1757   15   10-1733   41   10-1976   18   10-2600   22   10-2773   23   10-3143   31   10-3667   36   10-5213   4   10-5515   8   10-1571   15   10-1733   41   10-1980   18   10-2602   22   10-2773   23   10-3144   31   10-3677   36   10-5213   4   10-5516   8   10-1572   15   10-1735   41   10-1980   18   10-2602   22   10-2774   23   10-3146   31   10-3680   36   10-5214   4   10-5516   8   10-1577   15   10-1741   41   10-1982   18   10-2604   22   10-2778   23   10-3147   31   10-3680   36   10-5216   4   10-5516   8   10-1577   15   10-1741   41   10-1983   18   10-2605   22   10-2781   23   10-3147   31   10-3680   36   10-5216   4   10-5519   8   10-1578   15   10-1743   40   10-1984   18   10-2605   22   10-2781   23   10-3148   31   10-3697   36   10-5216   4   10-5520   8   10-1580   15   10-1743   40   10-1984   18   10-2600   22   10-2782   23   10-3145   30   10-3725   36   10-5216   4   10-5521   8   10-1580   15   10-1743   40   10-1985   18   10-2610   22   10-2783   23   10-3151   30   10-3725   36   10-5220   5   10-5522   8   10-1581   15   10-1750   41   10-1990   18   10-2613   22   10-2785   23   10-3153   30   10-3730   36   10-5220   5   10-5522   8   10-1583   15   10-1753   41   10-1991   18   10-2613   22   10-2795   23   10-3155   30   10-5042   14   10-5224   5   10-5526   8   10-1580   15   10-1753   41   10-1992   18   10-2616   22   10-2795   23   10-3155   30   10-5042   14   10-5225   5   10-5526   8   10-1580   15   10-1763   41   10-1994   18   10-2616   22   10-2797   23   10-3156   30   10-5042   14   10-5225   5   10-5526   8   10-1580   15   10-1763   41   10-1994   18   10-2616   22   10-2797   23   10-3156   30   10-5042   14   10-5225   5   10-5526   8   10-1580   15   10-1763   41   10-1994   18   10-2616   22   10-2797   23   10-3156   30   10-5042   14   10-5225   5   10-5526   8   10-1580   15   10-1765   41   10-1994   18   10-2																		
10-1570   15   10-1732   41   10-1975   18   10-2601   22   10-2772   23   10-3143   31   10-3670   36   10-5212   4   10-5515   8   10-1572   15   10-1733   41   10-1980   18   10-2601   22   10-2774   23   10-3145   31   10-3687   36   10-5214   4   10-5516   8   10-1573   15   10-1740   41   10-1981   18   10-2603   22   10-2775   23   10-3146   31   10-3686   36   10-5214   4   10-5516   8   10-1577   15   10-1740   41   10-1982   18   10-2604   22   10-2775   23   10-3146   31   10-3685   36   10-5214   4   10-5517   8   10-1577   15   10-1741   41   10-1983   18   10-2604   22   10-2780   23   10-3146   31   10-3685   36   10-5216   4   10-5519   8   10-1579   15   10-1743   40   10-1984   18   10-2606   22   10-2781   23   10-3148   31   10-3697   36   10-5216   4   10-5521   8   10-1580   15   10-1744   40   10-1985   18   10-2610   22   10-2782   23   10-3150   30   10-3720   36   10-5217   4   10-5521   8   10-1581   15   10-1745   42   10-1994   18   10-2610   22   10-2784   23   10-3150   30   10-3730   36   10-5221   5   10-5522   8   10-1581   15   10-1751   41   10-1991   18   10-2613   22   10-2796   23   10-3153   30   10-3737   36   10-5221   5   10-5524   8   10-1585   15   10-1752   41   10-1992   18   10-2613   22   10-2796   23   10-3155   30   10-5041   14   10-5224   5   10-5526   8   10-1584   15   10-1751   41   10-1993   18   10-2615   22   10-2796   23   10-3155   30   10-5041   14   10-5224   5   10-5526   8   10-1584   15   10-1751   41   10-1994   18   10-2615   22   10-2796   23   10-3155   30   10-5041   14   10-5224   5   10-5526   8   10-1584   15   10-1761   40   10-1994   18   10-2615   22   10-2797   23   10-3162   30   10-5041   14   10-5224   5   10-5526   8   10-1584   15   10-1761   40   10-1994   18   10-2615   22   10-2796   23   10-3162   30   10-5041   14   10-5224   5   10-5526   8   10-1584   15   10-1761   40   10-1994   18   10-2615   22   10-2816   23   10-3162   30   10-5041   14   10-5224   5   10-5526   8   10-1584   15   10-1766   40   10-1994   18   10-2					l .				-									
10-1572   15   10-1735   41   10-1980   18   10-2602   22   10-2774   23   10-3145   31   10-3680   36   10-5214   4   10-5516   8					l .			22										
10-1573   15   10-1740   41   10-1981   18   10-2603   22   10-2775   23   10-3146   31   10-3685   36   10-5215   4   10-5517   8   10-1577   15   10-1741   41   10-1982   18   10-2605   22   10-2780   23   10-3147   31   10-3680   36   10-5216   4   10-5520   8   10-1579   15   10-1742   41   10-1984   18   10-2605   22   10-2781   23   10-3150   30   10-3720   36   10-5218   4   10-5520   8   10-1579   15   10-1743   40   10-1985   18   10-2606   22   10-2783   23   10-3150   30   10-3720   36   10-5218   4   10-5521   8   10-1580   15   10-1744   40   10-1985   18   10-2610   22   10-2783   23   10-3151   30   10-3720   36   10-5218   4   10-5521   8   10-1581   15   10-1745   42   10-1990   18   10-2612   22   10-2785   23   10-3153   30   10-3730   36   10-5220   5   10-5522   8   10-1583   15   10-1751   41   10-1990   18   10-2613   22   10-2785   23   10-3153   30   10-3730   36   10-5222   5   10-5523   8   10-1583   15   10-1751   41   10-1991   18   10-2613   22   10-2795   23   10-3153   30   10-3730   36   10-5222   5   10-5524   8   10-1586   15   10-1752   41   10-1992   18   10-2613   22   10-2795   23   10-3153   30   10-3730   36   10-5222   5   10-5526   8   10-1586   15   10-1750   41   10-1993   18   10-2613   22   10-2795   23   10-3155   30   10-5040   14   10-5223   5   10-5526   8   10-1586   15   10-1760   40   10-1994   18   10-2616   22   10-2797   23   10-3165   30   10-5041   14   10-5224   5   10-5526   8   10-1586   15   10-1760   40   10-1995   18   10-2616   22   10-2810   23   10-3160   30   10-5041   14   10-5225   5   10-5526   8   10-1580   15   10-1763   41   10-1995   18   10-2620   22   10-2813   23   10-3160   30   10-5041   14   10-5226   5   10-5529   8   10-1580   15   10-1763   41   10-1995   18   10-2620   22   10-2813   23   10-3160   30   10-5041   14   10-5227   5   10-5540   9   10-1590   15   10-1764   40   10-1996   18   10-2620   22   10-2813   23   10-3160   30   10-5061   14   10-5227   5   10-5544   9   10-1590   15   10-1766   40   10-2271   20   10-2																		
10-1577   15					l .													
10-1579   15   10-1743   40   10-1984   18   10-2606   22   10-2782   23   10-3150   30   10-3720   36   10-5218   4   10-5521   8   10-1580   15   10-1744   40   10-1986   18   10-2610   22   10-2783   23   10-3152   30   10-3725   36   10-5220   5   10-5522   8   10-1582   15   10-1750   41   10-1990   18   10-2612   22   10-2785   23   10-3153   30   10-3737   36   10-5221   5   10-5523   8   10-1583   15   10-1751   41   10-1990   18   10-2612   22   10-2785   23   10-3153   30   10-3737   36   10-5222   5   10-5524   8   10-1584   15   10-1751   41   10-1991   18   10-2613   22   10-2795   23   10-3154   30   10-5040   14   10-5223   5   10-5525   8   10-1584   15   10-1752   41   10-1992   18   10-2614   22   10-2796   23   10-3156   30   10-5041   14   10-5223   5   10-5526   8   10-1586   15   10-1760   40   10-1994   18   10-2616   22   10-2810   23   10-3160   30   10-5041   14   10-5226   5   10-5527   8   10-1587   15   10-1760   40   10-1995   18   10-2612   22   10-2810   23   10-3160   30   10-5044   14   10-5226   5   10-5529   8   10-1587   15   10-1761   40   10-1995   18   10-2621   22   10-2812   23   10-3160   30   10-5044   14   10-5227   5   10-5540   9   10-1591   15   10-1763   40   10-2270   20   10-2622   22   10-2813   23   10-3162   30   10-5044   14   10-5227   5   10-5542   9   10-1592   15   10-1765   41   10-2272   20   10-2623   22   10-2813   23   10-3165   30   10-5064   14   10-5240   4   10-5243   4   10-5044   10-6060   28   10-1776   40   10-2273   20   10-2625   22   10-2813   23   10-3165   30   10-5062   14   10-5244   4   10-5244   4   10-5544   9   10-1602   28   10-1776   40   10-2274   20   10-2625   22   10-2813   23   10-3165   30   10-5068   14   10-5244   4   10-5544   9   10-1602   28   10-1776   40   10-2274   20   10-2625   22   10-2813   23   10-3166   30   10-5068   14   10-5244   4   10-5544   9   10-1602   28   10-1776   40   10-2274   20   10-2635   22   10-2812   23   10-3167   30   10-5068   14   10-5244   4   10-5555   9   10-1606   28   10-1786																		_
10-1580   15   10-1744   40   10-1985   18   10-2610   22   10-2783   23   10-3151   30   10-3725   36   10-5220   5   10-5522   8   10-1581   15   10-1745   42   10-1986   18   10-2611   22   10-2785   23   10-3152   30   10-3730   36   10-5221   5   10-5523   8   10-1583   15   10-1751   41   10-1991   18   10-2612   22   10-2795   23   10-3153   30   10-3730   36   10-5221   5   10-5522   8   10-1583   15   10-1751   41   10-1991   18   10-2612   22   10-2795   23   10-3153   30   10-3737   36   10-5222   5   10-5524   8   10-1584   15   10-1752   41   10-1992   18   10-2613   22   10-2795   23   10-3154   30   10-5040   14   10-5223   5   10-5525   8   10-1585   15   10-1753   41   10-1992   18   10-2615   22   10-2796   23   10-3156   30   10-5040   14   10-5224   5   10-5526   8   10-1586   15   10-1760   40   10-1994   18   10-2615   22   10-2797   23   10-3156   30   10-5042   14   10-5225   5   10-5527   8   10-1587   15   10-1761   40   10-1995   18   10-2620   22   10-2811   23   10-3160   30   10-5043   14   10-5225   5   10-5529   8   10-1588   15   10-1762   40   10-1996   18   10-2620   22   10-2812   23   10-3161   30   10-5044   14   10-5227   5   10-5540   9   10-1590   15   10-1763   40   10-2270   20   10-2622   22   10-2812   23   10-3163   30   10-5061   14   10-5228   5   10-5541   9   10-1592   15   10-1765   41   10-2272   20   10-2623   22   10-2814   23   10-3163   30   10-5061   14   10-5240   4   10-5543   9   10-1602   28   10-1766   40   10-2273   20   10-2625   22   10-2818   23   10-3168   30   10-5062   14   10-5244   4   10-5544   9   10-1602   28   10-1766   40   10-2275   20   10-2633   22   10-2812   23   10-3168   30   10-5064   14   10-5244   4   10-5544   9   10-1602   28   10-1766   40   10-2275   20   10-2633   22   10-2820   23   10-3168   30   10-5064   14   10-5246   4   10-5545   9   10-1603   28   10-1766   40   10-2281   20   10-2633   22   10-2820   23   10-3168   30   10-5008   14   10-5246   4   10-5555   9   10-1607   28   10-1766   40   10-2281   20   10-2															10-5217			
10-1581   15   10-1745   42   10-1986   18   10-2611   22   10-2784   23   10-3152   30   10-3730   36   10-5221   5   10-5523   8   10-1752   141   10-1990   18   10-2613   22   10-2785   23   10-3153   30   10-3737   36   10-5222   5   10-5524   8   10-1583   15   10-1751   41   10-1991   18   10-2613   22   10-2795   23   10-3153   30   10-5040   14   10-5223   5   10-5524   8   10-1751   41   10-1992   18   10-2614   22   10-2796   23   10-3155   30   10-5040   14   10-5223   5   10-5525   8   10-1585   15   10-1753   41   10-1993   18   10-2615   22   10-2796   23   10-3155   30   10-5041   14   10-5224   5   10-5526   8   10-1586   15   10-1760   40   10-1994   18   10-2616   22   10-2810   23   10-3160   30   10-5042   14   10-5225   5   10-5527   8   10-1586   15   10-1761   40   10-1995   18   10-2620   22   10-2811   23   10-3160   30   10-5043   14   10-5227   5   10-5540   9   10-1590   15   10-1763   40   10-1996   18   10-2621   22   10-2812   23   10-3163   30   10-5044   14   10-5227   5   10-5541   9   10-1591   15   10-1763   40   10-2270   20   10-2622   22   10-2813   23   10-3163   30   10-5060   14   10-5229   15   10-5542   9   10-1592   15   10-1765   41   10-2273   20   10-2623   22   10-2818   23   10-3166   30   10-5061   14   10-5240   4   10-5543   9   10-1601   28   10-1766   40   10-2273   20   10-2626   22   10-2818   23   10-3166   30   10-5061   14   10-5240   4   10-5544   9   10-1602   28   10-1766   40   10-2275   20   10-2626   22   10-2818   23   10-3166   30   10-5061   14   10-5243   4   10-5544   9   10-1602   28   10-1766   40   10-2275   20   10-2626   22   10-2818   23   10-3166   30   10-5068   14   10-5244   4   10-5544   9   10-1602   28   10-1766   40   10-2275   20   10-2626   22   10-2818   23   10-3166   30   10-5068   14   10-5244   4   10-5544   9   10-1606   28   10-1766   40   10-2275   20   10-2626   22   10-2812   23   10-3166   30   10-5068   14   10-5244   4   10-5545   9   10-1605   28   10-1766   40   10-2280   20   10-2633   22   10-2826   23   10																		
10-1583   15   10-1751   41   10-1991   18   10-2613   22   10-2795   23   10-3154   30   10-5040   14   10-5223   5   10-5525   8   10-1584   15   10-1752   41   10-1992   18   10-2614   22   10-2796   23   10-3155   30   10-5041   14   10-5224   5   10-5526   8   10-1586   15   10-1763   40   10-1994   18   10-2615   22   10-2810   23   10-3160   30   10-5042   14   10-5225   5   10-5527   8   10-1587   15   10-1761   40   10-1995   18   10-2620   22   10-2811   23   10-3160   30   10-5043   14   10-5225   5   10-5529   8   10-1588   15   10-1761   40   10-1995   18   10-2620   22   10-2811   23   10-3160   30   10-5044   14   10-5227   5   10-5540   9   10-1588   15   10-1763   40   10-1996   18   10-2620   22   10-2812   23   10-3162   30   10-5048   14   10-5228   5   10-5540   9   10-1590   15   10-1763   40   10-2270   20   10-2622   22   10-2813   23   10-3162   30   10-5060   14   10-5229   15   10-1564   9   10-1591   15   10-1764   40   10-2271   20   10-2623   22   10-2814   23   10-3165   30   10-5060   14   10-5240   4   10-5240   4   10-5240   4   10-10-1000   10-10000   10-10000   10-10000   10-10000   10-100000   10-10000000000																		
10-1584         15         10-1752         41         10-1992         18         10-2614         22         10-2796         23         10-3155         30         10-5041         14         10-5224         5         10-5526         8           10-1585         15         10-1760         40         10-1994         18         10-2616         22         10-2797         23         10-3156         30         10-5042         14         10-5225         5         10-5527         8           10-1586         15         10-1760         40         10-1994         18         10-2616         22         10-2810         23         10-3160         30         10-5043         14         10-5226         5         10-5529         8           10-1587         15         10-1761         40         10-1995         18         10-2621         22         10-2812         23         10-3160         30         10-5044         14         10-5227         5         10-5540         9           10-1590         15         10-1763         40         10-2270         20         10-2622         22         10-2813         23         10-3163         30         10-5060         14         10-5229																		
10-1585         15         10-1753         41         10-1993         18         10-2615         22         10-2797         23         10-3156         30         10-5042         14         10-5225         5         10-5527         8           10-1586         15         10-1760         40         10-1994         18         10-2616         22         10-2810         23         10-3160         30         10-5043         14         10-5226         5         10-5529         8           10-1587         15         10-1761         40         10-1996         18         10-2621         22         10-2812         23         10-3161         30         10-5044         14         10-5227         5         10-5540         9           10-1588         15         10-1762         40         10-1996         18         10-2621         22         10-2812         23         10-3163         30         10-5048         14         10-5227         5         10-5540         9           10-1591         15         10-1764         40         10-2271         20         10-2623         22         10-2814         23         10-3163         30         10-5061         14         10-5240																		
10-1586         15         10-1760         40         10-1994         18         10-2616         22         10-2810         23         10-3160         30         10-5043         14         10-5226         5         10-5529         8           10-1587         15         10-1761         40         10-1995         18         10-2620         22         10-2811         23         10-3161         30         10-5044         14         10-5227         5         10-5540         9           10-1588         15         10-1762         40         10-1996         18         10-2621         22         10-2812         23         10-3163         30         10-5060         14         10-5228         5         10-5541         9           10-1591         15         10-1764         40         10-2271         20         10-2623         22         10-2814         23         10-3163         30         10-5060         14         10-5229         15         10-1765         41         10-2272         20         10-2623         22         10-2814         23         10-3163         30         10-5061         14         10-5240         4         10-5543         9           10-1600																		
10-1588         15         10-1762         40         10-1996         18         10-2621         22         10-2812         23         10-3162         30         10-5048         14         10-5228         5         10-5541         9           10-1590         15         10-1763         40         10-2271         20         10-2622         22         10-2813         23         10-3163         30         10-5060         14         10-5229         15         10-5542         9           10-1591         15         10-1764         40         10-2271         20         10-2623         22         10-2814         23         10-3163         30         10-5060         14         10-5240         4         10-5543         9           10-1592         15         10-1765         41         10-2272         20         10-2624         22         10-2818         23         10-3165         30         10-5062         14         10-5241         4         10-5544         9           10-1600         28         10-1766         40         10-2273         20         10-2626         22         10-2819         23         10-3166         30         10-5063         14         10-5244															10-5226			
10-1590         15         10-1763         40         10-2270         20         10-2622         22         10-2813         23         10-3163         30         10-5060         14         10-5229         15         10-5542         9           10-1591         15         10-1764         40         10-2271         20         10-2623         22         10-2814         23         10-3164         30         10-5061         14         10-5240         4         10-5543         9           10-1600         28         10-1766         40         10-2273         20         10-2625         22         10-2819         23         10-3166         30         10-5062         14         10-5241         4         10-5544         9           10-1600         28         10-1766         40         10-2273         20         10-2625         22         10-2880         23         10-3167         30         10-5062         14         10-5242         4         10-5545         9           10-1601         28         10-1776         40         10-2275         20         10-2630         22         10-2820         23         10-3168         30         10-5068         14         10-5243																		
10-1592         15         10-1765         41         10-2272         20         10-2624         22         10-2818         23         10-3165         30         10-5062         14         10-5241         4         10-5544         9           10-1600         28         10-1766         40         10-2273         20         10-2625         22         10-2819         23         10-3166         30         10-5063         14         10-5242         4         10-5545         9           10-1601         28         10-1776         40         10-2275         20         10-2630         22         10-2820         23         10-3168         30         10-5068         14         10-5243         4         10-5546         9           10-1602         28         10-1776         40         10-2275         20         10-2630         22         10-2821         23         10-3168         30         10-5068         14         10-5243         4         10-5546         9           10-1603         28         10-1786         40         10-2280         20         10-2632         22         10-2822         23         10-3173         31         10-5071         14         10-5245																		
10-1600         28         10-1766         40         10-2273         20         10-2625         22         10-2819         23         10-3166         30         10-5063         14         10-5242         4         10-5545         9           10-1601         28         10-1770         41         10-2274         20         10-2626         22         10-2820         23         10-3167         30         10-5064         14         10-5243         4         10-5546         9           10-1602         28         10-1776         40         10-2275         20         10-2630         22         10-2821         23         10-3168         30         10-5068         14         10-5244         4         10-5549         9           10-1603         28         10-1786         40         10-2280         20         10-2821         23         10-3173         31         10-5070         14         10-5245         4         10-5550         9           10-1604         28         10-1796         40         10-2281         20         10-2632         22         10-2824         23         10-3171         31         10-5071         14         10-5246         4         10-5551																		
10-1601         28         10-1770         41         10-2274         20         10-2626         22         10-2820         23         10-3167         30         10-5064         14         10-5243         4         10-5546         9           10-1602         28         10-1776         40         10-2275         20         10-2630         22         10-2821         23         10-3168         30         10-5068         14         10-5244         4         10-5549         9           10-1603         28         10-1780         41         10-2276         20         10-2631         22         10-2822         23         10-3170         31         10-5070         14         10-5245         4         10-5550         9           10-1604         28         10-1796         40         10-2281         20         10-2632         22         10-2826         23         10-3171         31         10-5071         14         10-5246         4         10-5551         9           10-1605         28         10-1800         19         10-2282         20         10-2832         22         10-2826         23         10-3172         31         10-5088         14         10-5247																		
10-1602         28         10-1776         40         10-2275         20         10-2630         22         10-2821         23         10-3168         30         10-5068         14         10-5244         4         10-5549         9           10-1603         28         10-1780         41         10-2276         20         10-2631         22         10-2822         23         10-3170         31         10-5070         14         10-5245         4         10-5550         9           10-1604         28         10-1796         40         10-2281         20         10-2632         22         10-2824         23         10-3172         31         10-5071         14         10-5246         4         10-5551         9           10-1605         28         10-1796         40         10-2281         20         10-2826         23         10-3172         31         10-5088         14         10-5246         4         10-5551         9           10-1606         28         10-1800         19         10-2283         20         10-2634         22         10-2828         23         10-3173         31         10-5088         14         10-5246         4         10-5552													10-5064					
10-1604         28         10-1786         40         10-2280         20         10-2632         22         10-2824         23         10-3171         31         10-5071         14         10-5246         4         10-5551         9           10-1606         28         10-1800         19         10-2282         20         10-2633         22         10-2826         23         10-3172         31         10-5088         14         10-5246         4         10-5552         9           10-1606         28         10-1804         19         10-2283         20         10-2635         22         10-2830         23         10-3173         31         10-5089         14         10-5246         4         10-5553         9           10-1607         28         10-1804         19         10-2283         20         10-2635         22         10-2830         23         10-3173         31         10-5099         14         10-5251         7         10-5554         9           10-1609         28         10-1806         19         10-2284         20         10-2636         22         10-2832         23         10-3175         31         10-5091         14         10-5252	10-1602	28	10-1776	40	10-2275	20	10-2630	22	10-2821	23	10-3168	30	10-5068	14	10-5244	4	10-5549	9
10-1605         28         10-1796         40         10-2281         20         10-2633         22         10-2826         23         10-3172         31         10-5088         14         10-5247         4         10-5552         9           10-1606         28         10-1800         19         10-2282         20         10-2634         22         10-2828         23         10-3173         31         10-5089         14         10-5248         4         10-5553         9           10-1607         28         10-1804         19         10-2283         20         10-2635         22         10-2830         23         10-3174         31         10-5090         15         10-5251         7         10-5554         9           10-1608         28         10-1805         19         10-2284         20         10-2636         22         10-2832         23         10-3175         31         10-5091         14         10-5252         7         10-5555         9           10-1609         28         10-1806         19         10-2285         20         10-2640         22         10-2834         23         10-3176         31         10-5092         14         10-5253																		
10-1606     28     10-1800     19     10-2282     20     10-2634     22     10-2828     23     10-3173     31     10-5089     14     10-5248     4     10-5553     9       10-1607     28     10-1804     19     10-2283     20     10-2635     22     10-2830     23     10-3174     31     10-5090     15     10-5251     7     10-5554     9       10-1608     28     10-1805     19     10-2284     20     10-2636     22     10-2832     23     10-3175     31     10-5091     14     10-5252     7     10-5555     9       10-1609     28     10-1806     19     10-2285     20     10-2640     22     10-2834     23     10-3176     31     10-5092     14     10-5253     7     10-5556     9																		
10-1608         28         10-1805         19         10-2284         20         10-2636         22         10-2832         23         10-3175         31         10-5091         14         10-5252         7         10-5555         9           10-1609         28         10-1806         19         10-2285         20         10-2640         22         10-2834         23         10-3176         31         10-5092         14         10-5253         7         10-5556         9	10-1606	28	10-1800	19	10-2282	20	10-2634	22	10-2828	23	10-3173	31	10-5089	14	10-5248	4	10-5553	9
10-1609 28   10-1806 19   10-2285 20   10-2640 22   10-2834 23   10-3176 31   10-5092 14   10-5253 7   10-5556 9																		



					<b>•</b> •	sku-s	sor	ted	pac	ıe in	de	X					
item p	age	item	page	item p	oage		page	item			page		oage	item p	age	item	page
10-5560	9	10-5676	17	10-5845	7	11-0710	98	11-1149	97	11-1589		12-0240	63	12-0450	74	12-1107	86
10-5561 10-5562	9	10-5677 10-5678	17 17	10-5846 10-5847	7 7	11-0711 11-0712	98 98	11-1150 11-1151	97 97	11-1590 11-1591	106 106	12-0241 12-0243	63 63	12-0451 12-0470	74 71	12-1108 12-1109	86 86
10-5563	9	10-5679	17	10-5901	7	11-0712	98	11-1152	97	11-1592		12-0243	63	12-0470	79	12-1103	86
10-5564	9	10-5691	6	10-5902	7	11-0730	98	11-1153	97	11-1593		12-0247	63	12-0580	68	12-1112	86
10-5565	9	10-5692	6	10-5903	7 7	11-0740	98	11-1161	97 97	11-1594		12-0248	69 67	12-0583	68	12-1116	86
10-5566 10-5569	9	10-5693 10-5694	6	10-5904 10-5905	7	11-0780 11-0781	98 98	11-1170	97	11-1600 11-1601	96 96	12-0250 12-0251	67 67	12-0590 12-0591	68 68	12-1121 12-1122	86 86
10-5570	9	10-5695	6	10-5908	7	11-0782	98	11-1205	102	11-1603	96	12-0259	75	12-0592	68	12-1124	86
10-5571	9	10-5698	6	10-5911	7	11-0800	100	11-1226	102	11-1604	96	12-0261	67	12-0593	68	12-1125	86
10-5572 10-5573	9	10-5710 10-5711	8	10-5912 10-5913	7 7	11-0810 11-0811	100	11-1227 11-1240	102 107	11-1605 11-1606	96 96	12-0262 12-0263	67 67	12-0594 12-0595	68 68	12-1129 12-1140	86 86
10-5574	9	10-5712	8	10-5914	7	11-0830	100	11-1241	107	11-1607	96	12-0264	67	12-0600	62	12-1141	86
10-5575	9	10-5713	8	10-5915	7	11-0831	100	11-1242	107	11-1710	96	12-0266	67	12-0601	64	12-1149	80
10-5576 10-5579	9	10-5714 10-5715	8 8	10-5918 10-5921	7 7	11-0950 11-0951	100 100	11-1243 11-1250	107 104	11-1711 11-1712	96 96	12-0267 12-0281	69 62	12-0602 12-0603	66 66	12-1150 12-1151	80 80
10-5580	9	10-5716	8	10-5922	7	11-0952	100	11-1251	104	11-1713	96	12-0282	62	12-0604	62	12-1155	80
10-5581	9	10-5717	8	10-5923	7	11-0953	100	11-1252	104	11-1714	96	12-0286	62	12-0605	66	12-1156	80
10-5582 10-5583	9	10-5719 10-5720	8 8	10-5924 10-5925	7 7	11-0954 11-1000	100 104	11-1253 11-1257	104 104	11-1716 11-1718	96 96	12-0288 12-0289	75 74	12-0606 12-0607	66 64	12-1157 12-1180	80 80
10-5584	9	10-5720	8	10-5928	7	11-1000	104	11-1280	104	11-1710	96	12-0203	62	12-10007	76	12-1100	89
10-5585	9	10-5722	8	10-5941	7	11-1002	104	11-1281	106	11-1722	96	12-0291	62	12-1001	76	12-1191	89
10-5586 10-5589	9	10-5723 10-5724	8	10-5942 10-5943	7 7	11-1003 11-1004	104 104	11-1282 11-1283	106 106	11-1730 11-1744	96 96	12-0292 12-0293	62 62	12-1002 12-1005	76 76	12-1192 12-1193	88 88
10-5590	9	10-5724	8	10-5943	7	11-1004	104	11-1203	106	11-1744	96	12-0293	62	12-1005	76	12-1193	87
10-5591	9	10-5726	8	10-5945	7	11-1006	104	11-1285	106	11-1748	96	12-0298	74	12-1010	78	12-1205	87
10-5592	9	10-5727 10-5729	8	10-5948	7	11-1007 11-1009	104	11-1286	106	11-1750	96	12-0299	74	12-1011	78	12-1210	87
10-5593 10-5594	9	10-5729	10	10-5949 10-6000	15 37	11-1009	104 104	11-1287 11-1291	106 106	11-1752 11-1760	96 96	12-0300 12-0301	70 70	12-1012 12-1014	76 78	12-1220 12-1221	87 87
10-5595	9	10-5761	10	10-6001	37	11-1011	104	11-1292	106	11-1900		12-0302	70	12-1015	78	12-1230	87
10-5596	9	10-5762	10	10-6002	37	11-1012	104	11-1293	106	11-1901	109	12-0303	70	12-1016	78	12-1231	87
10-5599 10-5600	9	10-5763 10-5764	10 10	10-6003 10-6004	37 37	11-1013 11-1020	104 107	11-1294 11-1297	106 106	11-1911 11-1920	109 109	12-0304 12-0340	70 70	12-1020 12-1025	79 77	12-1320 12-1340	88 88
10-5601	6	10-5765	10	10-6020	37	11-1021	107	11-1298	106	12-0021	69	12-0341	70	12-1026	77	12-1341	88
10-5602	6	10-5766	10	10-6021	37	11-1029	99	11-1300	102	12-0023		12-0342	70	12-1027	77	12-1342	88
10-5603 10-5604	6	10-5767 10-5769	10 10	10-6022 10-6023	37 37	11-1030 11-1031	99	11-1301 11-1308	102 102	12-0025 12-0027	69 69	12-0343 12-0370	70 71	12-1028 12-1030	76 79	12-1343 12-1344	88 88
10-5605	6	10-5800	10	10-6023	37	11-1031	99	11-1310	102	12-0027		12-0370	71	12-1030	79	12-1345	88
10-5606	6	10-5801	10	10-6040	37	11-1033	99	11-1311	102	12-0052		12-0372	71	12-1032	79	12-1346	88
10-5607 10-5608	6	10-5802 10-5803	10 10	10-6050 11-0103	37 99	11-1034 11-1035	99	11-1312 11-1313	102 102	12-0053 12-0054		12-0376 12-0377	71 71	12-1035 12-1040	79 79	12-1347 12-1348	88 88
10-5610	6	10-5804	10	11-0105	99	11-1035	99	11-1313	102	12-0034		12-0377	71	12-1040	79	12-1349	88
10-5611	6	10-5805	10	11-0106	99	11-1037	99	11-1315	102	12-0101	66	12-0380	75	12-1042	79	12-1350	88
10-5612 10-5613	6 6	10-5806 10-5807	10 10	11-0107 11-0108	99 99	11-1038 11-1039	99 99	11-1316 11-1321	102 102	12-0102 12-0103		12-0381 12-0382	74 74	12-1043 12-1050	79 79	12-1351 12-1352	88 88
10-5614	6	10-5809	10	11-0100	99	11-1060	107	11-1360	102	12-0103	66	12-0302	75	12-1055	80	12-1352	88
10-5615	6	10-5810	10	11-0200	99	11-1061	107	11-1361	102	12-0111	66	12-0384	69	12-1056	80	12-1355	88
10-5616	6	10-5811	10	11-0201	99 99	11-1062	107	11-1362	102	12-0112		12-0385	71	12-1057	80 74	12-1390	82
10-5617 10-5618	6	10-5812 10-5813	10 10	11-0202 11-0220	99	11-1066 11-1067	107 107	11-1363 11-1364	102 102	12-0113 12-0114	66 67	12-0387 12-0388	70 70	12-1059 12-1060	78	12-1391 12-1392	82 82
10-5620	6	10-5814	10	11-0221	99	11-1068	107	11-1365	102	12-0115	67	12-0389	67	12-1061	75	12-1393	82
10-5621	6	10-5815	10	11-0222 11-0223	99	11-1069	107	11-1366	102	12-0116		12-0392	70	12-1062	75	12-1395 12-1396	82
10-5622 10-5623	6	10-5816 10-5817	10 10	11-0223	99	11-1088 11-1090	107 107	11-1367 11-1368	102 102	12-0117 12-0126		12-0393 12-0394	70 69	12-1063 12-1064	75 75	12-1396	82 82
10-5624	6	10-5819	10	11-0225	99	11-1092	107	11-1369	102	12-0127	66	12-0397	70	12-1065	75	12-1398	82
10-5625	6	10-5820	9	11-0226	99	11-1093	107	11-1376	102	12-0128		12-0398	70	12-1066	80	12-1406	82
10-5626 10-5627	6 6	10-5821 10-5822	9 9	11-0227 11-0240	99 100	11-1094 11-1095	107 107	11-1377 11-1390	102 103	12-0180 12-0181		12-0399 12-0400	69 71	12-1070 12-1071	80 80	12-1407 12-1408	82 82
10-5628	6	10-5823	9	11-0244	100	11-1096	107	11-1391	103	12-0182		12-0401	71	12-1076	81	12-1440	83
10-5629	15	10-5824	9	11-0245	100	11-1097	107	11-1397	103	12-0185		12-0402	71	12-1077	81	12-1441	83
10-5640 10-5641	6 6	10-5825 10-5826	9 9	11-0246 11-0400	100 100	11-1110 11-1119	97 97	11-1398 11-1399	103 103	12-0200 12-0201	64 64	12-0403 12-0406	71 71	12-1078 12-1080	81 81	12-1442 12-1443	83 83
10-5642	6	10-5829	9	11-0401	100	11-1120	101	11-1550	105	12-0202		12-0407	71	12-1082	81	12-1444	83
10-5643	6	10-5830	9	11-0402	100	11-1121	101	11-1552	105	12-0203	64	12-0410	71	12-1085	89	12-1450	83
10-5644 10-5645	6	10-5831 10-5832	9	11-0403 11-0404	100	11-1123 11-1130	101	11-1556 11-1557	105 105	12-0209 12-0211	64 75	12-0411 12-0412	71 71	12-1086 12-1090	89 81	12-1451 12-1455	83 83
10-5646	6	10-5833	9	11-0404	100	11-1130	101	11-1560	105	12-0211	63	12-0412	10	12-1090	81	12-1455	84
10-5647	6	10-5834	9	11-0406	100	11-1132	101	11-1561	105	12-0222	65	12-0422	10	12-1094	88	12-1461	84
10-5648	6 17	10-5835 10-5836	9	11-0500 11-0501	104 104	11-1133 11-1135	101	11-1562 11-1563	105	12-0226 12-0228		12-0424 12-0427	10 10	12-1095 12-1096	79 81	12-1462 12-1464	84 84
10-5670 10-5671	17	10-5839	9	11-0501	104	11-1136	101	11-1563	105 105	12-0220		12-0427	71	12-1096	80	12-1464	84
10-5672	17	10-5841	7	11-0506	104	11-1137	101	11-1565	105	12-0235	65	12-0429	10	12-1099	81	12-1466	84
10-5673 10-5674	17 17	10-5842 10-5843	7 7	11-0700 11-0701	98 98	11-1138 11-1140	101 101	11-1566 11-1567	105 105	12-0237 12-0238		12-0443 12-0445	71 71	12-1100 12-1101	86 86	12-1467 12-1468	84 84
10-5674	17		7	11-0701	98	11-1140	97			12-0236		12-0445		12-1101	86		84



			sku-sor	ted pag	ge inde	X		
item page	item page	item page	item page	item page		item page	item page	item page
12-1470 84	12-1920 85	12-4502 92	12-4600 95	13-3016 117	13-3320 116	15-1022 144	15-3262 151	15-4210 152
12-1471 84 12-1472 84	12-1921 85 12-1922 85	12-4503 92 12-4504 92	12-4601 95 12-4602 95	13-3019 116 13-3020 117	13-3350 118 13-3360 118	15-1023 144 15-1030 144	15-3263 151 15-3264 151	15-4220 152 15-4221 152
12-1472 84	12-1922 05	12-4505 92	12-4603 95	13-3020 117	13-5220 110	15-1030 144	15-3271 151	15-4266 131
12-1476 84	12-1924 85	12-4506 92	12-4604 95	13-3032 117	13-5221 110	15-1082 144	15-3272 151	15-4267 131
12-1477 84 12-1480 83	12-1926 85 12-1940 90	12-4507 92 12-4508 92	12-4605 95 12-4606 95	13-3061 108 13-3062 108	13-5223 110 13-5230 110	15-1083 144 15-1200 144	15-3273 151 15-3274 151	15-4268 131 16-1100 148
12-1481 83	12-1941 90	12-4509 94	12-4607 95	13-3063 108	13-5250 110	15-1201 144	15-3281 151	16-1102 148
12-1482 83	12-1942 90	12-4510 94	12-4608 95	13-3064 117	13-5251 110	15-1202 144	15-3282 151	16-1600 148
12-1485 75 12-1490 83	12-1950 90 12-1951 90	12-4511 94 12-4512 94	12-4609 95 12-4610 95	13-3065 117 13-3070 108	13-5252 110 13-5253 110	15-1203 144 15-1210 144	15-3283 151 15-3284 151	16-1601 148 16-1602 148
12-1492 83	12-1952 90	12-4513 94	12-4611 95	13-3071 108	14-1000 138	15-1211 144	15-3290 151	16-1603 148
12-1495 83 12-1500 84	12-1953 90 12-1960 85	12-4514 94 12-4515 94	12-4612 95 12-4613 95	13-3072 108 13-3073 108	14-1001 138 14-1002 138	15-1220 144 15-1221 144	15-3291 151 15-3292 151	16-1701 148 16-1710 148
12-1500 84	12-1900 03	12-4516 94	12-4614 95	13-3074 108	14-1002 138	15-1230 144	15-3293 151	16-1800 148
12-1502 84	12-2031 90	12-4517 94	12-4615 95	13-3075 108	14-1004 138	15-1231 144	15-3294 151	17-1101 148
12-1510 84 12-1511 84	12-2032 90 12-2040 90	12-4518 94 12-4519 94	12-4616 95 12-4617 95	13-3076 108 13-3077 108	14-1005 138 14-1006 138	15-1500 142 15-1502 142	15-3300 151 15-3301 151	17-1121 148 18-1135 101
12-1512 84	12-2041 90	12-4520 94	12-4618 95	13-3086 117	14-1007 138	15-1506 142	15-3302 151	18-1138 101
12-1514 84	12-2042 90	12-4521 94	12-4619 95	13-3087 117	14-1008 138	15-1515 142	15-3304 151	18-1165 101
12-1516 84 12-1517 84	12-2043 90 12-2044 90	12-4522 94 12-4523 94	12-4620 95 12-4621 95	13-3088 108 13-3089 108	14-1009 138 14-1020 138	15-1536 142 15-1537 142	15-3310 151 15-3311 151	18-1168 101 18-1181 101
12-1520 84	12-2045 90	12-4524 94	12-4622 95	13-3091 117	14-1030 138	15-1539 142	15-3313 151	18-1220 101
12-1521 84 12-1524 84	12-2046 90 12-2050 90	12-4525 94 12-4526 94	12-4623 95 12-4624 95	13-3099 117 13-3100 115	14-1031 138 14-1032 138	15-1543 142 15-2000 145	15-3314 151 15-3316 151	18-1225 101 18-1232 101
12-1524 84	12-2052 90	12-4527 94	13-1150 128	13-3101 115	14-1040 138	15-2001 145	15-3730 141	18-1233 101
12-1560 82	12-2100 90	12-4528 94	13-1151 128	13-3103 115	14-1050 138	15-2002 145	15-3731 141	18-1234 101
12-1562 82 12-1564 82	12-2151 85 12-2152 90	12-4529 93 12-4530 93	13-1152 128 13-1153 128	13-3105 115 13-3106 115	14-1051 138 14-1080 126	15-2003 145 15-2004 145	15-3732 141 15-3733 141	18-1235 101 18-1280 101
12-1565 82	12-2153 90	12-4531 93	13-1165 128	13-3107 115	14-1081 126	15-2016 145	15-3734 141	18-1320 101
12-1566 82 12-1570 84	12-2154 90	12-4532 93	13-1170 128	13-3120 115	14-1082 126	15-2017 145	15-3735 141 15-3736 141	18-1331 101
12-1570 84 12-1571 84	12-2155 90 12-2156 90	12-4536 93 12-4537 93	13-1171 128 13-1185 128	13-3121 115 13-3123 115	14-1083 126 14-1084 126	15-2018 145 15-2019 145	15-3736 141 15-3737 141	18-1370 101 18-1371 101
12-1572 84	12-2157 90	12-4539 93	13-1186 128	13-3124 115	14-1085 126	15-2026 145	15-4000 152	18-1372 101
12-1573 84 12-1574 84	12-2158 90 12-2159 90	12-4540 93 12-4541 93	13-1187 128 13-1188 128	13-3125 115 13-3126 115	14-1300 139 14-1301 139	15-2027 145 15-2028 145	15-4001 152 15-4004 152	18-1373 101 18-1380 101
12-1637 82	12-2200 85	12-4542 93	13-1287 111	13-3133 115	14-1310 139	15-2029 145	15-4005 152	18-1390 101
12-1638 82	12-2210 85	12-4543 93	13-1288 111	13-3144 115	14-1311 139	15-2030 145	15-4010 152	18-1391 101
12-1639 82 12-1640 82	12-2211 85 12-2220 85	12-4544 93 12-4545 93	13-1291 111 13-1293 111	13-3145 115 13-3146 115	14-1320 139 14-1330 139	15-2031 145 15-2032 145	15-4030 153 15-4031 153	18-1393 101 18-1395 101
12-1641 82	12-2240 85	12-4546 93	13-1294 111	13-3153 115	14-1335 139	15-2033 145	15-4037 153	19-1002 149
12-1642 82 12-1643 82	12-2241 85 12-2242 85	12-4547 92 12-4548 92	13-1295 111 13-1302 118	13-3154 115 13-3155 115	14-1336 139 14-1337 139	15-2034 145 15-2046 145	15-4038 153 15-4039 153	19-1003 149 19-1004 149
12-1644 82	12-2243 85	12-4549 92	13-1302 118	13-3156 115	14-1340 139	15-2047 145	15-4040 153	19-1004 149
12-1645 82	12-2244 85	12-4550 92	13-1304 118	13-3157 115	14-1341 139	15-2048 145	15-4041 153	19-1006 149
12-1646 82 12-1647 82	12-2250 85 12-2261 85	12-4551 92 12-4552 92	13-1305 118 13-1306 119	13-3158 115 13-3159 115	14-1350 139 14-1360 139	15-2049 145 15-2056 145	15-4042 153 15-4050 153	19-1007 149 19-1008 149
12-1648 82	12-2262 85	12-4553 92	13-1308 120	13-3160 115	14-1370 139	15-2057 145	15-4051 153	19-1009 149
12-1649 82 12-1650 82	12-2272 85 12-2273 85	12-4554 92 12-4555 94	13-1309 119 13-1313 120	13-3161 115 13-3171 118	14-1380 138 14-1381 138	15-2058 145 15-2059 145	15-4057 153 15-4058 153	19-1014 149 19-1015 149
12-1651 82	12-2273 85	12-4556 94	13-1314 119	13-3180 118	14-1390 138	15-2060 145	15-4060 153	19-1016 149
12-1652 82	12-3000 91	12-4557 94	13-1315 117	13-3181 118	14-1391 138	15-2062 145	15-4061 153	19-1017 149
12-1653 82 12-1654 82	12-3001 91 12-3002 91	12-4558 94 12-4560 95	13-1316 117 13-1317 119	13-3200 121 13-3201 121	14-1392 138 14-1393 138	15-2063 145 15-2072 145	15-4062 153 15-4063 153	19-1018 149 19-1019 149
12-1655 82	12-3005 91	12-4561 94	13-1318 119	13-3202 121	14-1400 138	15-2073 145	15-4064 153	19-1032 149
12-1656 82 12-1657 82	12-3010 91 12-3011 91	12-4562 95 12-4563 95	13-1322 120 13-1323 120	13-3210 121 13-3211 121	14-1500 126 14-1501 126	15-2074 145 15-2082 145	15-4065 153 15-4066 153	19-1033 149 19-1034 149
12-1657 82 12-1658 82	12-3020 91	12-4564 95	13-1348 120	13-3211 121	14-1501 126	15-2082 145	15-4067 153	19-1034 149
12-1659 82	12-3030 91	12-4565 95	13-1349 120	13-3252 108	14-1503 126	15-2084 145	15-4080 153	19-1036 149
12-1660 82 12-1662 82	12-3031 91 12-3050 91	12-4566 95 12-4567 95	13-1360 120 13-1361 120	13-3253 108 13-3254 108	14-1504 126 14-1505 126	15-2100 145 15-2101 145	15-4081 153 15-4082 153	19-1037 149 19-1038 149
12-1664 82	12-3051 91	12-4568 95	13-1362 120	13-3255 108	14-1506 126	15-2102 145	15-4090 153	19-1040 149
12-1666 82	12-3052 91	12-4569 95	13-1365 120	13-3256 108	14-1511 126	15-2103 145	15-4091 153	19-1041 149
12-1668 82 12-1670 82	12-3060 91 12-3080 91	12-4571 95 12-4572 95	13-1366 119 13-1367 119	13-3257 108 13-3258 108	14-1512 126 14-1513 126	15-2104 145 15-2123 144	15-4092 153 15-4140 153	19-1042 149 19-1043 149
12-1671 82	12-3081 91	12-4574 95	13-1368 119	13-3259 108	15-1000 144	15-3040 150	15-4142 153	19-1044 149
12-1710 89 12-1711 89	12-3082 91 12-3500 87	12-4575 95 12-4580 92	13-1369 119 13-1442 131	13-3260 108 13-3261 108	15-1001 144 15-1002 144	15-3041 150 15-3042 150	15-4150 153 15-4151 153	19-1045 149 19-1051 149
12-1711 89	12-3500 87	12-4581 92	13-3000 116	13-3300 116	15-1002 144	15-3042 150	15-4151 153	19-1051 149
12-1713 89	12-3502 87	12-4582 92	13-3001 116	13-3301 116	15-1010 144	15-3250 151	15-4200 152	19-1060 149
12-1714 89 12-1902 89	12-3504 87 12-3505 87	12-4583 92 12-4584 92	13-3009 116 13-3010 116	13-3310 116 13-3311 116	15-1011 144 15-1012 144	15-3251 151 15-3255 151	15-4201 152 15-4202 152	19-1061 149 19-1072 149
12-1910 85	12-3508 87	12-4585 92	13-3012 116	13-3312 116	15-1013 144	15-3256 151	15-4203 152	19-1100 149
12-1911 85 12-1912 85	12-4500 92 12-4501 92	12-4586 92 12-4587 92	13-3014 116 13-3015 117	13-3313 116 13-3314 116	15-1020 144 15-1021 144	15-3260 151 15-3261 151	15-4204 152 15-4205 152	19-1101 149 19-1102 149
12 1012 00	1 12 7001 02	1.2 4001 32	, 10 00 10 117	, 10 00 14 110	, 10 1021 174	, .0 0201 101	1 10 7200 102	1 10 1102 173



			sku-sor	ted pag	ge index	X		
item page	item page	item page	item page	item page	item page	item page	item page	item page
19-1103 149	20-4212 135	24-4573 164	24-7110 186	30-1054 49	30-1651 167	30-1833 49	30-2082 39	30-3000 170
19-1110 149 19-1111 149	20-4220 135 24-1600 165	24-4574 164 24-4575 164	24-7111 186 24-7245 186	30-1070 172 30-1071 172	30-1652 167 30-1670 167	30-1837 49 30-1840 49	30-2090 39 30-2091 39	30-3001 170 30-3002 170
19-1112 149	24-1601 165	24-4576 164	24-7246 186	30-1072 172	30-1671 167	30-1841 46	30-2094 39	30-3003 170
19-1113 149	24-1602 165	24-4577 164	24-7247 186	30-1074 172	30-1672 167	30-1842 49	30-2095 39	30-3004 170
20-1000 134 20-1001 134	24-1606 165 24-1607 165	24-4578 164 24-4579 164	24-7248 186 24-7255 186	30-1075 172 30-1076 172	30-1673 167 30-1684 167	30-1845 47 30-1846 47	30-2100 38 30-2101 38	30-3005 170 30-3006 170
20-1002 134	24-2600 164	24-4900 137	24-7256 186	30-1077 172	30-1686 167	30-1847 47	30-2102 38	30-3007 170
20-1003 134 20-1004 134	24-2601 164 24-2602 164	24-4910 137 24-4912 137	24-7257 186 24-7258 186	30-1078 172 30-1080 172	30-1687 167 30-1692 166	30-1848 49 30-1851 46	30-2103 38 30-2106 38	30-3008 170 30-3009 170
20-1005 134	24-2603 164	24-4914 137	24-7265 186	30-1081 172	30-1695 166	30-1852 46	30-2107 38	30-3010 168
20-1010 134 20-1011 134	24-2604 164 24-2634 164	24-4916 137 24-4920 137	24-7266 186 24-7267 186	30-1082 172	30-1696 166	30-1853 46 30-1854 46	30-2108 39 30-2109 39	30-3011 168 30-3012 168
20-1011 134	24-2635 164	24-4920 137	24-7268 186	30-1083 172 30-1085 172	30-1697 166 30-1698 166	30-1855 46	30-2109 39	30-3012 108
20-1013 134	24-2636 164	24-4960 137	24-7269 186	30-1086 172	30-1700 48	30-1867 45	30-2111 38	30-3014 168
20-1014 134 20-1015 134	24-2637 164 24-2638 164	24-4961 137 24-4962 137	24-7270 186 25-3410 137	30-1088 172 30-1089 172	30-1701 48 30-1702 48	30-1868 45 30-1869 45	30-2112 38 30-2113 38	30-3015 168 30-3016 168
20-1020 134	24-2640 164	24-4963 137	25-3412 137	30-1090 172	30-1703 48	30-1870 45	30-2116 38	30-3017 168
20-1021 134	24-2641 164	24-4964 137	25-3414 137	30-1091 172	30-1704 48	30-1876 48	30-2117 38	30-3018 168
20-1022 134 20-1030 134	24-2642 164 24-2643 164	24-4965 137 24-4966 137	25-3416 137 25-3420 137	30-1092 172 30-1093 172	30-1705 48 30-1709 48	30-1877 48 30-1878 48	30-2120 44 30-2130 44	30-3020 168 30-3021 168
20-1031 134	24-2644 164	24-7000 186	25-3422 137	30-1094 172	30-1713 48	30-1879 48	30-2131 44	30-3025 168
20-1032 134 20-1033 134	24-2650 165 24-2651 165	24-7001 186 24-7002 186	25-3424 137 25-3426 137	30-1095 172 30-1096 172	30-1720 48 30-1721 48	30-1880 48 30-1881 48	30-2132 44 30-2140 39	30-3026 168 30-3027 168
20-1033 134	24-2652 165	24-7002 186	25-3500 137	30-1090 172	30-1721 48	30-1882 48	30-2141 39	30-3027 108
20-1035 134	24-2653 165	24-7004 186	25-3510 137	30-1098 172	30-1723 48	30-1883 48	30-2142 39	30-3031 182
20-1040 134 20-1041 134	24-2656 165 24-2666 165	24-7005 186 24-7006 186	25-3520 137 25-3530 137	30-1170 172 30-1171 172	30-1724 48 30-1725 47	30-1884 48 30-1885 48	30-2143 39 30-2144 39	30-3032 182 30-3040 182
20-1042 134	24-2667 165	24-7007 186	25-3540 137	30-1172 172	30-1726 47	30-1891 47	30-2145 39	30-3050 175
20-1043 134	24-2668 165	24-7010 186	25-3550 137	30-1173 172	30-1727 47	30-1892 47	30-2150 38	30-3061 175
20-1044 134 20-1045 134	24-2669 165 24-2701 165	24-7011 186 24-7012 186	25-3560 137 25-3570 137	30-1174 172 30-1230 57	30-1728 47 30-1729 47	30-1893 47 30-1894 47	30-2151 38 30-2152 38	30-3062 175 30-3070 173
20-1050 134	24-2702 165	24-7013 186	25-3580 137	30-1231 57	30-1735 47	30-1900 45	30-2153 38	30-3071 173
20-1051 134 20-1052 134	24-2703 165 24-2710 165	24-7015 186 24-7016 186	25-3590 137 25-3600 137	30-1232 57 30-1233 57	30-1736 47 30-1737 47	30-1901 45 30-1902 45	30-2161 43 30-2162 43	30-3072 173 30-3074 173
20-1060 134	24-2711 165	24-7017 186	25-3610 137	30-1240 169	30-1740 47	30-1908 44	30-2163 43	30-3075 173
20-1061 134	24-2712 165	24-7018 186	25-3620 137	30-1241 169	30-1741 47	30-1910 44	30-2164 43	30-3076 173
20-1062 134 20-1063 134	24-2713 165 24-2714 165	24-7020 186 24-7021 186	25-3630 137 25-3640 137	30-1242 169 30-1243 169	30-1742 46 30-1743 46	30-1913 44 30-1915 44	30-2165 43 30-2166 43	30-3077 173 30-3078 173
20-1064 134	24-2721 165	24-7022 186	25-3650 137	30-1250 57	30-1744 46	30-1917 44	30-2170 42	30-3079 173
20-1065 134 20-1070 134	24-2722 165 24-2723 165	24-7023 186 24-7024 186	30-0990 169 30-0991 169	30-1251 57 30-1390 172	30-1745 46 30-1746 46	30-1925 45 30-1926 45	30-2171 42 30-2180 39	30-3080 173 30-3081 173
20-1070 134	24-4090 136	24-7025 186	30-0992 169	30-1395 172	30-1773 46	30-1935 33	30-2181 39	30-3082 173
20-1072 134	24-4091 136	24-7026 186	30-1000 57	30-1396 172	30-1776 46	30-1940 182	30-2190 43	30-3083 173
20-1080 134 20-1081 134	24-4100 136 24-4101 136	24-7027 186 24-7030 186	30-1001 57 30-1002 57	30-1410 159 30-1411 159	30-1777 46 30-1778 46	30-1961 46 30-1962 46	30-2191 43 30-2200 38	30-3085 173 30-3086 173
20-1082 134	24-4102 136	24-7031 186	30-1003 57	30-1580 167	30-1779 46	30-1963 46	30-2201 38	30-3089 173
20-1083 134 20-1084 134	24-4110 136 24-4111 136	24-7032 186 24-7033 186	30-1004 57 30-1005 57	30-1585 167 30-1590 166	30-1780 47 30-1781 47	30-1964 46 30-1965 46	30-2210 38 30-2211 38	30-3094 173 30-3097 173
20-1085 134	24-4112 136	24-7054 186	30-1006 57	30-1591 166	30-1782 47	30-1982 45	30-2220 39	30-3101 171
20-1090 134	24-4120 136	24-7055 186	30-1007 57	30-1592 166	30-1783 47	30-1983 45	30-2221 39	30-3105 171
20-1091 134 20-1092 134	24-4121 136 24-4122 136	24-7056 186 24-7057 186	30-1008 57 30-1009 57	30-1593 166 30-1594 166	30-1790 48 30-1791 48	30-1994 140 30-1995 140	30-2230 38 30-2231 38	30-3106 171 30-3107 171
20-1093 134	24-4130 136	24-7064 186	30-1010 57	30-1595 166	30-1792 48	30-1996 140	30-2235 38	30-3108 171
20-1094 134 20-1095 134	24-4131 136 24-4132 136	24-7065 186 24-7066 186	30-1011 57 30-1012 57	30-1596 166 30-1597 166	30-1795 48 30-1796 48	30-1997 140 30-1998 140	30-2236 38 30-2240 38	30-3111 171 30-3112 171
20-4001 135	24-4140 136	24-7067 186	30-1013 57	30-1598 166	30-1797 48	30-1999 140	30-2241 38	30-3140 180
20-4002 135	24-4141 136	24-7070 186	30-1014 57	30-1629 166	30-1800 46	30-2000 171	30-2245 38	30-3141 180
20-4003 135 20-4010 135	24-4142 136 24-4150 136	24-7071 186 24-7072 186	30-1015 57 30-1016 57	30-1630 166 30-1631 166	30-1801 46 30-1802 46	30-2001 171 30-2002 171	30-2246 38 30-2260 38	30-3142 180 30-3143 180
20-4011 135	24-4151 136	24-7073 186	30-1017 57	30-1632 166	30-1803 46	30-2003 171	30-2261 38	30-3152 182
20-4012 135 20-4021 135	24-4152 136 24-4540 164	24-7074 186 24-7075 186	30-1018 168 30-1020 168	30-1633 166 30-1634 166	30-1804 46 30-1805 46	30-2004 171 30-2020 171	30-2270 38 30-2271 38	30-3153 180 30-3160 180
20-4030 135	24-4541 164	24-7076 186	30-1020 100	30-1635 166	30-1806 46	30-2021 171	30-2280 38	30-3161 180
20-4040 135	24-4542 164	24-7077 186	30-1022 168	30-1636 166	30-1807 46	30-2022 171	30-2281 38	30-3171 173
20-4050 135 20-4051 135	24-4543 164 24-4545 164	24-7100 186 24-7101 186	30-1023 168 30-1024 168	30-1637 166 30-1640 166	30-1808 46 30-1809 46	30-2023 171 30-2024 171	30-2290 38 30-2291 38	30-3172 173 30-3174 173
20-4100 135	24-4550 164	24-7102 186	30-1025 168	30-1641 166	30-1825 47	30-2030 171	30-2310 53	30-3180 168
20-4101 135	24-4551 164 24-4552 164	24-7103 186 24-7104 186	30-1026 168 30-1027 168	30-1642 166 30-1643 166	30-1826 47 30-1827 47	30-2031 171 30-2032 171	30-2330 39 30-2331 39	30-3181 168 30-3182 168
20-4110 135 20-4200 135	24-4552 164	24-7104 186	30-1027 168	30-1644 166	30-1827 47	30-2032 171	30-2331 39	30-3182 168
20-4201 135	24-4554 164	24-7106 186	30-1039 49	30-1645 166	30-1829 49	30-2034 171	30-2333 39	30-3184 168
20-4202 135 20-4210 135	24-4570 164 24-4571 164	24-7107 186 24-7108 186	30-1047 49 30-1048 49	30-1646 166 30-1647 166	30-1830 49 30-1831 49	30-2040 171 30-2080 39	30-2340 39 30-2341 39	30-3185 168 30-3186 168
20-4211 135	24-4572 164	24-7109 186	30-1051 49	30-1650 167		30-2081 39	30-2401 53	



# sku-sorted page index

item page	item page	item page	item page	item page	item page	item page	item page	item page
30-3190 170	30-3640 181	31-3028 154	31-3832 178	32-2203 167	38-0314 55	38-2304 55	43-2054 156	43-2560 161
30-3191 170	30-3641 181	31-3200 176	31-3833 178	32-2204 167	38-0315 55	38-2305 55	43-2060 156	43-2565 161
30-3192 170	30-3642 181	31-3201 176	31-3834 178	32-2400 181	38-0316 55	38-2306 55	43-2061 155	43-2570 161
30-3193 170	30-3643 181	31-3210 176	31-3835 178	32-2401 181	38-0317 55	38-2307 55	43-2062 156	43-2571 161
30-3194 170	30-3644 181	31-3221 176	31-3836 178	32-2402 181	38-0400 53	38-2310 55	43-2065 155	43-2580 161
30-3195 170	31-1010 179	31-3224 176	31-3837 178	32-2403 181	38-0401 53	38-2311 55	43-2066 155	43-2581 161
30-3196 170	31-1011 179	31-3230 176	31-3838 178	32-2408 181	38-0402 53	38-2312 55	43-2067 155	45-1840 162
30-3197 170	31-1012 179	31-3250 176	31-3840 178	32-2410 181	38-0410 53	38-2313 55	43-2070 156	45-1841 162
30-3198 170 30-3210 169	31-1013 179 31-1014 179	31-3251 176 31-3255 176	31-3842 178 31-3843 178	32-2416 181 32-2432 181	38-0411 53 38-0412 53	38-2314 55 38-2315 55	43-2071 156 43-2072 156	45-1842 162 45-1850 162
30-3220 182	31-1016 179	31-3420 176	31-3844 178	32-2490 181	38-0500 53	38-2316 55	43-2075 156	45-1851 162
30-3221 182	31-1017 179	31-3430 176	31-3845 178	38-0000 54	38-0501 53	38-2317 55	43-2076 156	45-1852 162
30-3222 182	31-1018 179	31-3440 176	31-3846 178	38-0001 54	38-0502 53	38-2400 53	43-2100 157	45-2190 162
30-3230 170	31-1050 169	31-3441 176	31-3847 178	38-0002 54	38-0510 53	38-2401 53	43-2102 157	45-2191 162
30-3231 170	31-1120 178	31-3450 176	31-3848 178	38-0003 54	38-0511 53	38-2402 53	43-2104 157	45-2192 162
30-3232 170	31-1121 178	31-3460 176	31-3850 179	38-0004 54	38-0512 53	38-2410 53	43-2110 157	45-2195 162
30-3233 170	31-1122 178	31-3461 176	31-3851 179	38-0010 54	38-0550 55	38-2411 53	43-2112 157	45-2200 162
30-3240 169	31-1123 178 31-1124 178	31-3490 177 31-3497 177	31-3860 179 31-3861 179	38-0011 54 38-0012 54	38-0551 55 38-0552 55	38-2412 53 42-1051 133	43-2120 157 43-2122 157	45-2201 162 45-2202 162
30-3241 169 30-3242 169	31-1124 176	31-3498 177	31-3861 179 31-3870 179	38-0012 54	38-2000 54	42-1051 133	43-2122 157	45-2202 162 45-2205 162
30-3242 169	31-1131 178	31-3560 154	31-3880 178	38-0014 54	38-2001 54	42-1053 133	43-2135 157	45-2206 162
30-3250 169	31-1140 178	31-3561 154	31-3885 178	38-0022 54	38-2002 54	42-1054 133	43-2136 157	45-2207 162
30-3251 169	31-1141 178	31-3570 154	31-3890 178	38-0023 54	38-2003 54	42-1060 132	43-2140 157	45-2208 162
30-3252 169	31-1145 178	31-3580 154	32-1000 56	38-0024 54	38-2004 54	42-1061 132	43-2142 157	45-2210 162
30-3253 169	31-1146 178	31-3581 154	32-1001 56	38-0060 54	38-2010 54	42-1062 132	43-2150 157	45-2211 162
30-3280 174	31-1147 178	31-3600 155	32-1002 56	38-0061 54	38-2011 54	42-1070 132	43-2152 157	45-2212 162
30-3281 174	31-1148 178	31-3601 155	32-1003 56	38-0062 54	38-2012 54	42-1071 132	43-2160 157	45-2220 160
30-3282 174	31-1149 178	31-3602 155	32-1004 56	38-0063 54	38-2013 54	42-1072 132	43-2164 157	45-2221 160
30-3283 174 30-3290 174	31-2000 56 31-2001 56	31-3603 155 31-3610 155	32-1005 56 32-1006 56	38-0064 54 38-0070 54	38-2014 54 38-2022 54	42-1080 133 42-1211 132	43-2170 157 43-2171 157	45-2222 160 45-2223 160
30-3290 174	31-2001 56	31-3618 155	32-1000 50	38-0070 54	38-2023 54	42-1211 132	43-2171 157	45-2230 160
30-3292 174	31-2002 56	31-3630 155	32-1007 56	38-0072 54	38-2024 54	42-1214 132	43-2181 157	45-2231 160
30-3293 174	31-2004 56	31-3631 155	32-1009 56	38-0073 54	38-2060 54	42-1215 132	43-2200 158	45-2232 160
30-3305 174	31-2005 56	31-3632 155	32-1011 56	38-0074 54	38-2061 54	42-1216 132	43-2210 158	45-2240 160
30-3306 174	31-2006 56	31-3650 154	32-1012 56	38-0082 54	38-2062 54	42-1217 132	43-2221 158	45-2241 160
30-3307 174	31-2010 56	31-3651 154	32-1014 56	38-0083 54	38-2063 54	42-1223 132	43-2222 158	45-2242 160
30-3308 174	31-2011 56	31-3652 154	32-1015 56	38-0084 54	38-2064 54	42-1225 132	43-2230 158	45-2250 160
30-3310 174	31-2012 56	31-3653 154	32-1230 52	38-0090 54	38-2070 54	42-1228 132	43-2231 158	45-2251 160
30-3311 174 30-3312 174	31-2013 56 31-2014 56	31-3654 154 31-3655 154	32-1232 52 32-1234 52	38-0091 54 38-0093 54	38-2071 54 38-2072 54	42-1245 132 42-1251 132	43-2251 158 43-2252 158	45-2252 160 45-2257 160
30-3312 174	31-2015 56	31-3656 154	32-1234 52	38-0094 54	38-2073 54	42-1253 132	43-2260 158	45-2260 160
30-3320 174	31-2016 56	31-3657 154	32-1238 52	38-0100 54	38-2074 54	42-1255 132	43-2261 158	45-2290 162
30-3326 174	31-2017 56	31-3658 154	32-1240 52	38-0101 54	38-2082 54	42-1256 132	43-2280 158	45-2291 162
30-3327 174	31-2020 56	31-3660 155	32-1247 52	38-0102 54	38-2083 54	42-1258 132	43-2281 158	45-2300 163
30-3328 174	31-2021 56	31-3661 155	32-1248 52	38-0120 54	38-2084 54	42-1260 133	43-2285 158	45-2301 163
30-3329 174	31-2022 56	31-3662 155	32-1380 53	38-0121 54	38-2090 54	42-1400 133	43-2286 158	45-2302 163
30-3330 174	31-2023 56	31-3670 154	32-1383 53	38-0122 54	38-2091 54	42-1401 133	43-2290 158	45-2303 163
30-3331 174 30-3380 160	31-2024 56 31-2032 56	31-3671 154 31-3672 154	32-1386 53 32-1400 51	38-0200 55 38-0201 55	38-2093 54 38-2094 54	42-1405 133 42-1406 133	43-2291 158 43-2292 158	45-2310 163 45-2320 163
30-3381 160	31-2032 56	31-3673 154	32-1400 51	38-0202 55	38-2100 54	42-1400 133	43-2295 158	45-2330 163
30-3382 160	31-2034 56	31-3674 154	32-1410 51	38-0203 55	38-2101 54	42-1410 133	43-2330 160	45-2340 163
30-3390 174	31-2035 56	31-3680 155	32-1411 51	38-0204 55	38-2102 54	42-1411 133	43-2332 160	45-2350 163
30-3500 173	31-2050 56	31-3681 155	32-1420 51	38-0205 55	38-2120 54	42-1412 133	43-2336 160	45-2351 163
30-3501 174	31-2051 56	31-3682 155	32-1421 51	38-0206 55	38-2121 54	42-1420 133	43-2340 160	45-2352 163
30-3502 174	31-2052 56	31-3690 154	32-1430 51	38-0207 55	38-2122 54	42-1430 133	43-2342 160	45-2360 163
30-3503 174	31-2300 181	31-3691 154	32-1440 51	38-0210 55	38-2200 55	42-1431 133	43-2350 160	45-2361 163
30-3504 173 30-3505 173	31-2301 181 31-2600 181	31-3692 154 31-3693 154	32-1441 51 32-1450 51	38-0211 55 38-0212 55	38-2201 55 38-2202 55	42-1432 133 42-1440 133	43-2400 162 43-2402 162	45-2362 163 45-2370 163
30-3506 173	31-2670 181	31-3694 154	32-1450 51 32-1490 51	38-0213 55	38-2202 55	42-1440 133	43-2410 162	45-2371 163
30-3510 174	31-3000 155	31-3705 177	32-1491 51	38-0214 55	38-2204 55	42-1450 133	43-2412 162	45-2375 163
30-3511 174	31-3001 155	31-3706 177	32-1500 44	38-0215 55	38-2205 55	42-1460 133	43-2500 161	45-2380 163
30-3512 174	31-3002 155	31-3707 177	32-2011 180	38-0216 55	38-2206 55	43-1160 156	43-2501 161	45-2381 163
30-3513 174	31-3003 155	31-3710 177	32-2013 180	38-0217 55	38-2207 55	43-1161 156	43-2502 161	45-2390 163
30-3521 173	31-3004 155	31-3711 177	32-2015 180	38-0300 55	38-2210 55	43-1162 156	43-2510 161	45-2395 163
30-3522 173	31-3005 155	31-3720 177	32-2016 180	38-0301 55	38-2211 55	43-1163 156	43-2511 161	45-2400 163
30-3523 173 30-3601 167	31-3010 154 31-3012 154	31-3721 177 31-3730 154	32-2020 43 32-2021 43	38-0302 55 38-0303 55	38-2212 55 38-2213 55	43-1164 157 43-2000 156	43-2512 161 43-2520 161	50-0000 110 50-0002 110
30-3602 167	31-3012 154	31-3730 154	32-2021 43	38-0304 55	38-2214 55	43-2010 156	43-2521 161	50-0002 110
30-3603 167	31-3020 154	31-3732 154	32-2022 43	38-0305 55	38-2215 55	43-2010 156	43-2522 161	50-0003 110
30-3604 167	31-3021 154	31-3733 154	32-2050 180	38-0306 55	38-2216 55	43-2014 156	43-2530 161	50-0005 110
30-3605 167	31-3022 154	31-3800 178	32-2051 180	38-0307 55	38-2217 55	43-2020 156	43-2540 161	50-1000 122
30-3606 167	31-3023 154	31-3802 178	32-2111 182	38-0310 55	38-2300 55	43-2030 156	43-2541 161	50-1001 122
30-3607 167	31-3024 154	31-3810 178	32-2200 167	38-0311 55	38-2301 55	43-2040 156	43-2550 161	50-1002 122
30-3609 167	31-3025 154	31-3812 178	32-2201 167	38-0312 55	38-2302 55	43-2050 156	43-2551 161	50-1005 123
30-3610 167	31-3027 154	31-3830 178	32-2202 167	38-0313 55	38-2303 55	43-2052 156	43-2555 161	50-1007 123



#### sku-sorted page index

item page	item page	item page	item page	item page	item page	item page	item page	item page
50-1008 123	50-1201 147	50-1550 184	50-1851 159	55-1035 72	60-1011 192	61-0009 191	61-0071 191	86-0030 187
50-1009 123	50-1210 147	50-1560 185	50-1852 159	55-1036 72	60-1012 192	61-0014 191	61-0072 191	86-0090 187
50-1010 58	50-1211 147	50-1570 185	50-1880 159	55-1037 72	60-1030 192	61-0015 191	61-0073 191	86-0100 193
50-1023 58	50-1212 147	50-1580 185	50-3000 159	55-1038 72	60-1031 192	61-0020 191	61-0100 191	86-0101 193
50-1024 58	50-1213 147	50-1581 185	50-3001 159	55-1040 73	60-1032 192	61-0021 191	61-0110 191	86-0110 193
50-1035 58	50-1214 147	50-1590 184	50-3004 159	55-1041 73	60-1040 192	61-0022 191	61-0111 191	86-0120 193
50-1036 58	50-1215 147	50-1591 184	50-3005 159	55-1042 73	60-1041 192	61-0023 191	61-0112 191	86-0130 193
50-1040 122	50-1216 147	50-1592 184	50-3007 159	55-1050 73	60-1042 192	61-0030 191	61-0200 191	86-0200 193
50-1041 122	50-1217 147	50-1595 184	50-5120 195	55-1052 73	60-1043 192	61-0031 191	62-0100 190	86-0210 193
50-1045 122	50-1218 147	50-1596 184	50-5121 195	55-1053 73	60-1050 192	61-0032 191	62-0110 190	86-0220 193
50-1046 122	50-1300 159	50-1597 184	50-5130 195	55-1054 73	60-1051 192	61-0033 191	62-0111 190	86-0230 193
50-1050 123	50-1301 159	50-1598 184	50-5131 195	55-1055 73	60-1052 192	61-0036 191	62-0120 190	86-0300 187
50-1051 123	50-1302 159	50-1610 184	50-5132 195	59-0100 196	60-1060 192	61-0037 191	62-0130 190	86-0310 187
50-1052 123	50-1325 159	50-1620 184	50-5150 195	59-0101 196	60-1061 192	61-0038 191	62-0131 190	86-0311 187
50-1070 123	50-1326 159	50-1621 184	50-5500 165	59-0102 196	60-1070 192	61-0040 191	62-0140 190	86-0350 187
50-1075 123	50-1327 159	50-1630 184	50-5510 165	59-0109 196	60-1071 192	61-0041 191	62-0141 190	86-0360 187
50-1080 123	50-1360 159	50-1640 184	50-6055 140	59-0110 196	60-1072 192	61-0042 191	62-0150 190	86-0361 187
50-1081 123	50-1361 159	50-1641 184	50-6056 140	59-0111 196	60-1080 192	61-0043 191	62-0160 190	86-0370 187
50-1095 122	50-1362 159	50-1650 184	50-6057 140	59-0120 196	60-1085 192	61-0044 191	62-0161 190	86-1125 187
50-1096 122	50-1363 159	50-1651 184	50-6061 140	59-0121 196	60-1086 192	61-0046 191	62-0170 190	86-1126 187
50-1100 189	50-1500 185	50-1652 184	50-6067 140	59-0122 196	60-1090 192	61-0047 191	62-0171 190	86-1127 187
50-1101 189	50-1501 185	50-1660 184	50-6068 140	59-0200 196	60-1091 192	61-0048 191	62-0172 190	AD-P300 124
50-1102 189	50-1502 185	50-1800 159	50-6069 140	59-0210 196	60-1100 188	61-0050 191	62-0200 190	AD-P400 124
50-1104 189	50-1505 185	50-1801 159	55-1000 73	59-0220 196	60-1101 188	61-0051 191	62-0201 190	FL-U110 109
50-1120 189	50-1506 185	50-1802 159	55-1010 72	59-0230 196	60-1102 188	61-0052 191	85-0100 194	FL-U115 109
50-1121 189	50-1507 185	50-1810 159	55-1011 72	59-0240 196	60-1110 188	61-0053 191	85-0110 194	IN-T001 114
50-1122 189	50-1508 185	50-1811 159	55-1012 72	59-0250 196	60-1111 188	61-0055 191	85-0120 194	IN-T002 114
50-1130 189	50-1509 185	50-1812 159	55-1013 72	60-0000 188	60-1112 188	61-0056 191	85-0121 194	IN-T0CB 114
50-1131 189	50-1510 185	50-1820 159	55-1014 72	60-0001 188	60-1113 188	61-0057 191	85-0130 194	ME-D000109
50-1132 189	50-1511 185	50-1821 159	55-1015 72	60-0002 188	60-1200 192	61-0058 191	85-0140 194	
50-1133 189	50-1512 185	50-1822 159	55-1020 72	60-0020 188	60-1205 192	61-0060 191	85-0150 194	
50-1134 189	50-1515 185	50-1830 159	55-1021 72	60-0100 188	60-1206 192	61-0061 191	85-0151 194	
50-1139 189	50-1516 185	50-1831 159	55-1022 72	60-0101 188	60-1210 192	61-0062 191	85-0200 194	
50-1140 189	50-1517 185	50-1832 159	55-1030 72	60-0102 188	60-1250 192	61-0063 191	85-0201 194	
50-1141 189	50-1521 185	50-1840 159	55-1031 72	60-0103 188	61-0000 191	61-0066 191	86-0001 187	
50-1142 189	50-1522 185	50-1841 159	55-1032 72	60-1000 192	61-0001 191	61-0067 191	86-0002 187	
50-1150 189	50-1530 185	50-1842 159	55-1033 72	60-1001 192	61-0002 191	61-0068 191	86-0010 187	
50-1200 147	50-1540 185	50-1850 159	55-1034 72	60-1010 192	61-0003 191	61-0070 191	86-0020 187	

#### ordering

- \$50 minimum order
- · Order by phone, fax, email or mail.

#### payment

- For your convenience, FEI accepts Visa, Mastercard, Discover Card and American Express
- Your first purchase will be credit card, C.O.D. or cash with order unless prior credit arrangements are made.

#### shipping and handling

- Drop shipments made to your customer at no additional charge.
- All freight charges and special handling charges are the responsibility of the purchaser and will be added to your invoice.
- Freight charges can be charged directly to your account by using carrier's consignee billing option.
- Some products are shipped directly from the factory or made to order. These will normally be shipped within 10-15 work days.
- Rush shipments can be arranged when merchandise is available. Charges for this service will be added to your invoice.

#### damages

- You must initiate all claims with the delivering carrier.
- Damage to a shipping container must be reported to the carrier at time of delivery.
- If damage to the merchandise is obvious at time of delivery, note damage on the delivery bill and have driver sign it (acknowledging

#### the damage).

If damage to the merchandise is not obvious (hidden damage) at time of delivery, call the carrier within 10 days of delivery and request an inspection. Save the merchandise and all packaging materials for the inspection.

terms and conditions

#### returns

- Obtain a return authorization number (RMA) from FEI's customer service department.
- Return authorization number must appear on each returned package.
- Return the merchandise to the location given by FEI's customer service department.
- Keep a record covering the shipment of returned merchandise (this is for your protection in case you need to proof-of-delivery from the shipper). get a
- No unauthorized return will be accepted.
- No returns will be accepted after 30 days.
- FEI charges a 25% restocking fee for returned merchandise in resalable condition.
- The customer is responsible for return shipping expenses unless the merchandises is received as a result of FEI's error.
- FEI will credit your account for returned merchandise as soon as return is received and processed.
- Credits for return is issued only to original purchaser. Credit will be applied to your open balance or to future purchases. No cash refund will be issued.

· net 30 days

- U.S. funds drawn on U.S. bank
- returned checks subject to \$25.00 fee FEI reserves right to withhold shipment of pending orders or back-orders when account is overdue
- Overdue invoices shall be charged interest at the rate of 11/2% per month
- Customer shall be responsible for any cost of collection
- Terms are subject to change without notice pricing
  - Pricing is subject to change without notice
- Prices do not include freight charges product

· Equipment and supplies illustrated and specified in the FEI catalog are subject to modification consistent with product availability or other factors.

#### warrantv

- FEI offers a minimum 1-year warranty on all products. In cases where the manufacturer includes a warrantee card or an extended warranty, contact the manufacturer directly for complete warranty information.
- FEI reserves the right to repair or replace warranted products at its discretion.

#### sales tax

If your billing or shipping location is in New York, you will be charged sales tax. If you are exempt from paying New York State sales tax please send a copy of your resale certificate with your invoice payment.



	5	ku-sorted pag	ge index		
abductor wedge	168	body composition scale	86.88-89	crutches	155-156
Absolute+Axis <sup>™</sup> goniometer	77	body flexibility unit	89	CryoCup <sup>™</sup> cold therapy	107
accessories (band and tubing)	12-13	body level, measurement	81	CryoDerm <sup>™</sup> lotion	100
accordion fold mat	54-55	BodyLogic® bodyfat analyzer	86	CryoSpray <sup>™</sup> topical	99
AccuForce <sup>™</sup> exercise band	7	bolster, positioning	56-57,147,170	CryoStim <sup>™</sup> ice massager	107
AcuForce <sup>™</sup> massage tool	138	Bosu® balance exerciser	45	Cuff® weight/set	26-27
AcuAngle <sup>™</sup> inclinometer	80	Bow-Tie <sup>™</sup> resistive exercise	9	cushion, insole	125
ADA LOPS / LEAP monofilament	82	boxer shorts (disposable)	134	cushions, wheelchair	159
Adapta® table	124,142	BROM inclinometer	80	Cushy-Air <sup>™</sup> ball	47
adaptive positioning	166-183	bubble inclinometer	80	Cushy-Air™ mat	55
add-a-mat <sup>™</sup>	54	Buck hammer	84	D-ring handle/strap	12
adjustable exercise band system	10	Buhl spirometer	89	Detecto® professional scale	88
aesthesiometer	82-83	bungee cord (exercise)	10	dexterity test	91
AirCast® ankle brace	165	C-Trax <sup>™</sup> traction	122	diathermy	108
AirCast® CryoCuff® cold therapy	105	cable tensiometer	71	Disckson® paraffin	97
Airex® mats, and balance	44,52	calipers, breadth	87	Digi-Extend™ finger exerciser	20
alarm, patient	196	calipers, skinfold	86	Digi-Extend n' Squeeze™ exerc	
Amara® wheelchair cushion	159	Cambion™ insole	125	Digi-Flex® finger exerciser	18
ambulation training Amrex® electrotherapy	152-153 115,118	Cando® exercise product  • ball product	45-49	Dipsters® patientwear discriminator, sensory	134 83
anatomical model	92-95	ball product     band product	45-49 4-16	•	
anchor, exercise band accessory	13	hand exerciser	18-25	dispenser rack (exercise band)  DME	154-163
ankle	13	weight product	28-31	dolorimeter/algorimeter	75,83
brace	164-165	canes	156	Dolphin <sup>™</sup> chair	162
elastic cuff exerciser	104-103	carpal tunnel stretchers	20	dressing aid, ADL	187
Elgin® leg/ankle exerciser	32	Carrie <sup>™</sup> seats/accessories	160,174	drinking cup, ADL	192
stretching exerciser	59	Carpal-Trac™	106	drop-foot brace/AFO	165
Cuff® weight	26-27	car seat	179	dual handle medicine ball	30
aquatic therapy product	132-135	carts, electrotherapy	112-114,117	Dualer IQ, J-tech	74
arm elevator	159	Cats Claw <sup>™</sup> lotion	100	dumbbell, vinyl-coated	29
Armasport® mat	50-51	Cellex® medium for Fluidothe	erapy® 109	dumbbell, aquatic	135
arthrodial protractor	81	center-fold mats	53-55	Dura-Stick <sup>™</sup> electrode	11
Artromot® CPM	129-130	cervical pillows	147	Dycem® non-slip material	184-185
Auto-Therm® diathermy	108	cervical traction	122-125	dynamometer	
Babinski hammer	84	chair		<ul> <li>back-leg-chest</li> </ul>	70-71
Backnobber® massager	139	<ul> <li>adaptive positioning 15</li> </ul>	6-159,166-167	<ul> <li>cable tensiometer</li> </ul>	71
balance		• ball	49	<ul> <li>hand hydraulic</li> </ul>	61-63,66,69
• beam	44,182	• bath	149	<ul> <li>hand spring-type</li> </ul>	62
• pad/board	40-45	<ul><li>bean bag</li></ul>	158	<ul> <li>squeeze (pneumatic)</li> </ul>	62
• stone (inflatable)	45	• massage	129	manual muscle tester	68-70,74-75
ball pit	181	Chango® boards and paw	42	• push/pull	70-71
ball product	40	chest weight pulley	32	• wrist/forearm	67
• chair	49	chilling unit	104-105	Dynex® isometric exerciser	74
• inflatable	44-49	circumference tape	87	eating aid, ADL	192
medicine/plyometric     aguagable	30-31	cleanwheel pinwheel ColPac® cold pack	83	edema gauge Eggsercizer® hand therapy	87 21
<ul><li>squeezable</li><li>storage/inflation</li></ul>	21 49	Cold n'Hot compress	105 106	Ekho® Heart Rate watch	90
Bailey® treatment table	144	cold firriot compress cold pack, cover/accesory	104-107	electrode, electrotherapy	111,131
band (resistive exercise)	4-15	cold spray	11,16,98-100	elliptical bike	37
band (kinesiology taping)	137	Colombia® positioning produc		EMS	119
barrel roll	181	Collyer breadth caliper	87	Endorphin® parallel dip	32
Baseline® measurement	61-91	ComforTrac® traction	123	Endorphin® pedal exerciser	36
bathroom product	0.0.	Commander, J-Tech	75	English ice cap	107
• adult	160-163	commode	160	EnviroSafe <sup>™</sup> foam mat	54-55
pediatric	160,162	composite foam roller	38	ergogrip hand exerciser	19
• scales	88-89	compression sock	180	Ergobasic <sup>™</sup> treatment table	143
Be-better™ kit, band/spray	11	compression therapy	127	extremity strap	13
bedroom aid, ADL		conductive spray	111	EVA foam roller	38
Big Grip™ hand exerciser	21	Core-Tex <sup>™</sup> balance	43	face cradle	146
biofeedback	131	core-training dome	45	facilitation massager	138
BioFreeze® cold therapy	99	CPM units and accessories	128-130	fat caliper, skinfold	86
BP/pulse monitor	85	crater bowl	181	FCE	72-75
Bodyblade® exerciser	25	crawler	180	feeder seat/accessory	172-175



feet ball, inflatable   47   household gadqet, ADL   188   massage table/chair   140   1		sku-sorted pag	ge index	(	
					141
- Indidater (Indidater (Indidater (India ) open of the control open open open open open open open open	finger		61-63,66,69	massager, lotion/warmer	140
- ordinometer of month of the weight of the weig	• circumference gauge 87	Hudson® ergometer	35-37	mat	
· motion augue         76         Inchine board         118         · incline, various shapes         56.57           Fitter he balance         4.23         incline wedge         56.57         - thick (landing)         56.687           Fitter he balance         4.23         incline wedge         56.57         - thick (landing)         56.687           Floxt-Pac' hol/cold pack         106         - exercise ball         46.49         mechanical pinch gauge         108           6 hall         21         - exercise ball         46.49         medicine ball         21           1 cushion         159         rorl         3.70         40         metronome         9           4 handle         12.13         - seans ball         46         46         metronome         9           1 fluidher prival         3.83         7 roll         3.83         1.44         3.44         4         4           1 fluidher prival         1.55         4 seans ball         4.5         4.5         4.5         4.6         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7         4.7	• ladder (shoulder/finger) 58	Hydrocollator® hot/cold	103	<ul> <li>group exercise</li> </ul>	54-55
Filter   Palance   43-23   incline board   59   *lable   144-145			105,107	<ul> <li>individual rest/exercise</li> </ul>	50-53
Filter   Filtr   Filt	• motion gauge 78	IFC	118	<ul> <li>incline, various shapes</li> </ul>	56
Flexation   Flex	• weight 24	incline board	59	• table	144-145
IdeathBility lester (sit and reach)		incline wedge	56-57	<ul><li>thick (landing)</li></ul>	56
Flexis-Pack hot/coll pack   62   - sterping stone   46-49   medicaine ball   30-31	-		80		56-87
FlaxIP-Pac* hot/cold pack   16					
	0 0 0	•		•	
ball         21         - massage ball         140         metronome         90           - cushion         159         - roll         37,47         Mettler® electrotherapy         116-117,19           • handle         12-13         - sensi ball         46         MicroFet* muscle tester         74           • roller (foam roller)         30-39         - wodge         45         mirror (glassless, glass)         172-173,177           Fluidotherapy* dry heat         172-173,176-177         foot         instant cold pack         107         model anatomical         92-95           - brace (drop-foot)         151         stool         135         support/brace         150-151         instant ice, lopical         100         interferential         110         model foam rollers         9.97,101-103,106           G-5° percussor/massager galt belt galt variace electrotherapy         167         interferential         118         118         miderferential positioning system         41         114           G-5° percussor/massager galt belt galt variace electrotherapy         115,118,131         deschober* massager         130         deschober* massager         130         MeckPro* cervical traction         122         descknober* device electrotherapy         117         Mybro* rocker/wobble board         40         My	·				
- cushion				•	
+ nandle		<u> </u>			
- mal					
Fulidotherapy® dy heat					
Fluidotherapy® dry heat   109   inflator pump   49   model, anatomical   32-95		• •		, ,	
Floor sitter, positioning   172-173,176-177   Insole, shoe   125   Instant cold pack   107   Instant cold pack   107   Instant loe, topical   100   Instant loe	` ,	· ·			
foot         - brace (drop-foot)         151         Instant cold pack         107         moist heat         96-97,101-103,106           - stool         135         Instant ciot, topical         100         monifilament         82           - stool         135         Intelect <sup>®</sup> electrotherapy         111-114,118-121         multi axial positioning system         41           full body swing         167         intermittent compression therapy         118         multi axial positioning system         41           6-8° percussor/massager         126         intermittent compression therapy         110         MVP* rocker/wobble board         40           6-8° percussor/massager         126         Jacknobber* massager         139         neurologic hammer         40           gait trainer         155         Jamaré* measurement         62,84,66         Mcmorter* Tens         113           golivanic electrotherapy         115,118,131         job task analysis         67-72         Nimbo* Stander         154           golivaric electrotherapy         115,118,131         job task analysis         67-72         Nimbo* Stander         154           golivaric electrotherapy         115,118,131         job task analysis         67-72         Nimbo* Stander         154           plantin<					
• stool value         135         Intelect® electrotherapy         111-114,118-121         multi axial positioning system         41           full body swing         167         interferential         118         muscle strength testing         61-75           full body swing         167         intermittent compression therapy         127         Myssage® totion         111,140           G-8° percussor/massager         126         intermittent compression therapy         127         Myssage® totion         111,140           G-8° percussor/massager         126         intermittent compression therapy         120         MVP" rocker/wobble board         40           G-8° percussor/massager         129         Jacknobber® massager         139         neurologic hammer         40           galt belt         195         Jacknobber® massager         139         neurologic hammer         84           galvanic electrotherapy         1115,118,131         job task analysis         67-72         Nimb® Stander         154           genider         76-78         Jumb of sandsvis         66,67         NMES         118           gonimeter         - finger         76-78         kettlebell         30         Optifies® knee CPM         129-130           plastic         76-77 <t< td=""><td></td><td>·</td><td></td><td></td><td></td></t<>		·			
• support/brace         150-151         interferential         118         muscle strength testing         61-75           full body swing         167         intermittent compression therapy         127         Myossage® lotion         111,140           furniture (clinic, hospital)         141-153         iontophoresis         110         MYP" rocker/wobble board         40           G-5° percussor/massager         126         ion weight         32         Alexhoror         40           gait trainer         155         Jacknobber* massager         139         neuroteck*         184           gait trainer         155         Jacknobber* massager         139         neuroteck** TENS         119           gel balls         21         J-Tech** evaluation         66,67         Nimbo** Stander         154           goniometer         jumb polal, inflatable         47         Nistonde (pool)         135           + large joint         79         Kickboard, swim         135         Optifies* was expered         104           + x-ray visible         79         kickboard, swim         135         orthosis, AFO         165           gulk measurement lape         87         kinesiology tape         137         orthosis, AFO         165	` ' '	•			
full body swing         167         Intermittent compression therapy         127         Myossage® totion         111,140           furniture (clinic, hospital)         141-153         ionothoproresis         110         MCVP" rocker/wobble board         40           6-5° percussor/massager         126         iono weight         32         NeckPro" cervical traction         122           gait belt         195         Jacknobbere" massager         139         neurologic harmor         84           gait trainer         155         Jamare measurement         62,64,66         neurologic harmor         81           gel balls         21         J-Tech® evaluation         66,67         NMES         118           goniometer         jumbo floam roller         38         noodle (pool)         135           + large joint         79         Kaye® walker         154         Nylatex® wraps         104           + stainless steel         78-77         kettlebell         30         Orfti-spinting material         136           + stainless steel         78-79         kinesbology tape         137         orthosix, AFO         165           *-ray visible         79         kinesbology tape         137         orthosix, AFO         165           Grj-					
Furniture (clinic, hospital)				o o	
G-5º percussor/massager 126 iron weight 32 NeckPro" cervical traction 122 gait belt 195 Jacknobber® massager 139 neurologic hammer 84 pait trainer 155 Jamar® measurement 62,64,66 NeuroTech® TENS 119 galvanic electrotherapy 115,118,131 job task analysis 67-72 Nimbo "Stander 154 place balls 21 J-Tech® evaluation 66,67 NMES 1118 goniometer jumbo foam roller 38 noodle (pool) 135 inger 14 noodle (pool) 135 jump bali, inflatable 47 Nylatex® wraps 104 Nimbo "Stander 129-130 orthosix orthosix orthosix orthosix orthosix orthosix orthosix orthosix AFO 165 orthosix orthosix AFO 165 orthosix AFO 165 orthosix AFO 165 orthosix AFO 165 orthosix orthosix AFO 165 orthosix AFO 1	, 0	· ·	1 7		
gait belt         195         Jacknobber® massager         139         neurologic hammer         84           gait trainer         155         Jamar® measurement         62,64,66         NeuroTech® TENS         119           galvanic electrotherapy         115,118,131         job task analysis         67-72         NIMBO         154           gel balls         21         J-Tech® evaluation         66,67         NMES         118           goniometer         jumbo foam roller         38         noodle (pool)         135           + finger         76-78         jump ball, inflatable         47         Nylatex® wraps         104           + large joint         79         Kaye® walker         154         OptiFlex® knee CPM         129-130           + plastic         76-77         kettlebell         30         Orfit* splinting material         136           + x-ray visible         79         kickboard, swim         135         orthotix         76           gulick measurement tape         79         Kickboard, swim         132         orthotix         76           Grip-It Dycem® material         185         klips® band accessory         13         overdoor traction         122           Gymnic® ball/roll         48         <	, , ,	·			
gait trainer         155         Jamar® measurement         62,64,66         NeuroTech® TENS         119           galvanic electrotherapy         115,118,131         job task analysis         67-72         Nimbo® Stander         154           gel balls         21         J-Ech® evaluation         66,67         NMES         118           goniometer         jumb of fainger         76-78         jumb ball, inflatable         47         Nylatex® wraps         104           • large joint         79         Kaye® walker         154         Optifilex® knee CPM         129-130           • plastic         76-77         kettlebell         30         Orfit* splinting material         136           • stainless steel         78-79         kickboard, swim         135         orthosis, AFO         165           Grip-It Dycem® material         185         klips® band accessory         13         orthosis, AFO         162           Grip-It Dycem® material         185         klips® band accessory         13         overdoor pulley         58           gulick measurement tape         87         knee wedge         57         parallel dip         32           dymnic® ball/roll         48         knee wedge         57         parallel dip         32		_			
galvanic electrotherapy         115,118,131         job task analysis         67-72         Nimbo "Stander         154           gel balls         21         J-Tech® evaluation         66,67         NMES         118           goniometer         jump to foam roller         38         noodle (pool)         135           • Iarge joint         79         Kaye® walker         154         Optiflex® knee CPM         129-130           • Iarge joint         79         kaye® walker         154         Optiflex® knee CPM         129-130           • stainless steel         78-79         kickboard, swim         135         orthosis, AFO         165           • stainless steel         78-79         kickboard, swim         135         orthosis, AFO         165           graded pinch exerciser         24         Kinestek® CPM         128         Orthotix         76           Grip-It Dycem® material         185         klips® band accessory         13         overdoor pulley         58           gulkic measurement tape         87         knee supports         164         overdoor pulley         58           gulkic measurement tape         87         knee wedge         57         parallel dip         overdoor pulley         58           han		· ·		· ·	
gel balls         21         J-Tech® evaluation         66,67         NMES         118           goniometer         jumbo foam roller         38         noodle (pool)         135           • finger         76-78         jumb ball, inflatable         47         Nylatex® wraps         104           • large joint         79         kaye® walker         154         OptiFlex® knee CPM         129-130           • plastic         76-77         kettlebell         30         Orft" splinting material         136           • stainless steel         78-79         kickboard, swim         135         orthosix, AFO         165           • x-ray visible         79         kinesiology tape         137         orthotix         76           graded pinch exerciser         24         Kinetek® CPM         128         Otter" bathing chair         162           Grip-It Dycen® material         185         kips" band accessory         133         overdoor pulley         58           gulick measurement tape         87         knee supports         164         overdoor pulley         58           gulick measurement tape         87         knee supports         164         overdoor traction         122           Aymine* Dall/roll         48	9				
goniometer         jumbo foam roller         38         noodle (pool)         135           • finger         76-78         jump ball, inflatable         47         Nylatex® wraps         104           • large joint         79         Kaye® walker         154         OpfiFlex® kree CPM         129-130           • plastic         76-77         kettlebell         30         Orfit* splinting material         136           • x-ray visible         79         kickboard, swim         135         orthoxis, AFO         165           graded pinch exerciser         24         Kinetek® CPM         128         Otter* bathing chair         162           Grip-It Dycem® material         185         klips* band accessory         13         overdoor pulley         58           gulick measurement tape         87         knee supports         164         overdoor pulley         58           gulic ball/roll         48         knee wedge         57         pain relief, gel and spray         98-100           half roll         57,146,170         Knobble® massager         139         paraffin: wax and paraffin unit         96-97           hamd valuation set         66         lamp, exam (gooseneck)         148         patient room chair         148 <t< td=""><td></td><td>•</td><td></td><td></td><td>118</td></t<>		•			118
<ul> <li>• finger</li> <li>• large joint</li> <li>• large joint</li> <li>• large joint</li> <li>• Pastic</li> <li>• 76-77</li> <li>• kettlebell</li> <li>• Stainless steel</li> <li>•</li></ul>	· ·				
• plastic 76-77 kettlebell 30 Orfit™ splinting material 136 • stainless steel 78-79 kickboard, swim 135 orthosis, AFO 165 • x-ray visible 79 kinesiology tape 137 orthotix 76 graded pinch exerciser 24 Kinetek® CPM 128 Otter™ bathing chair 162 Grip-It Dycem® material 185 klips™ band accessory 13 overdoor pulley 58 gulick measurement tape 87 knee supports 164 overdoor traction 122 Gymnic® ball/roll 48 knee wedge 57 pain relief, gel and spray 98-100 half roll 57,146,170 Knobble® massager 139 paraffin: wax and paraffin unit 96-97 hammer, sensory evaluation 84 Lafayette® MMT/pegboard 75,91 parallel dip 32 hand dynamometer 61-63,66,69,75 lamp, exam (gooseneck) 148 handle, dynamometer 67-68,71 laser, therapeutic 121 patientwear/swimwear 134 handle, exercise 12 Legend®, Intelect® 112-114,121 pediatric positioning 168-180 hand weight 29 leg brace 165 hand/wrist exerciser 26-27 lifting boxes 72 Handly weight ball 29 loop, resistive exercise 7 head/neck support 172,174 loop stirrup 10ton, oulf*hoat therapy 98-100 head-neck support 172,174 loop stirrup 10ton, ould/hot therapy 98-100 heatr-rate monitor 90 lotion, cold/hot therapy 98-100 height measurement 88 low volt, electrotherapy 114-115 halder 12 manual muscles tester/MMT 68-70,74-75 high chair 148 hilder 172 manual muscles tester/MMT 68-70,74-75 hook and loop fastener 186 hread-neck suployed plototer 167 hook and loop fastener 186 halfies 186 halfies 186 horthoria 187 horthoria 187 horthoria 187 horthoria 187 horthoria 187 horthoria 188 horthoria 189 horthoria 189 horthoria 189 horthoria 189 horthoria 189 horthoria	_	jump ball, inflatable	47		104
* stainless steel 78-79 kickboard, swim 135 orthosis, AFO 165  * x-ray visible 79 kinesiology tape 137 orthotix 76 graded pinch exerciser 24 Kinetek® CPP 128 Otter bathing chair 162 Grip-lt Dycem® material 185 kijps® band accessory 13 overdoor pulley 58 gulick measurement tape 87 knee supports 164 overdoor traction 122 Gymnic® ball/roll 48 knee wedge 57 pain relief, gel and spray 98-100 half roll 57,146,170 Knobble® massager 139 paraffin: wax and paraffin unit 96-97 hammer, sensory evaluation 84 Lafayette® MMT/pegboard 75,91 parallel dip 32 hand dynamometer 61-63,66,69,75 lamp, exam (gooseneck) 148 patient room chair 148 hand evaluation set 66 lamp, UV, IR-UV, vitamin D 11 patientwear/swimwear 134 handle, dynamometer 67-68,71 laser, therapeutic 121 pediatric positioning 168-180 hand weight 29 leg brace 185 pedometer 90 HandMaster Plus® exerciser 20 lever handle, wrist dynamometer 67 peg board 60,91 hand/wrist exerciser 26-27 lifting boxes 72 PEP® pack 4,8 Handy weight ball 29 loop, resistive exercise 7 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100® exercise band 5-6 head protector 167 lotion, ultrasound 111-140 Purdue® peg board 91 heart-rate monitor 90 lotion, cold/hot therapy 98-100 heart-rate monitor 148 lumbar traction 104-106 pinwheels, sensory evaluation 83-84 hiRes® plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31	• large joint 79	Kaye <sup>®</sup> walker	154	OptiFlex® knee CPM	129-130
• x-ray visible 79 kinesiology tape 137 orthotix 76 graded pinch exerciser 24 Kinestek® CPM 128 Otter bathing chair 162 Orip-It Dycem® material 185 klips® band accessory 13 overdoor pulley 58 knee supports 164 overdoor traction 122 Gymnic® ball/roll 48 knee wedge 57 pain relief, gel and spray 98-100 half roll 57,146,170 Knobble® massager 139 paraffin: wax and paraffin unit 96-97 hammer, sensory evaluation 84 Lafayette® MMT/pegboard 75,91 parallel dip 32 hand dynamometer 61-63,66,69,75 lamp, exam (gooseneck) 148 patient room chair 148 hand evaluation set 66 lamp, UV, IR-UV, vitamin D 11 patientwear/swimwear 134 handle, dynamometer 67-68,71 laser, therapeutic 121 pedal/arm exercisers 33-37 handle, exercise 12 Legend®, Intelect® 112-114,121 pediatric positioning 168-180 hand weight 29 leg brace 165 pedometer 90 HandMaster Plus® exerciser 26-27 lifting boxes 72 PEP® pack 4,8 Handy weight ball 29 loop, resistive exercise 7 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100® exercise band 5-6 head protector 167 lotion, ultrasound 111-140 Purdue® peg board 91 heart-rate monitor 90 lotion, cold/hot therapy 98-100 pillow/support 147 heating pad 101 low powder, exercise band 4-5 pinch gauge 62-66,69,74-75 high chair 148 lumbar traction 104-106 pinwheels, sensory evaluation 83-84 HiRes® plastic goniometer 76-77 lumbar roll/support 122-125 pivotal therapy system 125 Hi Volt electrotherapy 114-115 Magnecizer® pedal, shoulder 46 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy/plyoball 30-31	• plastic 76-77	kettlebell	30	Orfit™ splinting material	136
graded pinch exerciser 24 Kinetek® CPM 128 Otter™ bathing chair 162 Grip-It Dycem® material 185 klips™ band accessory 13 overdoor pulley 58 gulick measurement tape 87 knee supports 164 overdoor traction 122 Gymnic® ball/roll 48 knee wedge 57 pain relief, gel and spray 98-100 half roll 57,146,170 Knobble® massager 139 paraffin: wax and paraffin unit 96-97 hammer, sensory evaluation 84 Lafayette® MMT/pegboard 75,91 parallel dip 32 hand dynamometer 61-63,66,69,75 lamp, exam (gooseneck) 148 hand evaluation set 66 lamp, UV, IR-UV, vitamin D 11 patientwear/swimwear 134 handle, dynamometer 67-68,71 laser, therapeutic 121 pedal/arm exercisers 33-37 handle, exercise 12 Legend®, Intelect® 112-114,121 pediatric positioning 168-180 hand/wrist exerciser 20 lever handle, wrist dynamometer 67 peg board 60,91 hand/wrist exerciser 26-27 lifting boxes 72 PEP™ pack 4,8 hand weight 129 loop, resistive exercise 7 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100™ exercise band 5-6 head protector 167 lotion, ultrasound 111-140 Purdue® peg board 91 heart-rate monitor 90 lotion, cold/hot therapy 98-100 pillow/support 147 heating pad 101 low powder, exercise band 4-5 pinch gauge 62-66,69,74-75 high chair 148 lumbar traction 104-106 pinwheels, sensory evaluation 83-84 HiRes™ plastic goniometer 76-77 lumbar roll/support 122-155 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31	• stainless steel 78-79	kickboard, swim	135	orthosis, AFO	165
Grip-It Dycem® material 185 klips™ band accessory 13 overdoor pulley 58 gulick measurement tape 87 knee supports 164 overdoor traction 122 Gymnic® ball/roll 48 knee wedge 57 pain relief, gel and spray 98-100 paraffin: wax and paraffin unit 96-97 hammer, sensory evaluation 84 Lafayette® MMT/pegboard 75,91 parallel dip 32 hand evaluation set 66 lamp, UV, IR-UV, vitamin D 11 patient wear/swimwear 134 handle, dynamometer 67-68,71 laser, therapeutic 121 padient exercises 12 Legend®, Intelect® 112-114,121 pedal/arm exercisers 33-37 handle, exercise 20 lever handle, wrist dynamometer 67 peg board 60,91 hand/wrist exerciser 26-27 lifting boxes 72 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100™ exercise band 5-6 head protector 167 lotion, ultrasound 111-140 pillow expresses 140 low powder, exercise band 4-5 heading pad 101 low powder, exercise band 4-5 high chair 148 lumbar traction 104-106 pinwheels, sensory evaluation 83-84 high chair 148 lumbar traction 104-106 plastic goniometer 76-77 hook and loop fastener 186 where ye poly paraffin: wax and paraffin unit 96-97 parallel dip 20 patient room chair 148 pedal/arm exercisers 33-37 pedal/arm exercisers 33-37 pedal/arm exerciser 32-62-77 pedometer 90 pedometer 90 pedometer 90 per postord 60,91 per pedometer 90 pedometer 90 per postord 90 per paraffin: wax and paraffin unit 96-97 pedal/arm exerciser 124 per pedal/arm exercise 149 pedal/arm exerciser 149 pedal/arm exerciser 140 pedal/arm	• x-ray visible 79	kinesiology tape	137	orthotix	76
gulick measurement tape Gymnic® ball/roll A8 knee wedge 57 pain relief, gel and spray 98-100 half roll 57,146,170 Knobble® massager 139 paraffin: wax and paraffin unit 96-97 hammer, sensory evaluation 84 Lafayette® MMT/pegboard 75,91 hand dynamometer 61-63,66,69,75 lamp, exam (gooseneck) 148 patient room chair 148 hand evaluation set 66 lamp, UV, IR-UV, vitamin D 11 patientwear/swimwear 134 handle, dynamometer 67-68,71 laser, therapeutic 121 pedal/arm exercisers 33-37 handle, exercise 12 Legend®, Intelect® 112-114,121 pediatric positioning 168-180 hand weight 29 leg brace 165 pedometer 167 peg board 169,91 hand/wrist exerciser 20 lever handle, wrist dynamometer 67 peg board 169,91 hand/wrist exerciser 20 loop, resistive exercise 7 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100™ exercise band 5-6 head protector 167 lotion, ultrasound 111-140 heating pad 101 low powder, exercise band 104-106 heating pad 101 low powder, exercise band 4-5 height measurement 88 low volt, electrotherapy 115 pinch exerciser 24 height measurement 176-77 high chair 148 lumbar traction 180 volt electrotherapy 115 pinch exerciser 126 pivotherapy ystem 125 pivotal therapy system 125 plastic goniometer 76-77 hook and loop fastener 186 wartin-type vigorimeter 62 plyotherapy/plyoball 30-31	graded pinch exerciser 24		128	Otter™ bathing chair	162
Gymnic® ball/roll48knee wedge57pain relief, gel and spray98-100half roll57,146,170Knobble® massager139paraffin: wax and paraffin unit96-97hammer, sensory evaluation84Lafayette® MMT/pegboard75,91parallel dip32hand dynamometer61-63,66,69,75lamp, exam (gooseneck)148patient room chair148handle, dynamometer67-68,71laser, therapeutic121pedal/arm exercisers33-37handle, exercise12Legend®, Intelect®112-114,121pediatric positioning168-180hand weight29leg brace165pedometer90HandMaster Plus® exerciser20lever handle, wrist dynamometer67peg board60,91hand/wrist exerciser26-27lifting boxes72PEP™ pack4,8Handy weight ball29loop, resistive exercise7percussor126head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch exerciser24high chair148lumbar traction104-106pinwheels, sensory evaluation83-84 <td></td> <td></td> <td>13</td> <td>. ,</td> <td>58</td>			13	. ,	58
half roll 57,146,170 knobble® massager 139 paraffin: wax and paraffin unit 96-97 hammer, sensory evaluation 84 Lafayette® MMT/pegboard 75,91 parallel dip 32 lamp, exam (gooseneck) 148 patient room chair 148 hand evaluation set 66 lamp, exam (gooseneck) 148 patient room chair 148 handle, dynamometer 67-68,71 laser, therapeutic 121 pedal/arm exercisers 33-37 handle, exercise 12 Legend®, Intelect® 112-114,121 pediatric positioning 168-180 hand weight 29 leg brace 165 pedometer 90 lever handle, wrist dynamometer 67 peg board 60,91 hand/wrist exerciser 26-27 lifting boxes 72 PEP™ pack 4,8 landy weight ball 29 loop, resistive exercise 70 loop stirrup 144 Perf-100™ exercise band 5-6 head/neck support 172,174 loop stirrup 144 Perf-100™ exercise band 5-6 head protector 167 lotion, ultrasound 111-140 Purdue® peg board 91 heart-rate monitor 90 lotion, cold/hot therapy 98-100 pillow/support 147 heating pad 101 low powder, exercise band 4-5 pinch exerciser 24 height measurement 88 low volt, electrotherapy 115 high chair 148 lumbar traction 104-106 pinchese, sensory evaluation 83-84 HoldRite™ handle 12 manual muscles tester/MMT 68-70,74-75 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31			164		122
hammer, sensory evaluation84Lafayette® MMT/pegboard75,91parallel dip32hand dynamometer61-63,66,69,75lamp, exam (gooseneck)148patient room chair148hand evaluation set66lamp, UV, IR-UV, vitamin D11patient room chair148handle, dynamometer67-68,71laser, therapeutic121pedal/arm exercisers33-37handle, exercise12Legend®, Intelect®112-114,121pediatric positioning168-180hand weight29leyer handle, wrist dynamometer67ped ometer90hand/wrist exerciser26-27lifting boxes72PEP™ pack4,8Handy weight ball29loop, resistive exercise7percussor126head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch exerciser24high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34 <td></td> <td></td> <td></td> <td></td> <td></td>					
hand dynamometer 61-63,66,69,75 lamp, exam (gooseneck) 148 patient room chair 148 hand evaluation set 66 lamp, UV, IR-UV, vitamin D 11 patientwear/swimwear 134 handle, dynamometer 67-68,71 laser, therapeutic 121 pedal/arm exercisers 33-37 handle, exercise 12 Legend®, Intelect® 112-114,121 pediatric positioning 168-180 hand weight 29 leg brace 165 pedometer 90 HandMaster Plus® exerciser 20 lever handle, wrist dynamometer 67 peg board 60,91 hand/wrist exerciser 26-27 lifting boxes 72 PEP™ pack 4,8 Handy weight ball 29 loop, resistive exercise 7 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100™ exercise band 5-6 head protector 167 lotion, ultrasound 111-140 Purdue® peg board 91 heart-rate monitor 90 lotion, cold/hot therapy 98-100 pillow/support 147 heating pad 101 low powder, exercise band 4-5 pinch exerciser 24 height measurement 88 low volt, electrotherapy 115 pinch gauge 62-66,69,74-75 high chair 148 lumbar traction 104-106 pinwheels, sensory evaluation 83-84 HoldRite™ handle 12 manual muscles tester/MMT 68-70,74-75 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31					
hand evaluation set66lamp, UV, IR-UV, vitamin D11patientwear/swimwear134handle, dynamometer67-68,71laser, therapeutic121pedal/arm exercisers33-37handle, exercise12Legend®, Intelect®112-114,121pediatric positioning168-180hand weight29leg brace165pedometer90HandMaster Plus® exerciser20lever handle, wrist dynamometer67peg board60,91hand/wrist exerciser26-27lifting boxes72PEP™ pack4,8Handy weight ball29loop, resistive exercise7percussor126head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch exerciser24high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goni		, ,			
handle, dynamometer67-68,71laser, therapeutic121pedal/arm exercisers33-37handle, exercise12Legend®, Intelect®112-114,121pediatric positioning168-180hand weight29leg brace165pedometer90HandMaster Plus® exerciser20lever handle, wrist dynamometer67peg board60,91hand/wrist exerciser26-27lifting boxes72PEP™ pack4,8Handy weight ball29loop, resistive exercise7percussor126head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter6262-				•	
handle, exercise12Legend®, Intelect®112-114,121pediatric positioning168-180hand weight29leg brace165pedometer90HandMaster Plus® exerciser20lever handle, wrist dynamometer67peg board60,91hand/wrist exerciser26-27lifting boxes72PEP™ pack4,8Handy weight ball29loop, resistive exercise7percussor126head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31				•	
hand weight29leg brace165pedometer90HandMaster Plus® exerciser20lever handle, wrist dynamometer67peg board60,91hand/wrist exerciser26-27lifting boxes72PEP™ pack4,8Handy weight ball29loop, resistive exercise7percussor126head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31				•	
HandMaster Plus® exerciser 20 lever handle, wrist dynamometer 67 peg board 60,91 hand/wrist exerciser 26-27 lifting boxes 72 PEP™ pack 4,8 Handy weight ball 29 loop, resistive exercise 7 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100™ exercise band 5-6 head protector 167 lotion, ultrasound 111-140 heating pad 101 low powder, exercise band 4-5 height measurement 188 low volt, electrotherapy 115 pinch gauge 62-66,69,74-75 high chair 148 lumbar traction 104-106 HiRes™ plastic goniometer 76-77 lumbar roll/support 122-125 Hi Volt electrotherapy 115 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 67 peg board 60,91 PeP™ pack 4,8 PeF™ pack 4,8 Perf-100™ exercise band 5-6 Perf-100™ exercise band 5-6 Porf-100™ exercise band 4-5 pinch exerciser 24 pinch gauge 62-66,69,74-75 pinch gauge 63-66,69,74-75 pinch gau		_			
hand/wrist exerciser26-27lifting boxes72PEP™ pack4,8Handy weight ball29loop, resistive exercise7percussor126head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31	<u> </u>	· ·		•	
Handy weight ball 29 loop, resistive exercise 7 percussor 126 head/neck support 172,174 loop stirrup 14 Perf-100™ exercise band 5-6 head protector 167 lotion, ultrasound 111-140 Purdue® peg board 91 heart-rate monitor 90 lotion, cold/hot therapy 98-100 pillow/support 147 heating pad 101 low powder, exercise band 4-5 pinch exerciser 24 height measurement 88 low volt, electrotherapy 115 pinch gauge 62-66,69,74-75 high chair 148 lumbar traction 104-106 pinwheels, sensory evaluation 83-84 HiRes™ plastic goniometer 76-77 lumbar roll/support 122-125 pivotal therapy system 125 Hi Volt electrotherapy 114-115 Magnecizer™ pedal, shoulder 34 plastic ball chair 48 HoldRite™ handle 12 manual muscles tester/MMT 68-70,74-75 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31		·		. •	
head/neck support172,174loop stirrup14Perf-100™ exercise band5-6head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31		· ·			
head protector167lotion, ultrasound111-140Purdue® peg board91heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31	, ,	1 7		•	
heart-rate monitor90lotion, cold/hot therapy98-100pillow/support147heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31		·			
heating pad101low powder, exercise band4-5pinch exerciser24height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31	•			. 0	
height measurement88low volt, electrotherapy115pinch gauge62-66,69,74-75high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31					
high chair148lumbar traction104-106pinwheels, sensory evaluation83-84HiRes™ plastic goniometer76-77lumbar roll/support122-125pivotal therapy system125Hi Volt electrotherapy114-115Magnecizer™ pedal, shoulder34plastic ball chair48HoldRite™ handle12manual muscles tester/MMT68-70,74-75plastic goniometer76-77hook and loop fastener186martin-type vigorimeter62plyotherapy / plyoball30-31	· .			•	
HiRes™ plastic goniometer 76-77 lumbar roll/support 122-125 pivotal therapy system 125 Hi Volt electrotherapy 114-115 Magnecizer™ pedal, shoulder 34 plastic ball chair 48 HoldRite™ handle 12 manual muscles tester/MMT 68-70,74-75 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31	· ·	•			
Hi Volt electrotherapy 114-115 Magnecizer™ pedal, shoulder 34 plastic ball chair 48 HoldRite™ handle 12 manual muscles tester/MMT 68-70,74-75 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31	S .			•	
HoldRite <sup>™</sup> handle 12 manual muscles tester/MMT 68-70,74-75 plastic goniometer 76-77 hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31	The state of the s				
hook and loop fastener 186 martin-type vigorimeter 62 plyotherapy / plyoball 30-31	1,7				
,, ,					
	·	,. °			
hot water bottle 101 MARV <sup>™</sup> band handle 12 PneuGel® wrist/arm support 164					164



		sku-sorted pag	e index		
pneumatic dynamometer	62	sitting/standing disc	45	ThermoTherapy® dry heat	109
pneumatic stool	148	skate, exercise	58	Thumper massagers®	162
Point-Relief <sup>™</sup> hot/cold pack	102-104,106	Skillbuilders®		Tufcoat™ foam roller	38
Point-Relief <sup>™</sup> hot/cold spot	11,98	<ul> <li>hi-low work table</li> </ul>	151	timer	90
Point-Relief <sup>™</sup> mini-massager	138	<ul> <li>mobility</li> </ul>	172	toilet seat, raised/support	161
Polar <sup>™</sup> heart and BP monitor	90	<ul> <li>positioning</li> </ul>	169-170,172	TouchTest <sup>™</sup> monofilament	82
Polar Frost® cold therapy	99	<ul> <li>seating/accessory</li> </ul>	172	traction	
pool, augatic therapy	135	skinfold caliper	86	<ul> <li>accessories</li> </ul>	125
positioning roll/wedge	56-57	Skyndex® caliper	86	<ul> <li>overdoor/cervical</li> </ul>	122
posture grid and evaluator	81	sled, weighted	72	<ul> <li>system/table</li> </ul>	123-124
potty chair, child	160	smedley spring dynamometer	62	transfer (sliding) board	159
Pouch™ adjustable weight	27	Sof-seam <sup>™</sup> knee support	164	transfer belt	195
powder board	58	SoftGrip <sup>™</sup> hand weight	29	Transport <sup>™</sup> , Intelect <sup>®</sup>	113,121
PresSsion® compression	127	Sonicator® electrotherapy	116-117	treadmill	37
Pressure Positive® massager	139	soundhead applicator	112,116	Tristander® pediatric stander	176
Pro®, finger exerciser	19	sphygmomanometer/cuff	85	TumbleForms 2 <sup>®</sup> 161,	176,180-182
Pronex® traction	123	Spider-Tech® taping	137	Turt-L-Shell® balance exerciser	182
PronPillo® face cradle	57,125	spirometer, lung capacity	89	tubing (resistive exercise)	8-10,17
Prossage® topical	110	splint bath and material	136	tuning fork	84
ProStretch® stretcher	59	spray, cold	98-100	Twist-n-Bend <sup>™</sup> bar	25
protective helmet	167	spring hand grip exerciser	21	TX® traction unit/table	106
protractor	81	stability trainer	42,44	ultrasound unit/accessories	111-116
Psytech™ exerciser/evaluation	24,64	stabilizer base, ball	49	universal inclinometer	80
pulley, various	14,32,58	Stabilizer <sup>™</sup> biofeedback	131	UV lamps	101
pulse oxymeter	85	stacker rings, inflatable ball	49	Valu-Form <sup>™</sup> roll/wedge	168
pump, electric/manual	49	stainless steel goniometer	79-79	Vari-Grip® finger exerciser	19
Purdue pegboard	91	stairs	152,181	velcro	186
putty, exercise	22-23	stander	176	Versa Form® pillows	171
PuttyCize <sup>™</sup> Theraputty <sup>™</sup> tool	23	stepper, inflatable	33	vestibular	
PVC inflatable ball rack	49	stethoscope	85	<ul> <li>boards, systems</li> </ul>	40-44
quad board	57	stimulator, muscle, EMS	114-119	<ul> <li>discs and wedges</li> </ul>	45
rack (band, weight, ball)		stirrup (for band and tubing)	13	• domes	181
• ball	49	Sroufe® support/brace	164	Vestibulator II®	167
<ul><li>band</li></ul>	15	stone (balance)	45	Via® finger exerciser	18
<ul><li>weight</li></ul>	26-32	stool, treatment/foot	148	vinyl-coated dumbbell	29
raised toilet seat	161	stopwatch	90	volumetric edema gauge	87
reaching aid	189	storage rack (exercise band)	15	wall-mounted goniometer	81
rebounder	31	stress ball	21	walker, pediatric/adult	154,157
refractometer	89	stretching product	59	wartenberg pinwheel	83
reider seifer tuning fork	84	StretchRite <sup>™</sup> strap	59	WalSlide <sup>™</sup> exercise station	14
Relief Pak™ hot/cold pack	102-104,106	stroller/accessories	173-174,177	WaxWel <sup>™</sup> paraffin	96
REP™ band/tubing	17	supine traction	122	web, hand exercise	24
rest mat	53	swing seat/accessories	166-167	wedge, positioning	56-57
robinson goniometer	79	Sys*Stim® electrotherapy	117,121	weight	
roll		tables		chest pulley	32
<ul> <li>positioning</li> </ul>	56-57,170	• massage	141	• Cuff®/Pouch™ wrist/ankle	26-27
<ul> <li>lumbar, cervical support</li> </ul>	147	• tilt	150	<ul> <li>dumbbell/kettlebell</li> </ul>	29
• inflatable	39,47-48	traction	124	• plate/rack	32
rocker board	40-44	<ul> <li>treatment, mat, and work</li> </ul>	142-145,151	<ul> <li>SoftGrip<sup>™</sup> hand weight</li> </ul>	29
rollator	157	tactile <sup>™</sup> monofilament	82	<ul> <li>weighted ball</li> </ul>	29-31
RolUp <sup>™</sup> bar	28	tape measure	87	West™ monofilament	82
rubber band	19	Taylor hammer	84	wheelchair/accessory	158-159
rulongmeter	76	tennis elbow support	164-165	whirlpool	132-133
saddle rolls (inflatable)	47	TENS therapy	119-120	wobble boards	40-43
Saunders® traction	123-125	Thera-Band® products		work hardening	72-73
scale, weight/height	88-89	<ul> <li>balls, hand exerciser</li> </ul>	21	wrap, hot/cold/electrotherapy	104
scoliosis meter	81	<ul> <li>band/accessories</li> </ul>	15-17	wrist	
scooter, adaptive mobility	180,182	• Flexbar®	25	<ul> <li>wrist dynamometer</li> </ul>	67
seat, feeder and floor sitter	172-175	<ul> <li>soft weights</li> </ul>	30	<ul> <li>wrist/hand exerciser</li> </ul>	18-25
sEMG unit/electrode	131	<ul> <li>stability trainer</li> </ul>	44	<ul> <li>wrist inclinometer</li> </ul>	79
Semmes-Weinstein monofilam		TheraBath® paraffin therapy	97	<ul> <li>wrist straps, accessory</li> </ul>	10,71
sensi ball	46	Theraputty <sup>™</sup> exercise material	22-23	<ul> <li>wrist weight</li> </ul>	26-27
	58,60	TheraTemp® mirowaveable pa	cks 106	X-ray goniometer	79



# Fabrication Enterprises order form

L DILL TO				0. 1457	1105	<b>25 5</b> 4		FNT						CVC c	ode
1: BILL TO:				3: MET								5)		1.00	
	· · · · · · · · · · · · · · · · · · ·			☐ Check			to Fab	rication	ı Enterpr	ises, Inc.)					
				☐ Bill our			ΠМ	lactor()	ard 🗇	Discover	ΠΔn	Eν			
Address:				Card Acco			<b>-</b> IVI	iastei O	alu 🗖		xpiration				
City:		State:	Zip:								П				
				Name (Plea	ase Print	:):									
our Name:				Credit Card											
				City:							St	ate:		Zip:	
				Signature:								Date:			
2: SHIP TO	(If different than bill to)	*Shipment cannot be made to	a post office box	4: SHIP	PING	СНА	RGI	=0.							
		. ,	•	• Drop shi					or at no	additional	charge				
				All freight								bility of the	purchas	er and wi	ill be add
				to your i	invoice.										
-					0-15 work			irectly i	ioiii tile	iactory or	made id	order. The	ese will i	normany	be snipp
City:		State:	Zip:					ged whe	en merch	andise is	available	. Charges	for this s	service wi	ll be add
				to your i	invoice.										
5: PRODUC	PRODUCT NUMBER	P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	R UNIT		тот	'AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	R UNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	R UNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	R UNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	R UNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	RUNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	RUNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	R UNIT		тот	AL
		P	PRODUCT DESCR	EIPTION					UNIT	PR	ICE PEI	RUNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	RUNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PR	ICE PEI	RUNIT		ТОТ	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PRI	ICE PEI	RUNIT		тот	AL
		P	PRODUCT DESCR	RIPTION					UNIT	PRI	ICE PEI	RUNIT		ТОТ	AL
		P	PRODUCT DESCR	RIPTION								E TOTAL		TOT	AL
QUANTITY			PRODUCT DESCR	RIPTION							ANDISI			TOT	AL

Please email, fax, mail, or phone the completed order form to the contact information below.

New York Sales Tax resale certificate must accompany orders that are "billed to" or "shipped to" New York state.





# industry leading brands from Fabrication Enterprises






















additional brands include: **ArmaSport® mats**, **Dipsters® patientware** and more visit www.Fab-Ent.com to view FEI products, brochures, and updates







Fabrication Enterprises Inc. PO Box 1500 White Plains, New York 10602 (USA)

tel: 914-345-9300 800-431-2830 fax: 914-345-9800 800-634-5370

www.Fab-Ent.com

PRESORT STANDARD
U.S. POSTAGE
PAID
FABRICATION
ENTERPRISES
INCORPORATED

© 2010 FEI, all rights reserved



















Relief Pak®

